## ENGLISH

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

### INTENDED USER PROFILE:

THE USER SHOULD BE ABLE TO:

- READ AND UNDERSTAND THE DIRECTIONS, WARNINGS AND CAUTIONS.

**INTENDED USE/INDICATIONS:** Restricts cervical spine flexion, extension and rotation.

**CONTRAINDICATIONS:** None

**WARNINGS AND CAUTIONS:** If you experience any pain, swelling, sensation changes, or if any unusual reactions while using this product, consult your medical professional immediately.

### SIZING INFORMATION:

1. For Adult sizing, use the SIZING GUIDE to determine the appropriate collar size by measuring the vertical distance from the tip of the chin to the top of the shoulder. For Pediatric sizing, reference the back of the SIZING GUIDE. (Figure 1)

### APPLICATION INFORMATION:

2. SUPINE APPLICATION: To fit the patient when laying down.
  - A. Slide the BACK PIECE under the patient's neck and center it. Ensure the black/grey foam padding extends beyond the edge of the plastic. (Figure A)
  - B. Flare the sides of the FRONT PIECE out, slide it up the chest wall and scoop it up under the chin. Sides of the collar FRONT PIECE should be oriented up, off the trapezius and toward the ears. (Figure B)
  - C. While holding the FRONT PIECE securely, curl the ends snugly against the patient's neck. Apply the Velcro® strap and secure the opposite side of the collar in the same fashion. Tighten straps alternately, one at a time, to an equal length on both sides. Make sure the FRONT PIECE of collar is inside the BACK PIECE. (Figure C)

**NOTE:** Collar application must be secure to ensure proper fit, maintaining prescribed alignment, and to prevent the chin from slipping inside. If the patient can slip his/her chin inside the collar, it is a clear indication that it is not snug enough.

3. UPRIGHT APPLICATION: To fit the collar to a sitting or standing patient.
  - A. With the patient in neutral head alignment, slide the FRONT PIECE of the collar up under the chin and properly position the chin support. Sides of the collar FRONT PIECE should be oriented up, off the trapezius and toward the ears. (Figure D)
  - E. Hold the FRONT PIECE firmly with one hand. Place BACK PIECE under brace of neck until panel is centered and attach both sides to FRONT PIECE. (Figure E)
  - F. Apply the Velcro strap and secure the opposite side of the collar in the same fashion. Tighten straps alternately, one at a time, to an equal length on both sides. Make sure the FRONT PIECE of collar is inside the BACK PIECE. (Figure F)

### OTHER APPLICATION INFORMATION:

**Angulation Buttons:** the ANGLULATION BUTTONS may be used to fine-tune the fit of the collar for anatomical conformity and pressure distribution. It is easier to adjust the ANGLULATION BUTTONS if you first remove the collar. (Figure G)

**Chin Support Positioning:** Slide up to help relieve chin tip pressure/discomfort. Slide down to help relieve jaw pressure/discomfort.

**Occiput Support Positioning Guidelines:** Slide up to help relieve occipital pressure/discomfort. Slide down to help enhance extension control. For a predominantly supine patient keep OCCIPITAL SUPPORT BUTTONS set to top position to help reduce pressure.

**CLEANING INSTRUCTIONS:** Separate the brace's removable pads before cleaning.

Hand wash in water (86°F/30°C) using soap. Air dry.

**MATERIAL CONTENTS: RIGID SHELL:** Low Linear Density Polyethylene; **SOFTGOODS:** UBL, Polyester, Polyurethane Foam, Interlock.

**WARRANTY:** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

**Rx ONLY. INTENDED FOR SINGLE PATIENT USE. NOT MADE WITH NATURAL RUBBER LATEX.**

NOTE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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