
Enovis™ Instruments and Instrument Cases  
 Instrumente und Instrumentenbehälter von Enovis™  
 Instruments chirurgicaux et boîtiers à instruments Enovis™  
 Instrumentos y estuches de instrumentos quirúrgicos de Enovis™  
 Strumenti chirurgici e custodie Enovis™  
 Χειρουργικά όργανα και θήκες οργάνων Enovis™  
 Enovis™ Aletleri ve Alet Kutuları

**Encore® Medical, L.P.**  
**9800 Metric Blvd.**  
**Austin, TX 78758-5445 USA**

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Distributors for countries where Enovis™ is not the distributor:

Country	Distributor Name	Distributor Address
South Korea	Biolmplant Technology	1029 Yeongdeok-dong, Giheung-gu, Yongin-si, Gyeonggi-do, Korea, 16950 #2101 U-Tower
Denmark	Innosurge Trauma	Alsikevej 16 DK-8920 Randers NV
Belgium	Orthogrow Distribution	DaVincilaan 1 1930 Zaventem Belgium
South Africa	MedHold	MSI Business Park 68 Rigger Road Spartan Kempton Park Gauteng 1619
Mexico	Daonsa	Miguel Hidalgo 2426 Pte, Obispado, 64010 Monterrey, N.L., Mexico
Columbia	Piemca-Venezuela / Piemca-Columbia	Cra. 15 No. 88 - 64 / Edificio Torre Zimma / Oficina 705 y 713 Bogotá D.C., Colombia
Saudi Arabia	Al-Ewan Medical Company	Prince Nasir Ibn Farhan Al Saud, Al Mursilat, Riyadh 12461, Saudi Arabia
Turkey	Armoni Medikal	Gazi Mahallesi Yavuz Kanat Sok. Pasifik Plaza No. 34 Yenimahalle 06560 Ankara, Turkey

EN

**1. Product Handling**

Devices not returned to Enovis™ should be treated as biohazardous material and disposed of in accordance with local laws and regulations.

**Recommendation for the Care and Handling for Enovis™ Instruments and Instrument Cases**

<b>REUSABLE INSTRUMENT DESCRIPTION</b>	Enovis™ instrumentation consists of devices and their accessories used in surgical procedures. Implantation of Enovis™ products should only be performed with Enovis™ instrumentation or instrumentation distributed by Enovis™. Enovis™ instruments and instrument cases are generally composed of titanium, stainless steel, aluminum, and/or polymeric materials. The cases may be multi-layered with various inserts to hold surgical instrumentation in place during handling and storage. The inserts may consist of trays, holders, and silicone mats. The instrument cases will allow sterilization of the contents to occur in a steam autoclave utilizing the cleaning, sterilization, and drying cycle that has been validated and detailed below. Instrument cases do not provide a sterile barrier and must be used in conjunction with FDA cleared sterilization wrap to maintain sterility. Instruments are provided non-sterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines listed below.
<b>WARNINGS</b>	Automated cleaning may not be thorough enough. Carefully inspect each instrument to ensure that all visible blood residue and other contaminants have been removed.
<b>CAUTION</b>	Federal Law (USA) restricts this device to sale by or on the order of a physician.
<b>REPROCESSING LIMITATIONS</b>	Enovis™ instruments can be steam sterilized and repeat sterilization will not adversely affect them. If problems related to instrument sets are identified when using our instruments or instrument cases, please bring it to the attention of Enovis™ for investigation. The lifetime of an instrument is typically limited by normal wear and damage due to use.
<b>DISCLAIMER</b>	Enovis™ instrument cases are intended to protect instrumentation and facilitate the sterilization process by allowing steam penetration and drying. Enovis™ has verified through laboratory testing that our instrument cases are suitable for the sterilization cycles listed in the sterilization section of the IFU. It is the user's responsibility to verify that equipment is performing as intended, and conditions are achieved.

**INSTRUCTIONS FOR USE**

<b>POINT OF USE PREPARATION</b>	Keep instruments moist and do not allow blood and/or bodily fluids to dry on the instruments. The decontamination process should begin immediately after the completion of the surgical procedure. If cleaning must be delayed, place instruments in a covered container with pH Neutral enzymatic detergent to delay drying. Instruments should be cleaned within 30 minutes of use to minimize the potential for drying prior to cleaning. Wash all instruments whether or not they were used or were inadvertently contacted with blood. Disassemble instruments with removable parts; loosen instruments with movable parts, as applicable.
<b>DECONTAMINATION</b>	Decontamination is for the purpose of microbial inactivation. Saturate the surface completely with full strength intermediate disinfectant/cleaner* (e.g. CaviCide) and allow to remain in contact with devices for 5 minutes.

<b>A. MANUAL CLEANING: ALL INSTRUMENTS</b>	<ol style="list-style-type: none"> <li><b>Pre-Cleaning:</b> Remove all visible soil by immersing the devices in room temperature neutral pH enzymatic cleaner* (e.g. MetriZyme) and disassemble/loosen instruments, if suitable. The majority of the surgical instruments and trial devices are simply constructed and will not require disassembly. However, some of the more complex instruments are made of several components and these should be disassembled into their individual parts prior to decontamination. Scrub with the appropriate soft bristle brush until visibly clean; actuate through the full range of motion.</li> <li><b>Washing:</b> Immerse devices in the ultrasonic washer/cleaner with room temperature neutral pH enzymatic cleaner* (e.g. MetriZyme) and sonicate for 10 minutes. Ultrasonic cleaners can be used with hot water per the manufacturers' recommended temperature; however, room temperature was qualified. Be aware that loading patterns, water temperature, and other external factors may change the effectiveness of the equipment.</li> <li><b>Rinsing:</b> Thoroughly rinse the devices with deionized or distilled water. For example, a minimum of 2 minutes three (3) times.</li> </ol> <p>* Do not use high acidic (pH &lt;4) or high alkaline (pH &gt;10) products for disinfection or cleaning, since these can corrode metal, cause discoloration or stress fractures. Enovis™ has qualified the above cleaning method with the provided solution examples, for a 3 Spore Log Reduction (SLR). Other cleaning/disinfection methods may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternate method using appropriate laboratory techniques.</p>
<b>B. MANUAL CLEANING: INSTRUMENTS WITH CANNULAS, LUMENS, OR HOLES</b>	<ol style="list-style-type: none"> <li><b>Pre-Cleaning:</b> Follow the "Pre-Cleaning" and "Washing" steps in Section A. Manual Cleaning – ALL INSTRUMENTS.</li> <li><b>Washing:</b> After ultrasonic cleaning, in a fresh enzymatic cleaning bath use a tight-fitting, soft, non-metallic cleaning brush or pipe cleaner to scrub any cannula, lumen, or hole(s). Push in and out, using a twisting motion to remove debris. Use a syringe filled with enzymatic neutral pH cleaning solution to flush hard to reach internal areas.</li> <li><b>Rinsing:</b> Flush the instrument paying special attention to the cannulations, lumens, and/or holes with deionized or distilled water. For example, a minimum of 2 minutes three (3) times.</li> </ol>
<b>C. MANUAL CLEANING: ARTICULATING INSTRUMENTS</b>	<ol style="list-style-type: none"> <li><b>Pre-Cleaning:</b> Follow the "Pre-Cleaning" and "Washing" steps in Section A. Manual Cleaning – ALL INSTRUMENTS.</li> <li><b>Washing:</b> After ultrasonic cleaning, immerse the instrument in fresh neutral pH enzymatic cleaning solution to avoid aerosol generation. Actuate moveable mechanisms through full range of motion, such as knobs, hinges, box locks, or spring-loaded/retractable features. For instruments with flexible shafts, bend or flex the instrument under the neutral pH cleaning solution while brushing the flexible areas. For instruments with internal cavities, after actuating components in the neutral PH cleaning solution, fully open components and use a tight-fitting, soft, non-metallic cleaning brush or pipe cleaner to scrub the internal cavities. Use a syringe filled with enzymatic neutral pH cleaning solution to flush hard to reach internal areas</li> <li><b>Rinsing:</b> Actuate and/or retract moveable parts while rinsing with deionized or distilled water. For example, a minimum of 2 minutes three (3) times. For instruments with flexible shafts, flex the instrument while rinsing.</li> </ol>

<b>AUTOMATED CLEANING</b>	<p>ALL INSTRUMENTS Manual Cleaning steps in Section A, and Section B for instruments with cannulas, lumens or holes and Section C for articulating instruments required before automated cleaning. The following minimum parameters required for automated cleaning of Enovis™ instruments. These parameters have been validated by Enovis™ under laboratory conditions.</p> <table border="1" data-bbox="529 191 1377 390"> <thead> <tr> <th>Treatment / Phase</th> <th>Minimum Time</th> <th>Minimum Temperature</th> <th>Cleaning Detergent Type</th> </tr> </thead> <tbody> <tr> <td>Pre-Wash</td> <td>1 Minute</td> <td>Cold Tap Water</td> <td>N/A</td> </tr> <tr> <td>Enzyme Wash</td> <td>3 Minute</td> <td>Warm Tap Water 38°C</td> <td>Neutral pH Enzymatic Detergent</td> </tr> <tr> <td>Wash</td> <td>5 Minute</td> <td>Warm Tap Water 38°C</td> <td>Neutral pH Detergent</td> </tr> <tr> <td>Rinse</td> <td>2 Minute</td> <td>Warm Tap Water 38°C</td> <td>N/A</td> </tr> <tr> <td>Final Rinse</td> <td>1 Minute</td> <td>Hot Deionized Water 82°C</td> <td>N/A</td> </tr> <tr> <td>Thermal Pure Water Rinse</td> <td>5 Minute</td> <td>Deionized Water 93°C</td> <td>N/A</td> </tr> </tbody> </table>	Treatment / Phase	Minimum Time	Minimum Temperature	Cleaning Detergent Type	Pre-Wash	1 Minute	Cold Tap Water	N/A	Enzyme Wash	3 Minute	Warm Tap Water 38°C	Neutral pH Enzymatic Detergent	Wash	5 Minute	Warm Tap Water 38°C	Neutral pH Detergent	Rinse	2 Minute	Warm Tap Water 38°C	N/A	Final Rinse	1 Minute	Hot Deionized Water 82°C	N/A	Thermal Pure Water Rinse	5 Minute	Deionized Water 93°C	N/A
Treatment / Phase	Minimum Time	Minimum Temperature	Cleaning Detergent Type																										
Pre-Wash	1 Minute	Cold Tap Water	N/A																										
Enzyme Wash	3 Minute	Warm Tap Water 38°C	Neutral pH Enzymatic Detergent																										
Wash	5 Minute	Warm Tap Water 38°C	Neutral pH Detergent																										
Rinse	2 Minute	Warm Tap Water 38°C	N/A																										
Final Rinse	1 Minute	Hot Deionized Water 82°C	N/A																										
Thermal Pure Water Rinse	5 Minute	Deionized Water 93°C	N/A																										
<b>DRYING</b>	<p>Ensure device is dry prior to inspection and sterilization preparation. Instruments must be thoroughly dried to remove residual moisture before they are stored. Filtered compressed air may be used prior to air drying if available.</p>																												
<b>MAINTENANCE INSPECTION AND TESTING</b>	<p>After cleaning, the instruments (disassembled, if applicable) should be visually inspected. Check for misalignment, burrs, bent, or fractured tips. Mechanically test the working parts (e.g. hinges) to verify that each instrument functions throughout its intended range of motion. Place instruments into appropriate configuration within instrument case and wrap with protective FDA cleared sterilization wrap according to AAMI / AORN guidelines.</p> <p>Surgical instruments and instrument cases are susceptible to damage from prolonged use, and through misuse or rough handling. Care must be taken to avoid compromising their performance. To minimize damage, conduct the following:</p> <ol style="list-style-type: none"> <li>1. Inspect instrument cases and instruments for damage when received and after each use and cleaning. Incompletely cleaned instruments should be re-cleaned, and those that need repair returned for servicing.</li> <li>2. Only use an instrument for its intended purpose.</li> <li>3. When handling sharp instruments use extreme caution to avoid injury. Consult with an infection control practitioner to develop safety procedures appropriate for all levels of direct instrument contact.</li> <li>4. If instruments appear to be damaged in such a way that may compromise the performance of the instrument, contact your Enovis™ representative for a replacement.</li> <li>5. Visually inspect the instrument and check for damage and wear, moveable parts should have smooth movement, locking mechanisms should fasten securely</li> </ol>																												
<b>TRANSPORT</b>	<p>Compliance with the general precautionary measures for handling contaminated/biologically hazardous materials is required.</p>																												
<b>STERILIZATION</b>	<p>Instruments supplied by Enovis™ have been thoroughly cleaned, inspected and tested for proper function prior to shipment. Unless otherwise indicated, these instruments are NOT STERILE and must be sterilized prior to use. Instruments provided outside of instrument sets should be fully loosened/disassembled and wrapped in FDA cleared sterilization wrap per AAMI ST:79/AORN Guidelines. Flash (immediate-use) steam sterilization by exposure at 132°C / 270°F should only be used as an emergency procedure. Instruments must be cleaned and disassembled prior to processing.</p> <p>The following are minimum cycles required for steam sterilization that has been validated by Enovis™ under laboratory conditions to achieve a SAL of 10<sup>-6</sup> with components loosened or disassembled. Enovis™ has data on file.</p> <p><b>Sterilization with a Pre-Vacuum Sterilizer (HI-VAC):</b> 270° F (132° C), 5-minute exposure time, and a 30-minute dry time.</p> <p><b>Sterilization with a Gravity Displacement Sterilizer:</b> 270° F (132° C), 30-minute exposure time, with a 30-minute dry time.</p>																												
<b>STORAGE/INSTRUMENT CARE</b>	<p>Instruments must be thoroughly dried to remove residual moisture before they are stored. Instruments or instrument cases that have been processed and wrapped to maintain sterility should be stored in a manner to avoid extremes in temperature and moisture. Care must be taken in handling wrapped instruments or instrument cases to prevent damage to the barrier. The user must be aware that maintenance of sterility is event-related and that the probability of occurrence of a contaminating event increases over time and with handling</p>																												
<b>CONTACT INFORMATION</b>	<p>Enovis™  <b>ATTN: Customer Service</b>  <b>9800 Metric Boulevard</b>  <b>Austin TX, 78758 USA</b>  <b>+ 1-800-456-8696</b></p>																												

## 2. Product Description

The devices covered by this IFU are the instruments used to implant Enovis™ implantable devices, as well as the instrument cases used to store said instruments for cleaning and transportation.

The following devices are covered by IFU 0400-0221 and not this IFU.

FA XALT TRIAL (803-99-039)  
FA FMP INSTRUMENT (803-99-018)  
FA RSP GLENOID (804-99-11)  
3D KNEE GAP BALANCER (S-200775)  
RSP HALF MOON REAMERS (804-06-012, 804-06-013, 804-06-014)  
FA K EXPRT REV 1 (800-99-092)  
FA K EXPRT REV 2 (800-99-093)  
FA DAA GENERAL (803-99-102)  
FA DAA RETRACTOR (803-99-103)  
FA TURON RETRACT (804-99-020)  
FA TURON HUMERAL STEM 2010 (804-99-117)  
FA TURON GLENOID 2010 (804-99-118)  
FA TURON HUMERAL HEAD 2010 (804-99-119)  
FA RSP HUMERAL (804-99-010)  
FA RSP GLENOID (804-99-011)  
FA RSP MONOBLOCK LTD RELEASE (804-99-025)  
FA RSP SIZE 44 INST (804-99-024)  
FA K EMP DRF FEM PREP (800-99-094)  
FA K EMP TIB PREP (800-99-095)  
FA K EMP PAT TOOL KIT (800-99-096)

FA K EMP BONUS KIT (800-99-097)  
 FA K EMP 3D TRL CORE LT (800-99-098)  
 FA K EMP 3D TRL CORE RT (800-99-099)  
 FA S ALTIVATE RSP HUM PREP (804-99-120)  
 FA S ALTIVATE RSP HUM TRLS (804-99-121)  
 FA K EMP 3D TRL PREP OUT SML (800-99-101)  
 FA K EMP 3D TRL PREP OUT LG (800-99-102)  
 FA K EMP INS TRL SPCRS (800-99-103)  
 FA LR FMP CUP INSERTER (803-99-098)  
 FA FMP ACET REAMER (803-99-003)  
 FA XALT INST (803-99-040)  
 FA MIS HIP (803-99-028)  
 FA TAPERFILL INSTRUMENTS (803-99-170)  
 FA TAPERFILL INSTRUMENTS – 7MM (803-99-170)  
 FA TAPERFILL BROACHES (803-99-171)  
 FA H TAPERFILL BROACH 7MM POST (803-99-172)  
 FA TAPERFILL BIG LUG FT BROACH (803-99-172)  
 FA K EMP 3D CR TRL PREP OUT LG  
 FA K EMP 3DCR TRL PREP OUT SML  
 FA K EMP CR TRL CORE LT  
 FA K EMP CR TRL CORE RT  
 FA K EMP PS FEM TRL CAP CORE (800-99-117)  
 FA K EMP PS FEM TRL CORE (800-99-117)  
 FA K EMP PS INS TRL CORE (800-99-118)  
 FA K EMP PS TRL PREP OUT L CAP (800-99-108)  
 FA K EMP PS TRL PREP OUT LG (800-99-108)  
 FA K EMP PS TRL PREP OUT S CAP (800-99-107)  
 FA K EMP PS TRL PREP OUT SML (800-99-107)  
 FA K EMP 3D + CR COMP OUT LG  
 FA K EMP 3D + CR COMP OUT SML  
 FA K EMP 3D COMP OUT LG  
 FA K EMP 3D COMP OUT SML  
 FA K EMP COMP PAT TOOL KIT  
 FA K EMP CR COMP OUT LG  
 FA K EMP CR COMP OUT SML  
 FA K EMP TIB PREP 2  
 FA S ALTIVATE RSP 44  
 FA S ALTIVATE RSP HUM TRLS SML  
 FA K EMP CEM STEM AUG PREP  
 FA K EMP TIB PREP COMPLETE  
 FA K EMP PS TRL PREP OUT SML 2  
 FA K EMP PS TRL PREP OUT LG 2  
 FA K EMP VVC INS TRL CORE  
 FA S ALTIVATE RSP REVISION  
 FA H EMPOWR ACET GENERAL INST  
 FA H EMPOWR ACET MIS HANDLES  
 FA H EMPOWR ACET TRL NEU 10DH  
 FA H EMPWR ACET OFFSET TR LNRS  
 FA H EMPOWR ACET ANCILLARY  
 FA H EMPOWR ACET GEN INST V2  
 FA H EMPOWR ACET ANCILLARY V2  
 FA S ALTIVATE RSP SHORT  
 FA K EMP PARTIAL  
 FA K EMP PARTIAL PREP 1 OF 2  
 FA K EMP PARTIAL TRIAL 2 OF 2  
 FA S ALT ANATOMIC CS EDGE  
 FA K EMP FEM PREP - ADJ SZR

The following devices are covered by IFU 0400-0248 and not this IFU.

FA K EMP PS FEM TRL CAP CORE (800-99-117)  
 FA K EMP PS FEM TRL CORE (800-99-117)  
 FA K EMP PS INS TRL CORE (800-99-118)  
 FA K EMP PS TRL PREP OUT L CAP (800-99-108)  
 FA K EMP PS TRL PREP OUT LG (800-99-108)  
 FA K EMP PS TRL PREP OUT S CAP (800-99-107)  
 FA K EMP PS TRL PREP OUT SML (800-99-107)  
 FA S ALT ANATOMIC HUM STEM (804-99-125)  
 FA S ALT ANATOMIC HUM HEAD (804-99-126)  
 FA S ALT ANATOMIC GLENOID (804-99-127)  
 FA S TURON KEELED GLENOID (804-99-128)  
 FA H EXPRT REV HIP REAMERS (803-99-095)  
 FA H EXPRT REV HIP GENERAL (803-99-096)  
 FA H FMP STR CUP INSERTER (803-99-101)  
 DOUBLE OFFSET BROACH HANDLE, LEFT (803-03-084)  
 DOUBLE OFFSET BROACH HANDLE, RIGHT (803-03-085)  
 FA H HXE NEU-10DH EXT D (803-99-111)

### 3. Indications

Reference the applicable implant IFU for Indications.

### 4. Intended Use

Surgical instruments are intended to aid the surgeon in performing a prosthesis implantation.

### 5. Clinical Benefit

Surgical instruments assist in the implantation of prostheses and do not have a direct therapeutic or diagnostic function.

**6. Contraindications**

Reference the applicable implant IFU for Contraindications.

**7. Precautions and Warnings**

Stainless steel surgical instrumentation may contain nickel which is a known sensitizer. Sound medical judgement should be used if nickel sensitivity or allergy is suspected. Stainless steel surgical instrumentation may contain cobalt at a concentration greater than 0.1%. Although the cobalt is classified as hazardous, when used in this application, there are no risks associated with the presence of cobalt in these instruments. Reference the applicable implant IFU for additional Precautions and Warnings.

Inspect all equipment before use and do not use damaged or worn instrumentation as this may lead to delay in surgery, revision surgery, pain/discomfort/disability to the patient.

**8. Preoperative Planning and Postoperative Care**

Reference the applicable implant IFU for Preoperative Planning and Postoperative Care.

**9. MRI Safety**

N/A

**10. Adverse Effects**

Refer to Precautions and Warnings section above and reference the applicable implant IFU for Adverse Effects.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the relevant Competent Authority as defined in EU 2017/745.

**11. Lifetime of Device**

Enovis™ does not define the maximum number of uses appropriate for re-usable instruments. While the expected lifetime of surgical instruments may be subject to a multitude of factors such as patient characteristics, surgeon experience, amount of use, and surgical technique, evaluating the time between the release of an instrument from production and the return of that instrument to the manufacturer from the user can give an indication of its expected lifetime.

Based on this information, Powered Instruments have shown to last as short as 5 days in the field or as long as 9.5 years in the field, with an average lifetime of 2.3 years. Non-Powered Impaction or Extraction Instruments have shown to last as short as 3.5 months in the field or as long as 18 years in the field, with an average lifetime of 2.8 years. Non-Powered Guide Instruments have shown to last as short as 35 days in the field or as long as 9.1 years in the field, with an average lifetime of 2.9 years. Non-Powered, Non-Impaction/Extraction, & Non-Guide Instruments have shown to last as short as 56 days in the field or as long as 9.3 years in the field, with an average lifetime of 3.5 years.

As product data continues to be collected, these lifetime estimates may be re-evaluated and adjusted if required. Users should note that careful inspection of the instrument before use is the best method of determining the end of serviceable life.

**12. Trademarks and Patents**

Reference the applicable implant IFU for Trademarks and Patents.

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