
DONJOY[®]

spinamic[®]

INSTRUCTIONS FOR USE

For translated Instructions for Use, please go to <https://ifu.enovis.com/>

enovis[™]

EN	ENGLISH	2
ES	ESPAÑOL.....	11
FR	FRANÇAIS.....	20
DE	DEUTSCH.....	29
IT	ITALIANO.....	38
NL	NEDERLANDS.....	47
CS	ČEŠTINA.....	56
SK	SLOVENSKÝ.....	65
DA	DANSK.....	74
FI	SUOMI.....	83
NO	NORSK.....	92
SV	SVENSKA.....	101
PT	PORTUGUÊS	110
EL	ΕΛΛΗΝΙΚΑ.....	119
TR	TÜRKÇE.....	128
PL	POLSKI.....	137
HI/BN	हिन्दी & বাংলা	146
AR	العربية.....	155
RO	ROMÂNĂ.....	164
SR	SRPSKI	173
HR	HRVATSKI.....	182
SL	SLOVENŠČINA.....	191
TH	ภาษาไทย.....	200
HU	MAGYAR.....	209

TABLE OF CONTENTS

INTENDED PURPOSE	3
INTENDED USER PROFILE.....	3
PRECAUTIONS FOR USE.....	3
PRECAUTIONS FOR HANDLING	3
WASH & CARE INSTRUCTIONS	3
WARRANTY & DISPOSAL	4
SYMBOLS USED ON THE LABEL.....	4
PRODUCT COMPONENTS.....	5
BODY MEASUREMENTS / PRODUCT SIZES	6

BEFORE APPLYING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED PURPOSE

The Spinamic hybrid scoliosis brace is a single-patient use, scoliosis orthosis intended to reduce or prevent scoliosis curve progression and to manage spinal curvature in patients with scoliosis.

INTENDED USER PROFILE

The intended user is the patient or patient's caregiver. The patient or caregiver should be able to understand, and be physically capable of performing the directions, warnings, and cautions in the instructions for use. Initial fitting and usage guidance should be provided by a doctor or healthcare professional.

PRECAUTIONS FOR USE

A. The Spinamic device should only be used under doctor's instructions and prescription.

B. If the patient experiences breathing difficulties or digestion problems while using the product, they should remove the product immediately and consult their doctor or healthcare professional.

C. Contraindications:

- Obesity, defined as BMI \geq 95th percentile for pediatric patients or BMI $>$ 30 for adult patients.
- Body shape or size that prevents proper fit of the brace.
- Pregnancy.
- Non-corrective scoliosis.
- Inability to communicate.
- Respiratory issues: Severe breathing problems or lung disease that could be worsened by chest compression.
- Skin problems: Open wounds, active infections, or severe skin sensitivities in the areas where the brace touches the body.

D. Side Effects:

- Mild soreness: Tight feeling in the back or ribs as the body adjusts.
- Skin irritation: Redness, itching, allergies, rashes, or blisters. Patients should check their skin daily for redness or irritation. If symptoms persist or worsen, patients should stop use and contact their healthcare provider.
- Digestive issues: Abdominal discomfort, early satiety, or mild indigestion after eating.
- Chest wall or breast discomfort.

PRECAUTIONS FOR HANDLING

- Keep the product away from flame or heat to avoid risk of product damage due to burning.
- Due to material characteristics, repetitive or intense friction may cause damage or shine of the product surface.
- If the product surface is touched or rubbed against rough objects, fluff or fraying may occur.

WASH & CARE INSTRUCTIONS

- Hand wash separately in cold water. Do not machine wash.
- If the product becomes soiled, please hand wash it immediately.
- Do not use bleach or immerse the product in detergent for more than 10 minutes.
- Please dry in the shade.
- Machine drying (heat drying) and ironing may damage the product.

F. Store the product at room temperature when not in use.

do not dry clean

do not iron

do not bleach

hand wash

dry in the shade

WARRANTY & DISPOSAL

WARRANTY: DJO, LLC will replace all or part of the unit and its accessories for material or workmanship defects for a period of twelve months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

DISPOSAL: Dispose of the device and package in accordance with local or national environmental regulations.


SYMBOLS USED ON THE LABEL

The following symbols may appear on the device labeling, packaging, or accompanying documentation.


	Conformity with applicable European Union legislation		Distributor of the device
	Medical device		Importer in the EU
	Batch/Lot number		Consult instructions for use
	Unique Device Identifier		Intended for use on a single patient
	Manufacturing date		Not made with natural rubber latex
	Legal manufacturer		Prescription device (US)
	Authorized representative in the EU		Manufacturer's catalogue number

PRODUCT COMPONENTS

SPINAMIC BRACE

PRESSURE REGULATOR

INDICATOR TAPES (2 EA)



BODY MEASUREMENTS / PRODUCT SIZES

All fitting, pad placement, and pressure calibration must be performed by a qualified healthcare professional.

HEALTHCARE PROVIDER (HCP) GUIDE



1. CHEST

Patient should breathe out. Measure by placing the measuring tape around the patient's chest, centered at the bust point.



2. HIP

Measure by placing the measuring tape around the patient's hip, centered at the most convex part of the hip.

3. BACK

Measure by placing the measuring tape from the back of neck (C7, the protruding point when the face is down) down to the most convex part of the hip while touching the natural waist in the middle.

A. VISUAL ASSESSMENT & INITIAL PAD SETTING

1. Adam's forward bend test


- Procedure: Have the patient stand with feet together and knees straight. Instruct the patient to bend forward at the waist with arms hanging naturally.
- Assessment: Identify the rib hump (the most prominent area of the back) through visual inspection and palpation.

2. Hump pad attachment

- Attachment: Attach the hump pad on the inside of the Spinamic vest, aligning it with the previously identified hump position. The hump pad should be placed at the base of the hump to aid in derotation.
- Orientation: Ensure the pad is oriented horizontally for the maximum coverage.

3. Flank pad attachment

- Attachment: Put down the flank strap on the flat table with the center hook and loop facing down. Attach flank pads on the flank strap so that each flank pad is placed between ribs and pelvis. The flank strap is worn to elongate the spine by increasing abdominal pressure.
- Orientation: The flank pads should be oriented horizontally or slightly at an angle. Make sure to check with the patient so that the pads are not causing any pain when the flank strap is worn.

B. DONNING THE SPINAMIC VEST

1. Putting on the flank strap


- Put on a shirt (thin compression wear, sold separately).
- Put on the flank strap with two flank pads by pulling it forward to the front.

2. Putting on Spinamic vest with hump pad attached

- Put on Spinamic vest and pelvic band, respectively by pulling it forward to the front. (Initial attachment of the hump pad should take place following the guidance of healthcare providers.)

3. Zipping-up of Spinamic vest

- Attach the hook and loop cover after zipping up the vest.
- Re-attach the shoulder end to create the best fit.

C. RADIOGRAPHIC ANALYSIS & PRECISE CALIBRATION

Radiographic assessment should be conducted as clinically indicated to evaluate alignment and brace effectiveness.

1. Post-fitting X-ray & alignment analysis




- **Imaging:** Take an X-ray while the patient is wearing the brace, with the indicator hook and loop tapes attached to the brace as shown in the image.
- **Pressure Regulator Localization:** Use the indicator tape eyelets as a coordinate system to calculate the necessary placement of the pressure regulator.

2. Regulator position fixation

- **Alignment:** Based on the X-ray analysis, slide the pressure regulators along the rail to match the patient's target apex.
 - **Thoracic curve:** Place the regulator so that it follows the ribcage to target the apex of the thoracic curve.
 - **Lumbar curve:** Place the regulator parallel to the apex of the lumbar curve. It is recommended to position it a little bit below the apex of the lumbar curve if the regulator on the thoracic curve is too close.
- **Fixation:** Secure the regulator directly by tightening the screws into the rail holes using the Phillips screwdriver provided.
- **Cap Application:** Once the screws are tightened, place the screw caps over the screw heads.

3. Wiring the dial of pressure regulator


- Push each dial of pressure regulators until "click" is heard.
 - Increase pressure strength by turning the dial clockwise.
 - Release pressure by pulling the dial until "click" is heard.
- **NOTE:** Do not turn the dials counter-clockwise.

D. ATTACHMENT OF SUPPLEMENTARY BANDS


Attach the supplementary bands by pulling them forward. Ensure that each supplementary band covers the center of the pressure regulator. The thoracic supplementary band is placed after donning the shoulder strap.

1. The thoracic supplementary band from the shoulder strap should be attached slightly upwards diagonally.
2. The lumbar supplementary band located in the back of the brace should be attached horizontally.


E. ATTACHMENT OF SHOULDER STRAP AND THE LUMBAR STRAP

1. Attach the shoulder strap by pulling it forward over the shoulder and through the arm. The shoulder strap should be placed on the same shoulder as the lower pressure regulator.
2. Attach the lumbar strap by pulling it forward. Make sure that the strap covers the supplementary band.

F. ATTACHMENT OF THE THORACIC STRAP


Attach the thoracic strap by pulling it forward. The thoracic strap should cover the supplementary thoracic band and be attached slightly upwards diagonally.

G. FINAL VALIDATION

Imaging: Take an X-ray while the patient is wearing the brace, with the indicator hook and loop tapes attached to the brace as shown in the image. Radiographic assessment should be conducted as clinically indicated to evaluate alignment and brace effectiveness.

Confirmation: Perform the final X-ray examination to ensure the pressure points are accurately targeting the spinal curvature.

NOTE: HOW TO REVERSE THE SHOULDER STRAP

1. Put the supplementary band through the hole. Located at the base of the supplementary band, where it attaches to the shoulder strap
2. Detach the rubber patches, and re-attach them in the reverse manner.

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©2026 DJO, LLC. All rights reserved.
VALUE & TRUST Co. Ltd
Javitech Consulting GmbH