ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS: The Procure One Piece cervical collar is designed to stabilize and relieve the cervical spine and may be suitable for the symptomatic relief of pain and stiffness in the cervical spine area following patient trauma or injury. It may help to minimize neck pain and discomfort. For further details please consult your prescribing physician. Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICATIONS: Do not use if you allergic to any of the materials contained within this product.

WARNINGS AND PRECAUTIONS:

This product is to be prescribed and fitted by a physician. Do not use over open wounds. Do not use this device if it was damaged and/or package has been opened. The collar should be snug but not tight enough to cause excessive pressure, pain or difficulty breathing or swallowing. If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact a medical professional immediately. There could be difficulties with Laryngoscopy when collars are used and the device may have to be removed. **NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INFORMATION:

1. **MEASURE HEIGHT:** With the patient in NEUTRAL ALIGNMENT, measure the vertical distance from the tip of the chin to the top of the shoulder. The measurement corresponds to the HEIGHT of the collar. If the patient's measurement falls between two consecutive heights, try fitting the shorter height first. **NOTE:** Cervical region (head, neck and back) of the patient must be in neutral alignment before a measurement is taken. (Fig.1)
2. **LOCK HEIGHT:** Determine the correct height for patient and adjust collar height to proper position. Press the locking mechanism into collar to secure collar. **NOTE:** The two locks on either side of collar must be set at the same heights. Verify collar is in the locked position by attempting to readjust. If movement occurs, repeat procedure until the collar is secure.
3. **TO ADJUST HEIGHT:** Release locks by pushing each lock from the inside of the collar until it is in unsecured position. Determine the proper height for the patient and press the lock back into the collar until secured. (Fig.3)
4. After adjusting to the proper height, apply the collar to the patient's neck. The chin should be secured by the molded chin support. If the patient is in a supine position, slide the BACK part of the collar behind the patient's neck until the hook-and-loop strap is visible. (Fig.4)
5. Bring the hook-and-loop strap around the back of the patient's neck and secure to the hook-and-loop fastener. Before attaching the hook-and-loop strap, make certain the patient's head and collar are in proper alignment. Tighten collar if chin is not sufficiently supported. Select shorter height if additional tightening results in patient extension. (Fig.5)
6. To remove the collar, make sure the patient's head is in neutral position and unfasten the hook-and-loop strap. Slowly slide the collar from behind the neck and remove. (Fig.2)

CLEANING INSTRUCTIONS: Hand wash with mild soap in lukewarm water. Rinse thoroughly. AIR DRY.

MATERIAL CONTENTS:

RIGID SHELL: Polyethylene.
SOFTGOODS: Polyurethane, Woven Nylon.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

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