VENAFLOW[®] Elite

Operator's Manual

Prophylaxis for Deep Vein Thrombosis

REF 30BI, 30BI-B

AIRCAS^T

VENAFLOW[®] *Elite*

ENGLISH	1
ESPAÑOL	31
DEUTSCH.	63
ITALIANO	95
FRANÇAIS	127
NEDERLANDS.	159
PORTUGUÊS.	191
简体中文	223
TÜRKÇE.	255
한국어	287
SVENSKA	319

Table of Contents

Design Philosophy	3
Function	3
Product Description	4
User Profile	4
Intended Use/Indications	4
Contraindications.....	4
Warnings	5
Cautions	6
Latex Information.....	6
Contents	6
General System Information and Application	
Pump Features	7
Symbol Definitions	8
Pump Set-up	9
Cuff Application	10
Pump Operation	11
Single Leg Operation	11
Patient Compliance Counter Reset	11
System Alarms	
Alarm Reset.....	12
Tube Alarm.....	12
Enabling Compliance Alarm.....	12
Single/Dual Leg Operation Alarm	12
Compliance Alarm	13
Call For Service Alarm	13
Battery Information and Operation	
Battery Operation	14
Battery Charging Instructions.....	14
Lithium Iron Phosphate Battery Maintenance Guidelines.....	14
Overview.....	14
Battery Maintenance	15
Storage	15
Handling Precautions.....	15
Low Battery Alarm.....	16
Replacement Kit Instructions	
Instructions for Installing/Removing/Replacing the Power Cord.....	17
Telescoping Bed Hanger Replacement.....	17
Instructions for Replacing Fuses	17

Maintenance

Cleaning Instructions	18
Maintenance	18
Factory Service	18

Warranty	19
----------------	----

Optional Tube Alarm Verification Test	19
---	----

Ordering Information	20
----------------------------	----

VenaFlow Elite International Plug Ordering Matrix	21
---	----

Compliance Statements

Electromagnetic Compatibility (EMC)	22
---	----

Electromagnetic Compatibility (EMC) Tables - RF Emissions Class A

Guidance and Manufacturer's Declaration - Electromagnetic Emissions	23
---	----

Guidance and Manufacturer's Declaration - Electromagnetic Immunity	24
--	----

Guidance and Manufacturer's Declaration - Electromagnetic Immunity	25
--	----

Recommended separation distances between portable and mobile RF communications equipment and the VenaFlow Elite	26
---	----

Specifications

Classifications	27
-----------------------	----

Storage/Transportation Conditions	27
---	----

Operating Conditions	27
----------------------------	----

Customer Care Contact Information	28
---	----

Design Philosophy

Pneumatic compression is a clinically proven modality for minimizing the risks associated with deep vein thrombosis.¹ The Aircast VenaFlow Elite System leverages the same clinically proven technology as the existing VenaFlow offering, but is presented in a low profile, light-weight design. VenaFlow Elite combines two proven technologies, rapid inflation and graduated sequential compression, via our new single tube Integrated Graduated Sequential Flow System (IGSF), to accelerate venous velocity and enhance fibrinolysis. This unique combination makes VenaFlow the only technology platform on the market that produces blood flows that mimic those achieved through ambulation.² Additionally, VenaFlow Elite incorporates asymmetric compression for the superior emptying of veins.³ The breathable and comfortable VenaFlow Elite calf, foot and thigh cuffs enhance the patient experience, assist in increasing compliance and are compatible with one universal pump.

Function

There are cuffs for the:

- 1) Foot
- 2) Thigh and
- 3) Calf

Either a single cuff on one leg or a cuff on each leg is allowed.

Once activated, the VenaFlow Elite system begins the compression cycle and inflates the cuffs one leg at a time, alternating between the two legs every 30 seconds. First the distal aircell inflates rapidly within less than 0.5 seconds, then the proximal aircell follows.

Foot Cuffs: The pressure peaks at approximately 130mmHg +/- 10%.

Calf Cuffs: The distal pressure peaks at approximately 73mmHg +/-15%, and the proximal pressure peaks at approximately 63mmHg +/- 15%.


Thigh Cuffs: The distal pressure peaks at approximately 73mmHg +/-15%, and the proximal pressure peaks at approximately 63 mmHg +/- 15%.

These pressures then settle at 45mmHg ±10% and after 6 seconds, the cuff deflates. In 54 seconds, the inflation cycle begins again. An icon will appear on the graphical display signaling when the inflation cycle is in progress (See 'Pump Features' section). Pressure and the inflate/deflate cycle are automatic and an alarm will be activated if the system is not functioning properly (See 'System Alarms' section).

1. Labropoulos N, OH D.S, Golts, E, et al: Improved Venous Return By Elliptical, Sequential and Seamless Air-cell Compression. Loyola University Medical Center, January 2003.

2. Whitelaw G, Oladipo O, Shah BP, et al: Evaluation of Pneumatic Compression Devices. Boston Orthopedics, March 2001.

3. Kamm R: Unsteady Venous Blood Flow Resulting From Different Modes of External Compression Cambridge, MIT, 1996

Product Description

The VenaFlow Elite System is a device that is for the Prophylaxis of Deep Vein Thrombosis. This system includes a pump, cuff (applied part), and connector tubing. The pump provides graduated, sequential compression to the aircells contained in the cuff via the connector tubing.

User Profile

The intended user should be a licensed medical professional. The user should be able to:

- read and understand the operator's manual, warnings and cautions
- manually apply the compression cuffs
- sense auditory and visual signals

Intended Use/Indications

The VenaFlow Elite System is intended as a prophylaxis for deep vein thrombosis (DVT). The System is for Adult Use only. The VenaFlow Elite System is designed to be used only with Aircast VenaFlow Elite Compression Cuffs. Battery operation use only with DJO, LLC battery packs.

Contraindications

The VenaFlow Elite System should not be used by persons with known or suspected deep vein thrombosis, severe congestive heart failure, pulmonary edema, thrombophlebitis, severe arteriosclerosis, or active infection. Do not use on extremities which are not sensitive to pain, where cuff will interfere with gangrene, on patients with vein ligation or recent skin grafts, or extreme deformity of the leg. Do not use the VenaFlow Elite System where increased venous or lymphatic return is undesirable.

Warnings

- The VenaFlow Elite device is intended to be used under the supervision of a licensed medical practitioner and administrated by a medical professional within a hospitable or clinical environment.
- The VenaFlow Elite cuffs are designed for single patient use only.
- Before applying the VenaFlow System, the operator should be educated in DVT Prevention and VenaFlow operation.
- Risk of explosion if used in the presence of flammable anesthetics.
- The equipment or VenaFlow device should not be used adjacent to or stacked with other equipment and, if adjacent or stacked use is necessary, the equipment or VenaFlow device should be observed to verify normal operation in the configuration in which it will be used.
- Device is designed to comply with electromagnetic safety standards. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. Harmful interference to other devices can be determined by turning this equipment on and off. Try to correct the interference using one or more of the following:
 - Reorient or relocate the receiving device
 - Increase the separation between the equipment
 - Connect the equipment to an outlet on a different circuit from that to which the other device(s) are connected and consult the factory field service technician for help
 - Consult your authorized DJO dealer for help
- Care must be taken when operating this equipment around other equipment to avoid reciprocal interference. Potential electromagnetic or other interference could occur to this or to the other equipment. Try to minimize this interference by not using other equipment in conjunction with this device.
- To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth (ensure electrical outlet is grounded).
- The mains plug is used as the primary means of disconnect. Ensure mains outlet is approachable while the pump is in operation.
- Periodically check VenaFlow devices to ensure battery is functioning properly. Replace battery as necessary.
- To avoid risk of battery leakage, remove battery prior to storage of the VenaFlow Elite System or if the system will not be used for two months.
- No modification of this equipment is allowed.
- If you experience pain, swelling, sensation changes or any unusual reactions while using this device, consult your medical professional immediately.

Warnings (Continued)

- Ingress of water or liquids could result in unpredictable failures
- Inspect equipment for obvious damage prior to use
- Do not walk with tubing connected to cuffs
- The device should never be opened by anyone other than a trained technician
- Use of device outside all recommended specifications could result in unpredictable failures.

Cautions

- Medical Electrical Equipment needs special precautions regarding EMC. Portable and mobile RF communication equipment can be affected by other medical electrical devices. If you believe interference is occurring, please consult Electromagnetic Compatibility (EMC) section.
- Cuffs used in combination with warming devices may cause skin irritation. Regularly check for patient comfort, compliance and skin irritation.
- Special attention should be given to those patients with neuropathies or tissue viability problems (i.e. diabetes, arterial or venous insufficiencies).
- To prevent extremity compartment syndrome, special attention should be given to patients who are positioned in the supine lithotomy position for extended lengths of time. This includes patients with or without compression cuffs.
- Contact customer service for information concerning the decommissioning of your equipment.

Latex Information

All components of the VenaFlow Elite System are not made with natural rubber latex. All VenaFlow Elite cuffs are not made with natural rubber latex and may be placed directly against the skin or over a light compression dressing.



Contents

Each package contains:

- System
- Power Cord
- Operator's Manual

VENAFLOW[®]Elite

Pump Features

- 1) Inflation cycle indicator - Indicates cuff inflation.
- 2) First 3 minutes (ramp up) indicator - Indicates 0-3 minutes allowed to select single leg operation.
- 3) Patient compliance counter - Displays treatment time in hours, minutes and seconds.
- 4) Patient compliance counter reset push button - Resets the patient compliance counter.
- 5) Single / Dual Leg operation indicator - Indicates single leg or dual leg operation is active.
- 6) Battery indicator* - Indicates battery charge.
- 7) Single / Dual Leg operation push button - Selects single leg or dual leg operation.
- 8) ON / OFF / RESET push button - Turns the device on or off and resets the alarms.



* Note: not all VenaFlow Elite Systems will contain a battery

Symbol Definitions


	Power on/off button and alarm reset		Manufacturer
	Dual or Single leg operation		Date of Manufacturer
	Compliance monitor reset		Dual calf/thigh
	Tube alarm		Single calf/thigh
	Call for service		Dual foot
	Inflation cycle indication (Bars indicate cuffs are inflating)		Single foot
	Type BF Equipment		Visual only compliance alarm enabled
	F2AL250V		Visual and audio compliance alarm enabled
	Battery		Warning or Caution
	Waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.		Non-ionizing electromagnetic radiation
	CE Mark of conformity with notified body number		Conforms to AAMI Std ES60601-1 Certified to CSA Std C22.2 No. 60601-1
	Catalogue number		
	Refer to Instruction Manual/Booklet		
	Keep dry		
	Fragile, handle with care		
	This end up		
	Temperature range		
	Humidity range		
	Atmospheric pressure range		
	Federal Law restricts this unit to sale by, or on the order of, a physician or licensed practitioner.		
	Authorized representative in the European Community		

Pump Set-up

- 1) Hang pump from bed frame (foot of bed), bed rail, or rest on floor or table.
To use the telescoping bed hanger, press release button on the back of the device and gently pull out bed hook to the desired width.

- 2) Connect tube assembly to the pump.
Make sure tubing connectors lock securely into pump connectors.


- 3) Plug in the power cord to electrical outlet.

Cuff Application

The VenaFlow Elite System will automatically detect the cuff that is attached to the device and apply appropriate pressures.


- 1) **Calf cuff** – Apply cuff with the tube pointed toward the foot. The aircell may be placed either on the back, side or front of the leg. It doesn't matter where the aircell sits on the leg. It will perform the same.


Foot cuff – Apply cuff with the aircell centered on bottom of foot and with the tube pointed to the left


Thigh cuff – Apply cuff with the distal and proximal aircells centered on either the back, side or front of the calf and thigh with the tube pointed toward the foot. It doesn't matter where the aircell sits on the leg. It will perform the same.

- 2) Connect cuffs to the tubing that is attached to the pump.

For varied patient sizes, trim calf and thigh cuffs as needed for proper fit.

Be sure to secure cuff straps snug but not tight. When using the calf cuff and thigh cuff, rotation of the aircell on the leg does not affect performance of the system.

Pump Operation

- 1) To turn the device on, press . The graphical display, green pump indicator lights and green light above  button will turn on.
- 2) Once device is powered on, system will immediately enter cuff detection mode which means the system will be detecting whether or not there are cuffs attached. Display will read "DETECTING CUFF" and will display percentage of detection complete. Once it reaches 100%, the display will then transition to the standard display screen.
- 3) To turn the device off, press . The graphical display, green pump indicator lights and green light above  button will turn off.
- 4) To shut down the power completely, disconnect the power plug by removing it from the AC socket.

Single Leg Operation

- 1) After system is turned on, within 3 minutes push the single leg operation button  to indicate that single leg use is desired. The single leg icon will appear at the top right corner of the graphical display.
- 2) If the single leg operation mode is not selected within 3 minutes of the system powering up and only one cuff is connected, the "LEG ALARM" will be activated.
- 3) If the single leg operation mode is selected within 3 minutes of powering on and two cuffs are connected instead of one, a "LEG ALARM" will be activated.
- 4) In single leg operation mode, the device inflates detected cuff every minute for 6 seconds.
- 5) Either port may be used for single leg operation. The system will automatically detect which one is in use.
- 6) To change from dual to single leg or single to dual leg operation mode after the 3 minute ramp up period, press the button to turn off the device, then again to power on again. The device will go into ramp up mode and user will have 3 minutes to change the dual leg operation to single leg mode.

Patient Compliance Counter Reset


To reset the patient compliance counter, press  for 1 second and release. The hours, minutes and seconds will reset.

System Alarms

Alarm Reset

To reset any alarm, press the  button and take steps to correct the alarm if necessary. Pressing the button does not fix the alarm but silences it until the problem is fixed. If you wish to turn the device completely off, you must press the power button once to silence alarm and a second time to power down the device.

Tube Alarm

- 1) If one or both tubes has a kink or leak, after 4-6 minutes this symbol  and the text "CHECK TUBES" will alternate on the graphical display, an audible alarm will sound and the pump indicator lights on the side of the system will flash red.
- 2) If a tube alarm occurs, reset alarm and then ensure that the tube connections are secure and that the tubing is kink-free. Check pump connectors and tubing/cuff connectors for cracks or damage. Check cuffs for damage as well.
- 3) If alarm occurs again, replace cuffs and tubing altogether.

Enabling Compliance Alarm

A compliance alarm can be enabled on the device as an option. To enable compliance feature, complete the following steps:

- 1) Start with the system powered off.
- 2) Press and hold the single/dual leg operation button  and while holding, press the power button. The graphical display screen will go blank with no backlight.
- 3) Press the single/dual leg operation button (insert button).
- 4) Select button underneath 'Change'.
- 5) Select button underneath 'Other'.
- 6) Select button under 'No alarm'.
- 7) Select the button under 'Done' to set the new configuration. The change is complete and the system is ready to operate. A "C" will appear on the display indicating that the compliance alarm is enabled.

Single/Dual Leg Operation Alarm

- 1) Once powered on, the VenaFlow Elite System will default to dual leg operation mode. In the event that the selected leg operation mode and the number of cuffs attached do not match, after 4 minutes, an audible alarm will sound and this symbol  will appear on the graphical display along with the text, "LEG

ALARM". Additionally, the pump indicator lights on the side of the system will flash red.

- 2) A leg alarm will occur if either of the following occur:
 - Single leg use is desired, one cuff is connected, but the single leg operation button is not pressed within 3 minutes after powering on the device.
 - Single leg operation button was pressed within 3 minutes of powering on the device, but two cuffs are connected.
- 3) When a single/dual operation "LEG ALARM" occurs, reset alarm, and then ensure that the icon at the top right corner of the graphical display is consistent with the number of cuffs that are connected to the tubing and adjust accordingly.

Compliance Alarm

- 1) If the cuffs are removed from the patient while the device is still functioning and the cuffs are not reapplied after 15 minutes, the system will go into alarm mode and the graphical display will read "COMPLIANCE ALARM". Additionally, the pump indicator lights on the side of the system will flash red.
- 2) Reset alarm first, and then reapply cuffs. You must reapply the cuffs within 9-15 minutes after resetting the alarm in order to prevent another compliance alarm from occurring.
- 3) If cuffs are reapplied without resetting the alarm, the compliance alarm will automatically reset after 9-15 minutes

Call For Service Alarm

- 1) If one of the system's critical components is not working properly, the Call For Service alarm symbol will appear on the graphical display and it will alternate with the text "CALL FOR SERVICE". Additionally, the pump indicator lights on the side of the system will flash red.
- 2) If Call For Service alarm occurs, unplug the device and call DJO Technical support. See Customer Care Contact Information section.

Battery Information and Operation

Battery Operation

The VenaFlow Elite can be installed with a battery upon request at the time of the initial system order. When the battery-installed system is powered on, a battery icon will appear on the right side of the graphical display beneath the single/dual leg mode icon. If the system does not have a battery installed, then there will be no battery icon on the graphical display at all. When the system with battery is unplugged and powered on, the system will automatically switch to battery power.

The battery must be charged before use. Once a system has been configured to run with a battery, it must always have a battery installed.

See 'Specifications' section for battery pack information.

Battery Charging Instructions

- 1) To charge the battery, simply plug in the system to an electrical outlet either while the system is on and functioning, or when it is powered off. The system will begin charging immediately after being plugged in.
- 2) The VenaFlow Elite System with a fully charged battery will last for approximately 2 - 4 hours. It takes an estimated 2 hours to fully charge the battery when the unit is off. When the unit is on, it takes an estimated 2 hours to charge the battery.

Lithium Iron Phosphate Battery Maintenance Guidelines

Overview

Do not leave batteries unused for extended periods of time because Lithium Iron Phosphate batteries continue to slowly discharge (self-discharge) when not in use or while in storage. Routinely check the battery's charge status. The product Operator's manual includes information on how to check battery status as well as battery charging instructions. When a battery has been unused for 6 months, check the charge status and charge or send the system back to DJO for battery replacement.

The typical estimated life of a Lithium Iron Phosphate battery is about three to five years or approximately 1,000 charge cycles, whichever occurs first. One charge cycle is a period of use from fully charged, to fully discharged, and fully recharged again. For batteries that do not run through complete charge cycles, there is a three to five year life expectancy. The intended operating temperature of the device is between 10°C and 40°C.

Rechargeable Lithium Iron Phosphate batteries have a limited life and will gradually lose their capacity to hold a charge. This loss of capacity (aging) is irreversible. As the battery loses capacity, the length of time it will power the product (run time) decreases.

Battery Maintenance

You should routinely check the battery's charge status. Carefully monitor batteries that are approaching the end of their estimated life. Consider shipping the device back to DJO to replace the battery if you note either of the following conditions:

1. The battery run time drops below about 80% of the original run time.
2. The battery charge time increases significantly.

If a battery is stored or otherwise unused for an extended period, be sure to follow the storage instructions in this document.

Storage

Charge or discharge the battery to approximately 50% of capacity before storage. Charge the battery to approximately 50% of capacity at least once every three months.

The battery is best stored at $23\pm 5^{\circ}\text{C}$. For periods of less than 1 month, the battery can be stored between -25°C to 70°C .

NOTE: The battery self-discharges during storage. Higher temperatures (above 20°C or 68°F) reduce the battery storage life.

Handling Precautions

- Do not disassemble, crush, or puncture a battery. Do not short the external contacts on a battery.
- Do not dispose of a battery in fire or water.
- Do not expose a battery to temperatures above 70°C (158°F). Keep the battery away from children.
- Avoid exposing the battery to excessive shock or vibration. Do not use a damaged battery.
- If a battery pack has leaking fluids, do not touch any fluids. Properly dispose of a leaking battery pack.
- In case of eye contact with fluid, do not rub eyes. Immediately flush eyes thoroughly with water for at least 15 minutes, lifting upper and lower lids, until no evidence of the fluid remains.
- Seek medical attention.

Low Battery Alarm

- 1) When there is approximately 10 minutes of charge left in the battery, an audible alarm will beep for 10 seconds per minute and a battery icon will appear on the graphical display and alternate with the text "LOW BATTERY". When there is approximately 5 minutes left, two 10 second beeps will sound per minute in addition to the graphical alert.
- 2) When either battery alarm occurs, immediately plug device in to electrical outlet to begin replenishing the charge.
- 3) If the device is not plugged in during the alarm period, the device will shut down and will fail to power on again until the device is plugged in again.


Replacement Kit Instructions

Replacement kits may be ordered through customer care.


Instructions for Installing/Removing/ Replacing the Power Cord

- 1) Make sure system is powered off and unplugged.
- 2) Remove the screw retaining the power cord cover using a 1/16" hex driver.
- 3) Remove the power cord cover.
- 4) Remove the power cord.
- 5) Install new power cord and reassemble power cord cover.


Telescoping Bed Hanger Replacement

- 1) The telescoping bed hanger piece can be removed and replaced if necessary.
- 2) To replace the bed hanger, hold down release button on the back of the system and firmly pull bed hanger out of the body of the system.
- 3) Discard the broken or ineffective hanger. To install the new bed hanger, slide the new hanger back into the body of the device and set at desired length.

Instructions for Replacing Fuses

- 1) Remove the power cord. Refer to the 'Instructions for installing/removing/replacing the power cord' section.
- 2) Remove fuse holder from power entry module.
- 3) Remove fuses from fuse holder.
- 4) Replace any blown fuses with 250V, 2A, IEC fast acting 5X20mm fuses.
- 5) Reassemble fuse holder and replace power cord per previous instructions.

Cleaning Instructions

The tube assembly is reusable. Both the pump case and the tube assembly can be cleaned with mild soap solution, antiseptic or disinfectant wipes. Apply using a dampened cloth or wipes. Do not submerge the device in any liquid.

No adverse effects were found using isopropyl alcohol and hydrogen peroxide (Oxivir Tb). It is advised to avoid using Dimethyl benzyl ammonium chloride or Dimethyl ethylbenzyl ammonium chloride (Tec-Surf II), or solutions containing chlorine.

Maintenance

Prior to each use, inspect the device, power cord, cuff and tubing for damage. If damage is detected, contact Customer Service. See Customer Care Contact Information section.

NOTE: The unit was tested during the manufacturing process and is ready to be placed into service upon delivery.

Factory Service

When the VenaFlow Elite System requires factory service, contact Customer Service. See Customer Care Contact Information section.

Warranty

For consigned systems: DJO, LLC will repair or replace all or part of the 30BI and 30BI-B systems for material or workmanship defects for the life of the customer contract. Warranty and maintenance terms are specified in each consignment contract.

For purchased systems: DJO, LLC will repair or replace all or part of the 30BI or 30BI-B systems for material or workmanship defects for 3 years from the date of sale.

For disposable cuffs, tubing and battery pack purchases: DJO, LLC will repair or replace all or part of the cuffs, tube units and battery packs for material or workmanship defects, not including normal depreciation of battery life or charge capacity, for a period of six months from the date of sale.

Optional Tube Alarm Verification Test

The Tube Alarm Verification Test is provided for facility verification purposes only and is not required prior to use.

Attach tube assembly to pump and connect power cord to electrical outlet. Simulate a kinked tube situation by folding the tube over onto itself once and holding it for 4-6 minutes. The second cycle after the simulated kink should activate an audible alarm and this symbol \triangle will appear with "CHECK TUBES" on the graphical display and the pump indicator lights on the side of the system will flash red.

Ordering Information

System:

Part #	Description
30BI	VenaFlow Elite System, International
30BI-B	VenaFlow Elite System, International with Battery

Accessories:

Part #	Description	Qty	Maximum Calf Circumference
3040	VenaFlow Elite Calf Cuff	Pair	48 cm (19 in.)
3042	VenaFlow Elite XL Calf Cuff	Pair	56 cm (22 in.)
3043	VenaFlow Elite Bariatric Calf Cuff	Pair	76 cm (30 in.)
3045	VenaFlow Elite Thigh Cuff	Pair	One size fits all
3046	VenaFlow Elite Foot Cuff	Pair	One size fits all
3050	VenaFlow Elite Foam Calf Cuff	Pair	48 cm (19 in.)
3008	Tube Assembly, 167.64 cm (66 in.)	Each	
3008XL	Tube Assembly, 259.08 cm (102 in.)	Each	
3008XXL	Tube Assembly, 320.04 cm (126 in.)	Each	
3008XXXL	Tube Assembly, 381 cm (150 in.)	Each	

System Replacement Parts

Part #	Description
3071	Bed Hanger
3072	Tube Attachment Tag
3073	Fuse

For powercords, see VenaFlow Elite International Plug Ordering Matrix.

VenaFlow Elite International Plug Ordering Matrix

Country	Language	Power		Plug	Part number
		Voltage	Freq.		
United Kingdom	English	230 V	50 Hz	BS1363	3048
Germany	German	230 V	50 Hz	CEE7/7	3049
Belgium	German French Dutch	230 V	50 Hz	CEE7/7	3049
Austria	German	230 V	50 Hz	CEE7/7	3049
Denmark	Danish	230 V	50 Hz	Danish 107-2-D1	3051
France	French	230 V	50 Hz	CEE7/7	3049
Netherlands	Dutch	230 V	50 Hz	CEE7/7	3049
Iberica	Portuguese Spanish	230 V	50 Hz	CEE7/7	3049
Spain	Spanish	230 V	50 Hz	CEE7/7	3049
Portugal	Portuguese	220 V	50 Hz	CEE7/7	3049
Italy	Italian	230 V	50 Hz	CEE7/7 CEI23-16/VII	3049
South Africa	English	220 V	50 Hz	BS546	3052
Malaysia	English Malay	240 V	50 Hz	BS1363 BS546	3048 3052
Thailand	English Thai	220 V	50 Hz	CEE7/7 UL-NEMA5-15 (250V)	3049 3047
Greece	Greek	230 V	50 Hz	CEE7/7	3049
Hong Kong	English	220 V	50 Hz	BS1363 BS546	3048 3052
India	English	230 V	50 Hz	BS546	3052
Jordan	English Arabic	230 V	50 Hz	UL-NEMA5-15 (250 V) BS1363 CEE7/7 BS546	3047 3048 3049 3052
Turkey	Turkish	230 V	50 Hz	CEE7/7	3049
Saudi Arabia	English Arabic	110 V 220 V	60 Hz	UL-NEMA5-15 UL-NEMA5-15 (250V) BS1363 CEE7/7	3070 3047 3048 3049
UAE	English Arabic	220 V	50 Hz	BS1363 BS546	3048 3052
Qatar	English Arabic	240 V	50 Hz	BS1363 BS546	3048 3052
Japan	Japanese	100 V	50 Hz 60 Hz	UL-NEMA5-15 UL-NEMA5-15 (250V)	3070 3047
China(Main land)	Mandarin (Chinese)	220 V	50 Hz	AS3112	3053
Canada	English	120 V	60 Hz	UL-NEMA5-15	3047

Compliance Statements

Electromagnetic Compatibility (EMC)

The VenaFlow Elite has been tested and found to comply with the electromagnetic compatibility (EMC) limits for medical devices to IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

Caution: Medical electrical equipment requires special precautions regarding EMC and must be installed and operated according to these instructions. It is possible that high levels of radiated or conducted radio-frequency electromagnetic interference (EMI) from portable and mobile RF communications equipment or other strong or nearby radio-frequency sources, could result in performance disruption of the system. Evidence of disruption may include image degradation or distortion, erratic readings, equipment ceasing to operate, or other incorrect functioning. If this occurs, survey the site of disruption, and take the following actions to eliminate the source(s).

- Turn equipment in the vicinity off and on to isolate disruptive equipment.
- Relocate or reorient interfering equipment.
- Increase distance between interfering equipment and your system.
- Manage use of frequencies close to the system frequencies.
- Remove devices that are highly susceptible to EMI.
- Lower power from internal sources within the facility control (such as paging systems).
- Label devices susceptible to EMI.
- Educate clinical staff to recognize potential EMI-related problems.
- Eliminate or reduce EMI with technical solutions (such as shielding).
- Restrict use of personal communicators (cell phones, computers) in areas with devices susceptible to EMI.
- Share relevant EMI information with others, particularly when evaluating new equipment purchases which may generate EMI.
- Purchase medical devices that comply with IEC 60601-1-2 EMC Standards (3V/

meter EMI immunity, limit interference level to 0.0014 V/meter).

Electromagnetic Compatibility (EMC) Tables - RF Emissions Class A

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The VenaFlow Elite is intended for use in the electromagnetic environment specified below. The customer or the user of the VenaFlow Elite should assure that it is used in such an environment.		
Emissions Tests	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1	The VenaFlow Elite uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The VenaFlow Elite is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The VenaFlow Elite is intended for use in the electromagnetic environment specified below. The customer or the user of the VenaFlow Elite should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/ Burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VenaFlow Elite requires continued operation during power mains interruptions, it is recommended that the VenaFlow Elite be powered from an uninterruptible power supply or a battery.
Power Frequency (50/60Hz) Magnetic Fields IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The VenaFlow Elite is intended for use in the electromagnetic environment specified below. The customer or the user of the VenaFlow Elite should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the VenaFlow Elite, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_T} \right] \sqrt{P}$ 150 KHz to 80 MHz $d = \left[\frac{3.5}{E_T} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_T} \right] \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VenaFlow Elite is used exceeds the applicable RF compliance level above, the VenaFlow Elite should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VenaFlow Elite.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the VenaFlow Elite

The VenaFlow Elite is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the VenaFlow Elite can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VenaFlow Elite as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80 MHz $d = \left[\frac{3,5}{f_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{f_1} \right] \sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{f_1} \right] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Specifications

Size: 19.3 cm (7.6 in.) x 21.6 cm (8.5 in.) x 11.7 cm (4.6 in.)

Weight: no battery: 1.8 kg (4 lbs.) with battery: 2.2 kg (4.75 lbs.)

Power input: 100 - 240~, 50/60 Hz, 75 VA

Cord: Hospital grade, 4.572 m (15 ft.)

Standards: IEC 60601-1:2005, IEC 60601-1-2:2007

Inflation mode: rapid inflation and graduated, sequential compression

Bed hanger: accommodates up to 9 cm (3.5 in.)

Battery: LiFePO₄ 12.8V 1500 mAh

Fuse: F2AL250V, Quantity 2

Classifications

Degree of Protection from Electrical Shock: Class 1

Type BF Applied parts (Cuffs)

IPX0 (Ordinary protection)

No sterilization is required

This device is not intended for use in an Oxygen Rich Environment

Mode of operation is Continuous Operation

Storage/Transportation Conditions

Temperatures: -13° F (-25° C) - 158° F (70° C)

Relative Humidity: 15% - 93%

Atmospheric Pressure: 500 hPa - 1060 hPa

Operating Conditions

Operating Temperatures: 10°C to 40°C

Relative Humidity: 30% to 75%

Atmospheric Pressure: 700 hPa to 1060 hPa

Altitude: below 3000 m

Customer Care Contact Information

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Tel: +32 (0) 14 24 80 90
Email: benelux.orders@DJOglobal.com

CANADA:

DJO Canada
6485 Kennedy Road
Mississauga

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Tel: +49 761 4566 251
Fax: +49 761 4566 55251
Email: infoservice@DJOglobal.com

INDIA:

DJO India
Regus RMZ Milenia Business Park 2
Campus 4B, 6th Floor, No 143,
Dr. MGR Road
Perungudi, Chennai-600 096

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Tel: +39 02 484 63386
Fax: +39 02 484 09217
Email: vendite@DJOglobal.com

SOUTH AFRICA:

DJO South Africa (Pty) Ltd
Unit 1, Brackengate Business Park
5 on London
Brackenfell, 7560
Cape Town

SOUTH AFRICA

Tel: +27 (0) 87 3102480
Fax: +27 (0) 86 6098891
Email: info.southafrica@DJOglobal.com

SPAIN:

DJO Ibérica
Carretera de Cornellá 144 1^o4^o
Esplugues de Llobregat
08950 Barcelona

SPAIN

Tel.: +34 93 480 32 02
Fax: +34 93 473 36 67
Email: es.info@DJOglobal.com

UK & IRELAND:

DJO UK Ltd.
1a Guildford Business Park
Broad Street Guildford Surrey
GU2 8XG

UK

Tel: +44 (0)1483 459 659
Fax: +44 (0)1483 459 470
Email: ukorders@DJOglobal.com

UNITED STATES:

DJO, LLC
5919 Sea Otter Place,
Suite 200
Carlsbad, CA 92010

U.S.A.

Tel: 1 800 336 6569
Fax: 1 800 936 6569

DJO Export

ASIA-PACIFIC:

DJO Asia-Pacific
Limited Unit 1905, 19/F
Tower II Grand Central Plaza
138 Shatin Rural Committee Road
Shatin

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
1. Labropoulos N, OH D.S, Golts, E, et al: Improved Venous Return By Elliptical, Sequential and Seamless Air-cell Compression. Loyola University Medical Center, enero de 2003.
2. Whitelaw G, Oladipo O, Shah BP, et al: Evaluation of Pneumatic Compression Devices. Boston Orthopedics, marzo de 2001.
3. Kamm R: Unsteady Venous Blood Flow Resulting From Different Modes of External Compression Cambridge, MIT, 1996.

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
1. Labropoulos N, OH D.S, Golts, E, et al: Improved Venous Return By Elliptical, Sequential and Seamless Air-cell Compression. Loyola University Medical Center, January 2003.
2. Whitelaw G, Oladipo O, Shah BP, et al: Evaluation of Pneumatic Compression Devices. Boston Orthopedics, March 2001.
3. Kamm R: Unsteady Venous Blood Flow Resulting From Different Modes of External Compression Cambridge, MIT, 1996.

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
1. Labropoulos N, OH D.S, Golts, E 외: Improved Venous Return By Elliptical, Sequential and Seamless Air-cell Compression. Loyola University Medical Center, January 2003.

2. Whitelaw G, Oladipo O, Shah BP 외: Evaluation of Pneumatic Compression Devices. Boston Orthopedics, March 2001.

3. Kamm R: Unsteady Venous Blood Flow Resulting From Different Modes of External Compression Cambridge, MIT, 1996.

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
1. Labropoulos N, OH D.S, Golts, E, et al: Improved Venous Return By Elliptical, Sequential and Seamless Air-cell Compression. Loyola University Medical Center, January 2003.
2. Whitelaw G, Oladipo O, Shah BP, et al: Evaluation of Pneumatic Compression Devices. Boston Orthopedics, Mars 2001.
3. Kamm R: Unsteady Venous Blood Flow Resulting From Different Modes of External Compression Cambridge, MIT, 1996.

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MDSS GmbH