CLEARCUT®

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RECOMMANDATIONS DE MISE AU REBUT

L'élimination du dispositif ou de ses composants doit s'effectuer en accord avec la réglementation locale des déchets en vigueur.
PRECAUTIONS OF USE
The surgeon should be trained in the surgical techniques used for these cement restrictors. The insertion of a foreign body in the tissue increases normal risks of infection associated with surgery during the postoperative period.
PRECAUTIONS OF USE
The surgeon should be trained in the surgical techniques used for these cement restrictors. The insertion of a foreign body in the tissue increases normal risks of infection associated with surgery during the postoperative period.
CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
Before using DJO Surgical products, the operator must thoroughly familiarise himself with the safety regulations figuring in the directions for use, as well as the information relating to each product (description, operating technique, leaflets, etc.). The corresponding information is available from DJO Surgical. The surgeon should also be informed of the residual risk associated with the product he considers using.
ALL-GENERAL INFORMATION
The implantation of DJO Surgical products should only be performed by qualified operators having a sound knowledge and full mastery of operative techniques specific to DJO Surgical products. Operative techniques may be acquired from the DJO Surgical distributors. The surgeon is responsible for any complications or harmful consequences which might result from an erroneous indication or operative technique, improper use of the equipment and failure to comply with the safety regulations provided in the directions for use. Neither the DJO Surgical manufacturer nor the authorized DJO Surgical representative can be held responsible for these complications.
Only combination of DJO Surgical products belonging to the system and implantation with corresponding DJO Surgical instruments are permitted. The use of DJO Surgical instruments for other purposes is not allowed.
INDICATIONS
CLEARCUT® cement restrictor, size 8 to 18, is a diaphyseal plug designed to occlude the medullary cavity before the introduction of acrylic cement during total hip arthroplasty.
CLEARCUT®, size 6, is a diaphyseal plug designed to occlude the medullary cavity before the introduction of acrylic cement during total elbow arthroplasty.*
CLEARCUT® prevents the cement from flowing through the diaphysis and therefore facilitates cement pressurization. Totally biocompatible, CLEARCUT® is completely reabsorbed within several days at 37°C and resists the temporary increase in temperature generated by the polymerization of acrylic cement.
***Size 6 is not covered by CE mark**
COMPOSITION
Gelatine (porcine-based)
Glycerol
Water
Methylparahydroxybenzoate
These components comply with European Pharmacopoeia standards.
INSTRUCTIONS FOR USE
The diameter of the CLEARCUT® can be predetermined by the surgeon by setting up the correspondence between the supplied templates of the prosthesis and diameters of the medullary canal.
Cement restrictor size 8 to 18 :
1 - Select the sizing trial device corresponding to the medullar canal diameter, screw it on the inserter and check that it is correctly secured.
2 - Lock the set screw.
3 - Introduce the sizing trial device into the medullar canal to check for a proper fit.
4 - Introduce the cement restrictor matching the sizing trial.
5 - Unlock the sizing trial from the inserter
Attach the restrictor to the inserter.
6 - Introduce the cement restrictor into the medullar cavity to the predetermined distance (distance corresponding to the distal tip of the stem and the greater trochanter. Add 1 cm to this distance).

CLEARCUT®

Bio-absorbable cement restrictor

Directions for use

Distributed by :

djosurgical
DJO Surgical
9800 Metric Blvd.
Austin, TX 78758

Manufactured in France by :

TEKNIMED S.A.S
8, rue du Corps Franc Pommies
65500 VIC EN BIGORRE - FRANCE
Tél. (33)5 62 96 88 38
Fax (33)5 62 96 28 72
www.teknimed.com

First EC marking date applied: 1998

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Before using DJO Surgical products, the operator must thoroughly familiarise himself with the safety regulations figuring in the directions for use, as well as the information relating to each product (description, operating technique, leaflets, etc.). The corresponding information is available from DJO Surgical. The surgeon should also be informed of the residual risk associated with the product he considers using.

GENERAL INFORMATION

The implantation of DJO Surgical products should only be performed by qualified operators having a sound knowledge and full mastery of operative techniques specific to DJO Surgical products. Operative techniques may be acquired from the DJO Surgical distributors. The surgeon is responsible for any complications or harmful consequences which might result from an erroneous indication or operative technique, improper use of the equipment and failure to comply with the safety regulations provided in the directions for use. Neither the DJO Surgical manufacturer nor the authorized DJO Surgical representative can be held responsible for these complications.

Only combination of DJO Surgical products belonging to the system and implantation with corresponding DJO Surgical instruments are permitted. The use of DJO Surgical instruments for other purposes is not allowed.

INDICATIONS

CLEARCUT® cement restrictor, size 8 to 18, is a diaphyseal plug designed to occlude the medullary cavity before the introduction of acrylic cement during total hip arthroplasty.
CLEARCUT®, size 6, is a diaphyseal plug designed to occlude the medullary cavity before the introduction of acrylic cement during total elbow arthroplasty.*
CLEARCUT® prevents the cement from flowing through the diaphysis and therefore facilitates cement pressurization. Totally biocompatible, CLEARCUT® is completely reabsorbed within several days at 37°C and resists the temporary increase in temperature generated by the polymerization of acrylic cement.

***Size 6 is not covered by CE mark**

COMPOSITION

Gelatine (porcine-based)
Glycerol

Water
Methylparahydroxybenzoate
These components comply with European Pharmacopoeia standards.

INSTRUCTIONS FOR USE

The diameter of the CLEARCUT® can be predetermined by the surgeon by setting up the correspondence between the supplied templates of the prosthesis and diameters of the medullary canal.
Cement restrictor size 8 to 18 :
1 - Select the sizing trial device corresponding to the medullar canal diameter, screw it on the inserter and check that it is correctly secured.
2 - Lock the set screw.
3 - Introduce the sizing trial device into the medullar canal to check for a proper fit.
4 - Introduce the cement restrictor matching the sizing trial.
5 - Unlock the sizing trial from the inserter
Attach the restrictor to the inserter.
6 - Introduce the cement restrictor into the medullar cavity to the predetermined distance (distance corresponding to the distal tip of the stem and the greater trochanter. Add 1 cm to this distance).

Note : a distance of 1cm is advised between the distal tip of the prosthesis and the restrictor.

4 - Unlock the set screw..

5 - Pull the handle upwards to leave the restrictor in situ.

The inserter may now be withdrawn, the cement restrictor is now positioned within the canal at the predetermined distance.

- Introduce the acrylic bone cement.

Cement restrictor size 6 :

- Introduce the sizing trial of the inserter into the medullar canal to check for a proper fit.

- Attach the cement restrictor at the other end of the inserter.

- Introduce the cement restrictor into the medullar cavity to the predetermined distance (distance corresponding to the tip of the stem and the prepared bone site. Add 1 cm to this distance).

Note : a distance of 1cm is advised between the distal tip of the prosthesis and the restrictor.

- Remove the inserter and leave the restrictor in situ.
- Introduce acrylic bone cement.

Caution : the instrument set must not be forced.

PRECAUTIONS FOR USE

The surgeon must be trained in the surgical techniques used for these cement restrictors. The insertion of a foreign body in the tissue increases normal risks of infection associated with surgery during the postoperative period.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

SIDE EFFECTS

No undesirable side-effects have ever been recorded.
Use during pregnancy and lactation : no contraindications.
Interaction with other medicaments : none known.

CONTRAINDICATIONS

The CLEARCUT® cement restrictor must not be used in patients allergic to the product compounds or in patients having known allergies.

STERILISATION

CLEARCUT® is sterilised by gamma radiation at a minimum dose of 25 kGy.
Prior to using, carefully check the protective packaging in order to ensure that it has not been damaged in a way that could affect its sterility.
When taking the product out of its packaging, be sure to follow the asepsis rules.
CLEARCUT® is delivered sterile and ready for use in the operating room.
Re-sterilization of the product is strictly prohibited.
Do not use after the use-by date.
Disposable. Do not use if the packaging is damaged.

PATIENT INFORMATION

The patient should be informed by the doctor of the potential consequences of the factors mentioned in the following paragraphs: contra-indications and side effects, that is, those liable to hinder the success of the operation, as well as possible complications which may arise. The patient should also be informed of the measures to be taken to diminish the possible consequences of these factors.

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CLEARCUT® is available in the following seven sizes: 6 - 8 - 10 - 12 - 14 - 16 - 18 mm.
Each size corresponds to the diameter of the medullary canal to occlude.
Store in its original packaging in a dry place below 30°C.

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