Orthopaedic Instruments

Instrument Care, Cleaning,
Disinfection and Sterilization

**INSTRUMENT CARE, CLEANING,
DISINFECTION AND STERILIZATION**
(in accordance with ISO 17664) – Rev. 02/2023

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INSTRUMENT CARE, CLEANING, DISINFECTION AND STERILIZATION (in accordance with ISO 17664) – Rev. 02/2023

Table of Contents

1.0	INTRODUCTION	pag.	4
2.0	CONSIDERATIONS	pag.	4
3.0	GENERAL WARNINGS, PRECAUTIONS AND RESTRICTIONS	pag.	6
4.0	INSPECTION	pag.	6
5.0	CLEANING AND DISINFECTION	pag.	7
	5.1 Restrictions	pag.	7
	5.2 Preparation at the point of use prior to processing	pag.	8
	5.3 Cleaning agents preparation	pag.	9
	5.4 Cleaning and Disinfection instructions	pag.	9
	5.5 Manual Cleaning Procedure	pag.	10
	5.6 Manual/Automated Cleaning/Disinfection Procedure	pag.	11
6.0	INSPECTION, MAINTENANCE AND LUBRIFICATION	pag.	13
7.0	STERILIZATION PROCEDURE	pag.	14
	7.1 Packaging	pag.	14
	7.1.1 Method of double wrapping ANSI/AAMI ST79:2017	pag.	15
	7.2 Sterilization procedure	pag.	16
8.0	STORAGE	pag.	18
9.0	HOSPITAL RESPONSIBILITIES FOR LIMACORPORATE LOANER INSTRUMENTS	pag.	18
10.0	FURTHER INFORMATION	pag.	19
11.0	SYMBOLGY	pag.	20
12.0	BIBLIOGRAPHY	pag.	21

1.0 INTRODUCTION

LimaCorporate instrumentation includes instruments according to the following classifications:

- class I and class IIa (according to the D.L.24/02/97 n.46 “Attuazione della direttiva 93/42/CEE concernente i dispositivi medici”) used in the implantation procedures of LimaCorporate prosthesis (Not applicable for U.S.);
- class I, Ir and class IIa (according to Regulation (EU) 2017/745, Not applicable for U.S.);
- class I and class II medical devices (according to the Federal Law - USA).

This document gives the reprocessing instructions, including care, maintenance, cleaning, disinfection and sterilization of orthopaedic re-usable surgical instruments provided by LimaCorporate. This series of activities are needed to prepare medical devices (new or used) to be used for a surgery.

This document offers a valid aid to users and distributors to perform an effective and safe reprocessing of surgical instruments provided by LimaCorporate.

The hospital staff, the sterilization centres operators and the surgical rooms staff can be involved in manipulation of re-usable surgical instruments provided by LimaCorporate.

The hospital directors and the leaders of the involved departments must know these instructions in order to ensure a safe and effective execution of the reprocessing procedures, and in order to prevent the possible damages or the uncorrected use/manipulation of the surgical instruments.

These instructions are applicable only to re-usable surgical instruments provided by LimaCorporate. Do not apply these instructions to single-use, sterile medical devices.

These “instrument care, cleaning, disinfection and sterilization” instructions and “instruments assembly/disassembly” instructions are available in LimaCorporate website.

2.0 CONSIDERATIONS

This document is referred to all the re-usable orthopaedic surgical instruments provided by LimaCorporate. Read carefully these instructions. This document replaces the previous ones.

In case of countries where the requirements are stricter than those reported in this document, equipments and parameters used for reprocessing procedure must be in accordance to the local laws.

The transportation and distribution of instruments (new or used) by LimaCorporate to the hospitals and to the users does not include the control and maintenance of sterilization or decontamination condition, therefore the users MUST perform again the entire reprocessing procedure.

Re-usable instruments (used or new) provided in non-sterile state must be reprocessed following this document, before use, in order to obtain a sterile state.

All LimaCorporate re-usable surgical instruments can be effectively and safely reprocessed applying manual or manual/automated instruction for cleaning and disinfection, described in this document.

During orthopaedic surgery, instruments are in direct contact with body fluids and tissues (blood, bone, bone marrow). The instruments can also be contaminated because of direct contact with body fluids and tissues infected by aetiological and pathogenic agents. All the users must pay attention to the management of cutting instruments during surgical operations and reprocessing activities.

The surgical technique, the use of instruments and characteristics of devices to be implanted must be well-known by surgeon before performing any surgical procedure.

It is forbidden modify or alter any instrument; furthermore, instruments must not be used for actions different than those indicated in the surgical technique (i.e. do not use the instruments as levers). If the surgical instruments are subjected to excessive loads, they can be damaged.

Any alteration of instruments performed by the user is forbidden. Any consequence resulting from an alteration of instruments by user is considered user's responsibility.

LimaCorporate designed a series of specific instruments sets for each type of implantable component: other instruments must not be used if not specified in the surgical technique. Do not use LimaCorporate instruments with implant devices of other competitors.

The orthopaedic instruments kits must be complete and in good conditions in order to be correctly used.

The instruments are provided assembled in kits, located on the specific trays, included into the boxes suitable for the transportation. A series of instruments with articulated components, removable components and trial components must be used to perform an orthopaedic surgery.

Additional instruments, or instruments provided in substitution of damaged instruments, can be provided by LimaCorporate in single packaging.

The instruments provided in single packaging must be reprocessed like the other instruments.

Re-processing by LimaCorporate for instruments received back by user: LimaCorporate checks the appropriate cleaning and disinfection of received instruments; LimaCorporate performs a cleaning and disinfection process in order to prepare the instruments for being used again. The re-processing performed by LimaCorporate does not include the sterilization phase.

The single-use instruments (provided by LimaCorporate in sterile state in single packaging) must not be re-used or re-processed. Check the labels to verify if the instrument is intended for single use.

Please follow these instructions during all the phases of transportation, cleaning, disinfection, sterilization, use and storage in order to maintain the instruments in good conditions.

LimaCorporate instruments are not intended to be implanted into the human body.

LimaCorporate surgical instruments are not normally used in surgical procedures where they contact high risk Transmissible Spongiform Encephalopathies (TSE) infective tissue, however if prion removal procedures may be required by local law and/or ordinances (e.g. following specific public health concerns), LimaCorporate provides some steps that may be beneficial in prion removal.

3.0 GENERAL WARNINGS, PRECAUTIONS AND RESTRICTIONS

Before to perform a surgery, the surgeon must know the surgical technique, the instruments, and the characteristics of devices.

The instruments must not be modified/altered; the instruments must not be used for applications not in compliance with their intended use.

Personal Protective Equipment (i.e. gown, mask, goggles or face shields, gloves latex free and shoe covers) shall be worn when handling or working with contaminated or potentially contaminated materials, devices and equipment. It is the hospital's responsibility to ensure the staff safety.

Pay attention during handling of instruments sharpened or intended to perform a cut.

Pay attention during handling and transportation of contaminated instruments to the cleaning and disinfection areas. The contaminated instruments shall be transported to the cleaning and disinfection areas using appropriate procedures in order to avoid staff infection and hospital contamination. It is the hospital's responsibility to ensure the staff safety.

The procedure described in this document has been developed with standard instrumentation, in respect of standard re-processing methods.

NOTE: the user must validate the equipment and the processes for cleaning/disinfection and sterilization according to the validated indications reported in this document.

4.0 INSPECTION

A set card showing the list of instruments is provided into the instruments box. The users must check the completeness of instruments according to the list reported in the set card.

Single-packaging instruments are supplied in a closed pouch with an identification label.

Check that instruments are the correct ones to be used for the surgery, then clean, disinfect and sterilize them. For major details about instrument inspection please refer to Section 6.0.

Check the instruments marking.

Most of instruments are intended for re-use, therefore a general wear is a normal consequence of their repeated use. If the instrument properties are not altered by wear (i.e. scratches) it is possible to re-use the instrument. The marking must be readable and the instruments must be uniquely identified. Any graduated scale marked on the instruments surface must be readable and clear.

Otherwise, **if the defects worsen the instruments performance** (i.g. edge loss, metals strain, corrosion), or if the marking is not readable, please contact LimaCorporate for instrument substitution. LimaCorporate recommends to **not re-sharpen instruments autonomously**.

Do not place heavy instruments over delicate devices.

5.0 CLEANING AND DISINFECTION

The process described into this document has been validated.

Alternative cleaning methods are not considered to be validated by the manufacturer.

Alternative cleaning methods must be validated by the users.

The user must ensure the cleaning and disinfection processes execution in compliance with the instructions provided in this document.

5.1 Restrictions

These instructions have been specifically implemented for re-usable instruments manufactured and or distributed by LimaCorporate. They are not applicable to trays of other manufacturers or containing devices that are not manufactured and/or distributed by LimaCorporate.

LimaCorporate instruments are provided in NON-STERILE state.

Cleaning, disinfection and sterilization instructions must be applied only to re-usable instruments (for first use and next re-uses). These instructions are not applicable to SINGLE-USE instruments. The organic residuals must be removed from instruments immediately after their use. **Do not allow drying of contaminant residuals on devices.** All the succeeding cleaning, disinfection and sterilization phases result to be more effective if the user avoids the drying of residuals (i.e. bone, body fluids, tissue debris) on used devices before performing manual or manual/automated cleaning and disinfection.

The user must perform, according to this document, a **manual cleaning** or, in alternative, a combination of **manual and automated cleaning**.

The only automated cleaning is not enough without a previous manual phase (the only automated cleaning is not effective in case of instruments having holes and complex parts).

Disinfection procedures shall be carried out using an automated washer-disinfector validated according to ISO 15883 series.

Metal brushes or scouring pads must not be used during manual cleaning procedures. These equipment could damage the surface or the external finishing of instruments. Soft-bristled, nylon brushes and pipe cleaners shall be used.

Use of hard water shall be avoided. Softened tap water (hardness <150mg/L) shall be used for initial rinsing. Critical water (obtained by ultra-filtering (UF), reverse-osmosis (RO) or deionization (DI) systems (or equivalent methods) according to AAMI TIR34) shall be used for final rinsing to eliminate mineral deposits on instruments. Please refer to the table below for the suggested critical water specification.

Do not use mineral oils or silicon lubricants.

Type of water		Critical water
Specification	Units	
Hardness	mg/L	<1
Conductivity	µS/cm	<10
pH		5-7
Chlorides	mg/L	<1
Bacteria	cfu/mL	<10
Endotoxin	EU/mL	<10

The repeated re-processing of the instruments according to the instructions included in this document has a minimal impact in terms of instruments life. Instruments life is mainly affected by wear and damages caused by surgical use.

LimaCorporate instruments are shipped clean but NOT sterile. Re-clean, disinfect and sterilize the instruments. LimaCorporate advises against the sterilization of instruments with ethylene oxide (EtO), gas plasma and dry-heat. Steam sterilization (autoclave) is suitable for the sterilization of LimaCorporate instruments.

5.2 Preparation at the point of use prior to processing

Instruments must be washed and disinfected as soon as possible after use, in order to minimize hazards of infections (for medical staff) and corrosion (for instruments). Immediately after surgery, excess body fluids and tissue must be removed with a suitable disposable, non-shedding wipe.

Instruments shall be cleaned within 30 minutes after their use to reduce the probability of drying. Submerge the instruments into the water at a temperature <43°C, in order to remove the organic material.

Used instruments must be transported to area suitable for the re-processing in closed or covered containers to prevent the contamination risk.

NOTE: In accordance with the current guidelines, disinfectant can be used in order to ensure the safety of the operators. Submerge the instruments into a disinfectant solution for surgical instruments according to the instructions provided by manufacturer about time, temperature and concentration.

NOTE: the instruments submersion into the enzymatic (having disinfectant activity) solution instead of submersion in the water reduces the risk for the operators and facilitates the cleaning/ disinfection of instruments with complex shapes. The enzymatic solution avoids adhesion of protein and haematic material on the instruments. To prepare the enzymatic solution, see the manufacturer's instructions.

NOTE: pre-soaking is indicated by current guidelines for reducing the risk of transmission of prion borne diseases (e.g. CJD, vCJD, TSE)

After pre-soaking, instruments are to be dried with a suitable non-shedding wipe.

5.3 Cleaning agents preparation

Enzymatic/disinfectant and cleaning agents with a neutral pH or included between 4.5 and 8.5 are recommended for cleaning of LimaCorporate re-usable devices. Do not use strong acids nor oxidizing agents or corrosive chemicals which can alter the instrument surface (i.e. chlorine, iodine, fluorine, organic solvents, ammonia, especially if they are plastic or other polymeric components).

Do **not** soak instruments in **Ringer solution or similar saline solutions**.

Alkaline agents outside the established limits (pH limits) may be used to clean stainless steel and some polymer instruments in countries where required by law or local ordinance; or where prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeldt–Jakob disease (CJD) are a concern. See the instructions provided by the manufacturer in order to verify the adverse effects of the agents on the different materials.

It is critical that alkaline cleaning agents be completely and thoroughly neutralized and rinsed from devices.

Strong basis must be avoided for the cleaning of trays made of aluminium and aluminium alloy.

For optimal performance of the cleaning agents, all cleaning agents shall be prepared at the use-dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare cleaning agents.

The powder cleaning agents must be completely dissolved before use, in order to avoid the instruments corrosion.

The choice of cleaning and disinfecting agents must be in compliance with this LimaCorporate instruction..

5.4 Cleaning and Disinfection instructions

If applicable, the instruments consisting of multiple components must be disassembled, cleaned and reassembled before the sterilization process. For some complex instruments, specific assembly/disassembly indications are available on website www.limacorporate.com

NOTE. fresh cleaning solutions shall be prepared when existing solutions become grossly contaminated (bloody and/or turbid) or when the solution is ineffective (as indicated by the chemical agent's manufacturer).

NOTE after their use, boxes and instruments are to be washed separately; the instrument boxes must be cleaned, disinfected, sterilized, according to the same processes performed for the instruments.

Place the instruments in the boxes only after the cleaning and disinfection phases.

Cleaning of instruments can be performed by **manual** or **manual/automated** process. Disinfection of instruments shall be performed solely with **automated** equipment.

5.5 Manual Cleaning Procedure

Cleaning

Phase M.1	Disassembly all instruments and check all threads are drawn back. It's recommended to keep together the parts of the disassembled instruments to facilitate the assembly. When in doubt follow the instruction for the disassembly/assembly made available by LimaCorporate.
Phase M.2	Completely submerge instruments for at least 15 minutes in a water solution with the proper cleaning agent (i.e.: detergent for surgical instruments*) following the specifications of the manufacturer of the agent. While the instruments are submerged, scrub them with a soft brush (i.e. soft nylon brush), paying special attention to remove all organic residuals and clean also the difficult to reach areas (threads, pivots, hinges). Do not use metal or abrasive brushes! NOTE: Use the syringe or water jet will improve flushing of difficult to reach areas and closely mated surfaces.
Phase M.3	Rinse the instruments in softened running water for at least 3 minutes, paying special attention to plenty flush lumens, cavities and all the difficult to reach areas.
Phase M.4	Place prepared cleaning agents in the sonication unit. Completely submerge devices in cleaning solution and sonicate ($\geq 40\text{kHz}$) for 10 minutes. It is recommended to use temperatures between 40°C and 50°C . Temperatures higher than 50°C (122°F) can cause encrustations when blood is still present.
Phase M.5	Rinse the instruments with running water for at least 5 minutes paying special attention to plenty sprinkle lumens, cavities and all the difficult to reach areas.
Phase M.6	Rinse the instruments with distilled/demineralized water for at least 5 minutes paying special attention to plenty sprinkle lumens, cavities and all the difficult to reach areas.
Phase M.7	Visually inspect the cleaned instruments for soil residuals. If residuals are present repeat the cleaning steps from Phase M.2. Otherwise proceed with disinfection steps.

*: the use of a disinfecting detergent, while not strictly requested, could be useful in this in order to ensure the safety of the operators and to enhance the following disinfection steps.

Disinfection

NOTE: only automated disinfection should be considered as a validated method. No manual method has been validated by LimaCorporate. Disinfection Phase shall be performed with an automated equipment as reported in Section 5.6 - Phase A.4.

Drying

Phase M.8	Immediately dry the instruments in a hot air oven (70°C / 158°F minimum) for at least 10 minutes or with a suitable disposable, non-shedding wipe paying special attention to slots, lumens, mated surfaces, connectors and difficult to dry areas (i.e. laser-markings). Long and thin brushes can be used to dry lumens. By doing so, the risk of corrosion caused by water drops on the instruments surface, can be significantly reduced.
------------------	---

5.6 Manual/Automated Cleaning/Disinfection Procedure

Cleaning

Before proceeding with automated cleaning/disinfection it is recommended to perform a careful manual cleaning according to the following phases:

Phase M.1	Disassembly all instruments and check all threads are drawn back. It's recommended to keep together the parts of the disassembled instruments to facilitate the assembly. When in doubt, follow the instruction for the disassembly/assembly made available by LimaCorporate.
Phase M.2	Completely submerge instruments for at least 10 minutes in a water solution with the proper cleaning agent (i.e.: detergent for surgical instruments) following the specifications of the manufacturer of the solution. While the instruments are submerged, scrub them with a soft brush, paying special attention to remove all organic residuals and clean also the difficult to reach areas (threads, pivots, hinges). Do not use metal or abrasive brushes! NOTE. the use of a ($\geq 24\text{kHz}$) sonicator allows a more accurate cleaning of devices. NOTE. Use the syringe or water jet will improve flushing of difficult to reach areas and closely mated surfaces.
Phase M.3	Rinse the instruments in running water for at least 5 minute, paying special attention to plenty flush lumens, cavities and all the difficult to reach areas.

After this phase, put the instruments in an automated cleaning/disinfection system (validated washer-disinfector in compliance to ISO 15883 series, or to an equivalent standard) and perform a standard cycle in a cleaning/disinfection system. Please follow washer-disinfector manufacturer's instruction about loading configuration, positioning, required accessories etc.. No specific requirements/limitation, in addition to the ones imposed by the manufacturer, are applicable to reprocess LimaCorporate surgical instruments.

NOTE: use of supports could be considered in the cleaning of components with blind lumens/holes. This could help in avoiding water accumulation after washing that would lead to an ineffective drying process.

Automated washing can be included as part of a validated washing, disinfection, and/or drying cycle in accordance to manufacturer's instructions (an example of standard cycle, used for the validation of the procedure, is following reported).

Phase A.1	Pre-cleaning, warm running water(1), softened, at least 2 minutes
Phase A.2	Cleaning with detergent for automatic washer/disinfector(2), hot running water(3), softened at least 360 seconds
Phase A.3	Rinse, hot running water(4) softened, at least 240 seconds

Disinfection and drying phases (see below)

(1) water temperature $\geq 35^{\circ}\text{C}/95^{\circ}\text{F}$

(2) for temperature and exposure times follow the manufacturer instructions

(3) water temperature $\geq 50^{\circ}\text{C}/120^{\circ}\text{F}$

(4) water temperature $\geq 45^{\circ}\text{C}/113^{\circ}\text{F}$

Disinfection

Phase A.4	Thermal disinfection in water ($90^{\circ}\text{C}/194^{\circ}\text{C}$) for at least 5 minutes
------------------	---

Drying

Phase A.5	Drying with hot air ($120^{\circ}\text{C}/248^{\circ}\text{F}$ to $140^{\circ}\text{C}/284^{\circ}\text{F}$) for at least 15 minutes
------------------	--

NB: LimaCorporate advises against the performance of the sole automated cleaning without previous manual cleaning steps as this procedure can be inefficacious for orthopaedic instruments.

NOTE: replace the detergent solution when considerably contaminated (blood-colored and/or turbid) or when the solution is ineffective (as indicated by the manufacturer of the chemical agent used).

After the procedure of cleaning/disinfection, assembly the instruments previously disassembled. For some complex instruments, specific assembly/disassembly indications are available on website www.limacorporate.com.

In case of instruments kits, place all the instruments into their dedicated box (as specified by set card)

6.0 INSPECTION, MAINTENANCE, AND LUBRICATION

All the re-usable instruments: carefully inspect each device to ensure that all visible contamination has been removed. Pay special attention to instruments with blind holes and/or cannulated instruments. If contamination is noted, repeat the cleaning/disinfection process. Verify that residual moisture and water drops have been removed by means of a visual inspection and/or (for instruments with holes/lumens) placing (and tapping) the devices onto a sheet of dry blotting paper, no moisture shall be noted on the sheet.

If residual moisture or water drops are noted, wipe the instruments with a disposable non-shedding wipe or use long and thin brushes (or compressed air) to dry lumens until all residual moisture and water drops are removed.

Visually inspect for damage or wear, including components in their disassembled state prior to re-assembly. Where instruments form part of a larger assembly, ensure a check is completed (e.g. that devices assemble readily with compatible components). Ensure components are re-assembled securely.

The marking must be readable and the instruments must be uniquely identified. Any graduated scale marked on the instruments surface must be readable and clear.

Instruments with moving parts: check the action of moving parts (e.g.: hinges, box-locks, connectors, sliding parts) to ensure smooth operation throughout the intended range of motion.

Sharp instruments and with cutting features: check that sharpened edges are not damaged, distorted and/or with large nicks. Edges should be continuous.

Reamer/drill bits: Inspect “chuck” end for burrs and distortion that might hinder insertion into a drill. Visual and functional inspect for defective coupling end, blunt and dull cutting flutes.

Beaters/Impactors: check that impact surfaces have not cracks or imperfections (burrs and large nicks).

Trial components: check that articular surfaces are without any irregularity (smooth and free of cracks and deep nicks).

Hinged devices: Check for smooth movement of the hinge without excessive “play”.

Locking mechanisms: Check for action.

Mating parts: Check to make sure that mating parts fit together without complications. Ensure components are re-assembled securely.

Metal surfaces: Inspect for corrosion and major deformation.

Instrument with long slender features (particularly rotating instruments): Check that they have not been deformed.

If an instrument malfunction or defect is observed, do not use the instrument and contact LimaCorporate for the replacement.

Hinged, rotating, or articulating instruments should be lubricated with a water soluble product (e.g.: Dr.Weigert neodisher IP spray or equivalent lubricant) **intended for surgical instruments that must be sterilized.** Some water-based instrument lubricants contain bacteriostatic agents which are beneficial, and could be useful in maintaining appropriate bioburden levels. Manufacturer's expiration dates shall be adhered to for both stock and use-dilution concentrations.

Follow the instructions provided by LimaCorporate for instruments disassembly/assembly.

After the procedure of cleaning and inspection it is necessary to perform the sterilization of instruments.

If the cleaned instruments are not immediately sterilized, check that components are perfectly dry and stored in a way suitable for ensure the bioburden limitation. Instruments shall be stored, inside their dedicated trays, in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes. An instruments storage before sterilization, with inadequate modalities, may cause a re-contamination and overcoming of bioburden limits, making the following sterilization ineffective.

7.0 STERILIZATION PROCEDURE

7.1 Packaging

Sterile packaging of single instruments

The single instruments can be packaged into wraps (double wrapping or sterilization boxes) or pouch intended for medical applications suitable for steam sterilization; the dimensions must be correct to ensure a double wrap of instruments.

The packaging must be arranged with the method of double wrapping AAMI or other equivalent method as explained in §7.1.1.

Commercially available, medical grade steam sterilization pouches or wrap shall be used to package individual instruments: check that the internal pouch is large enough to contain the instrument (without stretching the seals or damaging the packaging) and small enough to be put in a second pouch.

Sterile packaging of instruments kits

After inspection the instruments must be placed in the proper slots of their dedicated trays.

The instruments must be placed into their dedicated trays according to the indications of instrument set card, in order to avoid the contact each other and ensure the steam penetration to all instruments surfaces during the sterilization.

Put in the dedicated trays only their dedicated instruments according to the instruments set card.

Instruments shall not be stacked or placed in close contact.

For sterilization process, trays and cases with lids must be wrapped according to the following alternatives according ISO 11607-1:

- standard wraps for medical applications, applying the method of double wrapping AAMI or equivalent method as explained in §7.1.1;
- approved sterilization container with gasket lid for sterilization.

Follow the sterilization containers' manufacturer's instructions for inserting and replacing sterilization filters within sterilization containers.

The user must ensure that the instrument case is not tipped or the contents shifted once the devices are arranged in the case.

Safety instructions on weight:

The total weight of a wrapped instrument tray or box must not exceed 11.4Kg (according to AAMI ST77).

NOTE. slots designated for specific devices shall contain only LimaCorporate devices specifically intended for these slots.

NOTE. Use only cleared sterilization wrap or other cleared accessory (for U.S. use only, LimaCorporate recommends the use of FDA-cleared wrap or other FDA-cleared accessory (e.g., container, pouch)) that has been validated to allow sterilant penetration as well as and sterility maintenance.


7.1.1 Method of double wrapping ANSI/AAMI ST79:2017

The method used by LimaCorporate for the validation of sterilization with Sterisheet double wrapping following ANSI/AAMI ST79:2017 is the "envelope fold", explained in the following steps. Sequential wrapping refers to wrapping two layers of wrap material individually. Sequential wrapping using the envelope-fold technique shown in figures below.

The procedure is as follows:

- 1) Place the wrap on the table to form a diamond shape. Place the device(s) to be wrapped in the center of the wrap, parallel with the edge of the table.
- 2) Bring the lower corner up to cover the contents and fold the tip back on itself to form a tab or flap (which is used later to assist in opening the pack aseptically).
- 3) Fold the left corner over the contents and fold the tip back to form a tab.
- 4) Fold the right corner over the left fold and fold the tip back on itself to form a tab.
- 5) Bring the top corner down over the contents and tuck the corner under the right and left folds, leaving a small tab visible for easy opening.
- 6) Apply the second wrap by placing the single wrapped device in the center of the remaining wrap and repeating the wrap sequence to form a package within a package. Bring the lower corner up to cover the single wrapped item and fold the tip back on itself to form a tab or flap.
- 7) Fold the left corner over the single wrapped item and fold the tip back to form a tab.
- 8) Fold the right corner over the left fold and fold the tip back on itself to form a tab.
- 9) Bring the top corner down over the single wrapped item and tuck the corner under the right and left folds, leaving a small tab visible for easy opening. Secure with chemical indicator tape.

Other methods equivalent to this here described can be also used.

7.2 Sterilization procedure

The instruments provided by LimaCorporate are to be sterilized by means of steam sterilization (autoclave) validated according to the requirements of EN 285, European Pharmacopoeia or ISO 17665-1/-2/-3. In case the user chooses another sterilization method the individual and/or hospital location takes on the responsibility for the sterilization efficacy and possible damages of LimaCorporate instruments.

Make sure your sterilizing facility is properly serviced and the process is validated.

NOTE. Each sterilization equipment has its own operating parameters. The suitability of the parameters has to be validated by qualified and skilled personnel on sterilization procedures. Each change to the sterilization procedures is at the risk of the user.

Instrument sets shall be properly prepared and packaged in trays and/or cases that will allow steam to penetrate and make direct contact with all surfaces.

When sterilizing multiple instrument sets in one sterilization cycle, ensure that the manufacturer's maximum load is not exceeded.

Disinfection is only acceptable as a precursor to full sterilization for re-usable surgical instruments.

The following minimum sterilization parameters were validated to provide a 10^{-6} sterility assurance level (SAL):

Validated cycle (applicable worldwide)	
Cycle	Cycle 1
Cycle Type	Prevacuum / Pulsating Vacuum
Temperature	132° C / 270° F
Exposure Time	4 minutes
Drying Time	30 minutes

Validated cycles (not applicable for U.S.)		
Cycle	Cycle 2 (*)	Cycle 3 (*) (**)
Cycle Type	Prevacuum / Pulsating Vacuum	Prevacuum / Pulsating Vacuum
Temperature	134° C / 273° F	134° C / 273° F
Exposure Time	3 minutes	18 minutes
Drying Time	20 minutes	30 minutes

* This cycle is not intended for U.S.

** This cycle is to be used in scenarios where Transmissible Spongiform Encephalopathy (TSE) and Creutzfeldt–Jakob disease (CJD and vCJD) are a concern.

NOTE. The instructions of the sterilization equipment use, configuration and maximum load issued by the manufacturer must be strictly followed.

Once sterilized, the instruments must be handled with aseptic procedures.

8.0 STORAGE

Instruments are to be kept and carried in their boxes: this will provide suitable protection from impacts and damages, and at the same time will shield users from risks of cuts.

Instruments shall be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.

The sterile packaging of instruments must be carefully checked before opening to ensure that packaging integrity has not been compromised.

NOTE. Store instruments in their boxes. Protect components from involuntary contact with other instrument as it may lead to damaging their finishing. Before use check possible visible damages of each instrument.

NOTE. If damages are observed (lacerations, perforation, ...) on sterile packaging or if it has been exposed to moisture, the instrument set must be packaged and sterilized again.

NOTE. If the lid's gasket or the sterilization container's filters have been opened or damaged, it is needed to replace the sterile filter and re-sterilize the instrument set.

9.0 HOSPITAL RESPONSIBILITIES FOR LIMACORPORATE LOANER INSTRUMENTS

Instruments which no longer perform properly because of long use, mishandling, or improper care shall be returned to LimaCorporate who will provide of their maintenance or disposal.

Loaner sets shall undergo all steps of decontamination, cleaning, disinfection, inspection, and terminal sterilization before being returned to LimaCorporate.

Documentation of decontamination shall be provided with instruments being returned to LimaCorporate.

Instruments are to be returned in their boxes and the boxes are to be properly packaged inside a cardboard box suitable for preventing possible damages during transport.

10.0 FURTHER INFORMATION

The instructions supplied in this document are validated (according to ISO 17664 and FDA Guidance “Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling”) and they are intended to give the instructions for the preparation for use of re-usable orthopaedic instruments.

The users have the responsibility to ensure, according to these instructions for instrument care, cleaning, disinfection and sterilization, the correct re-processing of devices performed with the correct equipments and materials; the involved staff must be adequately trained.
The equipment must be validated and periodically monitored.

All the deviations from these standard instructions performed by the operators must be adequately evaluated, considering the possible consequences.

The agents used for the validation of cleaning, disinfection and sterilization processes are following reported:

Chemical Agent	Phase of manual cleaning/disinfection procedure	Phase of manual/automated cleaning/disinfection procedure
Prolystica 2x	Cleaning	Manual cleaning
Prolystica 2x	//	Manual/Automated cleaning
Dr. Weigert neodisher IP spray	Lubrication	Lubrication

The equipments used for the validation of cleaning and sterilization processes are following reported:

Equipment	Model
Automated washer/disinfector	DS 1000 Steelco
Sterilizing machine	STERIS AMSCO 600 SERIES FAHC7713-D

For further investigation and questions contact LimaCorporate S.p.A.
(<http://www.limacorporate.com/>).

11.0 SYMBOLOGY

Catalogue number

Medical Device¹

Batch code

Manufacturer

Date of manufacture

Non sterile

Unique device identifier - UDI

Contains hazardous substances

Distributed by

Imported in the European Community/European Union by

Model trade name

¹ Please note that in the ISO context, the 'MD' symbol is used to identify that the product in question is a medical device [ISO_15223-1]

12.0 BIBLIOGRAPHY

1	UNI EN 285:2006
2	ISO 17664:2017
3	ISO 17665-1: 2007
4	ISO 17665-2: 2009
5	ISO 17665-3: 2013
6	AAMI TIR 30:2011
7	AAMI TIR 12:2010
8	ANSI/AAMI ST77:2013
9	ANSI/AAMI ST79:2010/A4:2013
12	ASTM E 1766-15
13	Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff; 2017.
14	Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
15	European Pharmacopeia
16	AS/NZS 4187:2014
17	AAMI TIR 34:2014
18	ASTM E1837:96/(R)2014
19	EN 14885:2019
20	EN ISO 15883-4:2018
21	ISO TS 15883-5:2005

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1	UNI EN 285:2006
2	ISO 17664:2017
3	ISO 17665-1: 2007
4	ISO 17665-2: 2009
5	ISO 17665-3: 2013
6	AAMI TIR 30:2011
7	AAMI TIR 12:2010
8	ANSI/AAMI ST77:2013
9	ANSI/AAMI ST79:2010/A4:2013
12	ASTM E 1766-15
13	Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff; 2017.
14	Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
15	Европейска фармакопея
16	AS/NZS 4187:2014
17	AAMI TIR 34:2014
18	ASTM E1837:96/(R)2014
19	EN 14885:2019
20	EN ISO 15883-4:2018
21	ISO TS 15883-5:2005

BRIGA O INSTRUMENTU I NJEGOVO ČIŠĆENJE, DEZINFEKCIJA I STERILIZACIJA (u skladu s ISO 17664) – Rev. 02/2023

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12.0 BIBLIOGRAFIЈА

1	Uni EN 285:2006
2	ISO 17664:2017
3	ISO 17665-1: 2007
4	ISO 17665-2: 2009
5	ISO 17665-3: 2013
6	AAMI TIR 30:2011
7	AAMI TIR 12:2010
8	ANSI/AAMI ST77:2013
9	ANSI/AAMI ST79:2010/A4:2013
12	ASTM E 1766-15
13	Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff; 2017.
14	Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
15	European Pharmacopeia
16	AS/NZS 4187:2014
17	AAMI TIR 34:2014
18	ASTM E1837:96/(R)2014
19	EN 14885:2019
20	EN ISO 15883-4:2018
21	ISO TS 15883-5:2005

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12.0 BIBLIOGRAFIE

1	UNI EN 285:2006
2	ISO 17664:2017
3	ISO 17665-1: 2007
4	ISO 17665-2: 2009
5	ISO 17665-3: 2013
6	AAMI TIR 30:2011
7	AAMI TIR 12:2010
8	ANSI/AAMI ST77:2013
9	ANSI/AAMI ST79:2010/A4:2013
12	ASTM E 1766-15
13	Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry ss Food and Drug Administration Staff; 2017.
14	Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
15	European Pharmacopeia
16	AS/NZS 4187:2014
17	AAMI TIR 34:2014
18	ASTM E1837:96/(R)2014
19	EN 14885:2019
20	EN ISO 15883-4:2018
21	ISO TS 15883-5:2005

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12.0 KAYNAKÇA

1	UNI EN 285:2006
2	ISO 17664:2017
3	ISO 17665-1: 2007
4	ISO 17665-2: 2009
5	ISO 17665-3: 2013
6	AAMI TIR 30:2011
7	AAMI TIR 12:2010
8	ANSI/AAMI ST77:2013
9	ANSI/AAMI ST79:2010/A4:2013
12	ASTM E 1766-15
13	Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff; 2017.
14	Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
15	Avrupa Farmakopesi
16	AS/NZS 4187:2014
17	AAMI TIR 34:2014
18	ASTM E1837:96/(R)2014
19	EN 14885:2019
20	EN ISO 15883-4:2018
21	ISO TS 15883-5:2005

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Lima (Beijing) Medical Devices Co., Ltd.
Lima Orthopaedics Canada Inc.
Limacorporate S.p.A.