

ISOFORM[®] - TLSO PATIENT RE-APPLICATION INSTRUCTIONS:





Step 1. Re-Apply Brace

While standing, slide the TLSO on (like you might put on a jacket or backpack)-keeping head between the two shoulder straps. **(Figure A)**


Lower the TLSO into position. Refasten the quick release buckle at top. **(Figure B)**

Let the front/anterior panel and back/posterior panel hang over the torso. **(Figure C)**



Step 2. Close Panels

Close both of the compression side panels so that they are even and comfortably snug on the Anterior Panel. **(Figure D)**



Step 3. Activate Compression

Grasp pull-handle and pull straight out from the torso. After the desired level of compression has been achieved, fasten the pull handle to anterior portion of the compression panels. Adjust compression after five minutes, then re-adjusting as needed. **(Figure E & F)**

BRACE REMOVAL:

Step 1.

Release the compression handles **(Figure E)** then release one of the shoulder strap buckles. **(Figure B)**

Step 2.

Slide the orthosis off the other shoulder like removing a jacket or backpack. **(Figure A)**

HOT/COLD PACK (sold separately):

Follow instructions on hot/cold pack. Insert hot/cold pack into pocket and adhere to back panel.



CLEANING INSTRUCTIONS:

Separate the brace's removable REAR PANEL before cleaning. Hand wash in cold water, less than 86°F (30°C), using mild soap. Air dry. NOTE: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material. Do NOT iron / Do NOT tumble dry / Do NOT bleach.

INTENDED USER PROFILE:

The user should be able to: Read and understand the directions, warnings and cautions.

INTENDED USE/INDICATIONS:

The Isoform[®] TLSO brace is designed to support and stabilize the spine for acute pain relief through postural extension alignment and may be suitable for people with compression fractures of the thoracic spine, osteoporosis, thoracic mechanical pain and kyphosis.

CONTRAINDICATIONS:

Pregnancy and circulation, pulmonary, cardiovascular or skeletal conditions which have risk to be made worse as a result of compression and/or pressure

WARNINGS AND CAUTIONS:

If you experience any pain, swelling, sensation changes, excessive discoloration, irritation to the skin or any unusual reaction while using this product, consult your medical professional. If any distortion of the support panels, fraying or deterioration of lace cables is noticed, discontinue use. Brace is intended to be worn over clothing. Brace should be appropriately sized to the patient to ensure proper support and function. If the brace does not fit properly, it may be adjusted by a healthcare professional. Care should be taken with patients, especially the elderly, not to apply the brace too tightly, which MAY cause excessive pressure. Patients should be instructed on appropriate tensioning during ongoing use. Patient should see a healthcare professional for a replacement if a product quality issue occurs within the warranty timeframe. Brace should not be exposed to moisturizers, chemicals, or solvents, as they may affect durability.

MATERIAL CONTENTS:

Nylon, Polyester, Polyethylene, Polyurethane, Polypropylene, Acrylic, PVC, EVA, Acetal, Aluminum, Steel, Synthetic Rubber, Glass, Calcium Carbonate

NOT MADE WITH NATURAL RUBBER LATEX

WARRANTY:

DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

RX ONLY. INTENDED FOR SINGLE PATIENT USE.

Product Support:

Phone: 888.405.3251 / Email: product.specialists@djoglobal.com

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

