

# Economy Hinged Knee

## PATIENT APPLICATION



**BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

### INDICATIONS:

- Mild MCL and/or LCL sprains
- Mild arthritis
- Post arthroscopic meniscectomy
- Mild global instabilities.

### WARNINGS:

Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene by-products.

**WARNING:** If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

### CAUTION :

Please be advised that DJO, LLC, does not recommend simultaneous wearing of functional knee braces bilaterally (right and left legs), during activities (ie skiing) that may cause the braces to contact each other, resulting in potential loss of control and injury. Functional knee braces are not recommended to be worn over clothing. Consult your licensed healthcare professional for your activity limitations.

### WARRANTY:


DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

**NOT MADE WITH NATURAL RUBBER LATEX  
SINGLE PATIENT USE ONLY**

NOTICE : WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

## Application Information:

### WRAP AROUND





**A** Unfasten and loosen the straps at the thigh and calf.

**B** Grasp the edges of the sleeve and pull onto leg. If using wrap-around version, simply wrap around leg and fasten over the front of the leg. Position the sleeve so that the hinges are parallel to each other, and the top of the circular cut-out aligns with the top of the knee cap.

**C** Fasten the support straps at the thigh and calf.

### HINGED KNEE



**E REMOVAL:** To remove, unfasten straps and apply a downward pressure to the top of the sleeve at the hinges.

### USE AND CARE:

- Hand wash in cold water using mild soap, rinse thoroughly. AIR DRY. Hinges may be removed for washing. Note: If not rinsed thoroughly, residual soap may cause skin irritation and deterioration of the material.
- The brace may be used in fresh or salt water. Following use in water activities, drain thoroughly and rinse with clean, fresh water. AIR DRY.



DJO, LLC  
T 800.336.6569 D 760.727.1280 F 800.936.6569  
1430 Decision Street | Vista, CA 92081 | U.S.A.  
[DJOglobal.com](http://DJOglobal.com)

COPYRIGHT © DJO, LLC 2020

MKTIS00-10617 Rev A

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.