

DONJOY®

OA ADJUSTER™ 3
Mighty. Tough.



 **DJO**
GLOBAL

Mighty. Tough.

The improved OA Adjuster™ 3 is powered by nano MAG™ and reinforced with a tough aluminum frame making it the most durable and functional DonJoy OA knee brace.



New Contoured Design conforms specifically to the calf, allowing minimal migration and maximum comfort for even the toughest patients



Clinically-Proven*
DonJoy's Adjuster™ technology offloads a little—or a lot

nano MAG™ Technology on the offloading side provides lightweight function with the strength of aircraft-grade aluminum (19 oz)

Floating Hinges allow a contoured, personalized fit for all levels of activity, with double upright design

4 Points of Leverage™ system pioneered by DonJoy provides clinically proven ligament stability**

- 1 anterior thigh anchor
- 2 posterior tibia anchor
- 3 lower thigh strap pushes femur anteriorly
- 4 anterior tibia strap applies an active constant load to prevent anterior tibial translation

*Matthew Nadaud, Richard Komistek et al, 2006 (J Bone Joint Surg AM, 87:114-119, 2005)
**Fleming, Beynnon, et al 2000 (Am J Sports Med, 2000;28:815-824.)



90 Day Test Drive - Simply fit your patient in any Donjoy knee osteoarthritis brace. If they are not 100% satisfied with the brace you may return it for a full credit or refund within 90 days of brace fitting.



EASE OF USE - The Arthritis Foundation has honored every Donjoy knee osteoarthritis braces with its coveted "Ease-of-Use" Commendation after an independent research lab conducted rigorous testing to assure that they met the needs of people with functional limitations due to the effects of arthritis.

ORDERING INFORMATION PDAC ASSIGNED L1845

Part Number	Affected Compartment	Leg
11-1590-X	Medial	Right
11-1591-X	Medial	Left
11-1592-X	Lateral	Right
11-1593-X	Lateral	Left

Available in Graphite.

Lock and off-Load Technology™ | Floating Hinge | Swiveling D-Ring

PATIENT-READY

Brace Size	Thigh Circumference 6" Above Mid-Patella
x = 1 = XS	13 - 15.5" (33 - 39 cm)
x = 2 = S	15.5 - 18.5" (39 - 47 cm)
x = 3 = M	18.5 - 21" (47 - 53 cm)
x = 4 = L	21 - 23.5" (53 - 60 cm)
x = 5 = XL	23.5 - 26.5" (60 - 67 cm)
x = 6 = XXL	26.5 - 29.5" (67 - 75 cm)
x = 7 = XXXL	29.5 - 32" (75 - 81 cm)

ACCESSORIES

Part Number	Description	Size
25-1376-0-00000	Flexion/Extension Stop Kit	Universal
11-0075-X-00000	Cotton Lycra Suspension Sleeve	S - XXL

Call **800.336.6569** or visit www.djoglobal.com/donjoy



DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

DJGlobal.com/donjoy

Together in Motion™

DJO Global proudly sponsors:

