

DONJOY®

OA ADJUSTER™ 3  
**Mighty. Tough.**



 **DJO**  
GLOBAL

# Mighty. Tough.

The improved OA Adjuster™ 3 is powered by nano MAG™ and reinforced with a tough aluminum frame making it the most durable and functional DonJoy OA knee brace.



**New Contoured Design** conforms specifically to the calf, allowing minimal migration and maximum comfort for even the toughest patients



1

**Clinically-Proven\***  
DonJoy's Adjuster™ technology offloads a little—or a lot

**nano MAG™ Technology** on the offloading side provides lightweight function with the strength of aircraft-grade aluminum (19 oz)

3

**Floating Hinges** allow a contoured, personalized fit for all levels of activity, with double upright design

4

2

**4 Points of Leverage™** system pioneered by DonJoy provides clinically proven ligament stability\*\*

- 1 anterior thigh anchor
- 2 posterior tibia anchor
- 3 lower thigh strap pushes femur anteriorly
- 4 anterior tibia strap applies an active constant load to prevent anterior tibial translation

\*Matthew Nadaud, Richard Komistek et al, 2006 (J Bone Joint Surg AM, 87:114-119, 2005)  
\*\*Fleming, Beynnon, et al 2000 (Am J Sports Med, 2000;28:815-824.)



**90 Day Test Drive** - Simply fit your patient in any Donjoy knee osteoarthritis brace. If they are not 100% satisfied with the brace you may return it for a full credit or refund within 90 days of brace fitting.



**EASE OF USE** - The Arthritis Foundation has honored every Donjoy knee osteoarthritis braces with its coveted "Ease-of-Use" Commendation after an independent research lab conducted rigorous testing to assure that they met the needs of people with functional limitations due to the effects of arthritis.

## ORDERING INFORMATION PDAC ASSIGNED L1845

Part Number	Affected Compartment	Leg
11-1590-X	Medial	Right
11-1591-X	Medial	Left
11-1592-X	Lateral	Right
11-1593-X	Lateral	Left

Available in Graphite.

Lock and off-Load Technology™ | Floating Hinge | Swiveling D-Ring

## PATIENT-READY

Brace Size	Thigh Circumference 6" Above Mid-Patella
x = 1 = XS	13 - 15.5" (33 - 39 cm)
x = 2 = S	15.5 - 18.5" (39 - 47 cm)
x = 3 = M	18.5 - 21" (47 - 53 cm)
x = 4 = L	21 - 23.5" (53 - 60 cm)
x = 5 = XL	23.5 - 26.5" (60 - 67 cm)
x = 6 = XXL	26.5 - 29.5" (67 - 75 cm)
x = 7 = XXXL	29.5 - 32" (75 - 81 cm)

## ACCESSORIES

Part Number	Description	Size
25-1376-0-00000	Flexion/Extension Stop Kit	Universal
11-0075-X-00000	Cotton Lycra Suspension Sleeve	S - XXL

Call **800.336.6569** or visit [www.djoglobal.com/donjoy](http://www.djoglobal.com/donjoy)



**DJO, LLC | A DJO Global Company**

**T** 800.336.6569 **D** 760.727.1280 **F** 800.936.6569  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
[DJGlobal.com/donjoy](http://DJGlobal.com/donjoy)

*Together in Motion™*

DJO Global proudly sponsors:

