

ARMOR FOURCEPOINT

Indestructible Protection. Really?



4-POINTS-OF-LEVERAGE™

Indestructible means:

REDUCE ACL STRAIN BY 50% 4-POINTS-OF-LEVERAGE™

A unique cuff and strapping configuration that provides a posterior force to the tibia preventing anterior movement.



Indestructible means:

GUARANTEED PROTECTION

If a patient re-ruptures a reconstructed ACL while wearing a properly fitted DonJoy Defiance,™ Armor FourcePoint™ or FULLFORCE™ within one year of purchase, DJO Global will pay a surgical insurance deductible. Visit donjoy.com for details and eligibility.

FEATURING OUR LIGHTEST AND TOUGHEST HINGE

Indestructible means:

MINIMIZE RISK


FOURCEPOINT HINGE™

Knee flexion angles less than 30 degrees are considered "at risk" and contribute to ACL injuries. The FourcePoint™ hinge technology increases the knee flexion angles therefore reducing the time spent in the ACL "at risk" positions². The FourcePoint™ hinge enhances DonJoy's 4-Points-of-Leverage System™ by dampening knee extension, and reducing shear forces at the knee².

Indestructible means:

UNCOMPROMISED STRENGTH

Steel reinforced hinge plate and enhanced aluminum frame deliver a shield of proven protection.



*Independent test data from the Orthopedic Biomechanics Laboratory at Children's Hospital of San Diego.

Indestructible means:

PROVEN RESULTS

There are numerous publications to support the 4-Points-of-Leverage,™ FourcePoint™ hinge, and the uncompromised strength of Armor.







DONJOY - A BRACE FOR EVERY CONDITION

		Brace	Activity Level		
		Туре	High	Moderate	Mild
Brace Strength Brace	Defiance™	Custom			
	Armor FourcePoint™	Prefabricated			
	FULLFORCE™	Prefabricated			

CLINICAL EVIDENCE

STUDIES THAT SHOW THE 4-POINTS-OF-LEVERAGE™ CONFIGURATION SIGNIFICANTLY REDUCES ACL STRAIN:

Fleming BC, Renstrom PA, Beynnon BD, Rengstrom PA, Peura GD:
2000 The Influence of Functional Knee Bracing on the Anterior Cruciate Ligament Strain
Biomechanics in Weight Bearing and Non-Weight Bearing Knees,
The American Journal of Sports Medicine: Vol. 28 (6): 815-824

Beynnon BD, Fleming BC, Peura GD, Johnson RJ, Rengstrom PA, Nichols CE, Pope MH: An In Vivo Investigation of Anterior Cruciate Ligament Strain:

1995 The Effect of Functional Knee Bracing and Attachment Strap Tension, 41st Annual Meeting, Orthopaedic Research Society

Beynnon BD, Pope MH, Wertheimer CM, Johnson RJ, Fleming BC, Nichols CE, Howe JC:

1992 The Effect of Functional Knee Braces on Strain on the Anterior Cruciate Ligament In Vivo,
The Journal of Bone and Joint Surgery (Am.), Vol. 74-A (9): 1298-1312

² STUDIES THAT SHOW THE EFFECTIVENESS OF THE FOURCEPOINT™ HINGE IN REDUCING KNEE SHEAR FORCES:

Lin CF, Liu H, Garrett WE, Yu B:

2008 Effects of a Knee Extension Constraint Brace on Selected Lower Extremity Motion Patterns During a Stop-Jump Task, Journal of Applied Biomechanics, Vol. 24: 158-165

Yu B, Herman D, Preston J, Lu W, Kirkendall DT, Garrett WE:

2004 Immediate Effects of a Knee Brace with a Constraint to Knee Extension on Knee Kinematics and Ground Reaction Forces in a Stop-Jump Task,
The American Journal of Sports Medicine, Vol. 32 (5): 1136-1143

ORDERING INFORMATION:

Part Number	Description
11-144O-X	ACL, Short Calf Length, Right (15")
11-1441-X	ACL, Short Calf Length, Left (15")
11-1442-X	ACL, Standard Calf Length, Right (16")
11-1443-X	ACL, Standard Calf Length, Left (16")

Sizes Available: XS-XXXL Measurements taken 15cm above mid-patella

 $X = 1 = XS = 13^{\circ}-15 1/2^{\circ} (33-39 cm)$

x = 2 = S = 15 1/2"-18 1/2" (39-47 cm)

 $X = 3 = M = 18 \frac{1}{2} - 21'' (47-53 \text{ cm})$

x = 4 = L = 21"-23 1/2" (53-60 cm)

 $X = 5 = XL = 23 \frac{1}{2}$ - 26 $\frac{1}{2}$ (60-67 cm)

x = 6 = XXL = 26 1/2"-29 1/2" (67-75 cm)

x = 7 = XXXL = 29 1/2"-32" (75-82.5 cm)

ACCESSORIES

Part Number	Description	Size
11-0075-X-00000	Lycra® Suspension Sleeve	S-XXL
11-0122-X	Neoprene Suspension Sleeve	XS-XXL
11-1006-x-06000	Impact Guard	S-XXL
11-1525-X	Tru-Pull® Advance Attachment	XS-XXXL
11-0398-9-00000	Full Extension Stop	Universal





DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569 1430 Decision Street | Vista, CA 92081-8553 | U.S.A. DJOglobal.com

Together in Motion.

DJO Global proudly sponsors:









