


## CONCELLTRATE® 100

**ConCelltrate® 100** is a human bone derived, flowable inductive carrier matrix, with no additional carriers, indicated for homologous use in the treatment of surgically created or traumatic skeletal defects.

- May be hydrated with saline, blood, BMA, PRP, or other cellular components in accordance with a physician's well-informed medical judgement
- Extreme irrigation resistant properties intended to aid in maintaining final allograft placement/shape while sustaining its desired hydrating constituent
- In-vivo test results demonstrate all 5 bone forming elements present (chondrocytes, osteocytes, cartilage, bone marrow cells, and new bone)
- Every lot is verified for osteoinductivity post-sterilization
- Provided sterile via low dose, low temperature, gamma sterilization with a sterility assurance level (SAL) of 10<sup>-6</sup>
- Packaged with a shelf life of 5 years and may be stored in ambient temperatures
- Available in 1.0cc, 2.5cc, 5.0cc, and 10.0cc jars



Part Number	Volume
453010-TS	1.0cc jar
453025-TS	2.5cc jar
453050-TS	5.0cc jar
453100-TS	10.0cc jar