

AIRCAST®

AIRLIFT™ PTTD BRACE
Off the Shelf with a Customizable Fit





AIRLIFT™ PTTD BRACE

Off the Shelf with a Customizable Fit

The Aircast® AirLift™ PTTD Brace has a versatile stirrup design that aids in the treatment of posterior tibial tendon dysfunction (PTTD). The anatomically correct design and new intuitive Integrated Inflation System allow for customization of the air cell to help lift the foot arch for a more natural position. Along with Aircast® technology, durable cross strapping and easy step-in design, the AirLift provides an off-the-shelf customized fit for comfortable relief from PTTD symptoms.



Integrated Inflation System for Customizable Fit

- Allows patient to customize the air cell to help hold the arch in natural position



Medial Cross-strap

- Helps lift foot arch for additional support and stability



Stirrup Design

- Bilateral anatomically shaped shells help with medial/lateral support



Easy Step-in Design

- Simple, easy application with posterior step-in design

PART NUMBER		SHOE SIZE		
Left	Right	Size	Men	Women
02PSL-1	02PSR-1	Small	Up to 7	Up to 8.5
02PML-1	02PMR-1	Medium	7.5 – 11	9 – 12.5
02PLL-1	02PLR-1	Large	11.5 and up	13 and up

Recommended for those with:
Posterior tibial tendon dysfunction (PTTD) or adult acquired flatfoot (AAF)



T 800.336.6569 F 800.457.4221
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
djoglobal.com/aircast

Copyright © 2021 by DJO, LLC
MKT00-11411 Rev A

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.