

Preferred Products
CATALOG





As a DJO Global brand, DonJoy pioneered the concept of functional knee bracing over 30 years ago. We've led the profession of performance ever since by studying the body, listening to athletes, consulting physicians and pushing the envelope of innovation. Today, DonJoy braces are standard equipment for competitors in over 44 countries. They count on us to create world-class products, while medical professionals look to us for unmatched service and in-depth education. The result – complete confidence for those who wear DonJoy braces, and those who prescribe them.

- Knee Ligament Bracing..... 1
 - Custom 4
 - Prefabricated..... 7
- Osteoarthritis Bracing 11
 - Custom..... 12
 - Prefabricated..... 13
- Post-Op Bracing..... 15
- Patellofemoral Bracing..... 17
- Walking Braces..... 23
- Ankle Braces and Supports..... 27
- Lower Back Bracing..... 31
- Upper Extremity..... 33
- Cold Therapy 39
- Index
 - Numerical..... 45
 - HCPCS..... 47

In addition to fast and accurate delivery, you can expect a friendly and knowledgeable Customer Care Representative to take your order and answer any questions. We look forward to serving you between the hours of 4:30 a.m.- 6 p.m. PST.

How To Order

To order, contact Customer Service by:

- Phone:** 1-800-336-6569
(1-800-3-DONJOY)
- Email:** customer.care@djoglobal.com
- Fax** (24 hours a day): 1-800-936-6569
(1-800-9-DONJOY)
- Product Specialists:** 1-888-405-3251
Contact the DJO Product Specialists for assistance with brace recommendations, brace fitting and application, troubleshooting, and for product instruction.

Questions We Will Ask

Your Account Number:
If you do not know your account number, the representative will be happy to look it up for you.

Product:
Please state the product number, size and quantity (description if available).

Other Important Information

- Administration:** 1-800-321-9549
- Statement or Account Inquiries:** 1-800-236-3715
(7 a.m.- 5 p.m. PST)
- Customer Care:** 1-800-336-6569

Shipping

Standard Product Orders* received before 1 p.m. Pacific Standard Time will be shipped within 24 hours and can be shipped for overnight delivery. For order tracking, go to: www.djoglobal.com/tracking.

*Standard Product Orders include Rigid Bracing and other off the shelf items.

Special Requirements

The majority of orders will ship within 24 hours. Custom products and special orders (logos, replacement parts, special measurements) will ship within 3 - 10 business days. Please contact Customer Care for product details.

Warranty


Products have a limited warranty from the date of purchase against defects in materials and workmanship. Warranty varies with some products. Please contact Customer Care with specific questions.

Return Policy

Before returning or exchanging any product, please call a Customer Care Representative for a Return Authorization Number. They will assist you with product specific return policies. Non-warranty returns will be accepted within 180 days from date of purchase and will be subject to a 15% restocking fee.

Sizing

1=XS, 2=S, 3=M, 4=L, 5=XL, 6=XXL, 7=XXXL, 9=Universal



Off-The-Shelf Sizing

Measurements taken 6" above (thigh) or 6" below (calf) mid-patella.

Size	Thigh	Knee Center	Calf
XS	13"-15½" (33-39 cm)	12"-13" (30.5-33 cm)	10"-12" (25.5-30.5 cm)
S	15½"-18½" (39-47 cm)	13"-14" (33-35.5 cm)	12"-14" (30.5-35.5 cm)
M	18½"-21" (47-53.25 cm)	14"-15" (35.5-38 cm)	14"-16" (35.5-40.5 cm)
L	21"-23½" (53.25-59.5 cm)	15"-17" (38-43 cm)	16"-18" (40.5-47 cm)
XL	23½"-26½" (59.5-67.25 cm)	17"-19" (43-48.25 cm)	18"-20" (47-50.75 cm)
XXL	26½"-29½" (67.25-75 cm)	19"-21" (48.25-53.25 cm)	20"-22" (50.75-56 cm)
XXXL	29½"-32" (75-81.25 cm)	21"-23" (56-61 cm)	22"-24" (53.25-58.5 cm)

Knee Ligament Bracing



Clinically Proven Brace and Hinge Technology

4-POINTS-OF-LEVERAGE SYSTEM™

DonJoy pioneered **dynamic bracing** with the 4-Points-of-Leverage System, a unique cuff and strapping configuration that provides a net differential posterior force to the tibia, preventing anterior movement. The result is significantly reduced strain on the ACL.¹ Passive brace designs don't actively address instabilities.

POINT OF LEVERAGE 1

The femur is anchored by a cuff placed at the top of the anterior thigh.

Dynamic brace design protects at **all times**

POINT OF LEVERAGE 4

A strap placed on the anterior tibia applies an active constant load to prevent anterior tibial translation.



POINT OF LEVERAGE 3

A strap across the back of the lower thigh pushes the femur anteriorly.

POINT OF LEVERAGE 2

The tibia is anchored by a cuff placed at the bottom of the posterior calf.

FOURCEPOINT™ HINGE TECHNOLOGY

DonJoy has taken hinge design and ACL protection to another level with its patented FourcePoint technology; the only hinge design clinically proven to protect the ACL.²

FourcePoint Technology:

- Enhances DonJoy's 4-Points-of-Leverage System by dampening knee-joint extension and improving overall mechanical performance.
- Increases flexion angles (5°) while reducing shear forces (9%-13%)³ at the knee.
- Leaf spring mechanism applies gradual resistance, reducing time spent in the "at-risk" position (flexion angles < 60°).⁴
- Provides prophylactic benefits. Athletes who trained wearing a DonJoy 4-Points-of-Leverage brace equipped with a FourcePoint hinge were estimated to decrease their chance of non-contact ACL injury by over 50%.⁵

CLINICALLY-PROVEN PERFORMANCE

By combining 4-Points-of-Leverage with our FourcePoint hinge, DonJoy's bracing line offers the most widely accepted and clinically-proven solution for:

- Protecting a healing ACL graft
- Providing stability to an unstable ACL-deficient knee
- Prophylactic use

It's the powerful combination of performance and peace-of-mind that separates DonJoy from others and continues to make DonJoy the #1 brace company worldwide.

¹Beynon et al. 1992; Beynon et al. 1995; Beynon et al. 1997; Fleming et al. 2000. ²Garrett and Yu, AJSM, 2004. ³Nunley et al. 2003. ⁴Beynon; Markolf. ⁵Yu et al. 2008.

Introducing the World's First Knee Guarantee™

DJO is proud to stand behind the efficacy of the DonJoy bracing line by offering the world's first and only Knee Guarantee Program. If a new patient re-ruptures their reconstructed ACL while wearing a properly fitted custom Defiance/Defiance III within 1 year of purchase, DJO will pay a surgical insurance deductible up to \$1,000.

We're pleased to announce that two of our most popular prefabricated braces, Armor FourcePoint and FULLFORCE, have been added to the Knee Guarantee Program. DJO will pay a surgical insurance deductible up to \$750. Same patient requirements apply as stated above.

The brace must be purchased and registered within six months of surgery. Additional terms and conditions apply as detailed on www.donjoy.com/kneegarantee.



Defiance III® & Defiance III® Female Fource™



IDEAL APPLICATION: The Defiance III® represents DonJoy's flagship custom knee ligament brace, offering a perfect fit for extreme, water, contact, and non-contact sports. Indications include: hyperextension, moderate to severe ACL or PCL instabilities, ACL/PCL instabilities combined (CI), MCL/LCL instabilities, ACL or PCL reconstructions, prophylactic use. Female Fource is an advanced custom knee ligament brace designed specifically for females.

Ordering information

Each Defiance III is custom made.
Measurements taken with CCMI Mark III or casted in weight bearing.

Size	Thigh	Calf
Minimum	12½" (31.8 cm)	10¼" (26 cm)
Maximum	38¾" (97 cm)	26¼" (66.5 cm)

Accessories	Description	Size
11-0075-x-00000	Black Lycra® Undergarment	S-XXL
11-0122-x	Black Neoprene Undergarment	XS-XXL
11-1006-x-06000	Impact Guard	S-XXL
11-1367-9-06000	Calf Pinch Guard	Universal
11-0398-9-00000	Full Extension Stop	Universal
25-0422-0-00000	Flexion / Extension Stop Kit	Universal
11-0103-x-06000	Sports Cover (short)	S-XXL
11-0016-x-06000	Sports Cover (standard)	S-XXL
	Slide Shield	Universal
	Suede Strap Pad Kit	Universal
	Suede Liner Kit	Universal
	Chamois Strap Pad Kit	Universal
	Chamois Liner Kit	Universal
	Chamois Pneumatic	Universal
11-0400-9-00000	Pneumatic Condyle Pad	Universal
25-1923-0-06000	Air Condyle Kit	Universal
25-1480-0-00000	Brace Bag	Universal
11-1525-x	Tru-Pull Advance Attachment	XS-XXXL
11-1415	Sili-Grip Strap Pads	Universal
11-0543-x-06000	Neoprene Suspension Strap	XS-XXL

Proceeds from Female Fource sales will be donated to support breast cancer initiatives.

Features & Benefits

- Clinically-proven 4-Points-of-Leverage System
- Backed by DonJoy's Knee Guarantee Program
- Custom fabricated
- Strong, lightweight carbon composite frame
- Swooping Thigh Cuff for inner thigh clearance and bilateral brace use
- Custom Graphics Kit (FemaleFource only)
- FemaleFource only offered with FourcePoint hinge
- Lifetime warranty on frame and hinges
- 24 hour turnaround time

Options

- Patented clinically-proven FourcePoint hinge technology (stop kit not included)
- Over 40 custom frame colors and graphics
- Defiance III Extra includes reinforced carbon composite frame
- PowerCaps—metal strap tab caps (recommended for contact sports)
- High strength stainless steel gears (recommended for contact sports)
- ACL, PCL, or Combined Instability (CI) strap configurations
- Standard or short frames available





Defiance - Standard Shape Thigh

IDEAL APPLICATION: Defiance is the foundation of DonJoy's custom bracing line. Offering a true custom fit, DEFIANCE is recommended for extreme, water, contact and non-contact sports. Indications include: hyperextension, moderate to severe ACL or PCL instabilities, ACL/PCL combined instabilities (CI), MCL/LCL instabilities, ACL or PCL reconstruction and prophylactic use.

Features & Benefits

- Clinically-proven 4-Points-of-Leverage System
- Backed by DonJoy's Knee Guarantee Program
- Custom fabricated
- Strong, lightweight carbon composite frame
- Lifetime warranty on frame and hinges
- 24-hour turnaround time

Options

- Over 40 custom frame colors and graphics
- Patented clinically-proven FourcePoint hinge technology (stop kit not included)
- Defiance Extra includes reinforced carbon composite frame
- PowerCaps—metal strap tab caps (recommended for contact sports)
- High strength stainless steel gears
- ACL, PCL, or Combined Instability (CI) strap configurations
- Standard or short frames available

Ordering information

Measurements taken with CCMI Mark III or casted in weight bearing.

Size	Thigh	Calf
Minimum	12½" (31.8 cm)	10¼" (26 cm)
Maximum	38½" (97 cm)	26¼" (66.5 cm)

Accessories	Description	Size
11-0075-x-00000	Black Lycra® Undergarment	S-XXL
11-0122-x	Black Neoprene Undergarment	XS-XXL
11-1006-x-06000	Impact Guard	S-XXL
11-1367-9-06000	Calf Pinch Guard	Universal
11-0398-9-0000	Full Extension Stop	Universal
25-0422-0-00000	Flexion / Extension Stop Kit	Universal
11-0103-x-06000	Sports Cover (short)	S-XXL
11-0016-x-06000	Sports Cover (standard)	S-XXL
	Slide Shield	Universal
	Suede Strap Pad Kit	Universal
	Suede Liner Kit	Universal
	Chamois Strap Pad Kit	Universal
	Chamois Liner Kit	Universal
	Chamois Pneumatic	Universal
11-0400-9-00000	Pneumatic Condyle Pad	Universal
25-1923-0-06000	Air Condyle Kit	Universal
25-1480-0-00000	Brace Bag	Universal
11-1525-x	Tru-Pull Advance Attachment	XS-XXXL
11-1415	Sili-Grip Strap Pads	Universal
11-0543-x-06000	Neoprene Suspension Strap	XS-XXL

Defiance Color Choices

NON-METALLIC COLOR



METALLIC COLOR



CCMI Mark III - Custom Measurement Device



IDEAL APPLICATION: The Custom Contour Measuring Instrument (CCMI) is used to measure custom DonJoy Defiance and Defiance III braces. The device takes into account fourteen separate anatomical reference points from the patient's leg indicating thigh and calf circumference, thigh and calf contours, as well as knee width and varus/valgus.

Features & Benefits

- New sleek design to determine leg contour measurements for the construction of a DonJoy custom brace
- Retractable tape measures for quick and accurate measurements
- Standard and Short brace measuring options
- Easy to use and collapsible thigh and calf contour gauges

Ordering information

Part Number	Custom Brace Compatibility
11-0183	Defiance III, Defiance III Female Fource, Defiance, Adjustable OA Defiance

Accessories	Description
11-0025	Replacement Tape Module

GRADIENT COLORS



Combination of any 2 non-metallic colors only

GRAPHICS COLORS (extra charge)



Tiger Red, Butterfly, Stars & Stripes Non-metallic, Tiger Orange, Snake, Stars & Stripes Metallic, Camouflage, Leopard Blue



Armor

IDEAL APPLICATION: Represents DonJoy's strongest prefabricated knee ligament brace and is ideally suited for extreme and contact sports. Indications include mild to severe ACL, PCL, MCL and LCL instabilities.

Features & Benefits

- Clinically-proven 4-Points-of-Leverage System
- Backed by DonJoy's Knee Guarantee program, an industry first for prefabricated bracing (FourcePoint version only)
- Swooping Thigh Cuff for inner thigh clearance and bilateral brace use
- Strong aircraft grade aluminum frame (0.125" thick)
- 1 year warranty on frame and hinges

Options

- Patented clinically-proven FourcePoint hinge technology (stop kit not included)
- ACL, PCL, or Combined Instability (CI) strap configurations
- Standard or short frames available

Ordering information

Armor Action	Ligament	Hinge	Calf Length	Leg
11-1440-x	ACL	FourcePoint	Short	Right
11-1441-x	ACL	FourcePoint	Short	Left
11-1442-x	ACL	FourcePoint	Standard	Right
11-1443-x	ACL	FourcePoint	Standard	Left
11-1004-x-26000	ACL	Standard	Short	Right
11-1005-x-26000	ACL	Standard	Short	Left
11-1444-x	ACL	Standard	Standard	Right
11-1445-x	ACL	Standard	Standard	Left
11-1029-x-26000	CI	Standard	Short	Right
11-1030-x-26000	CI	Standard	Short	Left
11-1447-x	CI	Standard	Standard	Right
11-1448-x	CI	Standard	Standard	Left

Accessories	Description	Size
11-0075-x-00000	Black Lycra® Undergarment	S-XXL
11-0122-x	Black Neoprene Undergarment	XS-XXL
11-1006-x-06000	Patella / Shin Impact Guard	S-XXL
11-1095-9-06000	Calf Pinch Guard	Universal
11-0956-0-06000	Supra Condyle Pad	Universal
11-0398-9-00000	Full Extension Stop	Universal
25-0422-0-00000	Flexion / Extension Stop Kit	Universal
25-2428	Mesh Brace Bag	Universal
11-0016-x-06000	Sports Cover (standard)	S-XXL
11-1525-x	Tru-Pull Advance Attachment	XS-XXXL
11-1415	Sili-Grip Strap Pads	Universal
11-0543-x-06000	Neoprene Suspension Strap	XS-XXL

Sizes Available – page 10



FULLFORCE™

IDEAL APPLICATION: Combines DonJoy's proven 4-Points-of-Leverage System with our patented FourcePoint hinge technology, the most powerful and clinically-proven solutions available to protect the ACL^{1,2}. Created for all activity levels, FULLFORCE's sleek, low profile and exceptionally light weight have made it the go-to prefabricated brace for competitive and recreational athletes of all kinds. (¹Beynnon et al. 1992. Beynnon et al. 1995. Beynnon et al. 1997. Fleming et al. 2000. ²Garrett and Yu, AJSM, 2004)

Ordering information

Part Number	Ligament	Hinge	Calf Length	Leg
11-0258-x	ACL	FourcePoint	Standard	Right
11-0259-x	ACL	FourcePoint	Standard	Left
11-3220-x	ACL	FourcePoint	Short	Right
11-3221-x	ACL	FourcePoint	Short	Left
11-0264-x	CI	Fource Point	Standard	Right
11-0265-x	CI	FourcePoint	Standard	Left
11-3222-x	CI	FourcePoint	Short	Right
11-3223-x	CI	FourcePoint	Short	Left

Accessories	Description	Size
11-0075-x-00000	Black Lycra® Undergarment	S-XXL
11-0122-x	Black Neoprene Undergarment	XS-XXL
11-1006-x-06000	Patella / Shin Impact Guard	S-XXL
11-1095-9-06000	Calf Pinch Guard	Universal
11-0956-0-06000	Supra Condyle Pad	Universal
11-0398-9-00000	Full Extension Stop	Universal
25-0422-0-00000	Flexion / Extension Stop Kit	Universal
25-2428	Mesh Brace Bag	Universal
11-0016-x-06000	Sports Cover (standard)	S-XXL
11-1525-x	Tru-Pull Advance Attachment	XS-XXXL
11-1415	Sili-Grip Strap Pads	Universal
11-0543-x-06000	Neoprene Suspension Strap	XS-XXL

Sizes Available – page 10

Features & Benefits

- Clinically-proven 4-Points-of-Leverage System
- Patented clinically-proven FourcePoint hinge technology (stop kit not included)
- Backed by DonJoy's Knee Guarantee program, an industry first for prefabricated bracing
- Internally mounted swiveling strap tabs help accommodate leg movement
- Swooping Thigh Cuff for inner thigh clearance and bilateral brace use
- Strong, lightweight aircraft grade aluminum frame
- 1 year warranty on frame and hinges



	Brace Type	Activity Level		
		High	Moderate	Mild
Defiance™ 	Custom	[Blue bar spanning High to Moderate]		
Armor FourcePoint™ 	Prefabricated	[Red bar spanning High to Moderate]		
FULLFORCE™ 	Prefabricated		[Light blue bar spanning Moderate to Mild]	



Playmaker® – Drytex® & Neoprene



IDEAL APPLICATION: DonJoy's premier soft good knee ligament brace. Recommended for non-contact sports and activities of daily living. Suitable for mild to moderate ligament instabilities. Ideal intermediate-level product for sporting activities.

- Features & Benefits**
- Warm and compressive Neoprene or cool and breathable Drytex® construction
 - Sleeve or Wraparound styles for optimal fit and comfort
 - Popliteal cutout for improved fit and comfort
- Options**
- Patented clinically-proven FourcePoint hinge technology
 - Versatile strapping allows ACL, PCL, or Combined Instability (CI) configurations
 - IROM hinge ideal for intermediate post-op care

Ordering information

Part Number	Hinge	Style	Material
11-0557-x	Standard	Sleeve	Drytex®
11-0558-x	Standard	Sleeve w/popliteal cutout	Drytex®
11-1396-x	Standard	Sleeve w/patella donut	Drytex®
11-1395-x	Standard	Sleeve w/popliteal cutout & patella donut	Drytex®
11-0666-x	Standard	Wraparound	Drytex®
11-0759-x	Standard	Wraparound w/popliteal cutout	Drytex®
11-0860-x	Standard	Sleeve	Neoprene
11-0861-x	Standard	Sleeve w/popliteal cutout	Neoprene
11-0862-x	Standard	Sleeve w/patella donut	Neoprene
11-0863-x	Standard	Sleeve w/popliteal cutout & patella donut	Neoprene
11-0864-x	Standard	Wraparound	Neoprene
11-0865-x	Standard	Wraparound w/popliteal cutout	Neoprene

Accessories	Description	Size
25-0422-0-00000	Flexion / Extension Stop Kit	Universal
11-1524	FourcePoint Hinge Kit	Universal

Sizes Available – page 10

Sizes Available - Prefabricated Knee Bracing

Measurements taken 6" above (thigh) or 6" below (calf) mid-petella.

Size	Thigh	Knee Center	Calf
x = 1 = XS	13"-15½" (33-39 cm)	12"-13" (30.5-33 cm)	10"-12" (25.5-30.5 cm)
x = 2 = S	15½"-18½" (39-47 cm)	13"-14" (33-35.5 cm)	12"-14" (30.5-35.5 cm)
x = 3 = M	18½"-21" (47-53.25 cm)	14"-15" (35.5-38 cm)	14"-16" (35.5-40.5 cm)
x = 4 = L	21"-23½" (53.25-59.5 cm)	15"-17" (38-43 cm)	16"-18" (40.5-47 cm)
x = 5 = XL	23½"-26½" (59.5-67.25 cm)	17"-19" (43-48.25 cm)	18"-20" (47-50.75 cm)
x = 6 = XXL	26½"-29½" (67.25-75 cm)	19"-21" (48.25-53.25 cm)	20"-22" (50.75-56 cm)
x = 7 = XXXL	29½"-32" (75-81.25 cm)	21"-23" (56-61 cm)	22"-24" (53.25-58.5 cm)

Hinged Knee



IDEAL APPLICATION: The Drytex Hinged Knee, designed for mild to moderate medial lateral control, is a neoprene alternative constructed of a unique nylon core and polyester Lycra® fabric that allows for improved airflow. This fabric provides compression and support, yet allows for breathability. Ideal for patients who have allergies or irritations to neoprene or live in warm weather climates.

- Features & Benefits**
- Warm and compressive Neoprene or cool and breathable Drytex® construction
 - Sleeve or Wraparound styles for optimal fit and comfort

Ordering information

Part Number	Hinge	Style	Material
11-2006-x	Standard	Sleeve w/horseshoe buttress	Neoprene
11-2007-x	Standard	Sleeve w/popliteal cutout & horseshoe buttress	Neoprene
11-2008-x	Standard	Wraparound	Neoprene
11-2009-x	Standard	Wraparound w/popliteal cutout	Neoprene
11-0555-x*	Standard	Sleeve	Drytex®
11-0556-x*	Standard	Sleeve w/popliteal cutout	Drytex®
11-2010-x	Standard	Sleeve w/horseshoe buttress	Drytex®
11-2011-x	Standard	Sleeve w/popliteal cutout & horseshoe buttress	Drytex®
11-2012-x	Standard	Wraparound	Drytex®
11-2013-x	Standard	Wraparound w/popliteal cutout	Drytex®

*XS - XXL only

Accessories	Description	Size
25-0422-0-00000	Flexion / Extension Stop Kit	Universal

Sizes Available – page 10



Osteoarthritis Bracing



OA Assist

Custom Adjustable OA Defiance®

IDEAL APPLICATION: DonJoy's premier Custom OA product. It is designed for the active OA patient, is extremely lightweight, and provides durable support for moderate to severe levels of osteoarthritis. The bilateral hinge bars offer support for ligament instabilities and work in conjunction with a patented telescoping condyle pad which puts the load control in the hands of the patient.

Ordering information

Sizing

Every OA DEFIANCE is custom fabricated with a 24 hour turnaround time using either DonJoy's CCMI Mark III measuring device or a cast mold. See page 5 for further detail on CCMI Mark III.

Size	Thigh Circumference	Calf Circumference
Minimum:	12½" (31.8 cm)	10¼" (26 cm)
Maximum:	38⅝" (97 cm)	26¼" (66.7 cm)

Brace Length	Patient Height
Standard Calf:	5'4" or taller (163 cm or taller)
Short Calf:	5'4" or shorter (163 cm or shorter)

Accessories	Description	Size
11-0075-x-00000	Lycra® Undergarment/Oversleeve	S-XXL
11-0122-x	Neoprene Undergarment	XS-XXL
11-0016-x-06000	Brace Cover (standard)	S-XXL
11-0103-x-06000	Brace Cover (short)	S-XXL
11-1525-x	Tru-Pull® Advanced Patellofemoral Attachment	XS-XXL
11-0398-9-00000	Full Extension Stop	Universal
11-0400-9-00000	Pneumatic Condyle Pad	Universal
25-1923-0-06000	Air Condyle Pads	Universal
11-1107-0-00000	Floam Condyle Pads	Universal
11-1006-x-06000	Impact Guard	S-XXL
11-1095-9-06000	Calf Pinch Guard	Universal
11-0857-9-00000	Load Sensor Kit	Universal
25-1480-0-00000	Brace Bag	Universal
11-1415	Sili-Grip Strap Pads	Universal
11-0543-x-06000	Neoprene Suspension Strap	XS-XXL

Features & Benefits

- 4-Points-of-Leverage™ System pioneered by DonJoy provides ligament stability
- 1° to 7° of correction built into a strong, lightweight, carbon composite frame (3° is recommended)
- Patient controlled Telescoping Condyle Pad (up to 3° of additional correction) shifts compressive knee joint forces from the degenerative area to the healthy compartment
- 24 hour turnaround time

Options

- Over 40 custom frame colors and graphics
- Swooping Thigh cuff for inner thigh clearance and bilateral brace use
- Swiveling Strap Tabs to accommodate leg movement
- Reinforced carbon composite frame ACL, PCL, or Combined Instability (CI) strap configurations



Swooping Medial Thigh



Standard Thigh



Telescoping Pad



3-Point Load Off-Loading

Mechanically shift weight away from the damaged part of the knee



90 DAY TEST DRIVE

Braces that don't fit properly can be uncomfortable to wear, slip out of place, or produce excessive pressure that can lead to skin problems and soreness. DonJoy's custom and patient-ready OA braces accommodate most people and every brace comes with our no-risk TEST DRIVE offer.

If you are not 100% satisfied with your brace and it cannot be corrected by a DonJoy representative, return it for a full credit or refund within 90 days.



EASE OF USE

The Arthritis Foundation has honored every DonJoy knee osteoarthritis braces with its coveted "Ease-of-Use" Commendation after an independent research lab conducted rigorous testing to assure that they met the needs of people with functional limitations due to the effects of arthritis.

OA Adjuster™



IDEAL APPLICATION: The DonJoy OA Adjuster prefabricated knee osteoarthritis brace is clinically proven* to off-load the knee and is ideally suited for active patients with moderate to severe unicompartmental osteoarthritis, ligament instabilities or post-operative OCD rehabilitation. (* Matthew Nadaud et al. "In Vivo Three-Dimensional Determination of the Effectiveness of the Osteoarthritic Knee Brace: A Multiple Brace Analysis." Journal of Bone & Joint Surgery. 87:114-119, 2005.)

Features & Benefits

- 4-Points-of-Leverage™ System
- Lock and Off-Load Technology™ shifts compressive knee joint forces from the degenerative area to the healthy compartment
- Universal OA Key enables patient-controlled load adjustment to meet daily needs
- Floating Hinges allow a contoured fit with a double upright design

Ordering information

Part Number	Affected Compartment	Leg
11-0870-x	Medial	Right
11-0871-x	Medial	Left
11-0872-x	Lateral	Right
11-0873-x	Lateral	Left

Accessories	Description	Size
11-0075-x-00000	Lycra® Undergarment/Oversleeve	S-XXL
11-0857-9-00000	Load Sensor Kit	Universal
11-1100-x-06000	Strap Replacement Kit	XS-XXXL
11-1547-x	Suede Liner Replacement Kit (Medial, Right)	XS-XXXL
11-1548-x	Suede Liner Replacement Kit (Medial, Left)	XS-XXXL
11-1549-x	Suede Liner Replacement Kit (Lateral, Right)	XS-XXXL
11-1550-x	Suede Liner Replacement Kit (Lateral, Left)	XS-XXXL
11-1552-x	Suede Strap Pad Replacement Kit	XS-XXXL
11-1553	Suede Condyle Pad Replacement Kit	Universal
25-1376-0-00000	Flexion/Extension Stop Kit	Universal
11-1415	Sili-Grip Strap Pads	Universal
11-0543-x-06000	Neoprene Suspension Strap	XS-XXL

OA Assist



IDEAL APPLICATION: DonJoy's technologically advanced design combines a single upright frame and the 3DF hinge on the OA Assist. A patient-ready knee osteoarthritis brace that is ideally suited for activities of daily living.

Features & Benefits

- 3rd Degree of Freedom (3DF) hinge allows rotation on the calf cuff for an optimal fit. (+ / - 11° of rotation from the neutral position)
- Single upright frame with 3 point load that shifts compressive knee joint forces from degenerative area to healthy compartment
- Telescoping condyle pad (up to 6° of correction) enables patients to adjust load according to activity level
- Sili-Grip strap pads and liners promote brace suspension

Ordering information

Part Number	Affected Compartment	Leg
11-0947-x	Medial	Right
11-0948-x	Lateral	Right
11-0949-x	Medial	Left
11-0950-x	Lateral	Left

Accessories	Description	Size
11-0075-x-00000	Lycra® Undergarment/Oversleeve	S-XXL
11-0122-x	Neoprene Undergarment	S-XXL
11-0973-x	Strap Replacement Kit	XS-XXXL
11-1415-x	Sili-Grip Strap Pad Replacement Kit	XS-XXXL
11-0975-x	Liner Replacement Kit (Right Medial/Left Lateral)	XS-XXXL
11-0976-x	Liner Replacement Kit (Left Medial/Right Lateral)	XS-XXXL
11-1420	2 in 1 Strap Pad Kit	Universal
11-0543-x-06000	Neoprene Strap Kit	1=XS, 2=S/M, 3=L-XXL

OA Everyday™

IDEAL APPLICATION: The OA Everyday represents DonJoy's easiest to use, prefabricated knee osteoarthritis brace and is ideally suited for activities of daily living.

Ordering information

Part Number	Affected Compartment	Leg
11-1338-x	Medial	Right
11-1339-x	Lateral	Right
11-1340-x	Medial	Left
11-1341-x	Lateral	Left

Accessories	Description	Size
11-0075-x-00000	Lycra® Undergarment/Oversleeve	S-XXL
11-0122-x	Neoprene Undergarment	XS-XXL
11-1414-x	Strap Replacement Kit	XS-XXXL
11-1416-x	Liner Replacement Kit (Right Medial/Left Lateral)	XS-XXX
11-1417-x	Liner Replacement Kit (Left Medial/Right Lateral)	XS-XXXL
11-1415	Sili-Grip Strap Pad Replacement Kit	Universal

Features & Benefits

- 3 point load shifts compressive knee joint forces from the degenerative area to the healthy compartment
- Telescoping condyle pad (up to 6° of correction) enables patients to adjust load according to activity level
- Swiveling, quick release buckles facilitate ease of application/removal and accommodates leg movement
- Short brace length (12.5") accommodates short inseams



OA Lite™

IDEAL APPLICATION: The OA Lite is DonJoy's most comfortable, lightest weight, lowest profile, prefabricated knee osteoarthritis brace and is ideally suited for activities of daily living and sedentary lifestyles.

Ordering information

Part Number	Affected Compartment	Leg
11-1531-x	Medial	Right
11-1532-x	Medial	Left
11-1533-x	Lateral	Right
11-1534-x	Lateral	Left

Accessories	Description	Size
11-0644-9-00000	Goniometer/Bending Tool	Universal
25-0422-0-00000	Flexion/Extension Stop Kit	Universal
25-2436-2	Strap Replacement Kit	XS/S
25-2436-3	Strap Replacement Kit	M
25-2436-4	Strap Replacement Kit	L/XL
25-2436-5	Strap Replacement Kit	XXL/XXXL
25-2438	Padding Kit	Universal
11-1541	Flex Knee Therapy Kit	Universal

Features & Benefits

- Single upright is malleable to allow further load adjustments with thigh and calf paddles to mitigate hinge rotation
- Pre-loaded hinge (up to 6° of correction) incorporates flexion and extension control
- Arthritis friendly with reinforced strap ends, thigh grip, and pull loops
- Breathe-O-Prene® fabric eliminates skin irritation and provides comfort by wicking away moisture



Sizes Available - Prefabricated OA Knee Bracing

(Measurements taken 6" above mid-patella)

- x = 1 = XS...13" - 15½"...(33-39 cm)
- x = 2 = S...15½" - 18½"...(39-47 cm)
- x = 3 = M...18½" - 21"...(47-53 cm)
- x = 4 = L...21" - 23½"...(53-60 cm)
- x = 5 = XL...23½" - 26½"...(60-67 cm)
- x = 6 = XXL...26½" - 29½"...(67-75 cm)
- x = 7 = XXXL...29½" - 32"...(75-81 cm)

Post-Op Bracing



TROM Adjuster

Telescoping TROM™ Advance (Total Range of Motion)

IDEAL APPLICATION: DonJoy's premier telescoping technology to provide the ultimate in fit, comfort and adjustability while giving clinicians a broader range of motion options.

Ordering information

Part Number	Description	Size
11-9114-9	TROM Advance, Cool	< 26" (66 cm)
11-9116-9	TROM Advance, Cool, XL	> 26" (> 66 cm girth)
11-9115-9	TROM Advance, Full Foam	< 26" (66 cm)
11-9117-9	TROM Advance, Full Foam, XL	> 26" (> 66 cm girth)

Accessories	Description
11-0024-x	Ankle Stirrup, Right
11-0026-x	Ankle Stirrup, Left

Features & Benefits

- Tele-Fit allows independent strap movement proximal and distal of knee allowing an overall better fit
- Telescoping adjustability from 17 3/4" to 23 1/4"
- Quick-lock mechanism from (negative) -10° to 90°
- Quick release buckles
- Sili grip strap pads



Tele-Fit



TROM Adjuster

IDEAL APPLICATION: Create the optimal knee off-load with just a few turns. TROM Adjuster is the first post-op brace to offer controlled and consistent off-loading to create a better healing environment. This brace represents a comprehensive post operative solution for meniscal and chondral defect procedures.

Ordering information

Part Number	Description	Size
11-5115-9	Medial RT/Lateral LT	< 26" (66 cm)
11-5116-9	Medial LT/Lateral RT	< 26" (66 cm)
11-5117-9	Medial RT/Lateral LT, XL	26"-32" (66-80 cm) girth
11-5118-9	Medial LT/Lateral RT, XL	26"-32" (66-80 cm) girth


Accessories	Description
11-0024-x	Ankle Stirrup, Right
11-0026-x	Ankle Stirrup, Left

Features & Benefits

- Off Loading Telescoping hinge
- Tele-Fit allows independent strap movement proximal and distal of knee
- Quick-lock mechanism from (negative) -10° to 90°
- Sili Grip Strap pads



Off-loading adjustable hinge technology



+ adjustment



- adjustment



Patellofemoral Bracing



Tru-Pull Lite

REACTION Knee Brace

NEW!

IDEAL APPLICATION: Uses an innovative web approach for absorbing shock and relieving knee pain caused by Chondromalacia. Patellar Tendonitis/Tendinosis, general Patellofemoral tracking issues.



Features & Benefits

- Shock-absorbing web that shifts peak loads away from the painful area of the knee and stabilizes the patella on all sides to ensure proper tracking position.
- Flexible Dual-Axis Hinges that work in conjunction with the elastomeric web for optimal fit and support, providing energy dispersion to the knee
- Breathable Comfort from an open framework and mesh backing that makes this breathable and comfortable during use

Ordering information

Part Number	Knee Design	Size
11-0215-x	REACTION Knee Brace	XS – XXXL
11-0541-x	Replacement Under Sleeve	XS – XXXL

Sizes Available – page 21

Tru-Pull® Lite™

IDEAL APPLICATION: Latest patented technology from the Tru-Pull family of products. For acute mild to moderate patellar subluxation or dislocation and maltracking.



Features & Benefits

- Bifurcated strap dynamically pulls the patella into proper alignment reducing pain associated with patellofemoral dysfunction
- Dual Durometer Butress System interfaces directly with the patella applying a consistent corrective force to reduce pain
- Independent anchors at the thigh and calf prevent rotation of the brace
- New anti-microbial, wicking fabric is hypoallergenic, breathable, and light weight—thereby increases patient comfort and compliance
- Set of removable plastic hinges included with every brace to provide medial and lateral stability

Ordering information

Part Number	Knee Design	Size
11-0260-x	Right Sleeve	XS – XXXL
11-0261-x	Left Sleeve	XS – XXXL

Sizes Available – page 21



Featuring the Dual Durometer Butress System

Tru-Pull® Advanced System

IDEAL APPLICATION: Chronic patellar subluxation/dislocation.

Ordering information

Part Number	Knee Design	Size
11-1289-x	Right Sleeve	XS – XXXL
11-1290-x	Left Sleeve	XS – XXXL
11-1292-x	Right Sleeve, popliteal cutout	XS – XXXL
11-1293-x	Left Sleeve, popliteal cutout	XS – XXXL
11-2289-x	Right Sleeve, hinges	XS – XXXL
11-2290-x	Left Sleeve, hinges	XS – XXXL
11-2292-x	Right Sleeve, popliteal cutout, hinges	XS – XXXL
11-2293-x	Left Sleeve, popliteal cutout, hinges	XS – XXXL
11-1525-x	TruPull attachment to a rigid brace	XS – XXXL

Sizes Available – page 21

Features & Benefits

- Bifurcated, elastomeric strap dynamically pulls the patella as the knee extends
- Two Buttresses (¼" and ⅜") included, trimmable to provide superior, inferior, and lateral patellar stabilization
- Independent anchors prevent brace rotation
- Hinged version offers medial/lateral support and flexion/extension control



Tru-Pull Hinged Version

NEW! Hinged Lateral J

IDEAL APPLICATION: An improved design of a dependable patella support. Enhanced stability allows patients to return to function with more confidence and more comfort. For lateral patellar subluxation and dislocations, lateral patellofemoral malalignment, patellofemoral maltracking and pain.

Ordering information

Part Number	Knee Design	Size
11-0775-x	Neoprene, Right	XS – XXL
11-0776-x	Neoprene, Left	XS – XXL
11-0777-x	Drytex, Right	XS – XXL
11-0778-x	Drytex, Left	XS – XXL

Sizes Available – page 21

Features & Benefits

- Enhanced Knee Stability with malleable hinges
- Dependable Patella Support. Tubular lateral "J" buttress provides stabilization of the patella. Lateral to medial strapping allows an intimate fit of the "J" buttress along the lateral aspect of the patella for additional stabilization



Hinged Air DonJoy®

IDEAL APPLICATION: Chondromalacia. Patellar tendonitis. Hyperextension, medial/lateral support.

Ordering information

Part Number	Knee Design	Size
11-0855-x-06000	Right or Left Popliteal cutout, hinges	XS – XXL

Sizes Available – page 21

Features & Benefits

- Inflatable buttresses for medial, lateral, or full patellar stabilization
- Polycentric hinges offer medial/lateral support and flexion/extension control
- Wraparound thigh and sleeve calf design facilitates ease of application/removal
- Cool and breathable Drytex® construction
- Superior/Inferior straps for compression and support



Lateral "J"



IDEAL APPLICATION: Acute patellar subluxation/dislocation.

Features & Benefits

- Sewn-in, tubular "J" buttress provides inferior and lateral patellar stabilization
- Lateral to medial stabilization strap
- Dual spiral stays offer medial and lateral support
- Warm and compressive Neoprene or cool and breathable Drytex® construction

Ordering information

Part Number	Knee Design	Size
11-0659-x-06060	Right Drytex® sleeve	XS – XXL
11-0660-x-06060	Left Drytex® sleeve	XS – XXL
11-0320-x-06060	Right Neoprene sleeve	XS – XXL
11-0321-x-06060	Left Neoprene sleeve	XS – XXL

Sizes Available – page 21

Horseshoe Patella



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis.

Ordering information

Part Number	Knee Design	Size
11-0334-x-06060	Right or Left Sleeve, stays	XS-XXL
11-0344-x-06060	Right or Left Sleeve, popliteal cutout, stays	XS-XXL

Sizes Available – page 21

Features & Benefits

- Adjustable, horseshoe-shaped buttress for patellar stabilization
- Superior/Inferior straps for compression and support
- Dual spiral stays offer medial and lateral support
- Warm and compressive Neoprene construction

Hinged "H" Buttress



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis. Hyperextension, medial/lateral support.

Features & Benefits

- Sewn-in "H" buttress provides inferior, medial, and lateral patellar stabilization
- Polycentric hinges offer medial/lateral support and flexion/extension control
- Anterior cover for open patella and hinges
- Superior/Inferior straps for compression and support
- Warm and compressive Neoprene construction

Ordering information

Part Number	Knee Design	Size
11-0828-x-06060	Right or Left Sleeve, popliteal cutout, hinges	XS – XXXL

Sizes Available – page 21

Deluxe Donut



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis.

Ordering information

Part Number	Knee Design	Size
11-0327-x-06060	Right or Left Sleeve, stays	XS-XXL
11-0380-x-06060	Right or Left Sleeve, popliteal cutout, stays	XS-XXL

Sizes Available – page 21

Features & Benefits

- Sewn-in, trimable, tubular donut provides inferior, superior, medial, and lateral patellar stabilization
- Superior/Inferior straps for compression and support
- Dual spiral stays offer medial and lateral support
- Warm and compressive Neoprene construction

"H" Buttress



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis.

Features & Benefits

- Sewn-in "H" buttress provides inferior, medial, and lateral patellar stabilization
- Superior/Inferior straps for compression and support
- Dual spiral stays offer medial and lateral support
- Warm and compressive Neoprene construction

Ordering information

Part Number	Knee Design	Size
11-0032-x-06060	Right or Left Sleeve, stays	XS – XXL
11-0323-x-06060	Right or Left Sleeve, popliteal cutout, stays	XS – XXL

Sizes Available – page 21

Adjustable Patella Donut



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis.

Ordering information

Part Number	Knee Design	Size
11-0465-x-06000	Right or Left Neoprene sleeve	XS-XXL
11-0379-x-06000	Right or Left Neoprene sleeve, popliteal cutout	XS-XXL
11-0549-x-06000	Right or Left Drytex® sleeve	XS-XXL
11-0554-x-06000	Right or Left Drytex® sleeve, popliteal cutout	XS-XXL

Sizes Available – page 21

Features & Benefits

- Sewn-in, trimable, tubular donut provides inferior, superior, medial, and lateral patellar stabilization
- Warm and compressive Neoprene or cool and breathable Drytex® construction
- Available with or without popliteal covering

Deluxe Elastic Knee



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis.

- Features & Benefits**
- Sewn-in, silicone donut provides inferior, superior, medial, and lateral patellar stabilization
 - Single spiral stays offer medial and lateral support
 - Lightweight, breathable, elastic cotton construction
 - Available in beige or white

Ordering information

Part Number	Knee Design	Size
11-0345-x-04000	...Right or Left Beige, cotton sleeve, stays	...XS – XXL
11-0345-x-15000	...Right or Left White, cotton sleeve, stays	...XS – XXL

[Sizes Available – page 21](#)

Elastic Knee



IDEAL APPLICATION: Chondromalacia. Patellar tendonitis.

- Features & Benefits**
- Sewn-in, silicone donut provides inferior, superior, medial, and lateral patellar stabilization
 - Dual spiral stays offer medial and lateral support
 - Lightweight, breathable, elastic cotton construction
 - Available in white

Ordering information

Part Number	Knee Design	Size
11-0242-x Right or Left White, cotton sleeve, stays S – XL

[Sizes Available – page 21](#)

Cross Strap



IDEAL APPLICATION: Patellar tendonitis. Osgood-Schlatter disease.

- Features & Benefits**
- Provides inferior patellar stabilization
 - Alleviates pain by applying pressure to the patellar tendon
 - Warm and compressive Neoprene construction

Ordering information

Part Number	Knee	Straps/Pack
11-1327-x-06000Right or Left 1
11-1328-x-06000Right or Left 25

Sizes Available (Measurements around the knee at the base of the knee cap)
 x = 3 = M..... Up to 14" x = 4 = L..... Over 14"

Sizes Available - Patellofemoral Bracing

Thigh circumference taken 6" above mid-patella

- | | | | |
|---------------------------------------|--|--|-------------------------------------|
| x = 1 = XS...13" - 15½"...(33-39 cm) | x = 2 = S...15½" - 18½"...(39-47 cm) | x = 3 = M...18½" - 21"...(47-53 cm) | x = 4 = L...21" - 23½"...(53-60 cm) |
| x = 5 = XL...23½" - 26½"...(60-67 cm) | x = 6 = XXL...26½" - 29½"...(67-75 cm) | x = 7 = XXXL...29½" - 32"...(75-81 cm) | |



Ever since he heard his knee pop,
 Garrett McNamara has had difficulty
 getting out of the water.

Yes, life is good.

Garrett McNamara, World Champion Big Wave Surfer

Walking Braces

Aircast XP, FP and SP Walker™

IDEAL APPLICATION: Designed to provide pneumatic support for stable fractures of the foot and ankle, severe ankle sprains, soft tissue injury, hindfoot and heel protection, metatarsal fractures, forefoot and midfoot injury, bunionectomies and acute or post-operative use.

Ordering information

Part Number	Shoe Size		Size
	Men	Women	
XP Walker™ (extra pneumatic)			
01P-P	Up to 4	Up to 5	Pediatric
01P-S	4-7	5-8	Small
01P-M	7-10	8-11	Medium
01P-L	10-13	11-15	Large
01P-XL	13 +	15 +	X-Large
FP Walker™ (foam pneumatic)			
01F-P	Up to 4	Up to 5	Pediatric
01F-S	4-7	5-8	Small
01F-M	7-10	8-11	Medium
01F-L	10-13	11-15	Large
01F-XL	13 +	15 +	X-Large
SP Walker™ (short pneumatic)			
01A-P	Up to 4	Up to 5	Pediatric
01A-S	4-7	5-8	Small
01A-M	7-10	8-11	Medium
01A-L	10-13	11-15	Large
01A-XL	13 +	15 +	X-Large

Features & Benefits

XP Walker™ (extra pneumatic)

- Overlapping Duplex aircells provide intermittent pneumatic compression for edema reduction and callus formation
- Aircells can be custom inflated for a "total contact" fit
- Lightweight, durable, semi-rigid full-shell protection
- Includes two socks and a hand bulb
- Can be used with Brace-Lok™ system

FP Walker™ (foam pneumatic)

- Two adjustable distal aircells provide compression and support the malleoli
- Aircells can be individually inflated for a "total contact" fit
- Lightweight, durable, semi-rigid shell protection
- Includes a hand bulb

SP Walker™ (short pneumatic)

- Pre-inflated anterior aircell and two adjustable distal aircells provide compression and support the malleoli
- Aircells can be individually inflated for a "total contact" fit
- Lightweight, durable, semi-rigid shell protection



XP



FP



SP

MaxTrax® Air / MaxTrax® Air Ankle



IDEAL APPLICATION: Designed to increase stability and help decrease pain and edema following trauma or post-operative procedures. Cushioned inner/outer soles are designed to absorb shock upon heel strike and improve patient comfort. Low profile rocker bottom helps promote natural gait and reduce plantar pressures.

Features & Benefits
<ul style="list-style-type: none"> Adjustable pneumatic air chambers accommodate changes in edema during rehab and maximize patient compliance Integrated pump designed to make the liner simple to inflate and adjust Wide footbed increases stability and accommodates edema or bandages Swiveling D-rings allow for crisscross or parallel strapping Malleable uprights breakdown for shorter tibia/fibula on standard height walking brace

Ordering information		
Part Number	Description	Size
11-1371-x	MaxTrax Air	XS - XL
11-1380-x	MaxTrax Air Ankle	XS - XL

MaxTrax® ROM / MaxTrax® ROM Ankle

IDEAL APPLICATION: Provides protected range of motion following trauma or post-operative procedures. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5° and 30° plantar- and dorsi-flexion.

Ordering information		
Part Number	Description	Size
11-1383-x-00000	MaxTrax ROM	XS - XL
11-1385-x-00000	MaxTrax ROM Ankle	XS - XL

Features & Benefits
<ul style="list-style-type: none"> Same features and benefits as MaxTrax® ROM hinge promotes early range of motion



MaxTrax® ROM Air / MaxTrax® ROM Air Ankle



IDEAL APPLICATION: Designed to increase stability and help decrease pain and edema following trauma or post-operative procedures. Provides a protected range of motion following trauma or post-operative procedures. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5° and 30° plantar-flexion and dorsi-flexion.

Features & Benefits
<ul style="list-style-type: none"> Same features and benefits as MaxTrax Air ROM hinge promotes early Range of motion

Ordering information		
Part Number	Description	Size
11-1381-x	MaxTrax ROM Air	XS - XL
11-1390-x	MaxTrax ROM Air Ankle	XS - XL

Aircast XP Diabetic Walker™

IDEAL APPLICATION: Designed to provide pneumatic support for neuropathic ulcers, hindfoot and heel protection, soft tissue injury, stable foot and ankle fractures, and soft tissue injury. Includes two insoles, three stockings and one hand bulb with pressure gauge.

Ordering information			
Part Number	Shoe Size		Size
	Men	Women	
01PD-S	4-7	5-8	Small
01PD-M	7-10	8-11	Medium
01PD-L	10-13	11-15	Large
01PD-XL	13 +	15 +	X-Large

Features & Benefits
<ul style="list-style-type: none"> Optimized rocker sole for maximum off-loading of the diabetic foot Semi-rigid shells, full aircells and insole to help regulate shear stress Includes two insoles, three stockings and one hand bulb with pressure gauge Can be used with Brace-Lok™ system



MaxTrax® / MaxTrax® Ankle



IDEAL APPLICATION: Designed with a wide foot bed for comfort and stability. The low profile rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer soles provide shock absorption to help improve patient comfort during ambulation.

Features & Benefits
<ul style="list-style-type: none"> Wide footbed for increased stability and accommodation of edema or bandages Malleable uprights break down for shorter tibia/fibula on standard height walking brace One-piece liner design keeps seams away from post-surgical areas Anterior/posterior and medial/lateral shells available for compression and support

Ordering information		
Part Number	Description	Size
11-1382-x-00000	MaxTrax	XS - XL
11-1384-x-00000	MaxTrax Ankle	XS - XL

Accessories	Description
11-1372-x-00000	Liner Pad Kit
11-1373-x-00000	Replacement Liner, tall
11-1377-x-00000	Replacement Liner, short
11-1375-x-00000	Anterior/Posterior and Medial/Lateral Shells
79-95310	Walker strap extension set

MaxTrax® Diabetic Walker

IDEAL APPLICATION: Designed for the treatment of plantar ulcers with Impax™ Grid technology constructed of pre-cut layered Plastazote® foam. The bottom two layers of foam can be easily removed around the affected areas, redistributing load away from healing ulcers.

Ordering information		
Part Number	Description	Size
11-1546-x	MaxTrax Diabetic Walker	XS - XL

Accessories	Description	Size
25-1626-x	Foam Impax Grid	S-XL
25-1457	Toe Cover	Universal
25-1345-x	Replacement Boot Liner	S-XL
25-0241	Upright Padding Kit	Universal

Features & Benefits
<ul style="list-style-type: none"> Wide footbed for increased stability Impax Grid technology redistributes load away from healing ulcers while allowing for a continuous plantar surface Simple to remove or replace Impax foam blocks Plastazote foam has low shear and is simple to clean and maintain Soft liner with no seams in critical areas



Sizes Available - MaxTrax®

	Men Shoe	Women Shoe	Junior Shoe		Men Shoe	Women Shoe		Men Shoe	Women Shoe
x = 1 = XS	N/A	1 to 4	1 to 5	x = 2 = S	1 to 5	4½ - 6	x = 4 = L	10½ - 13½	11½ - 14½
				x = 3 = M	5½ - 10	6½ - 11	x = 5 = XL	14 - 17	N/A

Ankle Braces and Supports



Velocity ES

Velocity™ ES, MS and LS



ES

IDEAL APPLICATION: Velocity represents an evolution in DonJoy ankle brace technology, utilizing several new, proprietary design features, resulting in a new level of ankle brace performance and providing the foundation for successful ankle injury prevention and treatment.

- Features & Benefits**
- R³ (Rapid-Rigid-Ratcheted) Technology™ provides circumferential tibia/fibula compression
 - Prevents abnormal inversion and eversion
 - Rigid uprights protect during impact
 - SpeedWrap™ for secured ankle joint compression (ES & LS versions only)
 - Available in black or white (ES version only)

Ordering information

Part Number	Calf Cuff	Ankle	Size
Velocity™ ES (Extra Support)			
11-1497-x-xx000	Standard	Left	S - L
11-1498-x-xx000	Standard	Right	S - L
11-1499-x-xx000	Wide	Left	S - L
11-1500-x-xx000	Wide	Right	S - L
Velocity™ MS (Moderate Support)			
11-1493-x-06000	Standard	Left	S - L
11-1494-x-06000	Standard	Right	S - L
11-1495-x-06000	Wide	Left	S - L
11-1496-x-06000	Wide	Right	S - L
Velocity™ LS (Light Support)			
11-1489-x-06000	Standard	Left	S - L
11-1490-x-06000	Standard	Right	S - L
11-1491-x-06000	Wide	Left	S - L
11-1492-x-06000	Wide	Right	S - L

Calf Cuff (measure calf circumference 8" from ground):
 Standard Calf:.....7¼" – 10" (18 – 25 cm)
 Wide Calf:.....10" – 12½" (25 – 32 cm)
 Colors:xx = 06 = Black xx = 15 = White



MS

Sizes Available

	Men Shoe	Women Shoe		Men Shoe	Women Shoe		Men Shoe	Women Shoe
x = 2 = S	6 - 8	8 - 9½	x = 3 = M	8 - 12	9½ - 13½	x = 4 = L	12 +	13½ +

Aircast Air-Stirrup® Ankle Brace

IDEAL APPLICATION: Anatomically designed for acute ankle injury, post-operative use, chronic instability, and ankle sprains grade I, II, and III.

Ordering information

Part Number	Description	Height	Size
02JL	Left	6"	Pediatric
02JR	Right	6"	Pediatric
02CL	Left	8¾"	Small
02CR	Right	8¾"	Small
02BL	Left	9"	Medium
02BR	Right	9"	Medium
02AL	Left	10½"	Large
02AR	Right	10½"	Large

- Features & Benefits**
- Anatomically designed semi-rigid shells lined with the patented Duplex™ aircell technology for support and graduated compression during ambulation
 - The compression promotes efficient edema reduction in addition to helping accelerate rehabilitation from ankle sprains



Aircast Air-Stirrup Universe™

IDEAL APPLICATION: Anatomically designed for acute ankle injury, post-operative use, chronic instability, and ankle sprains grade I, II, and III.

Ordering information

Part Number	Description	Height	Size
02E		9½"	Universal
02E12	12 PK	9½"	Universal

- Features & Benefits**
- Universal size eliminates excess inventory
 - Designed with Duplex aircell technology to reduce swelling, encourage earlier ambulation for functional management



Aircast Air-Stirrup® Universe™ Care Kit

IDEAL APPLICATION: Anatomically designed as the complete system for functional management of an ankle sprain. Includes an Air-Stirrup Universe Ankle Brace, ankle wrap, cold pack, exercise band, and a DVD plus booklet showing instructions for treatment and rehabilitation.

Ordering information

Part Number	Height	Size
02EK	9½"	Universal

- Features & Benefits**
- Ankle wrap provides circumferential compression, cold pack to help reduce pain, exercise band to help strengthen the ankle, DVD and instruction booklet provides information for proper care
 - DVD and instruction booklet adhere to standards for patient education as identified by JCAHO



Stabilizing Ankle



IDEAL APPLICATION: The Stabilizing Ankle support features non-stretch nylon figure 8 straps which lock the calcaneus into place to help control abnormal eversion and inversion. Ideal for chronic ankle sprains and prophylactic use.

Features & Benefits
<ul style="list-style-type: none"> Fits right or left Sturdy cool mesh construction Removable medial/lateral plastic stays

Ordering information		
Part Number	Color	Size
11-0451-x-06000	Black	XXS - XXXL
11-0451-x-15000	White	XXS - XXXL

Sizes Available (Measurement is circumference, taken from the base of the heel around ankle, while wearing an athletic sock)

x = 0 = XXS 9" - 10" (23-25 cm)	x = 3 = M 12" - 13" (30-33 cm)	x = 6 = XXL 15" - 16" (38-41 cm)
x = 1 = XS 10" - 11" (25-28 cm)	x = 4 = L 13" - 14" (33-36 cm)	x = 7 = XXXL 16" - 17" (41-43 cm)
x = 2 = S 11" - 12" (28-30 cm)	x = 5 = XL 14" - 15" (36-38 cm)	

Aircast A60™ Ankle Support



IDEAL APPLICATION: Anatomically designed for prophylactic support, protection, and comfort with simplified application. May also be used for chronic ankle instabilities.

Features & Benefits
<ul style="list-style-type: none"> Incorporates a stabilizer molded at a 60 degree angle to help guard against sprains and prevent rollover Manufactured with breathable material, for cool comfort Light-weight anatomic design easily fits in athletic footwear Single strap replaces time-consuming lacing

Ordering information				
Part Number	Description	Men	Women	Size
02TSL	Left	Up to 7	Up to 8½	S
02TSR	Right	Up to 7	Up to 8½	S
02TML	Left	7½ - 11½	9 - 13	M
02TMR	Right	7½ - 11½	9 - 13	M
02TLL	Left	12 +	13½ +	L
02TLR	Right	12 +	13½ +	L

Aircast AirSport™



IDEAL APPLICATION: Anatomically designed for ankle sprains, chronic instability, prophylactic support, protection, and comfort for with simplified application.

Features & Benefits
<ul style="list-style-type: none"> Incorporates clinically proven semi-rigid shell and aircells to provide comfort and support Additional compression and stabilization is provided by anterior talofibular cross strap and integral forefoot and shin wraps "Step-in" design and automatic heel adjustment make the AirSport extremely fast and easy to put on

Ordering information				
Part Number	Description	Men	Women	Size
02MXSL	Left	Up to 5	Up to 5	XS
02MXSR	Right	Up to 5	Up to 5	XS
02MSL	Left	5½ - 7	5½ - 8½	S
02MSR	Right	5½ - 7	5½ - 8½	S
02MML	Left	7½ - 11	9 - 12½	M
02MMR	Right	7½ - 11	9 - 12½	M
02MLL	Left	11½ - 13	13 - 14½	L
02MLR	Right	11½ - 13	13 - 14½	L
02MXLL	Left	13½ +	15 +	XL
02MXLR	Right	13½ +	15 +	XL

AirHeel/DNS Care Kit - Round the Clock Relief

IDEAL APPLICATION: For daytime, the unique Aircast AirHeel is a clinically-proven solution that provides comfort and massaging relief for Plantar Fasciitis. For nighttime, the Aircast Dorsal Night Splint delivers both tailored flexion and optimal comfort for pain relief from Plantar Fasciitis.

Ordering information				
Part Number	Description	Men	Women	Size
09PKS	AirHeel/DNS Care Kit	Up to 7	Up to 8½	S
09PKM	AirHeel/DNS Care Kit	7½ - 10	9 - 11	M
09PKL	AirHeel/DNS Care Kit	10½+	11½+	L

Features & Benefits
<ul style="list-style-type: none"> Helps reduce swelling and discomfort, and enhance circulation Made of lightweight, breathable fabric Simple to apply and to adjust with one strap Anatomically designed to provide optimal comfort and ease-of-use



Aircast AirHeel™

IDEAL APPLICATION: This unique, clinically-proven solution provides comfort and massage relief for Plantar Fasciitis.

Ordering information				
Part Number	Description	Men	Women	Size
09AS	AirHeel Support	Up to 7	Up to 8½	S
09AM	AirHeel Support	7½ - 11	8½ - 12½	M
09AL	AirHeel Support	11½+	13+	L

Features & Benefits
<ul style="list-style-type: none"> Reduces swelling and discomfort while enhancing circulation Two interconnected aircells are located under the foot arch and in the back of the heel, applying pulsating compression throughout the day Fits easily into a shoe with an extremely low-profile design



Aircast Dorsal Night Splint

IDEAL APPLICATION: Delivers both tailored flexion and optimal comfort for pain relief from Plantar Fasciitis, allowing active patients to sleep and start their day on the right foot.

Ordering information				
Part Number	Description	Men	Women	Size
09DS	Dorsal Night Splint	5 - 9½	6 - 10½	S - M
09DL	Dorsal Night Splint	10 - 14	11 - 15	L - XL

Features & Benefits
<ul style="list-style-type: none"> Optimal dorsiflexion adjustment for custom stretch with single dorsiflexion adjustment strap, with quick release tab Maximum comfort for enhanced compliance with pliable softgoods, thermoformed shell liner and anatomically-designed forefoot pad cradle the foot and ankle Secure fit with the optional "assist" strap which secures the toes within the forefoot area Easy-to-use flexible hinge design assists in ease of application Three-point softgood adjustment allows the splint to slide on like a sock



Lower Back Bracing



Mechanical Advantage Pulley System

DonJoy Duel TLSO

IDEAL APPLICATION: Low Back Pain, Low Back Sprains/Strains, Lumbar Disc Displacement, Osteoporosis, Disc Herniation and Degeneration, Post-Operative Laminectomy, Post-Operative Disectomy, IDET Procedure, Spondylolisthesis, Post-Operative Fusion, Spinal Stenosis, Compression Fractures, and Scoliosis.

Features & Benefits

- Rigid anterior panel & posterior panels restrict forward flexion and trunk rotation. Panels are easily removed and can be thermoformed for customized fit
- Telescoping sternal Y-Bar can be adjusted to desired height or locked into position. The alloy can be bent and custom fitted to the desired angle
- Quick release shoulder buckles allow for easy donning and doffing
- Features the patented Mechanical Advantage Pulley System

Sizes Available (Waist Measurement)

- x = 2 = S...25" - 30" (60-75 cm)
- x = 3 = M...30" - 35" (75-89 cm)
- x = 4 = L...35" - 42" (89-107 cm)
- x = 5 = XL...42" - 48" (107-122 cm)
- x = 6 = XXL...48" - 57" (122-145 cm)



DonJoy Back Bracing

IDEAL APPLICATION: Low Back Pain, Low Back Sprains/Strains, Lumbar Disc Displacement, Osteoporosis, Disc Herniation and Degeneration, Post-Operative Laminectomy, Post-Operative Disectomy, IDET Procedure, Spondylolisthesis, Post-Operative Fusion, and Spinal Stenosis.

Ordering information

Part Number	Description	Size
DonJoy LSO with Chairback (8") 11-1562-x	Range of Spinal Support T9-S1	S - XXXL
DonJoy Regular Back (10") 11-1568-x	Range of Spinal Support L1-S1	S - XXXL
DonJoy Low Profile Back (8") 11-1569-x	Range of Spinal Support L1-S1	S - XXXL
DonJoy Lite Back (9") 11-1570-x	Range of Spinal Support L1-L5	S - XXL
Accessories	Description	Size
11-1563-x	Standard Chairback	S - XL
11-1564-x	Short Chairback*	M - XL
11-1565	7" Lordotic Insert (pair)	
11-1566	10" Lordotic Insert (pair)	
11-1567	12" Lordotic Insert (pair)	

*The Short Chairback is designed for shorter torsos, so the frame is 13" in height for all sizes.

Features & Benefits

Features the patented Mechanical Advantage Pulley System

DonJoy LSO with Chairback (8")

- Rigid anterior panel & posterior panels are constructed of heat gun moldable poly-blend plastic customizable to the patient for fluctuations in weight or distension
- Easy to use compression technology with simple hook loop closures
- Maximum comfort for patient compliance with the light weight and low profile construction
- Lordotic foam inserts available in varying degrees for comfort and increased lordosis

DonJoy Regular Back (10") and DonJoy Low Profile Back (8")

- Simple hook-loop closure for ease of application and removal
- Proprietary Mechanical Pulley System panels allow for a smooth and uniform compression
- Easy-to-use pull handle(s) enable patient to quickly adjust brace to desired level of comfort
- Breathable interior spacer fabric with outer mesh allows heat to escape without compromising brace performance

DonJoy Lite Back (9")

- The DonJoy Lite is a lighter version of the Low Profile, without the rigid anterior panel. The Lite provides support for mild lower back pain, from L1 - L5.



DonJoy LSO with Chairback (8")



DonJoy Low Profile Back



Standard Chairback Accessory



Lordotic Insert Accessory

Lordotic foam inserts fit easily into the posterior pockets of the DonJoy 8" and 10" braces.

Sizes Available - DonJoy LSO with Chairback (8"), DonJoy Regular Back (10") and DonJoy Low Profile Back (8") (Waist Measurement)

- x = 2 = S...25" - 30" (60-75 cm)
- x = 3 = M...30" - 35" (75-89 cm)
- x = 4 = L...35" - 40" (89-102 cm)
- x = 5 = XL...40" - 45" (102-114 cm)
- x = 6 = XXL...45" - 50" (114-127 cm)
- x = 7 = XXXL...50" - 55" (127-140 cm)

Sizes Available - DonJoy Lite Back (9") (Waist Measurement)

- x = 2 = S.....26" - 32" (66-81 cm)
- x = 3 = M.....32" - 38" (81-97 cm)
- x = 4 = L.....38" - 43" (97-109 cm)
- x = 5 = XL.....43" - 48" (109-122 cm)
- x = 6 = XXL.....48" - 53" (122-135 cm)

Upper Extremity



UltraSling III

UltraSling®

IDEAL APPLICATION: The mesh UltraSling® III Family provides immobilization for rotator cuff repairs, capsular shifts, Bankhart repairs, glenohumeral dislocations/subluxation and soft tissue repairs/strains.

Features & Benefits

- UltraSling® III / UltraSling® III AB**
- NEW! Moisture wicking mesh fabric for greater comfort, particularly in longer term use (15° of abduction with internal rotation)
 - Promotes auxiliary air exchange to reduce risk of secondary infections and prevents internal rotation contractures
 - AB Version: Large comfortable pillow gives greater abduction where required (45° or 60°)
 - The De-Rotation Strap hooks to the UltraSling® III and prevents internal rotation by securely holding the patient's arm in neutral position
 - Quick release buckles make it easy to take on and off for exercise/therapy

- UltraSling® III ER**
- Maintains the shoulder and arm in an externally rotated position
 - NEW! Moisture wicking mesh fabric for greater comfort, particularly in longer term use
 - 15° and 30° versions to suit your protocol
 - De-Rotation shoulder strap to keep product in place
 - Quick release buckles make it easy to take on and off for exercise/therapy

Ordering information

Part Number	Description	Size
UltraSling III 11-0449-x		S - XL
UltraSling III AB 11-0450-x		S - XL
Accessories	Description	Size
11-1530	De-Rotation Strap	S - XL
UltraSling® III ER		
11-1343-x	15°, Right	S - XL
11-1344-x	15°, Left	S - XL
11-1345-x	30°, Right	S - XL
11-1346-x	30°, Left	S - XL
11-9100	15° Wedge (each)	S - XL
11-9101	15° Wedge (5-pack)	S - XL

Sizes Available (Measure elbow crease to base of index finger, fits right or left arm)

x = 2 = S..Up to 11" (28 cm) x = 3 = M..11 - 13" (28 - 33 cm) x = 4 = L..13"-15" (33-38 cm) x = 5 = XL..15" + (38 cm +)

DonJoy® Shoulder Stabilizer

IDEAL APPLICATION: Provides shoulder immobilization and controlled range of motion for glenohumeral dislocations/ subluxations, rotator cuff tears and acromioclavicular separations. Designed to protect and stabilize the shoulder post-injury and post-operatively.

Ordering information

Part Number	Description	Size
11-1101-x-06000	Right	S - XXXL
11-1102-x-06000	Left	S - XXXL
Accessories	Description	Size
11-1106-0-06000	Replacement strap set (5 straps)	Universal

Features & Benefits

- Provides abduction control, external rotation control and support for acromioclavicular separations
- The range-of-motion control can be customized with easy strap adjustments
- Latex-free black polyester Lycra® material
- Available in bilateral version upon request



Sizes Available (Measurement taken chest circumference under arms and across nipples)

x = 2 = S..... 30" - 34" (75-85 cm) x = 4 = L..... 38" - 41"..... (95-102 cm) x = 6 = XXL.... 44" - 48"..... (110-120 cm)
 x = 3 = M..... 34" - 38"..... (85-95 cm) x = 5 = XL..... 41" - 44"..... (102-110 cm) x = 7 = XXXL . 48" - 52"..... (120-130 cm)

DonJoy® Shoulder Stabilizer SPA

IDEAL APPLICATION: The SPA (Shoulder Pad Attachment) is indicated for use by those with glenohumeral instabilities in contact sports where shoulder pads are worn.

Ordering information

Part Number	Size
11-1319-x-06000	M - XXXL

Features & Benefits

- Cool – just an arm band, no chest vest
- Easy to use – arm cuff and straps attach directly to the shoulder pads (must have holes drilled) which remain attached when not in use. Player can attach it themselves



Sizes Available (Measurement taken at center of flexed bicep muscle)

x = 3 = M..... Up to 14½" (37 cm) x = 5 = XL..... 16 - 18"..... (40-46 cm) x = 7 = XXXL . 20" - 22"..... (51-56 cm)
 x = 4 = L..... 14½ - 16"..... (37.5-40.5 cm) x = 6 = XXL.... 18" - 20"..... (46-51 cm)

Sully® Shoulder Stabilizer

IDEAL APPLICATION: Provides shoulder immobilization and controlled range of motion for anterior, multi-directional, inferior, and posterior instabilities, rotator cuff deceleration, shoulder separations, and muscle strains. The shoulder brace is designed to stabilize, assist or restrict movement of the shoulder post-injury and post-operatively.

Ordering information

Part Number	Size	Chest	Biceps
11-0525-1	XSmall	34" - 40"	9" - 10.5"
11-0525-2	Small	36" - 40"	11" - 13"
11-0525-3	Medium	38" - 42"	13" - 15"
11-0525-4	Large	42" - 46"	14" - 17"
11-0525-5	XLarge	46" - 51"	15" - 17"
11-0527	Bilateral		
11-0526	Custom		

(Measure bicep in non-flex position)

Accessories	Description	Size
11-0528	Sully AC Pad	Universal

Features & Benefits

- Pulls the humerus posteriorly as it restricts external rotation and abduction
- Elastic straps attach with Velcro at any point, in any direction and with any amount of force
- The smooth, controlled assistance or restraint follows the natural movement of the muscles and joints



DonJoy® Custom Edition Elbow

IDEAL APPLICATION: DonJoy's premium and strongest bilateral upright elbow brace. The carbon fiber frame in conjunction with bilateral polycentric hinges provide support for medial/lateral ligament instabilities and hyperextension control. The hinges have adjustable flexion/extension stops to address individual range of motion protocols and restrictions.

- Features & Benefits**
- Custom made from a cast for the ultimate in fit and function
 - Low profile design for an easy fit under uniforms and clothing
 - 24 hour turnaround time
 - Includes patented supra condylar pad medially and pneumatic air Condyle Pad laterally for comfort and to prevent migration

Ordering information

Each Custom Edition Elbow is custom made. Cast must be 16" in length (8" above & below the elbow) arm should be in 15° flexion.

Minimum Diameter Measure from elbow joint
 Biceps - 12½" (6½" above)
 Forearm - 10½" (6½" below)

Accessories	Description
11-0543-x-06000	Neoprene Suspension Strap
25-1923-0-06000	Air Condyle Pad
11-0956-0-06000	Supra Condylar Pad Kit
25-0422-0-00000	Flexion / Extension Stop Kit

Sizes Available (Measurement taken at mid biceps)

x = 2 = S..... 10" - 11½".....(25.4 - 29.21 cm)	x = 4 = L..... 13.1" - 14½".....(33.27 - 36.83 cm)	x = 6 = XXL..... 16.1" - 17½".....(40.89 - 44.45 cm)
x = 3 = M..... 11½" - 13".....(29.46 - 33.02 cm)	x = 5 = XL..... 14½" - 16".....(37.08 - 40.64 cm)	x = 7 = XXXL .. 17½" - 19".....(44.70 - 48.26 cm)

DonJoy® Elbow Guard

IDEAL APPLICATION: Designed with bilateral polycentric hinges to provide superior elbow support for medial/lateral ligament instabilities, while the cross strapping design at the elbow works to control hyperextension.

- Features & Benefits**
- Wraparound biceps/triceps portion with elastic strap for ease of application
 - Bi-lateral polycentric hinges for flexion/extension and medial/lateral control
 - "X" strap for great hyperextension control
 - Elbow pad for protection in contact sports

Ordering information

Part Number	Size
11-1003-x-06000	S - XXXL

Sizes Available (Measurement taken at mid biceps)

x = 2 = S..... 10" - 11½".....(25.4 - 29.21 cm)	x = 4 = L..... 13.1" - 14½".....(33.27 - 36.83 cm)	x = 6 = XXL..... 16.1" - 17½".....(40.89 - 44.45 cm)
x = 3 = M..... 11.6" - 13".....(29.46 - 33.02 cm)	x = 5 = XL..... 14.6" - 16".....(37.08 - 40.64 cm)	x = 7 = XXXL .. 17.6" - 19".....(44.70 - 48.26 cm)

IROM® Elbow

IDEAL APPLICATION: The IROM Elbow is ideal for fixed or controlled ROM in the treatment of stable supracondylar fractures, stable radial fractures, and lateral/medial epicondylitis.

- Features & Benefits**
- 10° incremental setting (0° - 120°)
 - Allows for controlled flexion/extension
 - Integrated straps for full circumferential compression
 - Allows for supination and pronation of forearm and wrist

Ordering information

Part Number	Description	Size
11-0180-x-13066	Right	S - L
11-0181-x-13066	Left	S - L

Sizes Available (Measurements taken wrist to elbow crease)

x = 2 = S..... < 8" (20 cm)	x = 4 = L..... > 8" (20 cm)
-----------------------------	-----------------------------

Aircast Mayo Clinic Elbow Brace

IDEAL APPLICATION: Fracture dislocation; bicep/tricep tendon rupture; ulnar nerve transposition; total elbow arthroplasty; ligament repair; radial fixation; osteoarthritis; epicondylitis.

Ordering information

Part Number	Description	Size
05EL.....	Left	Universal
05ER.....	Right	Universal

Accessory	Description	Size
05EW.....	ARC Forearm Rotational Brace	Universal

- Features & Benefits**
- Unique patented design provides static adjustable stretch of the elbow in flexion and extension.
 - Replaces the need for two braces for elbow contractures
 - Eliminates the need for serial casting for elbow reconstructions and acute fracture dislocations.
 - Arc™ Forearm Rotational Brace may be used for additional pronation and supination control.
- U.S. Patent # 7,517,329



Elbow Support

IDEAL APPLICATION: Provides warmth and compression for minor injuries such as tendonitis, bursitis and arthritis. The elbow support with pad provides added elbow protection.

Ordering information

Part Number	Description	Size
11-0001-x-06000		XS - XXL
11-0401-x-06000	w/Pad	XS - XXL

- Features & Benefits**
- 4 mm neoprene
 - Available with oval padding for added protection

Sizes Available (Measurement taken at circumference of elbow joint)

x = 1 = XS..... 7" - 8"..... (17.75 - 20.25 cm)	x = 3 = M..... 9" - 10"..... (22.75 - 25.25 cm)	x = 5 = XL..... 11" - 12"..... (28 - 30.5 cm)
x = 2 = S..... 8" - 9"..... (20.25 - 22.75 cm)	x = 4 = L..... 10" - 11"..... (25.25 - 28 cm)	x = 6 = XXL..... 12" - 13"..... (30.5 - 33 cm)



Spider® Pad Elbow

IDEAL APPLICATION: The unique honeycomb design contours for comfortable movement while providing cushioning protection to the elbow. Breathable Lycra® material is ideal for basketball, football and volleyball.

Ordering information

Part Number	Size
11-0392-x-06060	XS - XXXL

- Features & Benefits**
- High impact-resistant foam compared to conventional elbow pads
 - Low-profile and streamlined for protection during contact sports and job-related activities

Sizes Available (Measurement taken at circumference of elbow joint)

x = 1 = XS..... 7" - 8"..... (17.75 - 20.25 cm)	x = 3 = M..... 9" - 10"..... (22.75 - 25.25 cm)	x = 5 = XL..... 11" - 12"..... (28 - 30.5 cm)
x = 2 = S..... 8" - 9"..... (20.25 - 22.75 cm)	x = 4 = L..... 10" - 11"..... (25.25 - 28 cm)	x = 6 = XXL..... 12" - 13"..... (30.5 - 33 cm)
		x = 7 = XXXL .. 13" - 14"..... (33 - 35.75 cm)



DonJoy Functional Wrist Brace NEW!



IDEAL APPLICATION: The DonJoy Functional Wrist Brace helps restore function and structural characteristics to an injured wrist. The brace offers variable extension stops with the option to control radial/ulnar deviation. The open palm design allows for gripping motion for a return to work or sport. Available in Carbon Fiber (CXT) or a high strength composite material (SXT).

Features & Benefits

- Variable extension control (60°, 40°, 20°, 0°) to avoid excessive extension of the wrist
- Rigid stability, provides ideal support
- Open palm design, allows gripping motion
- Lateral motion control, optional radial ulnar stop
- Low profile fit can be worn with sporting gear
- Varying styles to support activity level. Lacer, Sport Wrap, or Speed Strap
- Cool, comfortable fit, easy to wear at work or sport

Ordering information

Part Number	Description	Size
11-2130-x	CXT Wrist Brace, Lacer	S - L
11-2131-x	CXT Wrist Brace, Speed Strap	S - L
11-2132-x	CXT Wrist Brace, Sport Strap	S - L
11-2135-x	SXT Wrist Brace, Lacer	S - L
11-2140-x	SXT Wrist Brace, Speed Strap	S - L
11-2141-x	SXT Wrist Brace, Sport Strap	S - L

Aircast A2™ Wrist Brace



IDEAL APPLICATION: Wrist ligament instability, sprain or muscle strain, control of Carpal Tunnel syndrome symptoms, post-op and post removal of casting.

Features & Benefits

- Durable, lightweight nylon foam construction for maximum patient comfort
- Malleable volar stay for customizable fit
- Contact closure straps for ease of application

Ordering information

Part Number	Description	Size
05WSL	Left	Small (5" - 6½")
05WSR	Right	Small (5" - 6½")
05WML	Left	Medium (6½" - 7¾")
05WMR	Right	Medium (6½" - 7¾")
05WLL	Left	Large (7½" - 9")
05WLR	Right	Large (7½" - 9")

Aircast A2™ Wrist Brace w/ Thumb Spica



IDEAL APPLICATION: Ligament instability, thumb injuries, sprains or muscle strain, GameKeep's Thumb, and deQuervain Syndrome symptoms.

Features & Benefits

- Durable, lightweight nylon foam construction for maximum patient comfort
- Malleable volar and thumb stay for customizable fit
- Contact closure straps for ease of application

Ordering information

Part Number	Description	Size
05WTSL	Left	Small (5" - 6½")
05WTSR	Right	Small (5" - 6½")
05WTML	Left	Medium (6½" - 7¾")
05WTMR	Right	Medium (6½" - 7¾")
05WTLL	Left	Large (7½" - 9")
05WTLR	Right	Large (7½" - 9")

Aircast Pneumatic Armband

IDEAL APPLICATION: Tennis elbow; Golfers elbow (epicondylitis)

Ordering information

Part Number	Description	Size
05A	Beige	8" - 14"
05A-B	Black	8" - 14"

Features & Benefits

- Single, pre-inflated aircell, concentrates compression on the extensor muscle, providing more support, less constriction
- Universal fit
- Available in beige and black



Aircast StabilAir™ Wrist Brace

IDEAL APPLICATION: Post-operative treatment for distal radius fractures.

Ordering information

Part Number	Description	Size
05FSL	Left	Small (5" - 6")
05FSR	Right	Small (5" - 6")
05FML	Left	Medium (6½" - 7¾")
05FMR	Right	Medium (6½" - 7¾")

Features & Benefits

- Adjustable Anatomic shells and adjustable aircells conform to arm for a custom fit
- Lightweight, water resistant brace
- Easy-to-use and can be applied quickly and easily
- Removable padded shells enable easy inspection or treatment.



Cold Therapy



IceMan

IceMan®



IDEAL APPLICATION: The IceMan cold therapy unit helps reduce pain and swelling and facilitates rehabilitation. It provides several hours of continuous cold therapy for a variety of conditions. Using a patented semi-closed loop system, the IceMan maintains a more constant and accurate temperature than competitor's systems.

Features & Benefits

- Two-position locking top
- Water tight seal
- Hose stirrup
- Inlet temperature control and monitoring
- Compatibility with a variety of pads

Ordering information

Part Number	Description
11-0494	IceMan Unit Model 1100
11-0657-9-00000	Bed Hanger
11-0674-9-02000	Dual Pad Connector
11-1421	IceMan w/ Universal Sterile, RH
11-1422	IceMan w/ Universal Loop, NS, RH
11-1423	IceMan w/ Universal, NS, RH
11-1424	IceMan w/ Universal Loop, NS, EH
11-1425	IceMan w/ Universal Sterile, EH
11-1426	IceMan w/ Rectangle, Sterile, EH, XL
11-9098	IceMan w/ Cold Pad, Shoulder
11-9099	IceMan w/ Cold Pad, McGuire
25-4882	IceMan Power Supply

Cold Pads



IDEAL APPLICATION: DonJoy offers a variety of pad sizes to provide cold therapy to various areas of the body.

Features & Benefits

- Even distribution throughout the pad
- Sterile and non-sterile pads

Ordering information

Part Number	Description (Sterile)
11-0813-9-00000	EH, Rectangle, XL (9.76" x 13.75")
11-0801-9-00000	RH, Universal (10.85" x 11.34")
11-0802-9-00000	EH, Universal (10.85" x 11.34")
11-0804-9-00000	EH, Universal, XL (12" x 11.34")

Part Number	Description (Non-Sterile)
11-0675-9-00000	RH, Universal (10.85" x 11.34")
11-0676-9-00000	EH, Universal (10.85" x 11.34")
11-0678-9-00000	EH, Universal, XL (12" x 11.34")



Large Rectangle

Small Rectangle

Hand/Wrist

Ankle

Wrap-On Pads

IDEAL APPLICATION: The wrap-on pads are designed with a hook engageable material to easily secure around the affected area. The universal size is ideal for the shoulder, knee and ankle. The wrap-on pads are not available in a sterile configuration.

Ordering information

Part Number	Description
11-0679-9-00000	RH Universal (10.85" x 11.34")
11-0680-9-00000	EH Universal (10.85" x 11.34")
11-0681-9-00000	RH, XL Universal (12" x 11.34")
11-0682-9-00000	EH, XL Universal (12" x 11.34")
11-1268-9-00000	McGuire Knee Wrap w/open patella (10.85" x 11.34")
11-1334-9-00000	Shoulder (14" x 23" x 12.5")
11-1397-9-00000	Shoulder, small (12.75" x 19" x 12.25")
11-1522	RH, Ankle
11-1521	RH, Forefoot
11-9111	Hip, RT
11-9112	Hip, LT

Note: RH = Regular Hose 6.5" EH = Extended Hose 15"

Features & Benefits

- Includes both cold pad and wrap



Universal



Shoulder

Hip

McGuire Knee

Forefoot

Ankle

Sterile Dressings

IDEAL APPLICATION: The sterile dressings provide a sterile moisture barrier between the skin and the cold therapy pad. The dressings are available in five sizes to use together with the various IceMan cold pads, including the wrap-on pads.

Ordering information

Part Number	Description
11-0820-9-00000	Rectangle, XLarge
11-0821-9-00000	Ankle
11-0822-9-00000	Universal
11-0823-9-00000	Universal, XLarge
11-1360-9-00000	Shoulder, XLarge

Features & Benefits

- Moisture barrier
- Five sizes
- Compatible with IceMan cold pads



Universal



Rectangle

Shoulder

Ankle

Aircast Cryo/Cuff™ IC



IDEAL APPLICATION: Intermittent cold and compression is clinically proven to help reduce post-operative swelling and regain range-of-motion. The Cryo/Cuff IC provides both gravity and motorized system for flexibility of care.

Features & Benefits

- An integrated pneumatic pump within cooler lid provides automated compression and cold therapy
- Simplicity of design and ease of operation makes it ideal for post-operative recovery, trauma, athletic training rooms and home use

Ordering information

Part Number	Description
51A	Cryo/Cuff IC Cooler
51A10A	Cryo/Cuff IC Cooler w/ Ankle Cryo/Cuff
51A11A	Cryo/Cuff IC Cooler w/ Knee Cryo/Cuff, Medium
51A11B	Cryo/Cuff IC Cooler w/ Knee Cryo/Cuff, Large
51A12A	Cryo/Cuff IC Cooler w/ Shoulder
51A14A	Cryo/Cuff IC Cooler w/ Back, Hip and Rib

Accessories

25-4882	Cryo/Cuff IC Cooler power supply
25-0238	Cryo/Cuff IC Lid

Aircast Cryo/Cuff™ Cooler



IDEAL APPLICATION: Cryo/Cuff utilizes gravity to deliver chilled water to anatomically specific cuffs.

Features & Benefits

- Cooler holds water and ice needed for 6 to 8 hours of cryotherapy
- Compatible with all Cryo/Cuff and Arctic Flow pads. Includes tube assembly and insulation disk.

Ordering information

Part Number	Description
2125	Cooler and Tube Assembly

Aircast Cryo/Cuff™

IDEAL APPLICATION: Anatomically designed cuff provides specific joint coverage for optimum treatment of swelling and pain.

Ordering information

Part Number	Description	Circumference*	Size
Shoulder Cryo/Cuff			
12A	Cuff and Cooler	32" – 48"	
12A01	Cuff Only	32" – 48"	
1220	Strap Only	32" – 48"	
12AXL	Cuff w/X-Long	42" – 54"	Strap and Cooler
12AXL01	Cuff Only	42" – 54"	w/X-Long Strap
1220XL	X-Long Strap Only	42" – 54"	

Back/Hip/Rib Cryo/Cuff

14A	Cuff and Cooler
14A01	Cuff Only

Knee Cryo/Cuff

11C	Cuff and Cooler	10" – 19"	S
11C01	Cuff Only	10" – 19"	S
11A	Cuff and Cooler	18" – 23"	M
11A01	Cuff Only	18" – 23"	M
11B	Cuff and Cooler	20" – 31"	L
11B01	Cuff Only	20" – 31"	L

* Measure 6" above patella

Ankle Cryo/Cuff

10A	Cuff & cooler
10A01	Cuff only

Foot Cryo/Cuff

10C	Cuff and Cooler	9" – 13"	M
10C01	Cuff Only	9" – 13"	M
10B	Cuff and Cooler	10" – 17"	L
10B01	Cuff Only	10" – 17"	L

* Measure forefoot circumference

Features & Benefits

- Available for Foot, Ankle, Knee, Shoulder and Back



Back/Hip/Rib Cryo/Cuff



Shoulder Cryo/Cuff



Knee Cryo/Cuff



Ankle Cryo/Cuff



Foot Cryo/Cuff

Aircast Knee Cryo/Cuff™ SC

IDEAL APPLICATION: A totally self-contained cuff that is filled directly with water and ice — no cooler required.

Ordering information

Part Number	Description
11ASC	Cuff (self contained)

Features & Benefits

- Cuff remains cold for about one hour.
- A hand bulb is included for cuff inflation and to add compression as needed
- Excellent fit for OA chronic pain treatment



Dura*Soft®



Shoulder Wrap

IDEAL APPLICATION: The universal design provides cold therapy to various areas of the body. The Dura*Soft products can effectively penetrate surgical dressings to decrease pain, edema and secondary hypoxic tissue injury.

Features & Benefits

- Inserts can be frozen and exchanged to provide extended cold therapy
- The nylon wraps and elastic straps provide a universal fit and effective compression
- Ultra-wick fabric with antimicrobial protection
- May be used repeatedly to provide years of long-term use for rehab and/or chronic conditions

Ordering information

Part Number	Description
Dura*Soft® Shoulder Wrap & Back Wrap	
11-0914-0-02000	Shoulder Wrap with 1 Ice Insert (12" x 13")
11-0915-0-02000	Shoulder Wrap with 2 Ice Inserts (12" x 13")
11-0916-0-02000	Back Wrap with 1 Ice Insert
11-2124-0-02000	Back Wrap with 2 Ice Inserts
Dura*Soft® Knee Sleeve & Knee Wrap	
11-0909-0-02000	Knee Sleeve with 2 Ice Inserts
11-0910-0-02000	Knee Sleeve with 4 Ice Inserts
11-0912-0-00000	Knee Wrap with 1 Ice Insert
11-0913-0-02000	Knee Wrap with 2 Ice Inserts (11" x 12")
Dura*Soft® Surgical Foot/Ankle	
11-0911-0-02000	Foot / Ankle Wrap with 2 Ice Inserts
11-2125-0-02000	Foot / Ankle Wrap with 4 Ice Inserts
Dura*Soft® Universal Wrap	
11-0918-0-00000	Universal Wrap with 1 Ice Insert
11-2123-0-02000	Universal Wrap with 2 Ice Inserts



Back Wrap



Knee Sleeve



Knee Wrap



Surgical Foot/Ankle



Universal Wrap



Nate Holland, US Snowboard Team

Dura*Kold®



Shoulder/Hip Wrap

IDEAL APPLICATION: The Dura*Kold products can effectively penetrate surgical dressings to decrease pain, edema and secondary hypoxic tissue injury.

Ordering information

Part Number	Description
Dura*Kold® Shoulder Wrap & Shoulder/Hip Wrap	
11-0885-9-02000	Shoulder Wrap, Regular
11-0886-9-02000	Shoulder Wrap w/Rotator Cuff
11-0898-2-02000	Shoulder / Hip Wrap, Small
Dura*Kold® Back Wrap	
11-0887-9-02000	Back Wrap (20" x 45" waist)
Dura*Kold® Wrist/Elbow Wrap	
11-0889-9-02000	Wrist / Elbow Wrap (12" x 15")
Dura*Kold® Leg & Arm Wrap	
11-0899-5-02000	Leg & Arm Wrap – XLarge (13" x 26")
11-1336-0-02000	Leg & Arm 4 Panel (13" x 18")
Dura*Kold® Surgical Knee Sleeve	
11-0904-0-02000	Surgical Knee Sleeve, 2 Ice Inserts (12" x 13")
11-0905-0-02000	Surgical Knee Sleeve, 4 Ice Inserts (12" x 13")
Dura*Kold® Arthroscopy Knee Wrap	
11-0896-9-02000	Arthroscopic Knee Wrap, Standard (13" x 15")
11-0897-4-02000	Arthroscopic Knee Wrap, Large (13" x 23")
11-0897-5-02000	Arthroscopic Knee Wrap, XLarge (13" x 28")
Dura*Kold® CPM Knee Wrap	
11-0895-9-02000	Universal
Dura*Kold® Foot & Ankle Wrap	
11-0900-9-02000	Standard
Dura*Kold® Surgical Foot Wrap	
11-0901-0-02000	Surgical Foot Wrap
Dura*Kold® Consumer Wrap	
11-0902-2-02000	Consumer Wrap (5" x 15")
11-0902-4-02000	Consumer Wrap (8" x 22")

Features & Benefits

- Provides 2 hours of cold therapy
- Removable inserts can be exchanged without removing the wrap to provide long-term cold therapy
- Keeps patients comfortable and dry while allowing mobility



Shoulder Wrap



Back Wrap



Wrist/Elbow Wrap



Leg & Arm Wrap



Surgical Knee Sleeve



Arthroscopy Knee Wrap



CPM Knee Wrap



Foot & Ankle Wrap



Surgical Foot Wrap



Consumer Wrap

Part No.	Description	Page No.	Part No.	Description	Page No.
Ankle Braces					
02AL	Aircast Air-Stirrup® Ankle Brace, Left, Large	28	11-0823-9-00000	Sterile Dressings - Universal, XLarge	40
02AR	Aircast Air-Stirrup® Ankle Brace, Right, Large	28	11-0885-9-02000	Dura*Kold® Shoulder Wrap & Shoulder/Hip Wrap, Shoulder Wrap, Regular	44
02BL	Aircast Air-Stirrup® Ankle Brace, Left, Medium	28	11-0886-9-02000	Dura*Kold® Shoulder Wrap & Shoulder/Hip Wrap, Shoulder Wrap w/Rotator Cuff	44
02BR	Aircast Air-Stirrup® Ankle Brace, Right, Medium	28	11-0887-9-02000	Dura*Kold® Back Wrap	44
02CL	Aircast Air-Stirrup® Ankle Brace, Left, Small	28	11-0889-9-02000	Dura*Kold® Wrist/Elbow Wrap	44
02CR	Aircast Air-Stirrup® Ankle Brace, Right, Small	28	11-0895-9-02000	Dura*Kold® CPM Knee Wrap	44
02E12	Aircast Air-Stirrup Universe™, 12 PK	28	11-0896-9-02000	Dura*Kold® Arthroscopy Knee Wrap, Arthroscopic Knee Wrap, Standard (13" x 15")	44
02E	Aircast Air-Stirrup Universe™	28	11-0897-4-02000	Dura*Kold® Arthroscopy Knee Wrap, Arthroscopic Knee Wrap, Large (13" x 23")	44
02EK	Aircast Air-Stirrup® Universe™ Care Kit	28	11-0897-5-02000	Dura*Kold® Arthroscopy Knee Wrap, Arthroscopic Knee Wrap, XLarge (13" x 28")	44
02JL	Aircast Air-Stirrup® Ankle Brace, Left, Pediatric	28	11-0898-2-02000	Dura*Kold® Shoulder Wrap & Shoulder/Hip Wrap, Shoulder / Hip Wrap, Small	44
02JR	Aircast Air-Stirrup® Ankle Brace, Right, Pediatric	28	11-0899-5-02000	Dura*Kold® Leg & Arm Wrap, Leg & Arm Wrap – XLarge (13" x 26")	44
02MLL	Aircast AirSport™, Left, L	29	11-0900-9-02000	Dura*Kold® Foot & Ankle Wrap	44
02MLR	Aircast AirSport™, Right, L	29	11-0901-0-02000	Dura*Kold® Surgical Foot Wrap	44
02MML	Aircast AirSport™, Left, M	29	11-0902-2-02000	Dura*Kold® Consumer Wrap, Consumer Wrap (5" x 15")	44
02MMR	Aircast AirSport™, Right, M	29	11-0902-4-02000	Dura*Kold® Consumer Wrap, Consumer Wrap (8" x 22")	44
02MSL	Aircast AirSport™, Left, S	29	11-0904-0-02000	Dura*Kold® Surgical Knee Sleeve, Surgical Knee Sleeve, 2 Ice Inserts (12" x 13")	44
02MSR	Aircast AirSport™, Right, S	29	11-0905-0-02000	Dura*Kold® Surgical Knee Sleeve, Surgical Knee Sleeve, 4 Ice Inserts (12" x 13")	44
02MXLL	Aircast AirSport™, Left, XL	29	11-0909-0-02000	Dura*Soft® Knee Sleeve & Knee Wrap, Knee Sleeve with 2 Ice Inserts	43
02MXLR	Aircast AirSport™, Right, XL	29	11-0910-0-02000	Dura*Soft® Knee Sleeve & Knee Wrap, Knee Sleeve with 4 Ice Inserts	43
02MXSL	Aircast AirSport™, Left, XS	29	11-0911-0-02000	Dura*Soft® Surgical Foot/Ankle, Foot / Ankle Wrap with 2 Ice Inserts	43
02MXSR	Aircast AirSport™, Right, XS	29	11-0912-0-00000	Dura*Soft® Knee Sleeve & Knee Wrap, Knee Wrap with 1 Ice Insert	43
02TLL	Aircast A60™ Ankle Support, Left, L	29	11-0913-0-02000	Dura*Soft® Knee Sleeve & Knee Wrap, Knee Wrap with 2 Ice Inserts (11" x 12")	43
02TLR	Aircast A60™ Ankle Support, Right, L	29	11-0914-0-02000	Dura*Soft® Shoulder Wrap & Back Wrap, Shoulder Wrap with 1 Ice Insert (12" x 13")	43
02TML	Aircast A60™ Ankle Support, Left, M	29	11-0915-0-02000	Dura*Soft® Shoulder Wrap & Back Wrap, Shoulder Wrap with 2 Ice Inserts (12" x 13")	43
02TMR	Aircast A60™ Ankle Support, Right, M	29	11-0916-0-02000	Dura*Soft® Shoulder Wrap & Back Wrap, Back Wrap with 1 Ice Insert	43
02TSL	Aircast A60™ Ankle Support, Left, S	29	11-0918-0-02000	Dura*Soft® Universal Wrap, Universal Wrap with 1 Ice Insert	43
02TSR	Aircast A60™ Ankle Support, Right, S	29	11-1268-9-00000	Wrap-On Pads - McGuire Knee Wrap w/open patella	40
09AS	AirHeel Support, S	30	11-1334-9-00000	Wrap-On Pads - Shoulder	40
09AM	AirHeel Support, M	30	11-1336-0-02000	Dura*Kold® Leg & Arm Wrap, Leg & Arm 4 Panel (13" x 18")	44
09AL	AirHeel Support, L	30	11-1360-9-00000	Sterile Dressings - Shoulder, XLarge	40
09DS	Dorsal Night Splint, S - L	30	11-1397-9-00000	Wrap-On Pads - Shoulder, small	39
09DL	Dorsal Night Splint, L - XL	30	11-1421	IceMan w/ Universal Sterile, RH	39
09PKS	AirHeel/DNS Care Kit, S	30	11-1422	IceMan w/ Universal Loop, NS, RH	39
09PKM	AirHeel/DNS Care Kit, M	30	11-1423	IceMan w/ Universal, NS, RH	39
09PKL	AirHeel/DNS Care Kit, L	30	11-1424	IceMan w/ Universal Loop, NS, EH	39
11-0451-x-06000	Stabilizing Ankle, Black	29	11-1425	IceMan w/ Universal Sterile, EH	39
11-0451-x-15000	Stabilizing Ankle, White	29	11-1426	IceMan w/ Rectangle, Sterile, EH, XL	39
11-1489-x-06000	Velocity™ LS - Standard, Left	27	11-1521	Wrap-On Pads - RH, Forefoot	40
11-1490-x-06000	Velocity™ LS - Standard, Right	27	11-1522	Wrap-On Pads - RH, Ankle	40
11-1491-x-06000	Velocity™ LS - Wide, Left	27	11-2123-0-02000	Dura*Soft® Universal Wrap, Universal Wrap with 2 Ice Inserts	43
11-1492-x-06000	Velocity™ LS - Wide, Right	27	11-2124-0-02000	Dura*Soft® Shoulder Wrap & Back Wrap, Back Wrap with 2 Ice Inserts	43
11-1493-x-06000	Velocity™ MS - Standard, Left	27	11-2125-0-02000	Dura*Soft® Surgical Foot/Ankle, Foot / Ankle Wrap with 4 Ice Inserts	43
11-1494-x-06000	Velocity™ MS - Standard, Right	27	11-9098	IceMan w/ Cold Pad, Shoulder	39
11-1495-x-06000	Velocity™ MS - Wide, Left	27	11-9099	IceMan w/ Cold Pad, McGuire	39
11-1496-x-06000	Velocity™ MS - Wide, Right	27	11-9111	Wrap-On Pads - Hip, RT	40
11-1497-x-xx000	Velocity™ ES - Standard, Left	27	11-9112	Wrap-On Pads - Hip, LT	40
11-1498-x-xx000	Velocity™ ES - Standard, Right	27	11A01	Knee Cryo/Cuff, Cuff Only, 18" – 23", M	42
11-1499-x-xx000	Velocity™ ES - Wide, Left	27	11A	Knee Cryo/Cuff, Cuff and Cooler, 18" – 23", M	42
11-1500-x-xx000	Velocity™ ES - Wide, Right	27	11ASC	Aircast Knee Cryo/Cuff SC	42
Cold Therapy					
10A01	Ankle Cryo/Cuff, Cuff only	42	11B01	Knee Cryo/Cuff, Cuff Only, 20" – 31", L	42
10A	Ankle Cryo/Cuff, Cuff & cooler	42	11B	Knee Cryo/Cuff, Cuff and Cooler, 20" – 31", L	42
10B01	Foot Cryo/Cuff, Cuff Only, 10" – 17", L	42	11C	Knee Cryo/Cuff, Cuff and Cooler, 10" – 19", S	42
10B	Foot Cryo/Cuff, Cuff and Cooler, 10" – 17", L	42	12A01	Shoulder Cryo/Cuff, Cuff Only, 32" – 48"	42
10C01	Foot Cryo/Cuff, Cuff Only, 9" – 13", M	42	12A	Shoulder Cryo/Cuff, Cuff and Cooler, 32" – 48"	42
10C	Foot Cryo/Cuff, Cuff and Cooler, 9" – 13", M	42	12AXL01	Shoulder Cryo/Cuff, Cuff Only, 42" – 54" w/X-Long Strap	42
11-0494	IceMan Unit Model 1100	39	12AXL	Shoulder Cryo/Cuff, Cuff w/X-Long, 42" – 54" Strap and Cooler	42
11-0657-9-00000	IceMan® Bed Hanger	39	14A01	Back/Hip/Rib Cryo/Cuff, Cuff Only	42
11-0674-9-02000	IceMan® Dual Pad Connector	39	14A	Back/Hip/Rib Cryo/Cuff, Cuff and Cooler	42
11-0675-9-00000	Cold Pads (Non-Sterile) - RH, Universal	39	25-0238	Cryo/Cuff IC Lid	41
11-0676-9-00000	Cold Pads (Non-Sterile) - EH, Universal	39	25-4882	Cryo/Cuff IC Cooler power supply	41
11-0678-9-00000	Cold Pads (Non-Sterile) - EH, Universal	39	25-4882	IceMan Power Supply	39
11-0679-9-00000	Wrap-On Pads - RH Universal	40	51A10A	Cryo/Cuff IC Cooler w/ Ankle Cryo/Cuff	41
11-0680-9-00000	Wrap-On Pads - EH Universal	40	51A11A	Cryo/Cuff IC Cooler w/ Knee Cryo/Cuff, Medium	41
11-0681-9-00000	Wrap-On Pads - RH, XL Universal	40	51A11B	Cryo/Cuff IC Cooler w/ Knee Cryo/Cuff, Large	41
11-0682-9-00000	Wrap-On Pads - EH, XL Universal	40	51A14A	Cryo/Cuff IC Cooler w/ Shoulder	41
11-0801-9-00000	RH, Universal (10.85" x 11.34")	39	51A	Cryo/Cuff IC Cooler	41
11-0802-9-00000	EH, Universal (10.85" x 11.34")	39			
11-0804-9-00000	EH, Universal, XL (12" x 11.34")	39			
11-0813-9-00000	EH, Rectangle, XL (9.76" x 13.75")	39			
11-0820-9-00000	Sterile Dressing, Rectangle	40			
11-0821-9-00000	Sterile Dressings - Ankle	40			
11-0822-9-00000	Sterile Dressings - Universal	40			

Part No.	Description	Page No.	Part No.	Description	Page No.
1220	Shoulder Cryo/Cuff , Strap Only, 32" – 48"	42	11-2011-x	Hinged Knee - Standard, Sleeve w/popliteal cutout & horseshoe buttress, Drytex®	9
1220XL	Shoulder Cryo/Cuff , X-Long Strap Only, 42" – 54"	42	11-2012-x	Hinged Knee - Standard, Wraparound, Drytex®	9
2125	Aircast Cryo/Cuff™ Cooler w/Tube Assembly	41	11-2013-x	Hinged Knee - Standard, Wraparound w/popliteal cutout, Drytex®	9
Knee Ligament Bracing-Custom					
11-0016-x-06000	Defiance Sports Cover	4, 5	11-3220-x	ACL FourcePoint Short Right	8
11-0025x	Replacement Tape Module	5	11-3221-x	ACL FourcePoint Short Left	8
11-0075-x-00000	Defiance Black Lycra Undergarment	4, 5	11-3222-x	CI FourcePoint Short Right	8
11-0103-x-06000	Defiance Sports Cover (short)	4, 5	11-3223-x	CI FourcePoint Short Left	8
11-0122-x	Defiance Neoprene Under-sleeve	4, 5	25-0422-0-00000	Flexion / Extension Stop Kit	7, 9
11-0183	CCMI Mark III - Custom Measurement Device	5	25-2428	Mesh Brace Bag	7, 8
11-0398-9-00000	Defiance Full Extension Stop	4, 5	Lower Back Bracing		
11-0400-9-00000	Defiance Pneumatic Condyle Pad	4, 5	11-1561-x	Range of Spinal Support T2/T3-S1	31
11-0543-x-06000	Defiance Neoprene Suspension Strap	4, 5	11-1562-x	DonJoy LSO with Chairback (8")	32
11-1006-x-06000	Defiance Impact Guard	4, 5	11-1563-x	Standard Chairback	32
11-1367-9-06000	Defiance Calf Pinch Guard	4, 5	11-1564-x	Short Chairback*	32
11-1415	Defiance Sili-Grip Strap Pads	4, 5	11-1565	7" Lordotic Insert (pair)	32
11-1525-x	Defiance Tru-Pull Advanced Attachment	4, 5	11-1566	10" Lordotic Insert (pair)	32
25-0422-0-00000	Defiance Flexion / Extension Stop Kit	4, 5	11-1567	12" Lordotic Insert (pair)	32
25-1480-0-00000	Defiance Brace Bag	4, 5	11-1568-x	DonJoy Regular Back (10")	32
25-1923-0-06000	Defiance Air Condyle Kit	4, 5	11-1569-x	DonJoy Low Profile Back (8")	32
			11-1570-x	DonJoy Lite Back (9")	32
Knee Ligament Bracing-Prefabricated					
5-0442-0-00000	Flexion / Extension Stop Kit	8	Osteoarthritis Bracing - Custom		
11-0016-x-06000	Sports Cover	7	11-0016-x-06000	Brace Cover (standard)	12
11-0016-x-06000	Sports Cover (standard)	8	11-0075-x-00000	Lycra Undergarment/Oversleeve	12
11-0075-x-00000	Black Lycra Undergarment	7, 8	11-0103-x-06000	Brace Cover (short)	12
11-0122-x	Black Neoprene Undergarment	7, 8	11-0122-x	Neoprene Undergarment	12
11-0258-x	ACL FourcePoint Standard Right	8	11-0398-9-00000	Full Extension Stop	12
11-0259-x	ACL FourcePoint Standard Left	8	11-0400-9-00000	Pneumatic Condyle Pad	12
11-0264-x	CI Fource Point Standard Right	8	11-0543-x-06000	Neoprene Suspension Strap	12
11-0265-x	CI FourcePoint Standard Left	8	11-0857-9-00000	Load Sensor Kit	12
11-0398-9-00000	Full Extension Stop	7, 8	11-1006-x-06000	Impact Guard	12
11-0543-x-06000	Neoprene Suspension Strap	7, 8	11-1095-9-06000	Calf Pinch Guard	12
11-0555-x	Hinged Knee - Standard, Sleeve, Drytex®	9	11-1107-0-00000	Floam Condyle Pads	12
11-0556-x	Hinged Knee - Standard, Sleeve w/popliteal cutout, Drytex®	9	11-1415	Sili-Grip Strap Pads	12
11-0557-x	Playmaker® - Standard, Sleeve, Drytex®	9	11-1525-x	Tru-Pull® Advanced Patellofemoral Attachment	12
11-0558-x	Playmaker® - Standard, Sleeve w/popliteal cutout, Drytex®	9	25-1480-0-00000	Brace Bag	12
11-0666-x	Playmaker® - Standard, Wraparound, Drytex®	9	25-1923-0-06000	Air Condyle Pads	12
11-0759-x	Playmaker® - Standard, Wraparound w/popliteal cutout, Drytex®	9	Osteoarthritis Bracing - Prefabricated		
11-0860-x	Playmaker® - Standard, Sleeve, Neoprene	9	11-0075-x-00000	Lycra® Undergarment/Oversleeve	13
11-0861-x	Playmaker® - Standard, Sleeve w/popliteal cutout, Neoprene	9	11-0122-x	Neoprene Undergarment	13, 14
11-0862-x	Playmaker® - Standard, Sleeve w/patella donut, Neoprene	9	11-0543-x-06000	Neoprene Strap Kit	13
11-0863-x	Playmaker® - Standard, Sleeve w/popliteal cutout & patella donut, Neoprene	9	11-0543-x-06000	Neoprene Suspension Strap	13
11-0864-x	Playmaker® - Standard, Wraparound, Neoprene	9	11-0857-9-00000	Load Sensor Kit	13
11-0865-x	Playmaker® - Standard, Wraparound w/popliteal cutout, Neoprene	9	11-0870-x	OA Adjuster™ (patient-ready) - Medial, Right	13
11-0956-0-06000	Supra Condyle Pad	8	11-0871-x	OA Adjuster™ (patient-ready) - Medial, Left	13
11-0956-0-06000	Supra Condylar Pad Kit	7	11-0872-x	OA Adjuster™ (patient-ready) - Lateral, Right	13
11-1004-x-26000	Armor - ACL, Standard, Short, Right	7	11-0873-x	OA Adjuster™ (patient-ready) - Lateral, Left	13
11-1005-x-26000	Armor - ACL, Standard, Short, Left	7	11-0947-x	OA Assist Medial Right	13
11-1006-x-06000	Patella / Shin Impact Guard	7, 8	11-0948-x	OA Assist Lateral Right	13
11-1025-x-26000	Armor - PCL, Standard, Short, Right	7	11-0949-x	OA Assist Medial Left	13
11-1029-x-26000	Armor - CI, Standard, Short, Right	7	11-0950-x	OA Assist Lateral Left	13
11-1030-x-26000	Armor - CI, Standard, Short, Left	7	11-0973-x	Strap Replacement Kit	13
11-1095-9-06000	Calf Pinch Guard	7, 8	11-0975-x	Liner Replacement Kit (Right Medial/Left Lateral)	13
11-1395-x	Playmaker® - Standard, Sleeve w/popliteal cutout & patella donut, Drytex®	9	11-0976-x	Liner Replacement Kit (Left Medial/Right Lateral)	13
11-1396-x	Playmaker® - Standard, Sleeve w/patella donut, Drytex®	9	11-1100-x-06000	Strap Replacement Kit	13
11-1415	Sili-Grip Strap Pads	7, 8	11-1338-x	OA Everyday™ (patient-ready) - Medial, Right	14
11-1440-x	Armor - ACL, FourcePoint, Short, Right	7	11-1339-x	OA Everyday™ (patient-ready) - Lateral, Right	14
11-1441-x	Armor - ACL, FourcePoint, Short, Left	7	11-1340-x	OA Everyday™ (patient-ready) - Medial, Left	14
11-1442-x	Armor - ACL, FourcePoint, Standard, Right	7	11-1341-x	OA Everyday™ (patient-ready) - Lateral, Left	14
11-1443-x	Armor - ACL, FourcePoint, Standard, Left	7	11-1414-x	Strap Replacement Kit	14
11-1444-x	Armor - ACL, Standard, Standard, Right	7	11-1415	Sili-Grip Strap Pads	12, 13, 14
11-1445-x	Armor - ACL, Standard, Standard, Left	7	11-1415-x	Sili-Grip Strap Pad Replacement Kit	13
11-1447-x	Armor - CI, Standard, Standard, Right	7	11-1416-x	Liner Replacement Kit (Right Medial/Left Lateral)	14
11-1448-x	Armor - CI, Standard, Standard, Left	7	11-1417-x	Liner Replacement Kit (Left Medial/Right Lateral)	14
11-1524	FourcePoint Hinge Kit	9	11-1420	2 in 1 Strap Pad Kit	13
11-1525-x	Tru-Pull Advance Attachment	8	11-1531-x	OA Lite™ (patient-ready) - Medial, Right	14
11-1525-x	Tru-Pull Advanced Attachment	7	11-1532-x	OA Lite™ (patient-ready) - Medial, Left	14

Part No.	Description	Page No.
11-1550-x	Suede Liner Replacement Kit (Lateral, Left)	13
11-1552-x	Suede Strap Pad Replacement Kit	13
11-1553	Suede Condyle Pad Replacement Kit	13
25-0422-0-00000	Flexion/Extension Stop Kit	14
25-1376-0-00000	Flexion/Extension Stop Kit	13
25-2436-x	Strap Replacement Kit	14
25-2438	Padding Kit	14

Patellofemoral Bracing

11-0032-x-06060	"H" Buttress - Right or Left Sleeve, stays	19
11-0215-x	REACTION Knee Brace	17
11-0541-x	Replacement Under Sleeve	17
11-0242-x	Elastic Knee - Right or Left White, cotton sleeve, stays	21
11-0260-x	Tru-Pull Lite - Right Sleeve	17
11-0261-x	Tru-Pull Lite - Left Sleeve	17
11-0320-x-06060	Lateral "J" - Right Neoprene sleeve	19
11-0321-x-06060	Lateral "J" - Left Neoprene sleeve	19
11-0323-x-06060	"H" Buttress - Right or Left Sleeve, popliteal cutout, stays	19
11-0327-x-06060	Deluxe Donut - Right or Left Beige, cotton sleeve, stays	20
11-0334-x-06060	Horseshoe Patella - Right or Left Sleeve, stays	20
11-0344-x-06060	Horseshoe Patella - Right or Left Sleeve, popliteal cutout, stays	20
11-0345-x-04000	Deluxe Elastic Knee - Right or Left Beige, cotton sleeve, stays	21
11-0345-x-15000	Deluxe Elastic Knee - Right or Left White, cotton sleeve, stays	21
11-0379-x-06000	Adjustable Patella Donut - Right or Left Neoprene sleeve, popliteal cutout	20
11-0380-x-06060	Deluxe Donut - Right or Left Sleeve, popliteal cutout, stays	20
11-0465-x-06000	Adjustable Patella Donut - Right or Left Neoprene sleeve	20
11-0549-x-06000	Adjustable Patella Donut - Right or Left Drytex® sleeve	20
11-0554-x-06000	Adjustable Patella Donut - Right or Left Drytex® sleeve, popliteal cutout	20
11-0775-x	Hinged Lateral "J" - Neoprene, Right	18
11-0776-x	Hinged Lateral "J" - Neoprene, Left	18
11-0777-x	Hinged Lateral "J" - Drytex, Right	18
11-0778-x	Hinged Lateral "J" - Drytex, Left	18
11-0828-x-06060	Hinged "H" Buttress - Right or Left Sleeve, popliteal cutout, hinges	19
11-0855-x-06000	Hinged Air DonJoy, XS - XXL	18
11-0260-x	Tru-Pull Lite™ - Right Sleeve	17
11-0261-x	Tru-Pull Lite™ - Left Sleeve	17
11-1289-x	Tru-Pull® Advanced System - Right Sleeve	18
11-1290-x	Tru-Pull® Advanced System - Left Sleeve	18
11-1292-x	Tru-Pull® Advanced System - Right Sleeve, popliteal cutout	18
11-1293-x	Tru-Pull® Advanced System - Left Sleeve, popliteal cutout	18
11-1327-x-06000	Cross Strap - Right or Left - Qty 1	21
11-1328-x-06000	Cross Strap - Right or Left - Qty 25	21
11-1525-x	Tru-Pull Advanced Attachment	18
11-2289-x	Tru-Pull® Advanced System - Right Sleeve, hinges	18
11-2290-x	Tru-Pull® Advanced System - Left Sleeve, hinges	18
11-2292-x	Tru-Pull® Advanced System - Right Sleeve, popliteal cutout, hinges	18
11-2293-x	Tru-Pull® Advanced System - Left Sleeve, popliteal cutout, hinges	18

Post-Op Bracing

11-0024-x	Ankle Stirrup, Right	16
11-0026-x	Ankle Stirrup, Left	16
11-5115-9	TROM Adjuster Medial RT/Lateral LT	16
11-5116-9	TROM Adjuster Medial LT/Lateral RT	16
11-5117-9	TROM Adjuster Medial RT/Lateral LT, XL	16
11-5118-9	TROM Adjuster Medial LT/Lateral RT, XL	16
11-9114-9	TROM Advance, Cool	16
11-9115-9	TROM Advance, Full Foam	16
11-9116-9	TROM Advance, Cool, XL	16
11-9117-9	TROM Advance, Full Foam, XL	16

Upper Extremity

05A	Aircast Pneumatic Armband	38
05EL	Aircast Mayo Clinic Elbow Brace, Left	36
05ER	Aircast Mayo Clinic Elbow Brace, Right	36
05EW	ARC Forearm Rotational Brace	36
05FML	Aircast StabilAir™ Wrist Brace, Left, Medium	38
05FMR	Aircast StabilAir™ Wrist Brace, Right, Medium	38
05FSL	Aircast StabilAir™ Wrist Brace, Left, Small	38
05FSR	Aircast StabilAir™ Wrist Brace, Right, Small	38
05WLL	Aircast A2™ Wrist Brace, Left, Large	37
05WLR	Aircast A2™ Wrist Brace, Right, Large	37
05WML	Aircast A2™ Wrist Brace, Left, Medium	37
05WMR	Aircast A2™ Wrist Brace, Right, Medium	37
05WSL	Aircast A2™ Wrist Brace, Left, Small	37
05WSR	Aircast A2™ Wrist Brace, Right, Small	37

Part No.	Description	Page No.
05WTL	Aircast A2™ Wrist Brace w/ Thumb Spica, Left, Large	37
05WTLR	Aircast A2™ Wrist Brace w/ Thumb Spica, Right, Large	37
05WTML	Aircast A2™ Wrist Brace w/ Thumb Spica, Left, Medium	37
05WTMR	Aircast A2™ Wrist Brace w/ Thumb Spica, Right, Medium	37
05WTS	Aircast A2™ Wrist Brace w/ Thumb Spica, Left, Small	37
05WTSR	Aircast A2™ Wrist Brace w/ Thumb Spica, Right, Small	37
11-0180-x-13066	IROM® Elbow - Right	35, 37
11-0181-x-13066	IROM® Elbow - Left	35
11-0449-x	UltraSling III	33
11-0450-x	UltraSling III AB	33
11-0525-x	Sully® Shoulder Stabilizer	34
11-0526	Sully® Shoulder Stabilizer Custom	34
11-0527	Sully® Shoulder Stabilizer Bilateral	34
11-0528	Sully® Shoulder Stabilizer Sully AC Pad	34
11-1101-x-06000	DonJoy® Shoulder Stabilizer - Right	34
11-1102-x-06000	DonJoy® Shoulder Stabilizer - Left	34
11-1106-0-06000	Replacement strap set (5 straps)	34
11-1319-x-06000	DonJoy® Shoulder Stabilizer SPA	34
11-1343-x-06000	Ultrasling® ER - 15°, Right	33
11-1344-x-06000	Ultrasling III ER - 15°, Left	33
11-1345-x-06000	Ultrasling® ER - 30°, Right	33
11-1346-x-06000	Ultrasling® ER - 30°, Left	33
11-1530	Derotation Strap	33
11-1530	Ultrasling III / Ultrasling III AB Derotation Strap	33
11-2135-x	sXt Functional Wrist Brace, S - L	37
11-XXXX-x	cXt Functional Wrist Brace, S - L	37
11-9100	Ultrasling® ER - 15° Wedge (each)	33
11-9101	Ultrasling® ER - 15° Wedge (5-pack)	33

Walking Braces

01A-L	SP Walker™ (short pneumatic), Large	24
01A-M	SP Walker™ (short pneumatic), Medium	24
01A-P	SP Walker™ (short pneumatic), Pediatric	24
01A-S	SP Walker™ (short pneumatic), Small	24
01A-XL	SP Walker™ (short pneumatic), X-Large	24
01F-L	FP Walker™ (foam pneumatic), Large	24
01F-M	FP Walker™ (foam pneumatic), Medium	24
01F-P	FP Walker™ (foam pneumatic), Pediatric	24
01F-S	FP Walker™ (foam pneumatic), Small	24
01F-XL	FP Walker™ (foam pneumatic), X-Large	24
01PD-L	XP Diabetic Walker™ System, Large	26
01PD-M	XP Diabetic Walker™ System, Medium	26
01PD-XL	XP Diabetic Walker™ System, X-Large	26
01P-L	XP Walker™ (extra pneumatic), Large	24
01P-M	XP Walker™ (extra pneumatic), Medium	24
01P-P	XP Walker™ (extra pneumatic), Pediatric	24
01P-S	XP Walker™ (extra pneumatic), Small	24
01P-XL	XP Walker™ (extra pneumatic), X-Large	24
11-1371-x	MaxTrax Air	25
11-1372-x-00000	MaxTrax® / MaxTrax® Ankle Liner Pad Kit	25
11-1373-x-00000	MaxTrax® / MaxTrax® Ankle Replacement Liner, tall	25
11-1375-x-00000	MaxTrax® / MaxTrax® Ankle Anterior/Posterior and Medial/Lateral Shells	25
11-1377-x-00000	MaxTrax® / MaxTrax® Ankle Replacement Liner, short	25
11-1380-x	MaxTrax Air Ankle	25
11-1381-x	MaxTrax ROM Air	25
11-1382-x-00000	MaxTrax	25
11-1383-x-00000	MaxTrax ROM	26
11-1384-x-00000	MaxTrax Ankle	25
11-1385-x-00000	MaxTrax ROM Ankle	26
11-1390-x	MaxTrax ROM Air Ankle	25
11-1546-x	MaxTrax® Diabetic Walker	26
25-0241	Upright Padding Kit	26
25-1345-x	Replacement Boot Liner	26
25-1457	MaxTrax® Diabetic Walker Toe Cover	26
25-1626-x	Foam Impax Grid	26
79-95310	MaxTrax® / MaxTrax® Ankle Walker strap extension set	25

Wrist & Elbow


11-0001-x-06000	Elbow Support	36
11-0392-x-06060	Spider® Pad Elbow	36
11-0401-x-06000	Elbow Support - w/Pad	36
11-0543-x-06000	DonJoy® Neoprene Suspension Strap	35
11-0956-0-06000	DonJoy® Supra Condylar Pad Kit	35
11-1003-x-06000	DonJoy® Elbow Guard	35
L-3915	DonJoy® CXT Wrist Brace, Lacer	37
L-3915	DonJoy® CXT Wrist Brace, Speed Strap	37

Part No.	Description	Page No.
L-3915	DonJoy® CXT Wrist Brace, Sport Strap	37
L-3915	DonJoy® SXT Wrist Brace, Lacer	37
L-3915	DonJoy® SXT Wrist Brace, Speed Strap	37
L-3915	DonJoy® SXT Wrist Brace, Sport Strap	37
11-3147	Palmer Stay Replacement Kit	37
11-3488-x	Hand Soft Goods Kit, Laser Sport	37
11-3489	SXT ROM Extension Kit	37
11-3522	CXT ROM Extension Kit	37
11-3490	SXT Ulnar/Radial Kit	37
11-3491-x	Speed Strap Kit	37
11-3492-x	Sport Strap Kit	37
11-3493-x	Lacer Strap Kit	37
11-3494-x	Lacer Strap Kit	37
11-3521	Over Padding Kit	37
25-0422-0-00000	DonJoy® Flexion / Extension Stop Kit	35
25-1923-0-06000	DonJoy® Air Condyle Kit	35

DJO provides possible coding suggestions based on publicly available information as a convenience to our customers. We do not make claims, promises or guarantees as to the availability of reimbursement for any DJO product. It is within the sole discretion of the customer to determine the appropriate billing code for a product, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payor.

Part No.	HCPCS	Part No.	HCPCS	Part No.	HCPCS
01A	L4360	11-0451	L1902	11-1268	A9270
01F	L4360	11-0494	E0218	11-1269	A9270
01P	L4360	11-0525	L3675	11-1289	A4466
01PD	L4360	11-0526	L3675	11-1290	A4466
02AL	L4350	11-0527	L3675	11-1292	A4466
02ALP	L4350	11-0528	L3675	11-1293	A4466
02AR	L4350	11-0541	A4466	11-1319	L3675
02ARP	L4350	11-0555	L1820	11-1326	A9900
02BL	L4350	11-0556	L1820	11-1327/1328	L2999
02BLP	L4350	11-0557	L1832	11-1334	A9270
02BR	L4350	11-0558	L1832	11-1336	A9270
02BRP	L4350	11-0666	L1832	11-1337	A9270
02CL	L4350	11-0675-9-00000	A9270	11-1338/1341	L1845
02CLP	L4350	11-0676-9-00000	A9270	11-1339	L1845
02CR	L4350	11-0678-9-00000	A9270	11-1340	L1845
02CRP	L4350	11-0679-9-00000	A9270	11-1343/1346	L3670
02DL	L4350	11-0680-9-00000	A9270	11-1344	L3670
02DR	L4350	11-0681-9-00000	A9270	11-1345	L3670
02E	L4350	11-0682-9-00000	A9270	11-1371	L4360
02F	L4350	11-0759	L1832	11-1380	L4360
02JL	L4350	11-0775	L1810	11-1381	L4360
02JR	L4350	11-0776	L1810	11-1382	L4386
02M	L1906	11-0777	L1810	11-1383	L4386
02T	L1902	11-0778	L1810	11-1384	L4386
03AL	L4370	11-0801-9-00000	A9270	11-1385	L4386
03AR	L4370	11-0802-9-00000	A9270	11-1390	L4360
03BL	L4370	11-0804-9-00000	A9270	11-1397	A9270
03BR	L4370	11-0813-9-00000	A9270	11-1398-9	E1399
03CL	L4370	11-0815-9-00000	A9270	11-1421	E0218
03CR	L4370	11-0817-9-00000	A9270	11-1422	E0218
03DL	L4370	11-0860	L1832	11-1423	E0218
03DR	L4370	11-0861	L1832	11-1424	E0218
05A	L3701	11-0864	L1832	11-1425	E0218
05E	L3760	11-0865	L1832	11-1426	E0218
05F	L3982	11-0870/0871	L1845	11-1440/1441	L1845
05W	L3908	11-0872/0873	L1845	11-1442	L1845
09AS	L2999	11-0883	A9270	11-1443	L1845
09DS	L4396	11-0884	A9270	11-1489/1492	L1971
09DL	L4396	11-0885	A9270	11-1490	L1971
09AM	L2999	11-0886	A9270	11-1491	L1971
09AL	L2999	11-0887	A9270	11-1493/1496	L1971
10A	A9270	11-0888	A9270	11-1494	L1971
10A01	A9270	11-0889	A9270	11-1495	L1971
10B	A9270	11-0895	A9270	11-1497/1500	L1971
10B01	A9270	11-0896	A9270	11-1498	L1971
10C	A9270	11-0897-4	A9270	11-1499	L1971
10C01	A9270	11-0897-5	A9270	11-1521	A9270
10D01	A9999	11-0898-2	A9270	11-1522	A9270
10P	A9270	11-0899-5	A9270	11-1531/1534	L1843
10P01	A9270	11-0900	A9270	11-1532	L1843
10W	A9270	11-0901	A9270	11-1533	L1843
10W01	A9999	11-0902-2	A9270	11-1546	L4386
11-0001	L3701	11-0902-4	A9270	11-1561	L0460
11-0032	A4466	11-0904	A9270	11-1562	L0637
11-0180/0181	L3760	11-0905	A9270	11-1568	L0627
11-0215	L1810 + L2795	11-0906	A9270	11-1569	L0627
11-0242	L1825	11-0909	A9270	11-1570	L0626
11-0258	L1845	11-0910	A9270	11-2006	L1820
11-0259	L1845	11-0911	A9270	11-2007	L1820
11-0260	L1810 + L2795	11-0912	A9270	11-2008	L1820
11-0261	L1810 + L2795	11-0913	A9270	11-2009	L1820
11-0264	L1845	11-0914	A9270	11-2012	L1820
11-0265	L1845	11-0915	A9270	11-2013	L1820
11-0320/0321	A4466	11-0916	A9270	11-2123	A9270
11-0323	A4466	11-0917	A9270	11-2124	A9270
11-0334	A4466	11-0918	A9270	11-2125	A9270
11-0344	A4466	11-0947	L1843	11-2126	A9270
11-0345	A4466	11-0948	L1843	11-2130	L3915
11-0392	L3701	11-0949	L1843	11-2131	L3915
11-0401	L3701	11-0950	L1843	11-2132	L3915
11-0449	L3670	11-1003	L3760	11-2135	L3915
11-0450	L3670	11-1101/1102	L3675	11-2140	L3915

Part No.	HCPCS
11-2141	L3915
11-2289	L1810 + L2795
11-2290	L1810 + L2795
11-2292	L1810 + L2795
11-2293	L1810 + L2795
11-3220	L1845
11-3221	L1845
11-3222	L1845
11-3223	L1845
11-5115-9	L1832
11-5116-9	L1832
11-5117-9	L1832
11-5118-9	L1832
11-9098	E0218
11-9099	E0218
11-9100	L3670
11-9101	L3670
11-9111	A9270
11-9112	A9270
11-9114/9116	L1832
11A	A9270
11A01	A9270
11A01EBI	A9999
11ASC	A9270
11B	A9270
11B01	A9270
11C	A9270
11C01	A9270
11D01	A9999
11P01	A9270
12A	A9270
12A01	A9270
12AXL	A9270
12AXL01	A9270
13A	A9270
13A01	A9270
13B01	A9270
13C	A9270
13C01	A9270
14A	A9270
14A01	A9270
15A	A9270
15A01	A9270
16A	A9270
16A01	A9270
51A	E0236
51A10A	E0218
51A11A	E0218
51A11B	E0218
51A12A	E0218
51A14A	E0218



DONJOY

Copyright © 2011 DJO, LLC. Air DonJoy, Air-Stirrup, Air-Stirrup Universe, A60, A2, AirSport, XP Walker, FP Walker, SP Walker, CI, Cryo/Cuff, Defiance, Drytex, Dura*Kold, Dura*Soft, FemaleFource, IceMan, Knee Guarantee, MaxTrax, Never Stop Getting Better, OA Lite, OAdjuster, Playmaker, StabilAir, Sully, TROM, TruPull, UltraSling, Velocity and Walkabout are trademarks of DJO, LLC.

Lycra is a registered trademark of E.I. du Pont de Nemours & Company, Velcro is a registered trademark of Velcro Industries B.V., FLOAM is a trademark of Mattel Inc., Breath-o-prene is a trademark of Accumed.

DONJOY®



DJO, LLC | A DJO Global Company
T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOglobal.com

Together in Motion™

DJO Global proudly sponsors:

