



# Exos™ Undersleeve

## Ideal Application

The soft, microfiber material in the Exos wrist undersleeve, enhances the comfort and protects the patients' skin from shear and moisture. Designed with a slight compression and a flat seam to eliminate skin irritation, the wrist undersleeve easily fits under the Exos wrist and fracture braces.

## Features and Benefits

- Provides protection from skin shear
- Soft and comfortable
- Wicks away moisture
- Keeps skin and inside of brace clean, dry and more hygienic
- Sold in pairs – one to wear and one to wash
- Reimbursable item – L-Code L3995

### Exos Wrist and Finger Undersleeve

The new wrist and finger undersleeve has added trim lines that provide the clinician the ability to easily trim to fit under Exos upper extremity fracture braces that include fingers.

*Added trim line to customize to brace type*



### Exos Wrist Undersleeve

The original wrist undersleeve design can be used under any of the Exos wrist and fracture braces.





# Exos™ Undersleeve

## Application Configurations

### Exos Wrist and Finger Undersleeve



Use with:

- Boxer's Fracture Brace
- Radial Gutter Fracture Brace
- Thumb Spica Fracture Brace
- Long Thumb Spica w/Boa I & II



### Exos Wrist Undersleeve



Use with:

- Short Arm Fracture Brace
- Short Arm Fracture Brace Open Thumb
- Pediatric Short Arm Fracture Brace
- Wrist Brace w/Boa



### Ordering Information

Part #	Size	Wrist Circumference	
<b>Exos Wrist Undersleeve</b>			
896-30-1111	PEDS	Min 3.5"(9 cm)	Max 5"(13 cm)
896-40-1111	XS/S	Min 5"(13 cm)	Max 6.5"(16 cm)
896-50-1111	M	Min 6.5"(16 cm)	Max 8"(20 cm)
896-60-1111	L/XL	Min 8"(20 cm)	Max 11"(28 cm)
<b>Exos Wrist and Finger Undersleeve</b>			
897-40-1111	XS/S	Min 5"(13 cm)	Max 6.5"(16 cm)
897-50-1111	M	Min 6.5"(16 cm)	Max 8"(20 cm)
897-60-1111	L/XL	Min 8"(20 cm)	Max 11"(28 cm)

For less compression, choose one size larger • For more compression, choose one size smaller



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

[DJOGlobal.com/exos](http://DJOGlobal.com/exos)