

ONE PIECE Cervical Collar

ONE PIECE injection molded collar provides proper support and helps restrict cervical spine flexion, extension and rotation to assist in patient recovery following trauma. Constructed with durable high density polyethylene lined with a soft, closed cell foam for maximum patient comfort. Both universal sizes have multiple heights to fit most adults and pediatrics. Oversized trachea opening provides quick access to carotid pulse monitoring and emergency tracheotomies. Stacks flat to preserve space and is X-ray, CT and MRI lucent.

PN'S	HEIGHT ADJUSTMENT*	NECK CIRCUM.	SIZE
ONE PIECE CERVICAL COLLAR			
79-83130	1.5", 1.75", 2"	8" – 18"	PEDIATRIC
79-83131	1.75", 2.5", 3", 3.5"	11" – 23"	ADULT

*Height = measurement from sternal notch to chin.

IDEAL FOR: Cervical spine immobilization and extraction.



ONE PIECE Cervical Collar
Adult



ONE PIECE Cervical Collar
Pediatric



T 800.526.8785 D 760.727.1280 F 800.457.4221
2900 Lake Vista Drive | Dallas, TX 75067 | U.S.A.
DJOGlobal.com

Copyright © 2018 by DJO, LLC
MKT00-8224 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.