

## PLAYMAKER™ II X-ROM®

Designed to help provide immobilization and protected range of motion to the knee associated with ligament and meniscal injuries and surgeries. The Playmaker II X-ROM combines the latest in DonJoy® technology post-operative care and softgoods, with easy range of motion settings built into the XROM hinge and a comfortable wrap around design to promote healing and recovery after surgery or injury.

Standard OTS sizing from XS-XXXL fits a wide range of adult, youth and pediatric patients.

#### INDICATIONS:

For immobilization and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs..



PRECISE ROM PROTECTION
Control ROM with settings in
10-degree increments from -10
degrees to 90 degrees extension and
-10 degrees to 120 degrees flexion.



**QUICK-LOCK DESIGN**Easily lock at 10-degree increments from -10 to 30 degrees.



COMFORTABLE FIT High-tenacity (Hi10™) nylon/spandex fabric offers unmatched comfort and fit.



BREATHABLE SPACER
MESH FABRIC
Lightweight spacer fabric on front and back panels promotes a "cool" environment support.



FAST, EASY APPLICATION
Molded strap tabs allow quick
and easy application and strap
adjustments to deliver a secure and
comfortable fit.



ANTI-MICROBIAL COATING
Neutralizes odor while increasing comfort and durability.



**CUSTOMIZABLE DESIGN**Trimmable, adjustable strap and buckle design allows off-the-shelf and customizable applications.



ACCOMMODATES SOME PEDIATRIC PATIENTS Refer to size chart for leg measurements

#### SIZED TO FIT MOST PATIENTS

PART NUMBER	DESCRIPTION	
11-0021-x	Playmaker™ II X-ROM®	
Y=1-7 YS-YYYI		

### PDAC Assigned Code: L1833

PART NUMBER	DESCRIPTION	
25-0021	Playmaker Standard Hinge Kit	

Call 800.336.6569 or visit djoglobal.com/donjoy

SIZE	THIGH	KNEE CENTER	CALF
XS (X=1)	13" – 15.5" (33 – 39 cm)	12" – 13" (30.5 – 33 cm)	10" – 12" (25.5 – 30.5 cm)
S (X=2)	15.5" – 18.5" (39 – 47 cm)	13" – 14" (33 – 35.5 cm)	12" – 14" (30.5 – 35.5 cm)
M (X=3)	18.5" – 21" (47 – 53.25 cm)	14" – 15" (35.5 – 38 cm)	14" – 16" (35.5 – 40.5 cm)
L (X=4)	21" - 23.5" (53.25 - 59.5 cm)	15" – 17" (38 – 43 cm)	16" – 18" (40.5 – 47 cm)
XL (X=5)	23.5" – 26.5" (59.5 – 67.25 cm)	17" – 19" (43 – 48.25 cm)	18" – 20" (47 – 50.75 cm)
XXL (X=6)	26.5" – 29.5" (67.25 – 75 cm)	19" - 21" (48.25 - 53.25 cm)	20" – 22" (50.75 – 56 cm)
XXXL (X=7)	29.5" – 32" (75 – 81.25 cm)	21" - 23" (56 - 61 cm)	22" - 24" (53.25 - 58.5 cm)

 $\text{Circumference measurements should be taken at knee center, 6" (15 \, \text{cm}) above knee center and 6" (15 \, \text{cm}) below knee center. }$ 

# enovis...

T 800.336.6569 F 800.457.4221

DJO, LLC 5919 Sea Otter Place, Suite 200 | I Carlsbad, CA 92010 | I U.S.A. enovis.com/donjoy

Copyright © 2023 DJO, LLC MKT00-12517 Rev B