

# PRE-SHAPED BONE WEDGES

Pre-shaped Bone Wedges offer surgeons convenient options for approaching opening wedge osteotomies in the foot and ankle, specifically shaped for Evans calcaneal osteotomies and Cotton osteotomies. The bicortical edge provides structural support and strength allowing impaction of the graft into its precise location while dense cancellous bone provides interconnected porosity to support bone incorporation.

## BIOLOGICS

- Convenient size-specific, pre-shaped design decreases the amount of trimming and shaping necessary, saving valuable O.R. time
- Sterilized through the BioCleanse® Tissue Sterilization Process and terminally sterilized using a validated gamma irradiation dose to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>
- Easy-to-use, all-in-one instrumentation set including size-specific trials, a distractor, osteotomes, and tamps



Cotton Wedge Part Number	Size Dimensions (LxWxH)
D18805-TS	22 x 16 x 5mm
D18806-TS	22 x 16 x 6mm
D18807-TS	22 x 16 x 7mm
D18808-TS	22 x 16 x 8mm

Evans Wedge Part Number	Size Dimensions (LxWxH)
D18006-TS	22 x 22 x 6mm
D18008-TS	22 x 22 x 8mm
D18010-TS	22 x 22 x 10mm
D18012-TS	22 x 22 x 12mm

Manufactured by RTI Surgical, Inc., 11621 Research Circle, Alachua, FL 31615. Trilliant is a distributor for RTI Surgical, Inc.  
 ® indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.



T 800.495.2919 F 877.778.3864  
 727 North Shepherd Drive, Suite 100 | Houston, TX 77007 | U.S.A.  
[djoglobal.com](http://djoglobal.com)

Copyright © 2021 by DJO, LLC  
 900-00-104 Rev C

Individual results may vary. DJO, LLC is a manufacturer of orthopedic implants and does not practice medicine. Only an orthopedic, or foot and ankle surgeon can determine what treatment is appropriate. The contents of this document do not constitute medical, legal, or any other type of professional advice. This material is intended for the sole use and benefit of the DJO, LLC sales force and physicians. It is not to be redistributed, duplicated, or disclosed without the express written consent of DJO, LLC. For more information on risks, warnings, and possible adverse side effects refer to the Instructions for Use provided with the device.