



# Exos™ Wrist Undersleeve

## Ideal Application

Low profile wrist undersleeve with integrated antimicrobial fabric is used to enhance comfort and protect skin from shear and moisture. Designed to fit easily under the EXOS wrist braces and other upper extremity products.

## Features and Benefits

- Provides protection from skin shear
- Soft and comfortable
- Wicks away moisture
- Keeps skin and inside of brace clean, dry and more hygienic
- Sold in pairs – one to wear and one to wash

**NEW  
Pediatric Size  
Now  
Available.**




Shown:  
Wrist Undersleeve  
alone



Shown:  
Wrist Undersleeve  
with EXOS brace



## How to Apply:



### Ordering Information

Part #	Size	
<b>Wrist Undersleeve</b>		
896-30-1111	PEDS	
896-40-1111	XS/S	
896-50-1111	M	
896-60-1111	L/XL	
<b>Sizing Chart</b>		
	<b>Wrist Circumference</b>	
Size	MIN	MAX
PEDS	3.5" (9 cm)	5" (13 cm)
XS/S	5" (13 cm)	6.5" (17 cm)
M	6.5" (17 cm)	8" (20 cm)
L/XL	8" (20 cm)	11" (28 cm)

For less compression, choose one size larger

For more compression, choose one size smaller



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569  
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
[DJOGlobal.com/exos](http://DJOGlobal.com/exos)

