

COME BACK STRONG

*DONJOY*<sup>®</sup>  
*X-ROM*<sup>™</sup>





## EASY-TO-USE POST-OP ROM PROTECTION WITH FAST APPLICATION AND COMFORTABLE STABILITY.

Our new easy-to-use post-op knee brace offers fast, simple application and adjustment through innovative strap, tab and buckle designs. The user-friendly X-ROM™ helps control and protect ROM with a quick-lock design and hinge pull-out stops. With four independently telescoping sliders and the ability to contour to the patient's leg to aid varus and valgus support, the X-ROM comfortably accommodates a wide range of patient shapes and sizes.

### Precise ROM Protection

- Updated hinge design: Precisely control ROM with settings in 10° increments from -10° to 90° extension and -10° to 120° flexion
- Quick-lock design: Easily lock at 10° increments from -10° to 30°

### Fast, Easy Application


- Easy-to-use design: Strap, tab and buckle design make post-op application and adjustment fast and simple
- Comfortably fits a wide range of patients: Easily accommodate patients from 4'11" to 6'5"
- Push button slider settings: Adjustments are easy and precise
- Ready to fit: Quickly fit patients right out of the box

### Comfortable Stability

- Tele-Fit™ technology: Four independently telescoping sliders with easy push button control settings help enhance brace placement, stability and comfort while avoiding the incision and peroneal nerve
- Optimal fit & stability: Bendable hinge bars help enhance varus/valgus contouring
- Secure contact points: Enhanced ankle pads and strategically placed contact points between the leg and brace help reduce slippage
- Lightweight strength: Strong yet lightweight design helps reduce brace migration
- Anti-microbial pads: Helps neutralize odor while increasing comfort and durability







### Ordering information

Part Number	Description	Size
11-2181-9	X-ROM™ Post-Op Knee Brace	Universal
25-7530	X-ROM Replacement Straps Kit	
25-7531	X-ROM Replacement Pad Kit	

Indications: To aid immobilization and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs



T 800.336.6569 D 760.727.1280 F 800.936.6569  
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.  
[DJOGlobal.com/donjoy](http://DJOGlobal.com/donjoy)

Copyright © 2020 by DJO, LLC  
MKT00-8796 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.