

enovis™

DONJOY®
X-ACT FIT SHOULDER™



DonJoy® X-ACT Fit™ Shoulder

The new DonJoy® X-ACT Fit™ Shoulder brace provides versatility, comfort, and functionality to help improve the postoperative or post-traumatic experience for you and your patients.

The unique design helps support the patient's arm and shoulder with a universal, customizable immobilization system that's easy to apply, adjust, and helps allow for a secure and functional recovery.

AB/AD
Abduction/Adduction

15° to 90°

EZ Axis Hinge and malleable armrest platform provides simple, customization for postoperative immobilization

Conforming Waist Band for improved fit and comfort

Easy-to-adjust, universal size helps reduce stocking units

Approach Recovery from Any Angle

The DonJoy® X-ACT Fit™ Shoulder offers a customizable, multi-angle system to help comfortably immobilize the patient and protect postoperative repairs during rehabilitation. The universal design fits right or left arm and features a waistband and arm envelope that can be adjusted and sized to accommodate a wide range of patients.

Patient Comfort and Compliance

- Lightweight, easy to don and doff while at home
- Soft, breathable mesh arm sling material helps improve airflow and comfort
- Precise-fit waist band and quick-release buckles make it easy to reapply brace at home or perform exercises without removing it
- Padded armrest platform to comfortably rest the forearm during rehabilitation
- Optional three-point off-loading shoulder ring to evenly distribute weight across the uninjured shoulder and can minimize stress on the neck
- Optional waist belt extension to accommodate larger patients



Adjustable waist belt with front adjustment strap provides a contoured, customized fit every time



Quick-release buckles and trimmable straps for ease of application and proper fit



Movable hip pad and shoulder ring allow for additional comfort and support during the healing process

IR/ER

Internal Rotation/External Rotation

0° to 45°+

DONJOY®

X-ACT FIT SHOULDER™

Brace Part Number	Description
11-0495	DonJoy® X-ACT Fit™ Shoulder Brace, Universal
Accessory Part Number	
25-0308	Arm Sling Shoulder Ring, Universal
25-0309	Waist Belt Extension
11-9098	Iceman Classic ³ ® with Shoulder Wrap-On Pad
11-1424	Iceman Classic ³ ® with Universal Wrap-On Pad



Indications:

Rotator cuff repairs; Bankart and capsular shift procedures; Soft-tissue strains; Shoulder dislocation/subluxation; Non or minimally displaced proximal humerus fractures; Pectoral tendon rupture and distal biceps rupture; AC-Joint instabilities; Proximal biceps tenodesis, SLAP-Repairs; ORIF clavicle fractures; Humeral shaft fractures; Scapular fractures; SC joint separation, UCL recon of the elbow and other elbow fracture indications; Non- and post-surgical rehabilitation.

PDAC approved L3960



enovis™

T 800.336.6569 F 800.4574221

DJO, LLC
5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A.
enovis.com/donjoy

Copyright © 2024 DJO, LLC
MKT00-13293 Rev A

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.