



enovis™

DONJOY®
OA DEFIANCE® PRO

Confidence Without Compromise



Defiance Pro
(shown in Black Gloss)



Black Gloss Black Matte Metallic Black Gloss Metallic Gunmetal Gloss Metallic Silver Gloss



Smoke Matte White Gloss Metallic White Gloss Bright Blue Gloss Metallic Bright Blue Gloss



Navy Gloss Metallic Navy Gloss Powder Gloss Seafoam Gloss Emerald Gloss Metallic Forest Gloss Metallic Gold Gloss



Volt Gloss Bumble Gloss Safety Gloss True Red Gloss Metallic True Red Gloss Fuchsia Gloss Cotton Candy Gloss



Crimson Gloss Grape Gloss Velvet Gloss Olive Drab Green Matte Desert Tan Matte Flat Dark Earth Matte Stars & Stripes Gloss

OA DEFIANCE[®] PRO

DonJoy's premiere custom knee brace for active patients with osteoarthritis just got better.

THE BEST CUSTOM OA CARBON FIBER BRACE IS NOW EVEN BETTER WITH:

- Lower profile design with internally mounted straps providing a more intimate fit.
- Swiveling straps create 100% contact and contour to the anatomy to grab the calf better which allows for removal of the anti-migration band and improved ease of application.
- C-6 Liner that is soft to the skin and provides moisture wicking and anti-microbial.
- Soft touch strap ends offer flexibility to form better to the leg while maintaining a secure closure.
- Memory foam condyle pads help provide even pressure on the OA offloading side and are great for all skin types.

COMBINED WITH THE ENHANCED LEGACY FEATURES:

- Strong, lightweight carbon composite frame is custom fabricated to each patient's specific measurements
- 0 to 7 degrees of correction can be built into the frame
- Additional offloading can be applied by the patient using a key as needed
- Telescoping condyle pad shifts the knee joint forces from the degenerative area to the healthy compartment
- Frame provides support for ligament instabilities
- Swooping thigh to provide better clearance for bilateral wear
- Medial or lateral compartment off-loading
- ACL, PCL, or combined instability (CI) strapping
- 30 custom frame colors and graphics
- Standard or short length frames can be combined to accommodate leg
- Flexion and extension control
- Option for legend strap ends
- Brace bag included

Extremely lightweight, the OA Defiance[®] PRO provides durable support for moderate to severe levels of osteoarthritis. The bilateral hinge bars offer support for ligament instabilities and work in conjunction with a patented telescoping condyle pad which puts the load control in the hands of the patient. Ideally suited for patients with active lifestyles or athletes.

4-POINTS-OF-LEVERAGE SYSTEM™

DonJoy[®] pioneered dynamic bracing with the 4-Points-of-Leverage System™, a unique cuff and strapping configuration that helps provide a net differential posterior force to the tibia, preventing anterior movement. The result is significantly reduced strain on the ACL.¹ Passive brace designs do not actively address instabilities.



WE BACK YOUR KNEES. EVERY DONJOY CUSTOM OA KNEE BRACE COMES WITH THE FOLLOWING WARRANTY:

Stainless steel	Lifetime
Non-stainless steel	5 years
Soft goods	6 months
Re-measure	4 months (no charge for first re-fit; additional re-fits at a charge)



OA DEFIANCE® PRO

Every OA Defiance PRO is custom fabricated based off measurements from the DONJOY Custom Contour Measuring Instrument (CCMI).

SIZE	THIGH CIRCUMFERENCE	CALF CIRCUMFERENCE
Minimum	12.5" (31.8 cm)	10.25" (26 cm)
Maximum	38.125" (97 cm)	26.25" (66.7 cm)
BRACE LENGTH	PATIENT HEIGHT RECOMMENDATION	
Standard	Calf 5'4" or taller (163 cm or taller)	
Short	Calf 5'4" or shorter (163 cm or shorter)	
RECOMMENDED HCPCS CODE:		
L1846		

ACCESSORIES

Ordering information

PART NUMBER	DESCRIPTION	SIZE
11-1525-x	Tru-Pull® Advanced Attachment	XS-XXXL
11-0398-9-00000	Full Extension Lockout Stop	Universal
11-0075-x-00000	Lycra Undergarment	S-XXL
11-0122-x	Neoprene Undergarment	XS-XXL
11-0016-x-06000	Sports Cover (standard)	S-XXL
11-0103-x-06000	Sports Cover (short)	S-XXL
25-1923-0-06000	Air Condyle Pads	Universal
11-1107-0-00000	Surround Condyle Pads	Universal
11-1006-x-06000	Impact Guard	S-XXL
11-1095-9-06000	Calf Pinch Guard	Universal
11-0059-9-00000	Atrophy Pad Kit	Universal

INDICATIONS:

- Moderate to severe unicompartmental osteoarthritis
- Moderate to severe ligament instabilities
- Post-operative rehabilitation
- Pre-operative rehabilitation

CCMI HOW TO FIT



T 800.336.6569 F 800.457.4221

DJO, LLC
5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A.
[enovis.com/donjoy](https://www.enovis.com/donjoy)

Copyright © 2025 DJO, LLC
MKT00-11146 Rev B

1. Palm HG, et al. Effects of Knee Bracing on Postural Control after ACL Ligament Rupture. Knee. 2011 Aug 24 [https://www.ncbi.nlm.nih.gov/pubmed/21871811]

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.