

AIRCAST®

PROCARE®



Product Catalog



AIRCAST/PROCARE/EXOS

DEDICATED TO HELPING PEOPLE REGAIN THE COMFORT OF NATURAL MOTION.

DJO Global is a leading designer, manufacturer and distributor of state-of-the-art orthopedic bracing and support products.

This catalog highlights our innovative Aircast®, Exos and ProCare® offerings. Aircast products feature our patented technology of graduated pneumatic compression that has set a standard of care for sprains and other injuries. Exos stands out as the only bracing system available for the treatment of fractures that's removable, adjustable, reformable and waterproof. For upper and lower extremity bracing and support, compression wear and patient care, ProCare has a complete line of rehabilitative products.

Count on DJO Global for scientifically based, technologically advanced, quality products. To learn more about our family of brands or our company, visit DJOglobal.com.



TABLE OF CONTENTS



SPECIALTY

Exos Bracing Systems	2
Cryo-Compression Therapy	11
Dynamic Compression	15
VenaFlow	16
ArterialFlow	19
Static Compression – Anti Embolism Stockings	21

UPPER EXTREMITY

Cervical Collars	23
Clavicle Supports	29
Shoulder Immobilizers	31
Arm Slings & Elevators	35
Elbow Braces & Supports	39
Wrist & Hand Supports	43
Wrist & Thumb Supports	50
Lumbosacral Supports	53
Abdominal & Rib Supports	57

LOWER EXTREMITY

Knee Immobilizers	60
Post-Op Knee Braces	63
Thigh, Knee & Calf Supports	67
Patellofemoral	71
Hinged Knee Braces	75
Ankle Braces & Support	79
Ankle Stirrups	79
Soft Ankle & Specialty	83
Plantar Fasciitis	87
Post-Op & Cast Shoes	89
Walking Braces	93

PATIENT CARE

Patient & Pressure Care	105
Traction & Patient Safety	114
Fracture Bracing & Splints	123
Compression Care	131

INDEX	136
-----------------	-----

EXOS BRACING SYSTEMS



Reformable, Removable Bracing System

Exos is the only bracing system offering you and your patients a removable, adjustable, reformable and waterproof solution for the treatment of fractures and other injuries requiring stabilization. Many Exos products utilize the proprietary Boa® Closure System for increased patient comfort and to accommodate swelling and atrophy.

- Completely reformable
- Easily removable
- Adjustable Boa® reel and lace system
- Waterproof and easy to clean
- Radiolucent

Short Arm Fracture Brace with Boa

Provides support and protection for injuries to the wrist and forearm. Used primarily for non-displaced fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as for post-operative support. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
310-31-XXXX	310-32-XXXX	XS
310-41-XXXX	310-42-XXXX	S
310-51-XXXX	310-52-XXXX	M
310-61-XXXX	310-62-XXXX	L
310-71-XXXX	310-72-XXXX	XL

See color options and sizing chart on page 8.





Short Arm Fracture Brace Open Thumb

The open thumb design allows for enhanced ease of fit and patient comfort while the lycra top cover allows for faster heating and cooling to speed the application process. The ability to heat form and custom fit ensures proper uniform support and comfortable fit. Provides support and protection for major injuries to the wrist. Used primarily for fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as post-operatively. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
312-31-XXXX	312-32-XXXX	XS
312-41-XXXX	312-42-XXXX	S
312-51-XXXX	312-52-XXXX	M
312-61-XXXX	312-62-XXXX	L
312-71-XXXX	312-72-XXXX	XL

See color options and sizing chart on page 8.



Pediatric Short Arm Fracture Brace Open Thumb

The smallest of patients can benefit from Exos technology. The open thumb design allows for enhanced ease of fit and patient comfort while the lycra top cover allows for faster heating and cooling to speed the application process. The ability to heat form and custom fit ensures proper uniform support and comfortable fit. Exos braces are radiolucent and do not need to be removed for x-ray imaging. Provides support and protection for major injuries to the wrist. Used primarily for fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as post-operatively. Does not contain latex.

Left	Right	Size
311-11-XXXX	311-12-XXXX	XXXS
311-21-XXXX	311-22-XXXX	XXS

See color options and sizing chart on page 8.



Wrist Brace with Boa

Intended for carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapho-lunate dislocations, triangular fibro-cartilage complex (TFCC) tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures. May be used to support stable wrist fractures. May be used to control wrist motion for other injuries that require stabilization. Does not contain latex.

Left	Right	Size
221-21-XXXX	221-22-XXXX	XXS
221-31-XXXX	221-32-XXXX	XS
221-41-XXXX	221-42-XXXX	S
221-51-XXXX	221-52-XXXX	M
221-61-XXXX	221-62-XXXX	L
221-71-XXXX	221-72-XXXX	XL

See color options and sizing chart on page 8.



Wrist Brace

Intended for carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapho-lunate dislocations, triangular fibro-cartilage complex (TFCC) tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures. May be used to support stable wrist fractures. May be used to control wrist motion for other injuries that require stabilization. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
220-31-1111	220-32-1111	XS
220-41-1111	220-42-1111	S
220-51-1111	220-52-1111	M
220-61-1111	220-62-1111	L
220-71-1111	220-72-1111	XL

Available in black only. See sizing chart on page 8.



Thumb Spica Fracture Brace

Designed to provide stabilization to the thumb, first metacarpophalangeal (MCP) joint, the first carpometacarpal (CMC) joint, wrist joint, distal radius and distal ulna. Useful for all injuries to the thumb, such as Bennett's or complex Rolando's fractures, suspected or confirmed navicular (scaphoid) fractures, and for other carpal bone injuries. Also useful for fractures to the distal radius and/or distal ulna, Colles fractures, Salter-Harris pediatric wrist and forearm injuries. Used for acute injuries as well as for post-op support. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
315-31-1111	315-32-1111	XS
315-41-1111	315-42-1111	S
315-51-1111	315-52-1111	M
315-61-1111	315-62-1111	L
315-71-1111	315-72-1111	XL

Available in black only. See sizing chart on page 8.



Long Thumb Spica with Boa

Provides primary stabilization for the first metacarpophalangeal (MCP) joint and first carpometacarpal (CMC) joints, as well as secondary support for the radiocarpal and radioulnar wrist joints. Provides stability for injuries to the base of the thumb at the CMC joint, such as Bennett's or Rolando's fractures, or to stabilize a painful CMC joint. Also useful for confirmed or suspected navicular (scaphoid) fractures, and for injuries or arthritis involving the trapezium or trapezoid bones of the wrist. Does not contain latex.

Left	Right	Size
231-21-XXXX	231-22-XXXX	XXS
231-31-XXXX	231-32-XXXX	XS
231-41-XXXX	231-42-XXXX	S
231-51-XXXX	231-52-XXXX	M
231-61-XXXX	231-62-XXXX	L
231-71-XXXX	231-72-XXXX	XL

See color options and sizing chart on page 8.



Long Thumb Spica

Provides primary stabilization for the first metacarpophalangeal (MCP) joint and first carpometacarpal (CMC) joints, as well as secondary support for the radiocarpal and radioulnar wrist joints. Provides stability for injuries to the base of the thumb at the CMC joint, such as Bennett's or Rolando's fractures, or to stabilize a painful CMC joint. Also useful for confirmed or suspected navicular (scaphoid) fractures, and for injuries or arthritis involving the trapezium or trapezoid bones of the wrist. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
230-31-1111	230-32-1111	XS
230-41-1111	230-42-1111	S
230-51-1111	230-52-1111	M
230-61-1111	230-62-1111	L
230-71-1111	230-72-1111	XL

Available in black only. See sizing chart on page 8.



Short Thumb Spica

Immobilizes the base of the thumb to stabilize injuries to the ulnar collateral ligament (UCL) commonly known as Gamekeeper's thumb, fractures or avulsion injuries to the proximal first phalanx or distal first metacarpal, severe sprains of the thumb, or for arthritic conditions involving the metacarpophalangeal (MCP) joint. Does not contain latex.

Left	Right	Size
210-21-XXXX	210-22-XXXX	XXS
210-31-XXXX	210-32-XXXX	XS
210-41-XXXX	210-42-XXXX	S
210-51-XXXX	210-52-XXXX	M
210-61-XXXX	210-62-XXXX	L
210-71-XXXX	210-72-XXXX	XL

See color options and sizing chart on page 8.



Extended Short Thumb Spica

Immobilizes the base of the thumb and the interphalangeal joint, to stabilize injuries to the ulnar collateral ligament (UCL) commonly known as Gamekeeper's thumb, fractures or avulsion injuries to the proximal first phalanx or distal first metacarpal, severe sprains of the thumb, or for arthritic conditions involving the metacarpophalangeal (MCP) joint, or the interphalangeal (I-P joint). May also be used to immobilize the thumb joints following surgical repair of soft tissue injuries. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
211-31-1111	211-32-1111	XS
211-41-1111	211-42-1111	S
211-51-1111	211-52-1111	M
211-61-1111	211-62-1111	L
211-71-1111	211-72-1111	XL

Available in black only. See sizing chart on page 8.



Boxer's Fracture Brace

Used to stabilize fractures to the 4th and/or 5th metacarpal (MC) bones, and also provides stabilization to the 4th and 5th interphalangeal and wrist joints. May be used for metacarpal (MC) head/neck fracture of the 4th and/or 5th metacarpal. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
325-31-1111	325-32-1111	XS
325-41-1111	325-42-1111	S
325-51-1111	325-52-1111	M
325-61-1111	325-62-1111	L
325-71-1111	325-72-1111	XL

Available in black only. See sizing chart on page 8.




Radial Gutter Fracture Brace

Used to stabilize fractures to the 2nd and/or 3rd metacarpal (MC) bones and the respective metacarpal-phalangeal (MCP) joints of the hand. This pre-fabricated, custom-molded orthosis is designed to provide superior fit and support compared to plaster, fiberglass, or pre-made braces. The radial gutter portion of the Radial Gutter Fracture Brace allows the caregiver to achieve greater than 70° of flexion at the MCP joints of the 2nd and 3rd digits. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
326-31-1111	326-32-1111	XS
326-41-1111	326-42-1111	S
326-51-1111	326-52-1111	M
326-61-1111	326-62-1111	L
326-71-1111	326-72-1111	XL

Available in black only. See sizing chart on page 8.



Hand Based Ulnar Gutter Brace

The Hand Based Ulnar Gutter Brace is used to stabilize fractures to the 4th and/or 5th metacarpals or phalanges. This molded orthosis is designed to provide superior fit and support compared to plaster, fiberglass, or pre-made braces. The ulnar gutter portion of the Hand Based Ulnar Gutter Brace allows the caregiver to achieve greater than 70° of flexion at the MCP joints of the 4th and 5th digits. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
125-31-1111	125-32-1111	XS
125-41-1111	125-42-1111	S
125-51-1111	125-52-1111	M
125-61-1111	125-62-1111	L
125-71-1111	125-72-1111	XL

Available in black only. See sizing chart on page 8.



Hand Based Radial Gutter Brace

The Hand Based Radial Gutter Brace is used to stabilize fractures to the 2nd and/or 3rd metacarpals or phalanges. This molded orthosis is designed to provide superior fit and support compared to plaster, fiberglass, or pre-made braces. The radial gutter portion of the Hand Based Radial Gutter Brace allows the caregiver to achieve greater than 70° of flexion at the MCP joints of the 2nd and 3rd digits. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
126-31-1111	126-32-1111	XS
126-41-1111	126-42-1111	S
126-51-1111	126-52-1111	M
126-61-1111	126-62-1111	L
126-71-1111	126-72-1111	XL

Available in black only. See sizing chart on page 8.



Humeral Fracture Brace



Used to stabilize fractures to the mid-shaft of the humerus. The Humeral Fracture Brace may be appropriate for some injuries to the more proximal humerus. Does not contain latex.

Left	Size
330-40-1111	S/Univ.
330-50-1111	M/Univ.
330-60-1111	L/Univ.

Available in black only. See sizing chart on page 8.

Use the H measurement to select the appropriate size using the size chart below.
H=Height

- SMALL 4¼"
- MEDIUM ... 5"
- LARGE 6"



Boa Locking Ring Kit

The Boa Locking Ring is intended to prevent unprescribed loosening of the Boa reel or removal of Exos braces.

	Description
800-LR	Boa Locking Ring Kit



Exos Top Loading Oven

Convenient sized oven used to heat the Exos Brace in order to gently mold the brace to the patient's limb.

	Description
800-OVEN3P	Exos Top Loading Oven
Accessories	
800-ER	Oven Extender Ring
800-OVENBAG	Oven Bag

Warm it, Fit it

The Exos brace is completely reformable.

REFORMABLE

The Exos brace allows:

- Modification to account for anatomical changes as healing progresses
- Sequential straightening
- One brace throughout duration of treatment
- Less material waste

ADJUSTABLE

Exos uses the exclusive Boa® reel and lace system for dynamic compression, allowing:

- Expansion to allow for initial swelling and to reduce the risk of Compartment Syndrome
- Tension modification for custom fit despite atrophy or changes in soft tissue
- Improved patient comfort and fit



1. Simply heat brace in our oven for 3-5 minutes.



2. Open and slide onto patient's limb.



3. Gently mold and adjust for a snug, comfortable fit.



4. Use the reel and lace to tighten as needed.



Wrist Undersleeve

Low profile wrist undersleeve with integrated antimicrobial fabric is used to enhance comfort and protect skin from shear and moisture. Designed to fit easily under the EXOS wrist braces and other upper extremity products.

- Provides protection from skin shear
- Soft and comfortable
- Wicks away moisture
- Keeps skin and inside of brace clean, dry and more hygienic
- Sold in pairs – one to wear and one to wash

Part Number	Size
896-40-1111	XS/S
896-50-1111	M
896-60-1111	L/XL

Wrist Circumference	MIN	MAX
XS/S	5" (13 cm)	6½" (16 cm)
M	6½" (16 cm)	8" (20 cm)
L/XL	8" (20 cm)	11" (28 cm)

For more compression, choose one size smaller.



Other Accessories

(Photos not shown)

	Description	Size	Qty.
800-SIZINGTAPE	Exos Sizing Tape		1
800-02-1111	Edge tape - Lycra	1" x 24"	10 pk.
800-1111	Edge tape - UBL - Black	1" x 24"	10 pk.
800-02-2265	Edge tape - Denim	1" x 24"	10 pk.
800-02-2293	Edge tape - Spider	1" x 24"	10 pk.
800-02-2388	Edge tape - Camo	1" x 24"	10 pk.
800-02-3285	Edge tape - Polka Dot	1" x 24"	10 pk.
800-02-4444	Edge tape - Pink	1" x 24"	10 pk.
800-40-1111	Sheet Materials .040" thick	18" x 24"	1
800-50-1111	Sheet Materials .050" thick	18" x 24"	1
800-65-1111	Sheet Materials .065" thick	18" x 24"	1
800-ATK	Boa Adj. Tensioning Kit		1

	Description	Size	Qty.
800-40-STRAP	Accessory Strap	4"	10 pk.
800-45-STRAP	Accessory Strap	4.5"	10 pk.
800-50-STRAP	Accessory Strap	5"	10 pk.
800-55-STRAP	Accessory Strap	5.5"	10 pk.
800-60-STRAP	Accessory Strap	6"	10 pk.



Bracing Sizing and Color Chart







EXOS BRACING SIZING CHART

Wrap tape around hand to measure size.

Brace Size	Measurement
3XSML	5¼" – <6½"
2XSML	6½" – <7¼"
XSML	7¼" – <8"
SML	8" – <9"
M	9" – <10"
LG	10" – <11"
XLG	11"

Description	Qty.
800-SIZINGTAPE Exos Sizing Tape	1

EXOS BRACING COLOR CHART

	
4444 - Pink	1111 - Black
	
2293 - Spider	2388 - Camo
	
3285 - Polka Dot	2265 - Denim



Exos Rapid Splint™ System

The Exos Rapid Splint System has a reinforced composite polymer that creates a lightweight, stable structure that quickly becomes moldable using dry heat with no water and no mess. Rapid Splint goes on smooth and dry, without sharp fiberglass edges, enhancing patient safety and comfort. It reaches stability quickly, is reformable and offers superior radiolucency.

Part Number	Description
Exos Rapid Splint	
100-001-210	Exos Rapid Splint, 2" x 10" 10 Pack
100-001-312	Exos Rapid Splint, 3" x 12" 10 Pack
100-001-336	Exos Rapid Splint, 3" x 36" 5 Pack
100-001-413	Exos Rapid Splint, 4" x 13" 10 Pack
100-001-436	Exos Rapid Splint, 4" x 36" 5 Pack
100-001-513	Exos Rapid Splint, 5" x 13" 10 Pack
100-002-326	Exos Rapid Splint Extra, 3" x 26" 5 Pack
100-002-426	Exos Rapid Splint Extra, 4" x 26" 5 Pack
100-002-526	Exos Rapid Splint Extra, 5" x 26" 5 Pack
Exos Rapid Splint with Padding	
100-001-P210	Exos Rapid Splint with Padding, 2" x 10" 5 Pack
100-001-P312	Exos Rapid Splint with Padding, 3" x 12" 5 Pack
100-001-P336	Exos Rapid Splint with Padding, 3" x 36" 5 Pack
100-001-P413	Exos Rapid Splint with Padding, 4" x 13" 5 Pack
100-001-P436	Exos Rapid Splint with Padding, 4" x 36" 5 Pack
100-001-P513	Exos Rapid Splint with Padding, 5" x 13" 5 Pack
100-002-P326	Exos Rapid Splint Extra with Padding, 3" x 26" 5 Pack
100-002-P426	Exos Rapid Splint Extra with Padding, 4" x 26" 5 Pack
100-002-P526	Exos Rapid Splint Extra with Padding, 5" x 26" 5 Pack
Exos Rapid Splint with Exothermix Heater Pack	
100-001-E301	Exos Rapid Splint + Exothermix, 3" x 12" 1 Pack
100-001-E305	Exos Rapid Splint + Exothermix, 3" x 12" 5 Pack
100-001-E401	Exos Rapid Splint + Exothermix, 4" x 13" 1 Pack
100-001-E405	Exos Rapid Splint + Exothermix, 4" x 13" 5 Pack
100-001-E501	Exos Rapid Splint + Exothermix, 5" x 13" 1 Pack
100-001-E505	Exos Rapid Splint + Exothermix, 5" x 13" 5 Pack
100-001-EKIT	Exos Rapid Splint + Exothermix Kit, 3"x12", 4"x13", 5" x 13" 3 Pack
Accessories (Photos not shown)	
800-HEATERBAG	Exos Rapid Splint Warmer
79-98233	Elastic Bandages - Self Closure, 5 yards x 2" 10 Pack
79-98235	Elastic Bandages - Self Closure, 5 yards x 3" 10 Pack
79-98237	Elastic Bandages - Self Closure, 5 yards x 4" 10 Pack
79-98238	Elastic Bandages - Self Closure, 5 yards x 6" 10 Pack
800-40-STRAP	Accessory Strap 4" (10.2 cm), 10 Pack



EXOS Form™ 621

The Exos FORM™ 621 provides superior stabilization of the sacroiliac (SI) joint through direct compression. Silicon grip strips prevents the brace from migrating.

- Adjustable (Angle) Belt Wings
- Semi-Universal Design
- Simple, Light Weight Low Profile Design
- Sacroiliac Pain.

Part Number	Description	Color	Size
200621-40			S/M
200621-60			L/XL
Accessories			
800-AP-40	Anterior Panel	Black	S/M
800-AP-60	Anterior Panel	Black	L/XL
800-1-WINGS-40	Belt Wing Set	Black	S/M
800-1-WINGS-60	Belt Wing Set	Black	L/XL
800-1-BWING-40	Bridge Wing	Black	S/M
800-1-BWING-60	Bridge Wing	Black	L/XL
800-1-INSERT-40	Lordotic Inserts	Black	S/M
800-1-INSERT-60	Lordotic Inserts	Black	L/XL
800-JITT	SI Belt Pad Set	Black	Univ.

Sizes Available: S – XL (Waist Measurement)

S/M: Less than 35" (89 cm)

L/XL: 36" - 55" (91-140 cm)

Recommended PDAC Code: L0621



EXOS Form™ 626

The Exos FORM™ 626 that provides relief from mild lower back pain. The Exos FORM™ 626 provides superior sagittal control and support from L-1 to L-5.

- Exos™ Thermoformable Panels
- Independent Superior and Inferior Compression
- Adjustable (Angle) Belt Wings
- Semi-Universal Design
- Simple, Light Weight Low Profile Design.
- Acute and Chronic Lower Back Pain
- Lower Back Strains / Strains
- Lumbar Disc Displacement
- Osteoporosis
- Disc Herniation and Degeneration
- Post-Operative Laminectomy
- Post-Operative Disectomy
- Spondylolisthesis

Part Number	Description	Color	Size
200626-40			S/M
200626-60			L/XL
Accessories			
800-AP-40	Anterior Panel	Black	S/M
800-AP-60	Anterior Panel	Black	L/XL
800-1-WINGS-40	Belt Wing Set	Black	S/M
800-1-WINGS-60	Belt Wing Set	Black	L/XL
800-1-BWING-40	Bridge Wing	Black	S/M
800-1-BWING-60	Bridge Wing	Black	L/XL
800-1-INSERT-40	Lordotic Inserts	Black	S/M
800-1-INSERT-60	Lordotic Inserts	Black	L/XL
800-JITT	SI Belt Pad Set	Black	Univ.

Sizes Available: S – XL (Waist Measurement)

S/M: 25" - 49" (91-125 cm)

L/XL: 49" - 60" (125-152 cm)

Recommended PDAC Code: L0626



EXOS Form™ 627

The Exos FORM™ 627 (with anterior panel) provides relief from a wide range of indications from mild to acute lower back pain, as well as post-operative support. The Exos FORM™ 627 provides superior sagittal control and support from L-1 to L-5.

- Exos™ Thermoformable Panels
- Independent Superior and Inferior Compression
- Adjustable (Angle) Belt Wings
- Semi-Universal Design
- Simple, Light Weight Low Profile Design
- Acute and Chronic Lower Back Pain
- Lower Back Strains / Strains
- Lumbar Disc Displacement
- Osteoporosis
- Disc Herniation and Degeneration
- Post-Operative Laminectomy
- Post-Operative Disectomy
- Spondylolisthesis

Part Number	Description	Color	Size
200627-40			S/M
200627-60			L/XL
Accessories			
800-AP-40	Anterior Panel	Black	S/M
800-AP-60	Anterior Panel	Black	L/XL
800-1-WINGS-40	Belt Wing Set	Black	S/M
800-1-WINGS-60	Belt Wing Set	Black	L/XL
800-1-BWING-40	Bridge Wing	Black	S/M
800-1-BWING-60	Bridge Wing	Black	L/XL
800-1-INSERT-40	Lordotic Inserts	Black	S/M
800-1-INSERT-60	Lordotic Inserts	Black	L/XL
800-JITT	SI Belt Pad Set	Black	Univ.

Sizes Available: S – XL (Waist Measurement)

S/M: 25" - 49" (63-125 cm)

L/XL: 36" - 55" (91-140 cm)

Recommended PDAC Code: L0627



EXOS Form™ 631

The Exos FORM™ 631 provides a higher degree of spinal support and relief from a wide range of indications from acute lower back pain to post-operative support. The Exos FORM™ 631 provides superior sagittal control and support from T9-S1.

- Exos™ Thermoformable Panels
- Independent Superior and Inferior Compression
- Adjustable (Angle) Belt Wings
- Semi-Universal Design
- Simple, Light Weight Low Profile Design
- Acute and Chronic Lower Back Pain
- Lower Back Strains / Strains
- Lumbar Disc Displacement
- Osteoporosis
- Disc Herniation and Degeneration
- Post-Operative Laminectomy
- Post-Operative Disectomy
- Spondylolisthesis

Part Number	Description	Color	Size
200631-40			S/M
200631-60			L/XL
Accessories			
800-AP-40	Anterior Panel	Black	S/M
800-AP-60	Anterior Panel	Black	L/XL
800-1-WINGS-40	Belt Wing Set	Black	S/M
800-1-WINGS-60	Belt Wing Set	Black	L/XL
800-1-BWING-40	Bridge Wing	Black	S/M
800-1-BWING-60	Bridge Wing	Black	L/XL
800-1-INSERT-40	Lordotic Inserts	Black	S/M
800-1-INSERT-60	Lordotic Inserts	Black	L/XL
800-JITT	SI Belt Pad Set	Black	Univ.

Sizes Available: S – XL (Waist Measurement)

S/M: 25" - 49" (63-125 cm)

L/XL: 49" - 60" (91-140 cm)

Recommended PDAC Code: L0631

CRYO-COMPRESSION THERAPY

AIRCAST®

Cryo/Cuff™ IC Cooler

Integrated cold and compression is clinically proven to help reduce post-operative swelling and regain range of motion. The Cryo/Cuff IC provides both gravity and motorized cold and compression therapy for flexibility of care. An integrated pneumatic pump within the cooler lid provides automated cold and compression therapy, with simplicity of design and ease-of-use operation.

Indication: Ideal for post-operative recovery, trauma, athletic training rooms, and home use.

	Description	Size
51A	IC Cooler Only	Univ.
51A10A	IC Cooler with Ankle Cryo/Cuff	Univ.
51A11A	IC Cooler with Knee Cryo/Cuff	M
51A11B	Cryo/Cuff IC Cooler with Knee Cryo/Cuff	L
51A12A	Cryo/Cuff IC Cooler with Shoulder	Univ.
51A12AXL	Cryo/Cuff IC Cooler with Shoulder, XLarge Strap	Univ.
51A14A	Cryo/Cuff IC Cooler with Back/Hip/Rib	Univ.
Accessories		
25-4882	Cryo/Cuff IC Cooler Power Supply	Univ.
25-0238	Cryo/Cuff IC Lid	Univ.
2130	Tube Assembly	Univ.
2002	Insulation Disk	Univ.





AIRCAST®

Cryo/Cuff™ Gravity Cooler

An integrated approach to cold therapy that combines cold and focal compression delivered using a simple, gravity-fed system. Designed with safety and ease-of-use in mind.

Indication: Ideal for post-operative recovery, trauma, athletic training rooms, and home-use.

	Description	Size
2125	Gravity Cooler Only	Univ.
10A	Gravity Cooler w/Ankle Cryo/Cuff	Univ.
10C	Gravity Cooler w/Foot Cryo/Cuff	S/M
10B	Gravity Cooler w/Foot Cryo/Cuff	L
10P	Gravity Cooler w/Pediatric Ankle Cryo/Cuff	Ped.
11C	Gravity Cooler w/Knee Cryo/Cuff	S
11A	Gravity Cooler w/Knee Cryo/Cuff	M
11B	Gravity Cooler w/Knee Cryo/Cuff	L
12A	Gravity Cooler w/Shoulder Cryo/Cuff	Univ.
12AXL	Gravity Cooler w/Shoulder XL Strap Cryo/Cuff	Univ.
14A	Gravity Cooler w/Back/Hip/Rib Cryo/Cuff	Univ.
15A	Gravity Cooler w/Elbow Cryo/Cuff	Univ.
16A	Gravity Cooler w/Hand/Wrist Cryo/Cuff	Univ.



AIRCAST®

Cooler Hanger

For use with any Aircast Cryo/Cuff coolers. Enables cooler to be hung from a bed rail.

	Description
2140	Cooler Hanger



AIRCAST®

Ankle Cryo/Cuff

Anatomically designed to completely fit the ankle. Minimizes hemarthrosis, swelling and pain. Universal adult size. Also available in pediatric size.

Indication: Ideal use for acute sprains, trauma, post-op, and rehabilitation.

	Description	Length x Height	Size
10A01	Ankle Cryo/Cuff Only	11.5" x 10"	Univ.
10P01	Ankle Cryo/Cuff Only	8" x 7"	Ped.



AIRCAST®

Foot Cryo/Cuff

Anatomically designed to provide complete foot coverage for optimal management of swelling and pain.

Indication: Ideal use for turf toe and post-op.

	Description	Circumference*	Size
10C01	Foot Cryo/Cuff Only	9" – 13"	M
10B01	Foot Cryo/Cuff Only	10" – 17"	L

*Measure forefoot circumference.



AIRCAST®

Knee Cryo/Cuff

Anatomically designed to completely fit the knee. Minimizes hemarthrosis, swelling, and pain. Available in three adult sizes. Also available in pediatric size.

Indication: Ideal use for trauma, post-op, rehabilitation, and sports injuries.

	Description	Circumference*	Size
11C01	Knee Cryo/Cuff Only	10" – 19"	S
11A01	Knee Cryo/Cuff Only	18" – 23"	M
11B01	Knee Cryo/Cuff Only	20" – 31"	L
11P01	Knee/Elbow Cryo/Cuff	Pediatric	P

*Measure 6" above the patella.



AIRCAST®

Self-Contained (SC) Knee Cryo/Cuff

A totally self-contained knee cuff that can be filled with ice and water - no cooler required. The SC Knee Cryo/Cuff remains cold for over one hour. A hand bulb is included to add compression as needed.

Indication: Ideal use for trauma, post-op, rehabilitation, and sports injuries.

	Description	Size
11ASC	Self-Contained Knee Cryo/Cuff	Univ.



AIRCAST®

Shoulder Cryo/Cuff

Anatomically designed to cover the shoulder area completely. Minimizes hemarthrosis, swelling, and pain. Can be ordered with an XL strap to accommodate larger torso/chest circumferences.

Indication: Ideal use for trauma, post-op, rehabilitation, and sports injuries.

	Description	Circumference*	Size
12B01	Shoulder Cryo/Cuff Only	Up to 34" and below	S
12A01	Shoulder Cryo/Cuff Only	34"+	Univ.
12AXL01	Shoulder Cryo/Cuff with X-Long Strap Only	42" – 54"	Univ.

Accessories

1200XL	Shoulder Cryo/Cuff XL Strap		Univ.
--------	-----------------------------	--	-------

*Measure torso circumference. (Chest circumference)



AIRCAST®

Calf Cryo/Cuff

Anatomically designed to completely fit the calf. Minimizes swelling and pain. 10" long.

Indication: Ideal use for acute injury, sports injuries, and soft tissue injury.

	Description	Circumference*	Size
13C01	Calf Cryo/Cuff Only	14" – 20"	Univ.

*Measure mid-calf circumference.



AIRCAST®

Thigh Cryo/Cuff

Designed to completely fit the thigh. Minimizes swelling and pain.

Indication: Ideal use for stress fractures, sports injuries, and soft tissue injury.

Description	Circumference*	Size
13A01 Thigh Cryo/Cuff Only	17" – 23"	S/M
13B01 Thigh Cryo/Cuff Only	20" – 27"	L/XL

*Measure thigh circumference.



AIRCAST®

Back/Hip/Rib Cryo/Cuff

Universal fit. The versatile design may be applied either to the back, hip, or rib area. Minimizes swelling and pain.

Indication: Ideal use for chronic/acute pain; post-op and sports injuries.

Description	Size
14A01 Back/Hip/Rib Cryo/Cuff Only	Univ. up to 45" waist



AIRCAST®

Elbow Cryo/Cuff

Anatomically designed to fit the elbow. Minimizes hemarthrosis, swelling, and pain. One size for adults. Also available in a pediatric size.

Indication: Ideal use for epicondylitis (tennis elbow), post-op, and throwing injuries.

Description	Size
15A01 Elbow Cryo/Cuff Only	Univ.
11P01 Elbow/Knee Cryo/Cuff Only	Ped.



AIRCAST®

Hand/Wrist Cryo/Cuff

Anatomically designed to cover the hand and wrist. Versatile design allows repositioning for desired hand coverage. Minimizes swelling and pain. For added support, aluminum support bars are incorporated and can be removed if increased range of motion is desired. 12" long.

Indication: Ideal use for carpal tunnel syndrome, tendonitis, post-op, and sports injuries.

Description	Circumference*	Size
16A01 Hand/Wrist Cryo/Cuff Only	12"	Univ.

*Measure forearm circumference.

DYNAMIC COMPRESSION - VENAFLOW

Vascular Systems

While intermittent pneumatic compression is clinically proven to reduce the risks of deep vein thrombosis, The VenaFlow system is the only device on the market proven to mimic ambulation. Through its unique technology, it has been shown to reduce DVT in comparison to other pneumatic compression systems.^{1, 2, 3, 4}

VENAFLOW®

VenaFlow Elite System

Designed as a prophylaxis for Deep Vein Thrombosis, the VenaFlow Elite System combines two proven technologies, rapid inflation and graduated sequential compression that work to mimic ambulation and accelerate venous velocity. For the superior emptying of veins, VenaFlow Elite also incorporates asymmetric compression. The breathable and comfortable calf, foot and thigh cuffs enhance the patient experience and are compatible with one universal pump.^{1, 2, 3, 4}

Indication: For the prevention of thrombus formation.

	Description
30B	Elite System
30B-B	Elite System with Battery (additional charge)

References

¹Eisele R, Kinzi L, Koelsch T: "Rapid-Inflation Pneumatic Compression For Prevention of Deep Venous Thrombosis." JBJS, 2008.

²Dai G, Gertler GP, Kamm RD: "The Effects of External Compression on Venous Blood Flow and Tissue Deformation in the Lower Leg." Journal of Biomechanical Engineering, 1999.

³Gardner AMN, Fox RH: "The Return of Blood to the Heart, Venous Pumps in Health and Disease." Second Edition, London, John Libbey & Company Ltd, 1993.

⁴Lachiewicz PF, Kelley SS, Haden LR: "Two Mechanical Devices for Prophylaxis of Thromboembolism After Total Knee Arthroplasty." JBJS, 2004.





VENAFLOW®

VenaFlow Elite Calf Cuff (30B, 30B-B Models)

Used exclusively with the VenaFlow Elite System, the VenaFlow Elite calf cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. Available in sterile variation. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3040	Calf Cuff	19"	11.5"	2
3041	Sterile Calf Cuff	19"	11.5"	1
3042	Calf Cuff XL	22"	11.5"	2
3043	Bariatric Calf Cuff	30"	11.5"	2
3044	Sterile Calf Cuff XL	22"	11.5"	1
3050	Foam Calf Cuff	19"	11.5"	2



VENAFLOW®

VenaFlow Elite System Tubing (30B, 30B-B Models)

Designed for use with the VenaFlow Elite System. The tubing connects the VenaFlow Elite System to a VenaFlow cuff and is available in a variety of lengths.

Indication: To be used in conjunction with the VenaFlow Elite System for the prevention of thrombus formation.

	Description	Length
3008	Regular	5.5'
3008XL	X-Long	8.5'
3008XXL	XX-Long	10.5'
3008XXXL	XXX-Long	12.5'



VENAFLOW®

VenaFlow Elite Foot Cuff (30B, 30B-B Models)

Used exclusively with the VenaFlow Elite System, the VenaFlow Elite foot cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3046	Foot Cuff	16"	9"	2



VENAFLOW®

VenaFlow Elite Thigh Cuff (30B, 30B-B Models)

Used exclusively with the VenaFlow Elite System, the VenaFlow Elite thigh cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circumference	Length	Qty.
3045	Thigh Cuff	29"	21"	2



VENAFLOW®

VenaFlow Elite System Repair Kits (30B, 30B-B Models)

Repair kits for use with the VenaFlow Elite tubing system.

Indication: To be used in conjunction with the VenaFlow Elite System for the prevention of thrombus formation.

	Description
3070	Power cord (U.S.)
3071	Bed hanger
3072	Tube attachment tag
3073	Fuse



VENAFLOW®

VenaFlow® System

Designed as a prophylaxis for Deep Vein Thrombosis, the VenaFlow System combines two proven technologies, rapid inflation and graduated sequential compression, to accelerate venous velocity and enhance fibrinolysis. For the superior emptying of veins, VenaFlow also incorporates asymmetric compression. The breathable and comfortable calf, foot and thigh cuffs enhance the patient experience and are compatible with one universal pump.

Indication: For the prevention of thrombus formation.

	Description
30A	VenaFlow System



VENAFLOW®

VenaFlow Foot Cuff (30A Model)

Used exclusively with the VenaFlow System, the VenaFlow cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3016-PL	Foot Cuff	16"	9"	2



VENAFLOW®

VenaFlow Calf Cuff (30A Model)

Used exclusively with the VenaFlow System, the VenaFlow cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. Available in sterile variation. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3010-PL	Calf Cuff	19"	11.5"	2
3011-PL	Sterile Calf Cuff	19"	11.5"	1
3012-PL	Calf Cuff XL	22"	11.5"	2
3013-PL	Bariatric Calf Cuff	30"	11.5"	2
3014-PL	Sterile Calf Cuff XL	22"	11.5"	1
3030-PL	Foam Calf Cuff	19"	11.5"	2



VENAFLOW®

VenaFlow Thigh Cuff (30A Model)

Used exclusively with the VenaFlow System, the VenaFlow cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3015-PL	Thigh Cuff	29"	21"	2



VENAFLOW®

VenaFlow System Repair Kits (30A Model)

Repair kits for use with the VenaFlow system.

Indication: To be used in conjunction with the VenaFlow System for the prevention of thrombus formation.

	Description
11-0201	Tube Storage Strap Repair Kit
11-0202	Power Cord Repair Kit
11-0203	Pump/Tubing Connector Repair Kit
11-0204	Pump Case (Back Half) Repair Kit
11-0205	Bed Hanger Foam Repair Kit
11-0206	Fuse Repair Kit



VENAFLOW®

VenaFlow System Tubing (30A Model)

Designed for use with the VenaFlow System. The tubing connects the VenaFlow System to a VenaFlow cuff and is available in a variety of lengths.

Indication: To be used in conjunction with the VenaFlow System for the prevention of thrombus formation.

	Description	Length
3007-PL	Regular	5.5'
3007XL-PL	X-Long	8.5'
3007XXL-PL	XX-Long	10.5'

DYNAMIC COMPRESSION - ARTERIALFLOW

ArterialFlow®

ArterialFlow® System

Designed to augment arterial blood flow and microcirculation, the Aircast ArterialFlow System is clinically proven to accelerate healing of ulcerative wounds and preserve limbs. ArterialFlow combines rapid inflation and graduated, sequential compression at high pressures and frequent inflation cycles, which works to empty the veins and reduce venous resistance. Decreasing venous resistance promotes unrestricted flow across capillaries, increasing blood to the extremities where healing can begin.^{5, 6}

Indication: The ArterialFlow System is intended for the improvement of blood circulation in the legs and arms to help prevent and reduce complications or poor circulation.



Description	
32A	ArterialFlow System®
*Rental program now available through DJO. For more information contact your local representative. ArterialFlow system and cuffs are PDAC approved.	

References

⁵Kavros SJ, Delis KT, Turner NS, et al: "Improving Limb Salvage in Critical Ischemia With Intermittent Pneumatic Compression: A Controlled Study With 18-Month Follow-up." *Journal of Vascular Surgery*, 2008.

⁶Montori VM, Kavros SJ, Walsh EE, et al: "Intermittent Compression Pump for Nonhealing Wounds in Patients With Limb Ischemia. The Mayo Clinic Experience." *International Angiology*, 2002.



ArterialFlow®

ArterialFlow Upper Extremity Cuff

The ArterialFlow cuffs are for use exclusively with the ArterialFlow System and are designed for multiple uses per patient. The cuffs are latex-free and can be placed directly against the skin.

Indication: For the improvement of blood circulation in the arms.

Description	
3225	ArterialFlow Upper Extremity Cuff
Accessories	
3204	ArterialFlow Upper Extremity Cover
3205	ArterialFlow Upper Extremity Aircell
3226	ArterialFlow Single Limb Operation Pad



ArterialFlow®

ArterialFlow Lower Extremity Cuff

The ArterialFlow cuffs are for use exclusively with the ArterialFlow System and are designed for multiple uses per patient. The cuffs are latex-free and can be placed directly against the skin.

Indication: For the improvement of blood circulation in the legs.

Description	
3227	ArterialFlow Lower Extremity Cuff
3229	ArterialFlow Lower Extremity Cuff XL
Accessories	
3217	ArterialFlow Lower Extremity Cover
3228	ArterialFlow Lower Extremity XL Cover
3218	ArterialFlow Lower Extremity Aircell
3226	ArterialFlow Single Limb Operation Pad



ArterialFlow®

ArterialFlow Tube Assembly

Designed for use with the ArterialFlow System. The tubing connects the ArterialFlow pump to the ArterialFlow cuffs and is available in three lengths.

Indication: To be used in conjunction with the ArterialFlow System for the improvement of blood flow and microcirculation.

	Description	Length
3007-PL	Regular	5.5'
3007XL-PL	X-Long	8.5'
3007XXL-PL	XX-Long	10.5'

STATIC COMPRESSION- ANTI EMBOLISM STOCKINGS

VENAFLOW[®]AES

VenaFlow Knee High Anti Embolism Stockings

Designed with quality, fit and safety in mind, the VenaFlow Anti Embolism stockings are key in preventing VTE. They enhance venous flow velocity while keeping patients comfortable during and after their hospital stay.

Indication: For the prevention of thrombus formation.

	Description	Qty.	Size
Regular [®]			
500R-2	Knee High Stocking	12 pair/pk.	S
500R-3	Knee High Stocking	12 pair/pk.	M
500R-4	Knee High Stocking	12 pair/pk.	L
500R-5	Knee High Stocking	12 pair/pk.	XL
500R-6	Knee High Stocking	12 pair/pk.	XXL
500R-7	Knee High Stocking	12 pair/pk.	XXXL
Long [®]			
500L-2	Knee High Stocking	12 pair/pk.	S
500L-3	Knee High Stocking	12 pair/pk.	M
500L-4	Knee High Stocking	12 pair/pk.	L
500L-5	Knee High Stocking	12 pair/pk.	XL
500L-6	Knee High Stocking	12 pair/pk.	XXL
500L-7	Knee High Stocking	12 pair/pk.	XXXL

[®]All stockings provide compression in mmHg: Compression at ankle 18 mmHg, Compression at calf 14 mmHg, Compression at knee 8 mmHg, Compression at lower thigh 10 mmHg, Compression at upper thigh 8 mmHg.





VENAFLOW[®]AES

VenaFlow Thigh High Anti Embolism Stockings

Designed with quality, fit and safety in mind, the VenaFlow Anti Embolism stockings are key in preventing VTE. They enhance venous flow velocity while keeping patients comfortable during and after their hospital stay.

Indication: For the prevention of thrombus formation.

Thigh High	Thigh High w/Belt	Qty.	Size
Regular[®]			
501R-2	502R-2	6 pair/pk.	S
501R-3	502R-3	6 pair/pk.	M
501R-4	502R-4	6 pair/pk.	L
501R-5	502R-5	6 pair/pk.	XL
501R-6	502R-6	6 pair/pk.	XXL
Long[®]			
501L-2	502L-2	6 pair/pk.	S
501L-3	502L-3	6 pair/pk.	M
501L-4	502L-4	6 pair/pk.	L
501L-5	502L-5	6 pair/pk.	XL
501L-6	502L-6	6 pair/pk.	XXL

[®]All stockings provide compression in mmHg: Compression at ankle 18 mmHg, Compression at calf 14 mmHg, Compression at knee 8 mmHg, Compression at lower thigh 10 mmHg, Compression at upper thigh 8 mmHg.

Thigh High	Thigh High w/Belt	Qty.	Size
Short[®]			
501S-2	502S-2	6 pair/pk.	S
501S-3	502S-3	6 pair/pk.	M
501S-4	502S-4	6 pair/pk.	L
501S-5	502S-5	6 pair/pk.	XL
501S-6	502S-6	6 pair/pk.	XXL

[®]All stockings provide compression in mmHg: Compression at ankle 18 mmHg, Compression at calf 14 mmHg, Compression at knee 8 mmHg, Compression at lower thigh 10 mmHg, Compression at upper thigh 8 mmHg.

VENAFLOW ANTI EMBOLISM KNEE LENGTH STOCKINGS SIZING CHART [®]					
	Size	1 - Calf Circumference	2 - Length		Color
			Regular	Long	
	S	< 30.5 cm / 12"	< 41 cm / < 16"	> 41 cm / < 16"	Yellow
M	30.5 – 38 cm / 12" – 15"	< 43 cm / < 17"	> 43 cm / < 17"	White	
L	38 – 44.5 cm / 15" – 17.5"	< 46 cm / < 18"	> 46 cm / < 18"	Blue	
XL	44.5 – 51 cm / 17.5" – 20"	< 46 cm / < 18"	> 46 cm / < 18"	Green	
XXL	51 – 58.4 cm / 20" – 23"	< 46 cm / < 18"	> 46 cm / < 18"	Purple	
XXXL	58.5 – 66 cm / 23" – 26"	< 46 cm / < 18"	> 46 cm / < 18"	Red	

[®]Knee Length Sizing: 1 - Measure calf circumference, 2 - Measure length.

VENAFLOW ANTI EMBOLISM THIGH LENGTH STOCKINGS SIZING CHART [®]							
	Size	1 - Thigh Circumference	2 - Calf Circumference	3 - Length			Color
				Short	Regular	Long	
	S	< 63.5 cm / 25"	< 30.5 cm / < 12"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Yellow
M	< 63.5 cm / 25"	30.5 – 38 cm / 12" – 15"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	White	
L	< 63.5 cm / 25"	38 – 44.5 cm / 15" – 17.5"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Blue	
XL	63.5 – 81.5 cm / 25" – 31.8"	44.5 – 55 cm / 17.5" – 21.6"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Green	
XXL	81 – 91 cm / 31.8" – 35.8"	55 – 66 cm / 21.6" – 26"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Purple	

[®]Thigh High Sizing: 1 - Measure Thigh Circumference, 2 - Measure Calf Circumference, 3 - Measure Length.

CERVICAL COLLARS

PROCARE®

ProCare® Transitional 172 Cervical Collar

The ProCare® Transitional 172 Cervical Collars were designed to provide optimal support and comfort. The collar provides substantial motion restriction while minimizing painful pressure points that can lead to skin maceration and/or breakdown. Enhanced comfort using soft, breathable foam padding helps promote and encourage patient compliance. Removable open cell foam pads wick moisture away from the skin to ensure proper patient hygiene and comfort while minimizing skin irritation. X-ray, CT and MRI lucent.

Indication: Restricts cervical spine flexion, extension and rotation to promote patient recovery.

Adult Collar Only	Adult Collar w/ Repl. Pads	Height	Neck Circ.	Size
Adult Collars				
79-83283	79-83363	2.25"	13" – 22"	Adult Short
79-83285	79-83365	3.0"	13" – 22"	Adult Regular
79-83287	79-83367	3.75"	13" – 22"	Adult Tall
79-83288	79-83368	4.5"	13" – 22"	Adult X-Tall
Adult Collar Accessories				
	Description			Size
79-83325	ACC Transitional 172 Collar Pad Set			Adult Universal
79-83383	ACC Transitional 172 Posterior Panel			Adult Small
79-83385	ACC Transitional 172 Posterior Panel			Adult Standard
79-83387	ACC Transitional 172 Posterior Panel			Adult Large
79-83393	ACC Transitional 172 Posterior Panel Pad			Adult Small
79-83395	ACC Transitional 172 Posterior Panel Pad			Adult Standard
79-83397	ACC Transitional 172 Posterior Panel Pad			Adult Large

Pediatric Collar Only	Pediatric Collar w/ Repl. Pads	Age	Height	Neck Circ.	Size
Pediatric Collars					
79-83272	79-83351	1-18 mo.	0.75"	8" – 12"	PED-1, Infant Short
79-83273	79-83352	9-24 mo.	1.25"	11" – 13"	PED-2, Infant
79-83275	79-83353	1-3 yrs.	2.0"	11.5" – 14"	PED-3, Child Toddler
79-83277	79-83354	2-5 yrs.	2.25"	12" – 14.5"	PED-4, Child Short
79-83278	79-83355	3-6 yrs.	2.5"	13" – 16"	PED-5, Child Regular
Pediatric Collar Accessories					
	Description				Size
	ACC Transitional 172 Collar Pad Set, PED-1/PED-2				PED-1/2, Infant
79-83295	ACC Transitional 172 Collar Pad Set, PED-3/PED-4/PED-5				PED-3/4/5, Child
79-83214	ACC Transitional 172 Posterior Panel, PED-1/PED-2				PED-1/2, Infant
79-83215	ACC Transitional 172 Posterior Panel, PED-3/PED-4/PED-5				PED-3/4/5, Child





PROCARE®

Aspen® Vista® TX Collar

Vista TX encompasses all of the features and benefits of the Vista collar, but with a thoracic extension. Removable, foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation. x-ray, CT and MRI lucent.

Indication: Motion control of cervical spine flexion, extension and rotation to promote patient recovery.

Collar Only	Collar w/Repl. Pads	Neck Circum.	Size
79-83390	79-83360	13" - 19"	Adult/Univ.
Accessories			
79-83372	Vista TX Repl. Pad Set	Adult	
79-83377	Back Only	Adult	



PROCARE®

Aspen® Vista® Colla

Innovative height adjustment technology offers six (6) collars in one. Two-piece polyethylene design enhances patient comfort and is easily adjusted with Velcro™ tabs. Anterior and posterior flex tabs on the front and back panels provide improved rigidity and durability. The large trachea opening provides functional access for the airway management and other essential procedures. Removable, foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation. x-ray, CT and MRI lucent.

Indication: Motion control of cervical spine flexion, extension and rotation to promote patient recovery.

Collar Only	Collar w/Repl. Pads	Neck Circum.	Size
79-83370	79-83380	13" - 19"	Adult/Univ.
Accessories			
79-83371	Vista Replacement Pad Set		Adult
79-83377	Back Only		Adult



PROCARE®

Universal Cervical Collar

Medium density foam collar with vinyl stabilizing panel helps support cervical spine in neutral position. Available in 2", 3" and 4" heights to fit a wide range of patient neck sizes. Cotton stockinette cover for better patient hygiene. Latex free.

Indication: Ideal for post-operative rehabilitation or cervical sprains and strains.

	Neck Circum.	Height*	Size
79-83030	10.5" - 22"	2"	Univ.
79-83000	10.5" - 22"	3"	Univ.
79-83010	10.5" - 22"	4"	Univ.

*Height = measurement from sternal notch to chin.



PROCARE®

Universal Clinic Collar

Medium density foam collar provides consistent fit and support. Fits a wide range of neck sizes. Additional stockinette cover provided for better patient hygiene. Latex free.

Indication: Ideal for supporting the cervical spine in neutral position.

	Height	Length	Size
79-83520	2.5"	24"	Univ.
79-83500	3"	24"	Univ.
79-83510	4"	24"	Univ.



PROCARE®

Form Fit Cervical Collar

Firm or medium density foam collar helps support cervical spine in neutral position. Contour design for use in either flexed or extended position. Cotton stockinette cover for better patient hygiene. Sold as each or in a pack of six (6). Latex free.

Indication: Ideal for post-operative rehabilitation; cervical sprains and strains.

Firm Density	Med. Density	Med. Dens. (6-pk.)	H [*] x L	Neck Circum.	Size
	79-83011		2" x 14"	7" – 12"	Ped./2XS
	79-83012		2" x 16.5"	9" – 13"	XS
79-83003	79-83013	79-83043	3" x 18.5"	11" – 16"	S
79-83004	79-83014	79-83044	3" x 22"	15" – 20"	S/Long
79-83005	79-83015	79-83045	4" x 20"	13" – 18"	M
79-83006	79-83016	79-83046	4" x 23"	16" – 21"	M/Long
79-83007	79-83017	79-83047	4.5" x 22.5"	15" – 20"	L
79-83008	79-83018	79-83048	4.5" x 24"	16" – 22"	XL
79-83009	79-83019		2" x 24"	17" – 22"	Long/Nar.

^{*}Height = measurement from sternal notch to chin.



PROCARE®

Serpentine Cervical Collar

Firmer support than Form Fit styles. Contour design helps maintain cervical spine in a neutral position and restricts motion. Cotton stockinette cover for better patient hygiene. Latex free.

Indication: Ideal for post-operative rehabilitation; cervical sprains and strains.

	Height [*]	Length	Size
Medium Density			
79-83073	3"	16"	S
79-83075	4"	16"	M
79-83078	4"	20"	L/Stout
Firm Density			
79-83093	3.25"	19"	S
79-83095	3.5"	21"	M
79-83097	4.5"	23.5"	L
79-83098	4"	20"	L/Stout

^{*}Height = measurement from sternal notch to chin.



PROCARE®

Clinic Low Density Collar

Special contour design collar with low density molded foam. Sold each or in a pack of six (6). Cotton stockinette cover for better patient hygiene. Latex free.

Indication: Ideal for immobilization of the neck.

6-Pack	Each	Height [*]	Length	Size
79-83503	79-83503S	3"	17"	S
79-83505	79-83505S	3.25"	19"	M
79-83507	79-83507S	3.75"	21.5"	L
79-83508	79-83508S	4.25"	23.5"	XL

^{*}Height = measurement from sternal notch to chin.



PROCARE®

Low Contour Cervical Collar

The special low, medium density contoured design offers a better fit for patients with short necks. Cotton stockinette cover for better patient hygiene. Latex free.

Indication: Ideal for post-operative rehabilitation; cervical sprains and strains.

	Height*	Length	Size
79-83303	2.5"	18.5"	S
79-83305	3"	20"	M
79-83307	3.25"	23"	L

*Height = measurement from sternal notch to chin.



PROCARE®

California Tracheotomy Collar

Pre-formed molded foam tracheotomy collar provides support for both extrication and rehabilitation. Two-piece polyethylene foam lamination. Collar maintains head in neutral position and supports occipital and mandibular areas. Available in four (4) heights. Collar can be incorporated into a plaster body jacket.

Indication: Ideal for cervical spine immobilization and extraction.

Height*				Size
2.25"	3.25"	4.25"	5.25"	Ped./XS
79-83122				
79-83123	79-83133	79-83143	79-83153	S
79-83125	79-83135	79-83145	79-83155	M
79-83127	79-83137	79-83147	79-83157	L

XS = 7" – 10"; S = 10" – 13"; M = 13" – 16"; L = 16" – 19".

*Height = measurement from sternal notch to chin.



PROCARE®

California Cervical Collar

Pre-formed, lightweight, two-piece polyethylene foam with rigid plastic supports front and rear. Collar maintains head in neutral position. May be used to provide support in combination with a cervical or scoliosis cast.

Indication: Ideal for cervical spine immobilization and extraction.

Height*				Size
2.25"	3.25"	4.25"	5.25"	
79-83161				Ped./XS
79-83163	79-83183	79-83193	79-83243	S
79-83165	79-83185	79-83195	79-83245	M
79-83167	79-83187	79-83197	79-83247	L

XS = 7" – 10"; S = 10" – 13"; M = 13" – 16"; L = 16" – 19".

*Height = measurement from sternal notch to chin.



PROCARE®

One-Piece Flat Collar

Contoured plastic/foam design provides immobilization for extrication and rehabilitation. One piece construction with large trachea opening stacks flat to preserve space. Color coded closure straps for quick size reference. CT and MRI compatible.

Indication: Ideal for cervical spine immobilization and extraction.

	Color	Height*	Size
79-83021	Pink	1.5"	Infant
79-83022	Aqua	1.75"	Pediatric
79-83028	Purple	1.75"	Low
79-83023	Blue	2.5"	Short
79-83025	Orange	3"	Regular
79-83027	Green	3.5"	Tall

*Height = measurement from sternal notch to chin.



PROCARE®

Patriot® Collar

One-piece injection molded collar provides proper support and helps restrict cervical spine flexion, extension and rotation to promote patient recovery following trauma. Constructed with durable high density polyethylene lined with a soft closed cell foam for maximum patient comfort. Both universal sizes have multiple heights to fit most adults and pediatrics. Oversized trachea opening provides quick access to carotid pulse monitoring and emergency tracheotomies. Stacks flat to preserve space and is x-ray, CT and MRI lucent. Latex free.

Indication: Ideal for cervical spine immobilization and extraction.

	Height Adjustment*	Neck Circum.	Size
79-83130	1.5", 1.75", 2"	8" – 18"	Pediatric
79-83131	1.75", 2.5", 3", 3.5"	11" – 23"	Adult

*Height = measurement from sternal notch to chin.

PROCARE®

Miami J® Collar

Designed to meet the needs of the adult ambulatory patient. Padding helps retard skin breakdown.

Indication: Ideal for motion control of cervical spine flexion, extension and rotation to promote patient recovery.

		Height*	Neck Circum.	Size
79-83201	MJ 200S	1"	10" – 20"	Super Short
79-83203	MJ 300	2"	10" – 20"	S
79-83205	MJ 400	2.5"	10" – 20"	M
79-83207	MJ 500	3"	10" – 20"	L
79-83208	MJ 200L	1.75"	14" – 24"	Stout
Standard Collar With Replacement Pad				
79-83221	MJR200S	1"	10" – 20"	Super Short
79-83232	MJR 250	1.25"	8" – 14"	XS
79-83233	MJR 300	2"	10" – 20"	S
79-83235	MJR 400	2.5"	0" – 20"	M
79-83237	MJR 500	3"	10" – 20"	L
79-83238	MJR 200L	1.75"	14" – 24"	Stout
Replacement Pads Only				
79-83224	MJP 100			Adult
79-83225	MJP 200L			Stout

*Height = measurement from sternal notch to chin.



PROCARE®

Miami Jr.® Collar

Designed to meet anatomical needs of children up to 12 years old. Padding helps retard skin breakdown.

Indication: Ideal for motion control of cervical spine flexion, extension and rotation to promote patient recovery.

	Age	Height*	Neck Circum.
Ped 0			
79-83281	MJ P0	0-6 mo.	
79-83206	MJR P0 w/repl. Pad	0-6 mo.	
79-83209	Repl. Pad Set Only		
Ped 1			
79-83251	MJ P1	6 mo. – 2 yrs.	1.125" ≤ 9"
79-83227	MJR P1 w/repl. pad	6 mo. – 2 yrs.	1.125" ≤ 9"
79-83253	Repl. Pad Set Only		
Ped 2			
79-83252	MJ P2	2 – 6 yrs.	1.25" 8" – 11"
79-83228	MJR P2 w/repl. pad	2 – 6 years	1.25" 8" – 11"
79-83254	Repl. Pad Set Only		
Ped 3			
79-83202	MJ P3	6 – 12 yrs.	1.75" 9" – 13"
79-83222	MJR P3 w/ repl. pad	6 – 12 yrs.	1.75" 9" – 13"

*Height = measurement from sternal notch to chin.



PROCARE®

911 First Response® Collar

Rigid radiolucent foam laminate controls flexion and extension of the neck and ensures support in lateral bend. Large tracheal hole allows for emergency patient monitoring. Features additional hook and loop strip that allows for spine board attachment. 10 per package. Disposable spine board straps are ideal for use by EMS personnel for spine board restraint. Constructed of nylon and feature easy release buckles.

Indication: Ideal for use by EMS personnel for spine board restraint.

	Neck Circum.	Height*	Qty.	Size
79-83911	20"	3.25"	10/pk.	XXS
79-83912	24.5"	3.25"	10/pk.	XS
79-83913	24.5"	4.5"	10/pk.	S
79-83915	24.5"	5.25"	10/pk.	M
79-83917	24.5"	6"	10/pk.	L

Accessories

Disposable Spine Board Straps

79-83750	5 foot, 2 piece	12/pk.
79-83760	9 foot, 1 piece	12/pk.

*Height = measurement from sternal notch to chin.

CLAVICLE SUPPORTS

PROCARE®**Clavicle Splint**

Adjustable prong buckles ensure correct fit and prevents slippage. Stockinette-covered felt shoulder straps. Washable.

Indication: Ideal for clavicular fractures and postural problems.

	Circum.*	Size
79-85001	18" – 20"	Ped./2XS
79-85002	20" – 24"	XS
79-85003	24" – 30"	S
79-85005	30" – 36"	M
79-85007	36" – 42"	L
79-85008	42" – 48"	XL
79-85009	48" – 54"	2XL

*Chest circumference measurement taken at bottom of sternum.





PROCARE®

Contoured Clavicle

Contoured, preformed straps for maximum comfort. Constructed of bonded polyurethane foam and covered with cotton stockinette. Figure 8 design with 2-way contact closure.

Indication: Ideal for clavicular fractures and postural problems.

Contact Closure	Circum.®	Size
79-98671	16" – 18"	Ped./3XS
79-98672	18" – 20"	2XS
79-98673	20" – 24"	XS
79-98674	24" – 30"	S
79-98675	30" – 36"	M
79-98677	36" – 42"	L
79-98678	42" – 48"	XL
79-98679	48" – 54"	2XL
79-85028	42" – 48"	XL

®Chest circumference measurement taken at bottom of sternum.



PROCARE®

Heavy Padded Clavicle

Combination felt and foam padded straps for added comfort. Contact closure for easier adjustment.

Indication: Ideal for clavicular fractures and postural problems.

	Circum.®	Size
79-85021	18" – 20"	Ped./2XS
79-85022	20" – 24"	XS
79-85023	24" – 30"	S
79-85025	30" – 36"	M
79-85027	36" – 42"	L
79-85028	42" – 48"	XL

®Chest circumference measurement taken at bottom of sternum.



PROCARE®

Clavicle Strap

Foam padded shoulder straps covered with stockinette for added patient comfort. 3-way adjustment buckle closure helps prevent strap slippage.

Indication: Ideal for clavicular fractures and postural problems.

	Circum.®	Size
79-85012	20" – 24"	XS
79-85013	24" – 30"	S
79-85015	30" – 36"	M
79-85017	36" – 42"	L
79-85018	42" – 48"	XL

®Chest circumference measurement taken at bottom of sternum.



PROCARE®

Universal Clavicle

Foam/nylon construction with contact closure for easy adjustment. Movable D-rings allow for proper adjustment.

Indication: Ideal for clavicular fractures and postural problems.

	Circum.®	Size
79-85100	Up to 40"	Univ.

®Chest circumference measurement taken at bottom of sternum.

SHOULDER IMMOBILIZERS

AIRCAST®**Quick-Fit® Shoulder Immobilizer**

Universal design fits a variety of patients providing a custom fit, easy to use option for shoulder protection post-trauma. Pillow offers 15° of abduction to help prevent internal rotation contractures. Arm sling is constructed with cool mesh material allowing for proper air flow and fully padded waist and shoulder straps offer greater patient comfort.

Indication: Provides Support for Post-operative rotator cuff repairs, Postoperative Bankart procedures, Capsular shifts, Glenohumeral, dislocation/subluxation, Soft tissue strains/repairs.

		Arm Length*	Depth	Size
06AB		10.5" - 17"	7"	Univ.
06A	Sling Only	10.5" - 17"	7"	Univ.

*Measure from elbow crease to base of index finger.





PROCARE®

Shoulder Immobilizer w/Foam Straps

Sling is constructed with soft cotton/poly for patient comfort. Designed to exert upward pressure and immobilize arm and shoulder. Waist and shoulder straps are a plush foam laminate with contact closure. Fits right or left arm.

Indication: Ideal for immobilization and support of the shoulder and elbow joints.

	Envelope Size Depth x Length	Size
79-84162	6" x 11"	XS
79-84163	7" x 13"	S
79-84165	8" x 15"	M
79-84167	9" x 18"	L
79-84168	10" x 20"	XL



PROCARE®

Clinic Shoulder Immobilizer

High quality, durable poly cotton sling with contact closure adjustment on shoulder and body straps. Fits right or left arm and helps prevent shoulder rotation. Shoulder pad provides additional comfort.

Indication: Ideal for immobilization and support of the shoulder and elbow joints.

	Envelope Size Depth x Length	Size
79-84011	5" x 9"	2XS
79-84012	6" x 11"	XS
79-84013	7" x 13"	S
79-84015	8" x 15"	M
79-84017	9" x 18"	L
79-84018	10" x 20"	XL



PROCARE®

Deluxe Shoulder Immobilizer

Constructed of fiber laminate for patient comfort. Adjustable wrist cuff and shoulder strap may be positioned as needed for proper control and custom fit. Adjustable humeral cuff also available on one style. Available in child or adult models. Fits right or left and may be worn without shoulder strap.

Indication: Ideal for treatment of traumatic or post-surgical immobilization.

		Waist Measurement	Size
79-84100		≤ 52"	Univ.
79-84108		52" – 62"	Univ./XL
79-96820	w/detachable humeral strap	≤ 52"	Univ.
79-84100C	Contoured	≤ 52"	Univ.
79-84101	Child	≤ 38.5"	XS



AIRCAST®

Arm Immobilizer

Engineered with innovative features designed to provide secure and comfortable support. Enhanced design allows two options to secure underarm strap to maximize patient comfort. Available with or without abduction pillow. Pillow provides 15° of abduction to help prevent internal rotation and immobilize shoulder and arm as required.

Indication: Arm Support for a variety of pre- and post-op indications involving the shoulder, elbow, proximal or mid-humerus and wrist.

		Arm Length*	Size
06GMA		Up to 10"	M
06GLA		9.5" +	L
06GM	Sling Only	Up to 10.5"	M
06GL	Sling Only	10" +	L

*Measure from elbow to ulnar styloid.



PROCARE®

Elastic Shoulder Immobilizer

Elastic waist band construction for compressive support. Durable yet comfortable foam wrist and nylon humeral cuff with contact closure for ease of application. Fits right or left arm.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

		Rib Measurement	Size
79-84032	Male	24" – 32"	XS
79-84033	Male	32" – 36"	S
79-84035	Male	36" – 40"	M
79-84037	Male	40" – 44"	L
79-84038	Male	44" – 48"	XL
79-84039	Male	48" – 52"	2XL
79-84043	Female	24" – 28"	S
79-84045	Female	28" – 32"	M
79-84047	Female	32" – 36"	L
79-84048	Female	36" – 42"	XL
79-84049	Female	42" – 48"	2XL



PROCARE®

Deluxe Sling and Swathe

Cool, lightweight pile/foam construction provides patient comfort during rehabilitation following surgery or dislocation. Two over-the-shoulder criss-cross straps with a slide buckle for proper strap adjustment. Shoulder straps attach to the sling with contact closure to ensure proper weight distribution and ease of use. One size fits up to 54" waist. Fits right and left.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

		Size
79-84230		Univ.
79-99600	Swathe only	54" length



PROCARE®

Quick Release Economy Shoulder Immobilizer

Constructed of breathable Tietex mesh material. Buckle closure. Thumb loop prevents wrist drop. Sized for precise fit.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

	Envelope Length	Waist Strap Length	Size
79-84342	12"	32"	XS
79-84343	14"	34.5"	S
79-84345	15.5"	37.5"	M
79-84347	18"	41"	L
79-84348	20.5"	44.5"	XL



PROCARE®

Sling and Swathe

Sling and swathe constructed with durable cotton/duck material. Slide buckle closure on soft cotton woven shoulder strap. Available in four (4) sizes.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

	Swathe Size	Envelope	Size
79-84243	4" x 44"	9" x 14"	S
79-84245	4" x 49"	10" x 16"	M
79-84247	4" x 54"	11" x 17"	L
79-84248	4" x 59"	12" x 19"	XL



PROCARE®

Universal Mesh Shoulder Immobilizer

Constructed of open mesh knit with Tec-Clips™ for proper one-time size adjustment. Universal sizing to fit most patients. Hook and loop closures on shoulder and waist straps. Maintains glenohumeral joint, forearm and wrist in neutral position.

Indication: Ideal for immobilization; S/P glenohumeral acromioplasty; rotator cuff repair.

	Qty.	Size
79-84350	14/pk.	Univ.



PROCARE®

Shoulder Abduction Kit

Allows the shoulder and arm to remain in a neutral position. Quality foam/nylon construction, universal design to fit most patients. Fits right or left.

Indication: Provides abduction where required (45° or 60°) to provide post-injury or post-surgical support.

	Size
79-84500	Univ.



PROCARE®

Hot/Cold Therapy Wrap

Safe, non-toxic pack is microwave and freezer safe. Ideal for use on shoulder, wrist, elbow, back, knee and ankle.

Indication: Therapeutic wrap for muscle strains and sprains.

	Size
79-97210	Univ.
Replacement Hot/Cold Pack Only	
79-97220	

ARM SLINGS AND ELEVATORS

PROCARE®**Blue Vogue Arm Sling**

Generous size envelope constructed of durable poly/cotton provides patient with comfortable, custom fit yet helps ensure immobility. Fully padded 2" wide foam strap with contact closure and D-rings for easy adjustment and application. Thumb loop helps prevent sling from migrating.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-84153	7.5"	12.5"	S
79-84155	8"	15"	M
79-84157	8.5"	17.5"	L
79-84158	9"	20"	XL

*Measure length from elbow crease to base of index finger.





PROCARE®

Deluxe Arm Sling w/Pad

Generous size envelope to support the arm. Cotton/poly construction with web strap and contact closure. Thumb loop helps prevent sling from migrating. Shoulder pad included.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-84001	6.5"	8.25"	Pediatric/2XS
79-84002	7"	11"	Pediatric/XS
79-84003	7.5"	12.5"	S
79-84005	8"	15"	M
79-84007	8.5"	17.5"	L
79-84008	9"	20"	XL

*Measure length from elbow crease to base of index finger.



PROCARE®

Specialty Arm Sling

Constructed with cool mesh material and double D-rings. Hook and loop closure. Thumb loop prevents wrist drop.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	With Padded Strap	Depth	Length*	Size
79-84272	79-84262	6.5"	13"	Pediatric/XS
79-84273	79-84263	7"	16"	S
79-84275	79-84265	7.5"	17"	M
79-84277	79-84267	8"	19.5"	L
79-84278	79-84268	8.5"	21.5"	XL

*Measure length from elbow crease to base of index finger.



PROCARE®

Super Arm Sling

Durable cotton/poly sling with double D-rings allows for an open envelope and ease of application. Comfortable thumb loop helps prevent sling from migrating. Web shoulder strap with contact closure.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-84112	7"	11"	Pediatric/XS
79-84113	7.5"	12.5"	S
79-84115	8"	15"	M
79-84117	8.5"	17.5"	L
79-84118	9"	19"	XL

*Measure length from elbow crease to base of index finger.



PROCARE®

Clinic Arm Sling

Comfortable cotton/poly envelope with web shoulder strap and slide buckle closure. Packaged 6 per box.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Qty.	Size
79-84021	4.25"	10"	6/pk.	Pediatric/2XS
79-84022	5.25"	12.5"	6/pk.	Pediatric/XS
79-84023	5.5"	14.5"	6/pk.	S
79-84025	6.5"	16.5"	6/pk.	M
79-84027	7.5"	18.5"	6/pk.	L
79-84028	8.5"	22.5"	6/pk.	XL

*Measure length from elbow crease to base of index finger.



PROCARE®

Universal Arm Sling

Contact closure adjustment on envelope accommodates sizing for most patients. Durable cotton/poly material, double O-rings, and contact closure shoulder strap for proper adjustment. Available with or without padded strap.

Indication: Ideal for immobilization and support of the shoulder, arm and elbow joints.

	Depth	Length	Size
79-84300	7"	18"	Univ.
79-92070	w/padded strap 7"	18"	Univ.

**Measure length from elbow crease to base of index finger.*



PROCARE®

Quick Release Arm Sling

Breathable Tietex material. Deep pocket design to accommodate larger casts. Buckle closure strap for ease of application and removal.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-84292	6.5"	13"	XS
79-84293	7"	16"	S
79-84295	7.5"	17"	M
79-84297	8"	19.5"	L
79-84298	8.5"	21.5"	XL

**Measure length from elbow crease to base of index finger.*



PROCARE®

Chieftain™ Arm Sling

Breathable Tietex mesh material with hook and loop closure. Thumb loop prevents wrist drop.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Qty.	Size
79-84172	6.5"	13"	12/pk.	XS
79-84173	7"	16"	12/pk.	S
79-84175	7.5"	17"	12/pk.	M
79-84177	8"	19.5"	12/pk.	L
79-84178	8.5"	21.5"	12/pk.	XL

**Measure length from elbow crease to base of index finger.*



PROCARE®

Deep Pocket Economy Arm Sling

Generous size poly/cotton envelope supports entire area and provides maximum patient comfort. Adjustable web shoulder strap with slide buckle for ease of application and removal.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-99152	7"	11"	XS
79-99153	7.5"	12.5"	S
79-99155	8"	15"	M
79-99157	8.5"	17.5"	L
79-99158	9"	20"	XL

**Measure length from elbow crease to base of index finger.*



PROCARE®

Healthcare Bear Arm Sling

Cotton/poly envelope with Healthcare Bear design. Cotton web strap with quick-release buckle for easy application.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-99132	4.5"	11"	Pediatric
79-99133	6.5"	13"	Child

*Measure length from elbow crease to base of index finger.



PROCARE®

Arm Elevator Sling with Pockets

Easily converts to arm sling. Manufactured in sturdy, moisture-proof, mildew-resistant Vivatex materials with sewn-in ice pockets. Available with or without ice bags.

Indication: Ideal for patients requiring arm elevation following surgery.

		Size
79-84052	Pockets w/Ice Bags	Univ.

ELBOW BRACES & SUPPORTS



X-Act™ ROM Elbow

The X-Act ROM Elbow brace is a versatile, easy to apply, comfortable telescoping brace, allowing for ideal fixed or protected range of motion to help improve patient outcomes.

Indication: Designed for Ulnar Collateral Ligament (Tommy John), Bicep Tendon Repair, Supracondylar Fractures, Radial Elbow Fractures, and Lateral/Medial Epicondylitis. Single, lateral aluminum hinge and arm design maximize protection and minimize brace weight

		Elbow
11-9121	X-Act ROM Elbow	Right
11-9122	X-Act ROM Elbow	Left
Accessories		
25-6753	Shoulder Sling	
25-6757	ComfortForm Replacement Pad Kit	
25-6758	Elbow Strap Replacement Kit	
25-6755	Brace Lock Kit	
25-6756	Pronation/Supination Kit, RT	
25-6759	Pronation/Supination Kit, LT	





PROCARE®

Universal Surround® Elbow w/FLOAM™

Comfortable neoprene strap features FLOAM bladder and wraparound universal design with loop-lock closure. FLOAM is intended to reduce impact stress to the epicondyle and helps dissipate muscle force over a greater area. Durable to withstand wear during exercise, sporting activities, and repetitive work conditions.

Indication: Ideal for bursitis, sprains, strains, tendonitis, Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.°	Size
79-82570	≤ 18"	Univ.

*Forearm circumference.



PROCARE®

ElbowRANGER® Motion Control Splint

For non-operative application, collateral ligament strains and ruptures, and ulnar nerve transposition. May be used following early cast removal to stabilize fractures, and for post-operative immobilization. Removable hand grip assists in the prevention of pronation/supination. Aluminum R.O.M. bilateral hinges allow adjustments in 30° increments with formed Kydex® shells. Universal application fits right or left arm.

Indication: Ideal for post-operative immobilization and protected range-of-motion.

	Circum.°	Size
79-94200	8" – 11.5"	S
79-94210	11.5" – 14.5"	M
79-94220	14.5" – 19"	L

*Measurement at mid-bicep.



AIRCAST®

Pneumatic Armband

A single, pre-inflated aircell concentrates compression on the extensor muscle, not around the arm—providing more support, less constriction. Available in black and beige. Universal fit.

Indication: Ideal for Tennis elbow; Golfers elbow (epicondylitis).

	Circum.°	Color	Size
05A-B	8" – 14"	Black	Univ.
05A	8" – 14"	Beige	Univ.

*Forearm circumference.

0501 Replacement Aircell



AIRCAST®

Mayo Clinic Elbow Brace

Provides static adjustable stretch of the elbow in flexion and extension. Replaces the need for two braces for elbow contractures and eliminates the need for serial casting for elbow reconstructions and acute fracture dislocations. ARC™ Forearm Rotational Brace may be used for additional pronation and supination control. Max bicep cir. 23.5". Max forearm cir. 13". US Patent #7,517,329.

Indication: Ideal for fracture dislocation; bicep/tricep tendon rupture; ulnar nerve transposition; total elbow arthroplasty; ligament repair; radial fixation; osteoarthritis; epicondylitis.

Left	Right	Size
05EL	05ER	Univ.
ARC Forearm Rotational Brace		
05EW		Univ.
05E142	Replacement Knob	



PROCARE®

Arm Band w/Compression Pad

Quality nylon construction with reinforced neoprene pad to alleviate stress over lateral/medial epicondyle to address tendonitis, bursitis, and strains. D-ring with contact closure for adjustable compression. One size fits all.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.°	Size
79-97000	≤ 18"	Univ.

°Forearm circumference.



PROCARE®

Elbow Sleeve

Neoprene sleeve provides maximum warmth and comfort without roll or ride-up. Offers support and relief to elbow joint, muscles and tendons. 8.5" length.

Indication: Ideal for mild support and compression of the elbow.

	Circum.°	Size
79-82312	8" – 9"	XS
79-82313	9" – 10"	S
79-82315	10" – 12"	M
79-82317	12" – 14"	L
79-82318	14" – 16"	XL
79-82319	16" – 18"	2XL

°Forearm circumference.



PROCARE®

Tennis Elbow Support w/FLOAM™

Combines neoprene band with FLOAM bladder for compression over tendon insert point. FLOAM has no memory, which prohibits pressure from being displaced to bony prominences.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.°	Size
79-81183	6" – 8"	S
79-81185	8" – 11"	M
79-81187	11" – 14"	L
79-81188	14" – 16"	XL

°Forearm circumference.



PROCARE®

Clinic Tennis Elbow

Two-inch wide nylon strap with D-ring closure to allow comfortable fit and compression to address lateral/ medial epicondylitis and tendonitis. One size fits most patients.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Qty.	Size
79-81031	6/pk.	Univ.



PRO CARE®

Deluxe Tennis Elbow

Perforated vinyl construction with counter-force strap and contact closure provides a comfortable fit to address medial and lateral epicondylitis. Universal size fits most patients.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Size
79-81030	Univ.



PRO CARE®

Tennis Elbow Support

Combines soothing warmth of neoprene elbow sleeve with a counterforce strap that restricts full extension of muscles and tendons. Contact closure strap with D-ring provides exact fit and adjustment.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.*	Length	Size
79-82322	8" – 9"	8"	XS
79-82323	9" – 10"	8"	S
79-82325	10" – 12"	8"	M
79-82327	12" – 14"	8"	L
79-82328	14" +	8"	XL

*Forearm circumference.



PRO CARE®

IMAK Elbow Wrap

Provides comfort and pain relief while restricting painful movements during rest and sleep. Removable splint allows maximum immobilization with the stay inserted, or gentle support with the stay removed. Two adjustable straps for easy change of size and compression. Constructed with breathable foam covered with soft poly flannel for support and warmth. Latex free.

Indication: Ideal for relief of cubital tunnel syndrome, Medial epicondylitis (Golfer's elbow), Lateral epicondylitis (Tennis elbow), and other elbow conditions.

	Circum.*	Size
79-81190	≤ 11"	Univ.
79-81210	11" +	XL/Univ.

*Forearm circumference.



PRO CARE®

Elastic Elbow Support

Pull-on cotton/elastic elbow support for mild compression to address minor elbow strains, sprains, or edema.

Indication: Ideal for mild compression and support of the elbow.

	Circum.*	Size
79-81213	8" – 9"	S
79-81215	9" – 10"	M
79-81217	10" – 11.5"	L
79-81218	11.5" – 13"	XL

*Forearm circumference.

WRIST AND HAND SUPPORTS

PROCARE®**Quick-Fit® Wrist II**

The simple, single-pull lace closure design of the Quick-Fit Wrist II and the adjustable, malleable dorsal stay pod provide variability and a contoured, custom fit to accommodate the changing wrist anatomy during rehabilitation, post-trauma or injury. Universal size for both left and right help limit inventory to improve efficiencies.

US Patent #:8147438

Indication: Ideal for the treatment of sprains, strains and control of Carpal Tunnel Syndrome.

Right	Left	Wrist Circum.	Length	Size
79-87560	79-87570	≤ 10"	8"	Univ.
79-87561	79-87571	10" +	8"	Univ./XL





PROCARE®

ComfortFORM™ Wrist

Durable lightweight foam laminate is Lycra® lined for breathability and patient comfort. Preformed aluminum stay and loop/lock closure helps provide anatomically correct fit and proper support.

Indication: Ideal for sprains, strains and control of Carpal Tunnel syndrome symptoms.

Right	Left	Wrist Circum.	Length	Size
79-87281	79-87291	≤ 4.5"	7.0"	Ped./2XS
79-87282	79-87292	4.5" – 5.5"	7.5"	Ped./XS
79-87283	79-87293	5.5" – 6.5"	7.5"	S
79-87285	79-87295	6.5" – 7.5"	8"	M
79-87287	79-87297	7.5" – 8.5"	8.5"	L
79-87288	79-87298	8.5"+	8.5"	XL



PROCARE®

Quick-Fit® Wrist & Forearm

Comfortable stabilization is provided with cool, soft construction and an easy single-pull lace closure for convenient application. An adjustable dorsal stay pod can be moved proximal or distal for desired controlled range-of-motion and the secondary dorsal stay may be removed to accommodate rehabilitation needs. US Patent #: 8147438

Indication: Ideal for post surgical immobilization, post fracture casting, wrist/forearm sprains and strains.

Right	Left	Wrist Circum.	Length	Size
79-87500	79-87510	≤ 13"	10"	Univ.
79-87501	79-87511	13" +	10"	Univ./XL



AIRCAST®

StabilAir™ Wrist Brace

The StabilAir Wrist Brace is a lightweight, water resistant brace designed for the post-operative treatment of wrist fractures and for a variety of wrist injuries. The easy-to-use brace can be applied quickly and easily in the office, hospital or clinic. Anatomic shells and adjustable aircells conform to arm for a custom fit. Removable padded shells enable inspection or treatment. Available in 2 sizes.

Indication: Ideal for post-operative-use for distal radius fractures.

Right	Left	Wrist Circum.	Size
05FSR	05FSL	5" – 6"	S
05FMR	05FML	6" – 7.75"	M
01PG	Replacement inflation bulb		



AIRCAST®

A2™ Wrist Brace

The Aircast A2 Wrist Brace is designed to restrict motion and provide support for a variety of wrist injuries. Easy to apply, the A2 can be worn during most daily activities. The A2 Wrist Brace with Thumb Spica is available for thumb injuries.

Indication: Ideal for ligament instability, sprain or muscle strain, Carpal Tunnel Syndrome, post-op, post removal of casting or sprain, thumb injuries.

Right	Left	Wrist Circum.	Size
05WSR	05WSL	5" – 6.5"	S
05WMR	05WML	6.5" – 7.5"	M
05WLR	05WLL	7.5" – 9"	L



PRO-CARE®

Quick-Fit® Wrist

Cool, soft felt/brush nylon construction with a single one-pull design fits on a wide range of wrist sizes. Dorsal stay within adjustable stay pod can be moved proximal or distal for desired controlled range-of-motion. Secondary dorsal stay may be removed to accommodate rehabilitation needs. Contoured palmar stay may be adjusted for specific patient needs. Brace is designed to fit below palmar crease for full finger dexterity.

Indication: Ideal for sprains, strains and Carpal Tunnel Syndrome.

Right	Left	Wrist Circum.	Length	Size
79-87460	79-87470	≤ 10"	7.5"	Univ.
79-87461	79-87471	10" +	7.5"	Univ./XL



PRO-CARE®

Contoured Wrist Support

Perforated suede construction with soft flannel lining and cotton stockinette for patient comfort. D-rings with loop/lock closure provides secure fit. Aluminum palmar stay can be contoured to fit each patient. Designed to prevent brace from extending beyond palmar crease.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right	Left	Wrist Circum.	Length	Size
79-87002	79-87012	≤ 5.5"	9"	XS
79-87003	79-87013	5.5" – 6.5"	9"	S
79-87005	79-87015	6.5" – 7.5"	9"	M
79-87007	79-87017	7.5" – 8.5"	9"	L
79-87008	79-87018	8.5"+	9"	XL



PRO-CARE®

CTS Wrist Support

Contoured cotton/elastic construction provides comfortable compression for proper support following mild sprains, strains and cast removal. Circumferential strap for improved wrist stabilization. Preformed aluminum stay provides anatomically correct fit. Contact closure.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right	Left	Wrist Circum.	Size
79-87152	79-87162	≤ 5.5"	Ped./XS
79-87153	79-87163	5.5" – 6.5"	S
79-87155	79-87165	6.5" – 7.5"	M
79-87157	79-87167	7.5" – 8.5"	L
79-87158	79-87168	8.5"+	XL



PRO-CARE®

Elastic Wrist Brace

Soft, durable cotton/elastic construction provides compressive support. Contoured elastic body and preformed aluminum stay for anatomically correct fit and proper support following mild sprain, strain or cast removal. Contact closure straps for ease of application. Designed to prevent brace from extending beyond palmar crease allowing full finger range-of-motion.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right	Left	Circum.	Length	Size
79-87072	79-87082	≤ 5.5"	7"	Ped./XS
79-87073	79-87083	5.5" – 6.5"	7"	S
79-87075	79-87085	6.5" – 7.5"	7"	M
79-87077	79-87087	7.5" – 8.5"	7"	L
79-87078	79-87088	8.5"+	7"	XL



PROCARE®

Lace-Up Wrist Support

Perforated suede/flannel construction. Single pull laced closure for easy application. 7" and 10" lengths available.

Indication: Ideal for wrist sprains, strains, and Carpal Tunnel Syndrome.

Right 7"	Left 7"	Right 10"	Left 10"	Size*
79-87202	79-87212	79-87222	79-87232	Ped./XS
79-87203	79-87213	79-87223	79-87233	S
79-87205	79-87215	79-87225	79-87235	M
79-87207	79-87217	79-87227	79-87237	L
79-87208	79-87218	79-87228	79-87238	XL

*Wrist Circumference: XS = ≤ 5.5", S = 5.5" – 6.5", M = 6.5" – 7.5, L = 7.5" – 8.5", XL = 8.5"+.



PROCARE®

Universal CTS Wrist Brace

Panel elastic construction. Converts from right to left by moving aluminum palmar stay. Circumferential wrist strap provides improved stabilization. Contact closure.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

	Wrist Circum.	Size
79-97012	≤ 5.5"	Ped./XS
79-97013	5.5" – 6.5"	S
79-97015	6.5" – 7.5"	M
79-97017	7.5" – 8.5"	L
79-97018	8.5" – 9.5"	XL
79-97019	9.5" – 10.5"	2XL



PROCARE®

Universal Elastic Wrist Brace

Cotton/elastic construction provides compressive support. Contoured elastic body and aluminum stay intended to provide anatomically correct fit. Converts from right to left by changing location of aluminum stay. Contact closure straps for easy adjustment and application.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

	Wrist Circum.	Size
79-87092	≤ 5.5"	Ped./XS
79-87093	5.5" – 6.5"	S
79-87095	6.5" – 7.5"	M
79-87097	7.5" – 8.5"	L
79-87098	8.5"+	XL



PROCARE®

Deluxe Wrist & Forearm Supports

Contoured design and durable perforated material provides improved patient fit and comfort. Removable, preformed aluminum palmar stay for proper support. Loop/lock closure. Available in 7" and 10" length.

Indication: Ideal for post surgical immobilization, post fracture casting, wrist sprains and strains.

Right 7"	Left 7"	Right 10"	Left 10"	Size*
79-87032	79-87042	79-87052	79-87062	Ped./XS
79-87033	79-87043	79-87053	79-87063	S
79-87035	79-87045	79-87055	79-87065	M
79-87037	79-87047	79-87057	79-87067	L
79-87038	79-87048	79-87058	79-87068	XL

*Wrist Circumference: XS = ≤ 5.5", S = 5.5" – 6.5", M = 6.5" – 7.5, L = 7.5" – 8.5", XL = 8.5"+.



PROCARE®

IMAK Pii-O-Splint®

Provides exceptional comfort and maximum pain relief during rest and sleep. Soft cotton lining extends to fingertips to keep fingers from moving while sleeping. Two rigid plastic dorsal stays immobilize the hand, wrist, and fingers in a neutral position. Constructed with a breathable foam covered with soft poly flannel promotes warmth and good circulation. Palm rests on a tiny "pillow" for added comfort. Latex free. Adjustable Pii-O-Splint allows for easy change of size.

Indication: Ideal for Carpal Tunnel Syndrome, forearm tendinitis, arthritis, and other wrist and hand pain.

	Circum.	Length	Size
79-87540	≤ 12"	12"	Standard/Univ.
79-87550	12"+	12"	Adjustable/Univ.

*Forearm measurement.



PROCARE®

Canvas B.A.T.H.™ Wrist Splint

Constructed of heavy duty canvas for durability. Splint features a double contour stay for precise fit. Slip on design for easy application. Does not extend past palmar crease which does not limit finger dexterity.

Indication: Ideal for wrist sprains or strains; Carpal Tunnel Syndrome and Tendonitis.

Right	Left	Wrist Circum.	Length	Size
79-87441	79-87481	3.5" – 4.5"	5.25"	Ped.
79-87444	79-87484	4" – 5"	5.25"	Child
79-87442	79-87482	4.5" – 5.5"	7.5"	XS
79-87443	79-87483	5.5" – 6.5"	8"	S
79-87445	79-87485	6.5" – 7.5"	8"	M
79-87447	79-87487	7.5" – 8.5"	8"	L
79-87448	79-87488	8.5"+	8"	XL



PROCARE®

Universal Wrist & Forearm Supports

Breathable nylon/fiber with flannelette liner. Loop/lock closure. Adjustable aluminum palmar and dorsal stays for support. Available in 7" and 10" length.

Indication: Ideal for post surgical immobilization, post fracture casting, wrist sprains and strains.

Right	Left	Wrist Circum.	Length	Size
79-87050	79-87060	≤ 9"	10"	Univ.
79-87070	79-87080	≤ 9"	7"	Univ.



PROCARE®

Safety Wrist

Ergonomically designed brace allows for full hand dexterity during repetitive daily activities. Durable, cotton/elastic construction, contact closure straps and plastic palmer stay provide support to help address wrist sprains, tendonitis, and Carpal Tunnel Syndrome.

Indication: Ideal for wrist sprains, tendonitis and Carpal Tunnel Syndrome.

Right	Left	Wrist Circum.	Length	Size
79-82223	79-82233	5.5" – 6.5"	6"	S
79-82225	79-82235	6.5" – 7.5"	6"	M
79-82227	79-82237	7.5" – 8.5"	6"	L
79-82228	79-82238	8.5" – 9.5"	6"	XL
79-82229	79-82239	9.5"+	6"	2XL



PROCARE®

Universal Foam Wrist and Forearm Splint

Constructed of open cell foam that can be trimmed to fit. Universal size fits most patients and helps to reduce stocking inventory. Removable dorsal stay to ensure appropriate fit. Hook and loop closure straps for ease of application and removal.

Indication: Ideal for wrist sprains, strains, Carpal Tunnel Syndrome and post-operative rehabilitation for wrist.

Right	Left	Size
79-87030	79-87040	Univ.



PROCARE®

Universal Wrist-O-Prene™

Wraparound neoprene brace provides compression and support for treatment of Carpal Tunnel Syndrome (CTS) and wrist sprains and strains. Removable contoured aluminum palmar stay. Inside contact closure allows for ease of application and proper fit. Durable construction withstands wear during exercise, sporting activities and repetitive work conditions.

Indication: Ideal for compression and support for treatment of Carpal Tunnel Syndrome (CTS) and wrist sprains and strains.

Right	Left	Wrist Circum.	Size
79-82470	79-82471	≤11"	Univ.



PROCARE®

Cock-Up Splint

Durable, canvas construction with malleable aluminum palmar stay. Lace closure for ease of application.

Indication: Ideal for post-fracture casting, wrist sprains, strains and Carpal Tunnel Syndrome.

Right	Left	Length	Size*
79-87321	79-87331	4.5"	Ped./3XS
79-87341	79-87351	6"	Ped./2XS
79-87342	79-87352	6"	Ped./XS
79-87343	79-87353	6"	S
79-87345	79-87355	6"	M
79-87347	79-87357	6"	L
79-87348	79-87358	6"	XL
79-87381	79-87391	8"	Ped./2XS
79-87382	79-87392	8"	Ped./XS
79-87383	79-87393	8"	S
79-87385	79-87395	8"	M
79-87387	79-87397	8"	L
79-87388	79-87398	8"	XL

*Wrist Circumference: 3XS = ≤ 3.5", 2XS = 3.5" – 4.5", XS = 4.5" – 5.5", S = 5.5" – 6.5", M = 6.5" – 7.5", L = 7.5" – 8.5", XL = 8.5"+.



PROCARE®

Foam Wrist Splint

Foam tricort construction for breathable comfort. Dorsal and palmar malleable aluminum stays may be removed and contoured to proper angulation.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right	Left	Wrist Circum.	Length	Size
79-87422	79-87432	4.5" – 5.5"	10"	Ped./XS
79-87423	79-87433	5.5" – 6.5"	10"	S
79-87425	79-87435	6.5" – 7.5"	10"	M
79-87427	79-87437	7.5" – 8.5"	10"	L
79-87428	79-87438	8.5" +	10"	XL



PROCARE®

IMAK Smart Glove®

Bead-filled wrist pad protects and cushions the underside of the wrist while encouraging neutral positioning of the hand and wrist. Constructed with breathable cotton Lycra® removable dorsal stay for flexible support and ergonomically correct wrist position. Reversible to fit left or right hand. Latex free.

Indication: Ideal for the prevention and relief of repetitive motion injuries. Indications include the treatment and early prevention of Carpal Tunnel Syndrome, arthritis, and other wrist and hand pain.

	Palm Width	Size
79-87542	up to 2.75"	XS
79-87543	2.75" – 3.25"	S
79-87545	3.25" – 3.75"	M
79-87547	3.75" – 4.25"	L
79-87548	4.25" +	XL



PROCARE®

Loop-Lock Cock-Up Splint

Malleable aluminum palmar stay for custom anatomical position. Sturdy canvas construction with loop-lock closure. Splint allows full finger function above palmar crease.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right 6"	Left 6"	Right 8"	Left 8"	Size*
79-87362	79-87372	79-87402	79-87412	Ped./XS
79-87363	79-87373	79-87403	79-87413	S
79-87365	79-87375	79-87405	79-87415	M
79-87367	79-87377	79-87407	79-87417	L
79-87368	79-87378	79-87408	79-87418	XL

*Wrist Circumference: XS = ≤ 5.5", S = 5.5" – 6.5", M = 6.5" – 7.5", L = .5" – 8.5", XL = 8.5"+.



PROCARE®

Universal Cock-Up Splint

Two piece design constructed of durable perforated vinyl, foam, and tricot. Available in 7" and 10" length. May be worn on either hand. Dorsal and palmar aluminum stays provide proper patient support.

Indication: Ideal for immobilization of the wrist and post-fracture casting.

7"	10"	Size
79-87010	79-87020	Univ.
79-87021		Univ./XL



PROCARE®

Wrist Wrap

Adjustable neoprene wraparound design provides a secure compressive support. Fits right or left. Two styles.

Indication: Ideal for compressive support for weak or injured wrist.

		Size
79-82050	Wistlet	Univ.
79-82300	Nylon 2	Univ.

WRIST AND THUMB SUPPORTS

PROCARE®**Quick-Fit® W.T.O.**

Universal wrist, thumb orthosis with a single pull lace-closure design allows for ease of application and a secure, customizable fit delivering ease of use and functional treatment. Lightweight foam laminate maximizes patient comfort and the pre-contoured aluminum palmer and thumb stay may be adjusted for proper angulation to address wrist and thumb pathologies.

U.S. Patent 8147438.

Indication: Ideal for DeQuervain's Syndrome, Gamekeeper's thumb, scaphoid injuries, sprains and strains.

Right	Left	Wrist Circum.	Size
79-87480	79-87490	≤ 12"	Univ
79-87581	79-87591	12+"	Univ./XL





PRO-CARE®

ComfortFORM™ Wrist w/Abducted Thumb

Durable lightweight fabric is Lycra® lined for breathability and patient comfort. Dual aluminum stays support palmar surface of wrist and extensor surface of thumb. Stays can be adjusted for proper angulation. Contact closure straps.

Indication: Ideal for wrist sprains, strains, scaphoid injuries, Carpal Tunnel, Gamekeeper's thumb, and deQuervain's Syndrome symptoms.

Right	Left	Wrist Circum.	Length	Size
79-87302	79-87312	≤ 5.5"	7.5"	Ped./XS
79-87303	79-87313	5.5" – 6.5"	7.5"	S
79-87305	79-87315	6.5" – 7.5"	8"	M
79-87307	79-87317	7.5" – 8.5"	8.5"	L
79-87308	79-87318	8.5"+	8.5"	XL



PRO-CARE®

Universal Thumb-O-Prene™

Neoprene, wrap-around design offers two products in one with a low profile, custom fit. Two stays, rigid and spiral provide needed immobilization for various thumb pathologies or may be removed for step-down rehabilitation. The addition of unique, movable elastic strap provides compression to the base of the thumb to help stabilize the CMC joint.

Indication: Instability of the CMC (Carpometacarpal) joint, Basal Joint Arthritis and thumb tendonitis.

	Size
79-82700	Univ.
79-82701	Univ./XL



PRO-CARE®

ProCare® Thumb Splint

Unique perforated neoprene design with dorsal, radial and palmar stay pockets for a custom fit and maximum stabilization of the thumb. Each pocket contains a malleable aluminum stay and spiral stay that may be removed or shaped for desired support and custom fit.

Indication: Ideal for Gamekeeper's/Skier's thumb, Basal Joint Arthritis and thumb tendonitis.

	Size
79-82710	Univ.



PRO-CARE®

Wrist/Thumb Wrap

Neoprene construction provides compressive wrist support with thumb abduction. One size fits most patients. Right or left.

Indication: Ideal for soft tissue injuries to the thumb and wrist.

	Size
79-82600	Univ.



PROCARE®

ThumbSPICA™

Lightweight, perforated outer foam shell with a cool, cotton-terry liner. Elastic, circumferential, contact closure straps for proper positioning of splint.

Indication: Ideal for soft tissue injury to the thumb, arthritis, tendonitis, deQuervain's Syndrome.

Right	Left	Length	Size*
79-87113	79-87114	9"	S/M
79-87117	79-87118	9"	L/XL

*S/M - Fits most women and smaller men, L/XL - Fits most men and larger women.



PROCARE®

Thumb Splint

Perforated suede with flannel liner. Malleable aluminum stay helps abduct thumb in proper position. Elastic contact closure straps. One size, fits right or left.

Indication: Ideal for soft tissue injuries to the thumb, arthritis, tendonitis, deQuervain's Syndrome.

	Length	Size
79-92171	5"	XS
79-92170	7"	Univ.
79-96100	9"	Univ.



PROCARE®

ThumbGUARD™

Lightweight protection and stabilization for collateral ligament injuries to the thumb. Constructed with preformed plastic outer shell and thin foam inner liner. May be trimmed with cast scissors for anatomically correct fit.

Indication: Ideal for dislocations, collateral ligament injuries (GamesKeeper) and arthritis of the thumb.

Right	Left	Width*	Size
79-87262	79-87272	1.75" – 2.25"	XS
79-87263	79-87273	2.25" – 2.75"	S
79-87265	79-87275	2.75" – 3.25"	M
79-87267	79-87277	3.25" – 3.75"	L
79-87268	79-87278	3.75" – 4.25"	XL

*Width across palm of hand.



PROCARE®

Abducted ThumbSPICA™

Volara™ foam laminate construction with elastic contact closure straps.

Indication: Ideal for soft tissue injury to the thumb, arthritis, tendonitis, deQuervain's Syndrome.

Right	Left	Size
79-87100	79-87101	Univ.

LUMBROSACRAL SUPPORTS

PROCARE®**ComfortFORM™ Back**

Overlapping double pull straps help support the lumbar and abdominal area to promote proper spinal alignment and help prevent and protect injuries including strains and sprains. Silicone gripper strips on either side of back pad help prevent migration. Optional shoulder straps available. Latex free.

Indication: Helps prevent and protect injuries including strains and sprains; Aids in proper alignment while sitting and standing for prolonged periods. Support may be worn under or over clothing.

	Waist Measurement	Size
79-89352	20" – 25"	XS
79-89353	25" – 30"	S
79-89355	30" – 34"	M
79-89357	34" – 38"	L
79-89358	38" – 42"	XL
79-89359	42" – 56"	2XL
79-89359-10	56" – 62"	3XL
79-89350	Shoulder Straps	Univ.





PROCARE®

Sacro-Lumbar Support w/Compression Straps

Low profile 9" elastic support with double side pull compression straps. Neoprene posterior panel accommodates ComfortFORM™ Moldable Insert #79-89120. Foam lumbar pad is included with 79-8250x style. No stays.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

Black	White	Waist Meas.	Size
79-82502		18" – 24"	XS
79-82503	79-89263	24" – 33"	S
79-82505	79-89265	33" – 39"	M
79-82507	79-89267	39" – 45"	L
79-82508	79-89268	45" – 53"	XL
79-82509	79-89269	53" – 59"	2XL
79-82509-10		59" – 65"	3XL
79-82509-11		65" – 72"	4XL



PROCARE®

Clinic Retention Support with Compression Straps

Foam/pile, with 9" panel elastic construction and double side pull elastic compression straps. Two sizes. Accommodates ComfortFORM™ moldable insert #79-89120.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89240	24" – 40"	S/M
79-89241	40" – 54"	L/XL



PROCARE®

10" Criss-Cross Support

Universal style criss-cross elastic support with two plastic posterior stays.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89040	28" – 50"	Univ.



PROCARE®

Criss-Cross Support with Compression Straps

10" elastic criss-cross support with compression straps. Four 9" stays—two aluminum that can be easily formed to help ensure proper support and two plastic to maintain position of support.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89183	24" – 30"	S
79-89185	30" – 36"	M
79-89187	36" – 42"	L
79-89188	42" – 48"	XL
79-89189	48" – 52"	2XL



PROCARE®

10" Double-Pull Sacro-Lumbar Support

Premium all-elastic construction with compression straps. Four 9" posterior plastic stays for stability and support.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89003	31" – 34"	S
79-89005	35" – 38"	M
79-89006	36" – 42"	M/L
79-89007	39" – 42"	L
79-89008	43" – 46"	XL
79-89009	47" – 54"	2XL
79-89009-10	55" – 60"	3XL



PROCARE®

Sacro-Lumbar Support

Elastic construction with felt back pad. Four posterior preformed aluminum stays. Three adjustable web compression straps with front clip buckle closure for custom fit.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

13"	Waist Measurement	Size
79-89023	28" – 32"	S
79-89024	32" – 36"	M
79-89026	36" – 40"	M/L
79-89027	40" – 44"	L
79-89028	44" – 48"	XL
79-89029	48" – 52"	2XL



PROCARE®

Sacro-Lumbar Support with Mesh Back

10" elastic construction with ventilated mesh back panel. Four 9" plastic posterior stays. Side pull elastic compression straps with contact closure.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89203	24" – 30"	S
79-89205	30" – 36"	M
79-89207	36" – 42"	L
79-89208	42" – 48"	XL
79-89209	48" – 54"	2XL



PROCARE®

Neoprene/Elastic Sacro-Lumbar Support w/Compression Straps

Low profile 10" elastic construction with neoprene back panel and side pulls for warmth, compression and adjustability. Two plastic side stays prevent rolling and 3 spiral posterior stays provide added support. Available in sized and universal.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89013	28" – 32"	S
79-89015	34" – 38"	M
79-89016	36" – 40"	M/L
79-89017	40" – 44"	L
79-89018	44" – 48"	XL
79-89019	48" – 52"	2XL
79-89010	32" – 52"	Univ.



PROCARE®

Industrial Back Support

8" black, mesh back support with contact closures. Double-pull elastic compression suspenders. Choice between sewn-on or optional attachable shoulder suspension straps.

Indication: Ideal for helping to prevent and protect back injuries.

Black w/o Suspenders	Blk. w/Suspenders	Waist	Size
	79-89142	22" – 24"	XS
79-89083	79-89143	24" – 30"	S
79-89085	79-89145	30" – 36"	M
79-89087	79-89147	36" – 42"	L
79-89088	79-89148	42" – 50"	XL
79-89089	79-89149	50" – 56"	2XL
	79-89149L	56" – 64"	3XL
	79-89149X	64" – 72"	4XL
79-89190	Optional Attachable Suspenders		Univ.



PROCARE®

Industrial Back Support w/Compression Pad

Elastic construction with 8" molded back compression pad. Two 6" plastic stays to ensure proper support. Contact closure.

Indication: Ideal for helping to prevent and protect back injuries.

	Waist Measurement	Size
79-97033	24" – 30"	S
79-97035	30" – 36"	M
79-97037	36" – 42"	L
79-97038	42" – 50"	XL
79-89350	Suspenders	Univ.



PROCARE®

Clinic Industrial Back Support

Lightweight, low profile 8" back support with double-pull elastic compression straps and four sewn-in 7" plastic stays. Contact closure. Sewn-on adjustable suspenders.

Indication: Ideal for helping to prevent and protect back injuries.

Right	Left	Size
79-89172	18" – 24"	XS
79-89173	24" – 30"	S
79-89175	30" – 36"	M
79-89177	36" – 42"	L
79-89178	42" – 50"	XL
79-89179	50" – 56"	2XL



PROCARE®

Maternity Belt

8" plush elastic belt helps relieve discomfort of lower abdomen and back during all stages of pregnancy. Abdominal strap supports stomach while elastic back provides compression and support for lumbar area. Optional 8" ComfortFORM™ moldable back insert #79-89120 provides additional support.

Indication: Ideal for helping relieve lower back discomfort during pregnancy.

	Hip Measurement*	Size
79-89293	32" – 42"	S
79-89295	42" – 52"	M
79-89297	52" – 62"	L
79-89298	62" – 72"	XL

*Measurement taken during pregnancy.

ABDOMINAL & RIB SUPPORTS

PROCARE®

Premium Panel Elastic Binder

Constructed with panel elastic and soft cotton type lining for patient comfort. Universal style, flannel lined panel elastic support with contact closure. Multi-panel design helps prevent panel rolling. Available in each or in a pack of six (6). 9", 12" or 15" height.

Indication: Ideal for providing compression and support for post-natal or post-abdominal surgery, abdominal strains and weakness.

9" 3-Panel	12" 4-Panel	15" 5-Panel	Waist Meas.	Qty.	Size
79-89070	79-89090		30" – 45"		Univ.
79-89071	79-89091		45" – 62"		Univ.
79-89210	79-89220	79-92149	62" – 74"		XL
79-89270	79-89280		72" – 84"		2XL
79-89250	79-89260		82" – 100"		3XL
79-99430	79-99440		30" – 45"	6/pk.	Univ.
79-99431	79-99441		45" – 62"	6/pk.	Univ.





PROCARE®

Abdominal Binder

Elastic binder with four sewn-in plastic stays provides abdominal and lumbar support. Abdominal, post-natal or post-surgical application. Easy adjustment with contact closure. Latex free.

Indication: Ideal for providing compression and support for strains and weakness of the abdominal area.

10"	12"	14"	Waist Meas.	Size
79-89043	79-89323	79-89333	24" – 30"	S
79-89045	79-89325	79-89335	30" – 36"	M
79-89047	79-89327	79-89337	36" – 42"	L
79-89048	79-89328	79-89338	42" – 48"	XL
79-89049	79-89329	79-89339	48" – 54"	2XL



PROCARE®

Personal Abdominal Binder

Constructed of multi-panel, vented elastic. Multi-panel helps binder conform to patient's body. 10" height is an all elastic panel. Latex free.

Indication: Ideal for providing compression and support for strains and weakness of the abdominal area.

	Height	Waist Meas.	Size
79-89364	9"	20" – 42"	S/M
79-89366	9"	36" – 65"	M/L
79-89376	10"	32" – 65"	M/L
79-89384	12"	20" – 42"	S/M
79-89386	12"	36" – 65"	M/L



PROCARE®

10" Universal Abdominal Support

10" elastic and pile/foam laminate with three 9" removable plastic stays. Contact closure to ensure ease of application and removal. Latex free.

Indication: Ideal for providing compression and support for strains and weakness of the abdominal area.

	Waist Meas.*	Size
79-89080	28" – 50"	Univ.

*Measurements taken from male waist or female hip.



PROCARE®

Universal Deluxe Rib Belt

6" two-panel elastic belt with flannel lining and contact closure. Choice of male or female styles. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

Male	Female	Rib Meas.	Qty.	Size
79-89150	79-89160	28" – 50"		Univ.
79-99040	79-99030	28" – 50"	6/pk.	Univ.
79-89151		50" – 65"		XL



PROCARE®

Universal Rib Belt

6" elastic and pile/foam construction. Contact closure. Choice of male or female styles. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

		Rib Meas.	Size
79-89050	Male	24" – 50"	Univ.
79-89060	Female	24" – 50"	Univ.



PROCARE®

Deluxe Rib Belt

6" elastic and foam/flannel construction contoured for male or female styles. Contact closure for ease of application. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

Male	Rib Meas.	Female	Rib Meas.	Size
79-89063	32" – 36"	79-89073	24" – 27"	S
79-89065	36" – 40"	79-89075	27" – 33"	M
79-89067	40" – 46"	79-89077	33" – 37"	L
79-89068	46" – 50"	79-89078	37" – 44"	XL
79-89069	50" – 54"	79-89079	44" – 50"	2XL



PROCARE®

8" Universal Rib Belt

Panel elastic with foam/flannel lining. Contact closure for proper adjustment. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

		Rib Meas.	Size
79-89100	Male	28" – 50"	Univ.
79-89110	Female	28" – 50"	Univ.

KNEE IMMOBILIZERS

PRO CARE®

3-Panel Knee Splint

Universal size, 3-piece nylon/fiber laminate construction provides post-operative knee support and immobilization. Adjustable side panels with I-beam stays and patella strap allow for correct fit and easy application. Three posterior aluminum stays for proper immobilization. Loop lock closure. Standard sizes fits up to 29" and XL up to 36" maximum thigh circumference. Adjustable foam pop pad included. 20" and 24" length are available with cotton/terry liner and removable stays in 10-pack configuration only.

Indication: Ideal for post-operative immobilization of the knee joint; grade 1 and grade 2 collateral ligament sprains; patellar subluxations.

Standard	10 Pack	Max Length	Thigh Circum.*	Size
79-80000		12"	29"	Univ.
79-80010		14"	29"	Univ.
79-80110		16"	29"	Univ.
79-80020	79-80720	20"	29"	Univ.
79-80021		20"	36"	XL
79-80040		22"	29"	Univ.
79-80041		22"	36"	XL
79-80030	79-80730	24"	29"	Univ.
79-80031		24"	36"	XL
79-80220 with 35° bend		20"	29"	Univ.

*Circumference 4" above patella.





PROCARE®

Super Knee Splint

Comfortable, sturdy, foam/mesh construction with alternating web straps and loop/lock contact closure. Medial/lateral I-beam stays and three double bend posterior stays for proper immobilization. Adjustable foam pop pad included.

Indication: Ideal for immobilization of the knee joint.

12"	16"	20"	Circum.°	Size
79-80002			12" – 14"	Ped./XS
79-80003	79-80013	79-80023	14" – 16"	S
79-80005	79-80015	79-80025	16" – 19"	M
79-80007	79-80017	79-80027	19" – 22"	L
79-80008	79-80018	79-80028	22" – 26"	XL
79-80019	79-80029		26" – 30"	2XL

22"	24"	26"	Circum.°	Size
79-80053	79-80033	79-80043	14" – 16"	S
79-80055	79-80035	79-80045	16" – 19"	M
79-80057	79-80037	79-80047	19" – 22"	L
79-80058	79-80038	79-80048	22" – 26"	XL
79-80059	79-80039		26" – 30"	2XL

°Circumference 4" above patella.



PROCARE®

Quick-Fit® Basic Knee Splint

Foam laminate construction with 3" elastic contact closure straps provides compressive support and immobilization. Dual contoured, I-beam posterior stays and moveable medial/lateral double stay casings for correct fit.

Indication: Ideal for immobilization of the knee joint.

	Length	Max Thigh Circum.°	Size
79-96011	10"	29"	Univ.
79-96012	12"	29"	Univ.
79-96014	14"	29"	Univ.
79-96016	16"	29"	Univ.
79-96018	18"	36"	Univ.
79-96019	20"	36"	Univ.
79-96021	22"	36"	Univ.
79-96024	24"	36"	Univ.
79-96026	26"	36"	Univ.

°Circumference 4" above patella.



PROCARE®

Clinic 3-Panel Splint

Universal three-piece design for ease of adjustment and correct size adjustment. Foam posterior panel may be trimmed for a customized fit. Foam/mesh side panels with I-beam stays and encircling loop-lock straps for maximum stability. Straight stay version features popliteal pad.

Indication: Ideal for post-operative immobilization; grade 1 and grade 2 collateral ligament sprain; patellar subluxations.

	Length	Max Thigh Circum.	Size
Straight Stays			
79-80810	16"		
79-80820	20"		
79-80830	24"		
Contoured Stays			
79-80410	12"	32"	Univ.
79-80420	16"	32"	Univ.
79-80170	20"	32"	Univ.
79-80180	24"	32"	Univ.
79-80290	26"	32"	Univ.



PROCARE®

Universal Foam Knee Splint

Constructed of open cell foam and detachable medial/lateral aluminum I-beam stay panels for easy adjustment and custom fit. Three double bend posterior aluminum stays. Alternating web straps provide dual adjustability, using a D-ring buckle and a hook tab.

Indication: Ideal for post-operative immobilization, grade 1 and grade 2 collateral ligament sprain, and patella subluxations.

	Length	Size
79-80600	8"	Univ.
79-96870	10"	Univ.
79-96871	12"	Univ.
79-96914	14"	Univ.
79-96880	16"	Univ.
79-96881	18"	Univ.
79-96882	20"	Univ.
79-96890	22"	Univ.
79-96891	24"	Univ.
79-96900	26"	Univ.



PROCARE®

Universal Foam Knee with Loop Lock Closure

Constructed of durable open cell foam with moveable medial/lateral stay casings for proper positioning. Available with contoured or straight posterior stays. Immobilizer can be trimmed for custom fit. 4" elastic patella strap for compressive support. Straight stay version features open patella cut-out and popliteal pad.

Indication: Ideal for post-operative immobilization; grade 1 and grade 2 collateral ligament sprain; patellar subluxations.

Contour Stays	Straight Stays	Length	Max Thigh Circ.	Size	Strap Color
79-80611	79-80610	12"	22.5"	Univ.	White
79-80360		14"	29"	Univ.	Black
79-80621	79-80620	16"	24"	Univ.	White
79-80370		19"	29"	Univ.	Black
79-80631	79-80630	20"	29"	Univ.	White
79-80380		22"	29"	Univ.	Black
79-80641	79-80640	24"	29"	Univ.	White
79-80390		26"	29"	Univ.	Black
79-80651	79-80650	28"	29"	Univ.	White



PROCARE®

Compression Knee Dressing

Trimable foam laminate body with 6" wide multi-panel plush elastic straps for full compression to aid in prevention of edema. Moveable medial/lateral stay casings for proper positioning. Adjustable foam pop pad included.

Indication: Ideal for providing compression and immobilization of the knee following arthroscopic and reconstructive knee surgery; grade 1 and grade 2 collateral ligament sprain; patellar subluxations.

Contoured Stays	Straight Stays	Without Stays	Length	Max Thigh Circ.	Size
79-80090	79-80320	79-80430	16"	29"	Univ.
79-80190	79-80330	79-80440	20"	30"	Univ.
79-80260	79-80340	79-80450	24"	32"	Univ.
79-80280			28"	32"	Univ.



AIRCAST®

Aircast Knee Immobilizer

Stabilizes the knee in full extension. Shaped to provide space behind the knee, even at full extension. May be used with the Knee Cryo/Cuff™ to add the benefits of cold and compression.

Indication: Ideal for post-op stabilization.

	Length	Inseam	Size
07C	23"	25" – 28"	S
07A	26.5"	27" – 33"	M
07B	29.5"	32" – 36"	L

POST-OP KNEE BRACES



X-Act™ ROM Knee

The X-Act ROM is remarkably quick and easy to apply, allowing efficient use of OR time. The anti-migration foam and dynamic suspension frame naturally contour this lightweight brace to the patient's leg, providing a precise fit after surgery.

Indication: DonJoy X-Act™ ROM Knee orthosis is designed to provide precise protection and support post-operatively for patients following knee procedures.

		Size
11-2151-9	X-Act ROM Knee Brace	Universal
Accessories		
25-6755	X-Act ROM Knee Zip Ties	





X-Act™ ROM Hip

The X-Act ROM Hip Brace includes a Compression Belt with a dual pull-tab system that is simple for the patient to tighten and provides support to the lower back. Our Three Point Design helps push the femoral head into the acetabulum to minimize dislocation.

Indication: DonJoy X-Act™ ROM Hip orthosis is designed to provide precise protection and support post-operatively for patients following various hip procedures.

	Description
11-3233	X-Act™ ROM HIP
25-7525	Kit, Replacement Pad, X-Act™ ROM Hip
25-7528	Replacement Belt, X-Act™ ROM Hip
25-7529	Assembly, Hip Component, X-Act™ ROM Hip



KneeRANGER® II - Universal

Post-operative immobilization with protected range-of-motion. Dual axis, polycentric hinge allows for convenient dial-in adjustment with flexion and extension stops at 0°, 15°, 30°, 60° and 90°. Available in two (2) lengths. Trimmable foam leg wraps for customized fit. Medial/lateral plastic shells riveted to hinge bars provide additional support following post-operative procedures. Accommodates Ankle/Shoe Insert for added rotational control.

Indication: Ideal for providing protected range-of-motion for post-operative immobilization.

	Length	Max Thigh Circum.	Qty.	Size
79-94230	17"	36"		Univ./Short
79-94240	26"	36"		Univ./Long
Accessories				
79-94260	Replacement wraps, Short		Set	Univ.
79-94261	Replacement wraps, Long		Set	Univ.
79-94270	Replacement Strap		Set (6)	Univ.
79-96170	Foam Condyle Pads		Set	Univ.



KneeRANGER® - Universal

Post-operative immobilization and protected range-of-motion with 7.5° incremental settings between 0° and 120° of flexion and 0° and 70° of extension. Hinge bar breakdown allows for up to 5" reduction in length to convert from long leg to short for long-term rehabilitation. Trimmable medial/lateral plastic shells. Accommodates Ankle/Shoe Insert for added rotational control. Replaceable foam leg wraps.

Indication: Ideal for patients recovering from injury or surgery providing full and limited range-of-motion settings.

	Circum.*	Size
79-94130	16" overall length	XS
79-94000 ^o	Short, 32" inseam or less	S
79-94001 ^o	Long, 32" plus	L
Accessories		
79-94020	Replacement wraps, short	Set
79-94021	Replacement wraps, long	Set
79-96170	Foam Condyle Pads	Set

^oBreakdown bars. Maximum circumference: thigh 31", calf 22"



KneeRANGER® - Sized

24" length of offers post-operative immobilization and protected range-of-motion with settings between 0° - 120° of flexion and 0° - 70° of extension in 7.5° increments. One-piece easy on/off design. Breathable nylon/fiber construction with soft inner liner. Accommodates Ankle/Shoe Insert for added rotational control.

Indication: Ideal for immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs.

	Circum.*	Size
79-94003	14" - 16"	S
79-94005	16" - 18"	M
79-94007	18" - 20"	L
79-94008	20" - 22"	XL
79-94009	22" - 24"	2XL

^oCircumference measurement taken 4" above mid-patella.



PROCARE®

KneeRANGER® II - Sized

One-piece, 24" length design allows for easy application and adjustment. Dual axis hinge provides protected, full range-of-motion control with flexion and extension stops at 0°, 15°, 30°, 60°, and 90°. Breathable nylon/fiber/flannelette construction. Integrated strap/hinge bar design accommodates Ankle/Shoe Insert for added rotational control.

Indication: Ideal for immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs.

	Circum.°	Size
79-94253	14" – 16"	S
79-94255	16" – 18"	M
79-94257	18" – 20"	L
79-94258	22" – 24"	XL

[°]Circumference measurement taken 4" above mid-patella.



PROCARE®

KneeRANGER® Lite

Economical, lightweight, universal leg brace for warmer climates. Narrow band circumferential cuffs replace full cut leg wraps. Integrated strap/hinge bar design. Provides immobilization and protected range-of-motion in 7.5° increments from 0° – 120° of flexion and 0° – 70° of extension.

Indication: Ideal for immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs.

	Length*	Qty.	Size
79-94040	24"		Univ.
Accessories			
79-94041	Replacement wraps	Set	
79-96170	Foam Condyle Pads	Set	

*Maximum circumference: thigh 31", calf 22".



PROCARE®

KneeRANGER® Range-of-Motion

Provides protected range-of-motion for postoperative immobilization. Incremental settings between 0° and 120° of flexion and 0° and 70° of extension in 7.5° increments. Trimmable nylon/foam leg wraps. Adjustable contact closure straps. Accommodates Ankle/Shoe Insert for added rotational control.

Indication: Ideal for providing protected range-of-motion for post-operative immobilization.

	Circum.°	Qty.	Size
79-94060	26"		Univ.
Accessories			
79-94021	Replacement wraps, Long	Set	Univ.
79-96170	Foam Condyle Pads	Set	Univ.

*Maximum circumference: thigh 31", calf 22".



PROCARE®

KneeRANGER® Rehab Splint

Full foam wrap may be trimmed for customized fit. Anterior plastic shells provide additional compression and support. Protected range-of-motion in 7.5° increments with settings between 0° – 120° of flexion and 0° – 70° extension. Universal size. Fits right or left.

Indication: Ideal for immobilization and protected range-of-motion of the knee.

	Length*	Qty.	Size
79-94120	16"		Univ.
Accessories			
79-94150	Replacement wraps	Set	

*Maximum circumference: thigh 31", calf 22"



PROCARE®

KneeRANGER® – Neoprene

Designed for wear when compressive support, motion control, and heat retention are indicated. Contoured 3/16" neoprene sleeve. Range-of-motion settings between 0° – 120° of flexion and 0° – 70° extension in 7.5° increments for protected range-of-motion.

Indication: Ideal for mild to moderate MCL and/or LCL sprains, post arthroscopic meniscectomy, and mild global instabilities.

	Circum.*	Size
79-94132	13.5" – 15.5"	XS
79-94133	15.5" – 18"	S
79-94135	18" – 20.5"	M
79-94137	20.5" – 23"	L
79-94138	23" – 25.5"	XL
79-94139	25.5" – 28"	2XL
79-94139-10	28" – 31"	3XL

*Measurement taken 6" above mid-patella.



PROCARE®

Dual Axis Hinge

3/16" of set polycentric design provides stability and full range-of-motion control with convenient dial-in flexion/extension stops at 0°, 15°, 30°, 60° and 90°. Patent No. 5, 443,444. Sold in pairs.

Indication: Ideal for casting, fabrication and range-of-motion control.

		Length
79-94400	Dual Axis Hinge	10.75"
79-94410	Dual Axis Hinge	14.5"
79-94550	w/rubber condyle pad	10.75"
79-94560	w/rubber condyle pad and cap cover	10.75"
79-94580	w/rubber condyle pad	14.5"
79-94590	w/rubber condyle pad and cap cover	14.5"
79-93820	w/rubber condyle pad and cap cover	21.5"



PROCARE®

KneeRANGER® Ankle/Shoe Insert

Designed to be used with long leg post-op braces. Attaches inside hinge bars to extend brace leverage for added rotational control. Contoured heel cup. Unrestricted dorsi/plantar flexion.

Indication: To be used in conjunction with any of the IROM/ROM long leg splints for added rotational control.

Right	Left	Male Shoe Size	Fem. Shoe Size	Length	Size
79-94153	79-94163	3 – 4	4 – 5	12.5"	S
79-94155	79-94165	5 – 9	6 – 10	12.5"	M
79-94157	79-94167	10 – 14	11 – 15	12.5"	L



PROCARE®

KneeRANGER® Fracture Hinges

A high-performance, low-cost alternative to traditional cast hinges. The ultimate in hinged motion control. Range-of-motion settings between 0° – 120° flexion and 0° – 70° extension in 7.5° increments.

Indication: To be used in conjunction with standard casting materials to combine an adjustable hinge mechanism at the knee joint with circumferential casting above and below the joint.

	Qty.	Size
79-94140	1 pair	Univ.
79-96035	1 pair	Univ./without paddles

THIGH, KNEE & CALF SUPPORTS

PROCARE®**Sport Knee Sleeve**

1/8" neoprene sleeve designed for durability and greater comfort. Extra length pad at the knee of offers localized heat retention, compressive support, and improved protection. Choice between open or closed patella. 13" Length.

Indication: Ideal for general compression and support.

Closed Patella	Open Patella	Open Pop.	Circum.*	Size
79-82002	79-82012	79-82122	13.5" – 15.5"	XS
79-82003	79-82013	79-82123	15.5" – 18"	S
79-82005	79-82015	79-82125	18" – 20.5"	M
79-82007	79-82017	79-82127	20.5" – 23"	L
79-82008	79-82018	79-82128	23" – 25.5"	XL
79-82009	79-82019	79-82129	25.5" – 28"	2XL
79-82009-10	79-82019-10	79-82129-10	28" – 31"	3XL
79-82019-11	79-82129-11		31" – 34"	4XL

*Measurement taken 6" above mid-patella.





PROCARE®

Clinic Knee Sleeve

1/8" single-sided nylon neoprene with contoured design for controlled, uniform compression. 10" length.

Indication: Ideal for general compression and support.

Open Patella	Circum.*	Size
79-82612	13.5" – 15.5"	XS
79-82613	15.5" – 18"	S
79-82615	18" – 20.5"	M
79-82617	20.5" – 23"	L
79-82618	23" – 25.5"	XL
79-82619	25.5" – 28"	2XL



PROCARE®

Knee Support with Closed Patella

Constructed of high quality, reversible 1/8" neoprene for comfortable compression and warmth. 10" length.

Indication: Ideal for general compression and support.

	Circum.*	Size
79-82623	15.5" – 18"	S
79-82625	18" – 20.5"	M
79-82627	20.5" – 23"	L
79-82628	23" – 25.5"	XL
79-82629	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Knee Support with Reinforced Patella

1/8" neoprene provides compressive support and warmth. Reversible. Blue/black neoprene. 10" length.

Indication: Ideal for Patellofemoral (knee cap) disorders.

	Circum.*	Size
79-82633	15.5" – 18"	S
79-82635	18" – 20.5"	M
79-82637	20.5" – 23"	L
79-82638	23" – 25.5"	XL
79-82639	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Universal Knee Wrap with Closed Patella

1/8" neoprene wraparound style with medial and lateral spiral stays for compressive support following strains, sprains or other minor indications.

Indication: Ideal for mild to moderate compression and patellar control.

		Circum.*	Size
79-82431	Nylon one side	21"	U/XL
79-82441	Nylon two sides	21"	U/XL

*Measurement taken 6" above mid-patella.



PROCARE®

Arthroscopy Knee Support

1/8" neoprene wraparound design provides compression and support with easy adjustment and application. Bilateral design fits right and left. Optional 1/8" neoprene Lateral Release Pad offers additional support following arthroscopic procedures.

Indication: Ideal for post-arthroscopic compression and support to the knee.

	Circum.®	Size
79-82402	13.5" – 15.5"	XS
79-82403	15.5" – 18"	S
79-82405	18" – 20.5"	M
79-82407	20.5" – 23"	L
79-82408	23" – 25.5"	XL
79-82400	Lateral Release Pad	Univ.

®Measurement taken 6" above mid-patella.



PROCARE®

Elastic Knee Support

Sturdy, heavy-duty elastic provides cool, compressive support. Closed patella.

Indication: Ideal for mild support and compression; mild patella bursitis.

	Circum.®	Size
79-80193	15.5" – 18"	S
79-80195	18" – 20.5"	M
79-80197	20.5" – 23"	L
79-80198	23" – 25.5"	XL

®Measurement taken 6" above mid-patella.



PROCARE®

Elastic Cartilage Knee Support

Heavy knit cotton/elastic construction with open patella design and sewn-in medial/lateral cartilage pads and stays provide added compression and support following arthroscopic procedures.

Indication: Ideal for mild chondromalacia; osteoarthritis; grade 1 collateral ligament sprain; patella femoral dysfunction support and compression.

	Circum.®	Size
79-80243	15.5" – 18"	S
79-80245	18" – 20.5"	M
79-80247	20.5" – 23"	L
79-80248	23" – 25.5"	XL
79-80249	25.5" – 28"	2XL

®Measurement taken 6" above mid-patella.



PROCARE®

Elastic Spiral Knee Support

Cotton/elastic construction provides cool, comfortable compressive support. Six (6) sewn-in medial/lateral spiral stays for additional support.

Indication: Ideal for osteoarthritis; bursitis; grade 1 collateral ligament sprain.

	Circum.®	Size
79-80233	15.5" – 18"	S
79-80235	18" – 20.5"	M
79-80237	20.5" – 23"	L
79-80238	23" – 25.5"	XL
79-80239	25.5" – 28"	2XL

®Measurement taken 6" above mid-patella.



PROCARE®

Thigh Sleeve

1/8" neoprene thigh sleeve provides firm compression and soothing warmth for hamstring and quadriceps strain. 12" length.

Indication: Ideal for strains or contusions to the quadriceps, hamstring, or groin muscles.

	Circum.*	Size
79-82332	18" – 19"	XS
79-82333	19" – 20"	S
79-82335	20" – 22"	M
79-82337	22" – 25"	L
79-82338	25" – 28"	XL
79-82339	28" – 31"	2XL

*Circumference measurement taken at mid-thigh.



PROCARE®

Calf Sleeve

1/8" neoprene sleeve contours to assure proper fit and correct degree of compression and support for strains and sprains. 12" Length.

Indication: Ideal for strains; shin splints.

	Circum.*	Size
79-82342	12" – 13"	XS
79-82343	13" – 14"	S
79-82345	14" – 15"	M
79-82347	15" – 17"	L
79-82348	17" – 19"	XL

*Circumference measurement taken at mid-calf.

PATELLOREMORAL

PROCARE®**Stabilized Knee Support**

1/8" neoprene with medial/lateral spiral stays for added support. Superior and inferior compression straps help prevent migration and provide additional support. Universal neoprene buttress may be adjusted to allow for proper patella management and stabilization.

Indication: Ideal for mild patella tendonitis, chondromalacia or subluxations.

Closed Pop.	Open Pop.	Circum.*	Size
79-82722	79-82752	13.5" – 15.5"	XS
79-82723	79-82753	15.5" – 18"	S
79-82725	79-82755	18" – 20.5"	M
79-82727	79-82757	20.5" – 23"	L
79-82728	79-82758	23" – 25.5"	XL
79-82729	79-82759	25.5" – 28"	2XL
79-82729-10		28" – 31"	3XL
79-82729-11		31" – 34"	4XL
79-82491	Horseshoe Buttress		Univ.

*Measurement taken 6" above mid-patella.





PROCARE®

Patella Stabilizer with Buttress

1/8" neoprene support with choice of adjustable felt horseshoe-buttress or sewn-in donut buttress to accommodate mild patella tendonitis, chondromalacia, or subluxation. Donut buttress may be trimmed to allow for medial, lateral, inferior, or superior buttress. Medial/lateral spiral stays and superior/inferior compression straps provide additional stability. Open popliteal available with donut buttress only.

Indication: For Patello Femoral (knee cap) disorders.

Horseshoe Buttress	Done Buttress w/Open Pop.	Circum.*	Size
79-92852	79-92842	13.5" – 15.5"	XS
79-92853	79-92843	15.5" – 18"	S
79-92855	79-92845	18" – 20.5"	M
79-92857	79-92847	20.5" – 23"	L
79-92858	79-92848	23" – 25.5"	XL
79-92859	79-92849	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Sport Knee Reinforced

1/8" neoprene sleeve with universal neoprene buttress offers full patellar control for chondromalacia, tendonitis, and patella tracking abnormalities. Neoprene buttress may be adjusted to allow for medial, lateral, inferior or superior patella stabilization. Choice between closed or open popliteal. Anterior oval compression pad with patella cutout provides additional warmth and compression.

Indication: Ideal for patella femoral tracking dysfunction; chondromalacia.

Closed Pop.	Open Pop.	Circum.*	Size
79-82702	79-82052	13.5" – 15.5"	XS
79-82703	79-82053	15.5" – 18"	S
79-82705	79-82055	18" – 20.5"	M
79-82707	79-82057	20.5" – 23"	L
79-82708	79-82058	23" – 25.5"	XL
79-82709	79-82059	25.5" – 28"	2XL
79-82709-10		28" – 31"	3XL
79-82491	Horseshoe Buttress		Univ.

*Measurement taken 6" above mid-patella.



PROCARE®

Reinforced Patella Stabilizer

1/8" neoprene with sewn-in horseshoe tubular buttress helps stabilize and improve patella tracking. Adjustable superior/inferior straps with contact closure for additional compression and stability.

Indication: Ideal for patella femoral tracking dysfunction; chondromalacia, patella tendonitis.

	Circum.*	Size
79-94433	15.5" – 18"	S
79-94435	18" – 20.5"	M
79-94437	20.5" – 23"	L
79-94438	23" – 25.5"	XL
79-94439	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Hinged Patella Stabilizer

1/8" neoprene construction with sewn-in tubular patella buttress and removable dual-axis polycentric hinges for additional medial/lateral support. Adjustable superior/inferior straps for compression and stability.

Indication: Ideal for medial/lateral instabilities with patella tendonitis, chondromalacia, subluxation, grade 1 collateral ligament sprains.

	Circum.*	Size
79-94483	15.5" – 18"	S
79-94485	18" – 20.5"	M
79-94487	20.5" – 23"	L
79-94488	23" – 25.5"	XL
79-94489	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Lateral Patella Stabilizer

1/8" neoprene with tubular sewn-in "J" buttress and dynamic side-pull compression strap supports patella malalignment, dislocation, and subluxation. Sold in right and left configuration.

Indication: Ideal for supporting lateral patella subluxation, dislocation or lateral Patellofemoral malalignment.

Right	Left	Circum.*	Size
79-94463	79-94473	15.5" – 18"	S
79-94465	79-94475	18" – 20.5"	M
79-94467	79-94477	20.5" – 23"	L
79-94468	79-94478	23" – 25.5"	XL
79-94469	79-94479	25.5" – 28"	2XL
79-94469-10	79-94479-10	28" – 31"	3XL
79-94469-11	79-94479-11	31" – 34"	4XL

**Measurement taken 6" above mid-patella.*



PROCARE®

Universal Patella Knee with Buttress

1/8" neoprene wrap with open patella and spiral stays. Sleeve can accommodate buttress pad for additional patellar support.

Indication: Ideal for mild support and compression; mild patella bursitis.

	Circum.*	Size
79-82460	≤ 21"	Univ.
79-82461	≤ 25"	Univ./XL

**Circumference of knee fully extended.*



PROCARE®

Lites Visco Knee Support

Constructed of lightweight four-way stretch elastic, and two medial/lateral stays for compression and support. Silicone buttress provides patella stabilization.

Indication: Ideal for chondromalacia; osteoarthritis; patellofemoral tracking dysfunction; grade I collateral ligament sprain.

	Circum.*	Size
79-80162	14.5" – 16"	XS
79-80163	16" – 17"	S
79-80165	17" – 18"	M
79-80167	18" – 19.25"	L
79-80168	19.25" – 20.5"	XL
79-80169	20.5" – 22"	2XL
79-80169-10	22" – 23"	3XL

**Measurement taken 6" above mid-patella.*



PROCARE®

Surround® Patella Strap

The primary body of the strap is a soft, latex free nylon/foam laminate with contact closures. A unique cylindrical FLOAM™ bladder is encased inside the body of the primary patella strap and may be adjusted for added compression. Floam provides uniform pressure and a customized fit while prohibiting pressure from being displaced to bony prominences. Patent pending.

Indication: Helps relieve anterior knee pain caused by patellar tendonitis, chondromalacia and Osgood Schlatter by applying pressure to the patella tendon.

	Circum. Just Below Kneecap	Size
79-80322	≤10"	XS
79-80323	10.5" – 12.5"	S
79-80325	13" – 15"	M
79-80327	15.5" – 17.5"	L
79-80328	18"+	XL



PRO-CARE®

Dynamic Pull Universal Knee Wrap

1/8" neoprene universal knee wrap provides compressive support and warmth. Adjustable medial/lateral patella pull strap to help improve patella tracking. Choice of nylon one side or nylon both sides.

Indication: Ideal for patellofemoral tracking dysfunction and chondromalacia.

		Circum.*	Size
79-82150	Nylon one side	≤24"	Univ.
79-82160	Nylon two sides	≤24"	Univ.

*Circumference of knee fully extended.



AIRCAST®

Infrapatellar Band

An aircell focuses compression on the patellar tendon, helping to decrease stress at the tibial tubercle. Available in black only. Universal fit.

Indication: Ideal for Osgood Schlatter; Patella tendonitis.

	Circum.	Color
08A-B	10" – 17"	Black

HINGED KNEE BRACES



Reaction Knee® Brace

The Reaction Web™ knee brace is a responsive, webbed approach to anterior knee pain that gives a distinct alternative to the basic knee sleeve. The elastomeric web is an innovative solution that provides energy dispersion to the painful knee, allowing comfortable motion in a brace that actually stays in place.

Indication: The elastomeric web design helps to reduce pain by dispersing energy across the knee. The web absorbs shock and shifts the peak loads away from the painful area of the knee. The elastomeric web acts to dynamically stabilize the patella on all sides, bringing the patella into proper tracking position to reduce pain for patients with general patellofemoral instabilities.

	Measurement*	Size
11-0215-2-XXXXX	13"-18.5" (33-47 cm)	XS/S
11-0215-3-XXXXX	18.5"-23.5" (47-60 cm)	M/L
11-0215-4-XXXXX	23.5"-29.5" (60-75 cm)	XL/XXL
11-0215-5-XXXXX	29.5"-31" (75-79 cm)	XXXL

Replacement Under Sleeve

11-0541-2	13" -18.5" (33-47 cm)	XS/S
11-0541-3	18.5"-21" (47-53.25 cm)	M
11-0541-4	21"-23.5" (53.25-59.5 cm)	L
11-0541-5	23.5"-26.5" (59.5-67.25 cm)	XL
11-0541-6	26.5"-29.5" (67.25-75 cm)	XXL
11-0541-7	29.5"-32" (75-81.25 cm)	XXXL

Color Information = XXXXX

Orange	14000	Red	01000	Green	11000
Pink	09000	Yellow	08000	Blue	02000
Grey	n/a (If color is not specified web will be grey)				

*Measurements listed below are the circumference taken around the thigh, 6" above mid-patella.





PROCARE®

Hinged Knee Support

Compressive neoprene knee support is available with and without buttress. Universal buttress may be adjusted to allow for medial, lateral, inferior, or superior patella stabilization. Sturdy removable dual axis hinges and superior/inferior contact closure straps provide medial/lateral ligament support.

Indication: Ideal for moderate medial/lateral knee instability, mild ACL/PCL sprain, osteoarthritis, patellofemoral dysfunction.

3/16" w/Open Pop.	1/8" w/Reinf. Pad	Circum.*	Size
79-82152	79-82162	13.5" – 15.5"	XS
79-82153	79-82163	15.5" – 18"	S
79-82155	79-82165	18" – 20.5"	M
79-82157	79-82167	20.5" – 23"	L
79-82158	79-82168	23" – 25.5"	XL
79-82159	79-82169	25.5" – 28"	2XL
79-82159-10	79-82739-10	28" – 30.5"	3XL
	79-82739-11	30.5" – 33"	4XL

*Measurement taken 6" above mid-patella.

1/8" w/Univ. Butt.	1/8" w/Open Pop.	Circum.*	Size
79-82732	79-82742	13.5" – 15.5"	XS
79-82733	79-82743	15.5" – 18"	S
79-82735	79-82745	18" – 20.5"	M
79-82737	79-82747	20.5" – 23"	L
79-82739	79-82749	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Reddie® Brace

1/8" neoprene wraparound configuration and strap tabs provide ease of application. Posterior strap adjustment allows for one-time adjustment and proper positioning of hinges. Removable dual axis polycentric hinges provide medial/lateral support. Available with open popliteal only.

Indication: Ideal for mild MCL and/or LCL sprains, mild arthritis, post arthroscopic meniscectomy, or mild global instabilities.

	Circum.*	Size
79-82392	13.5" – 15.5"	XS
79-82393	15.5" – 18"	S
79-82395	18" – 20.5"	M
79-82397	20.5" – 23"	L
79-82398	23" – 25.5"	XL
79-82399	25.5" – 28"	2XL
79-82399-10	28" – 30.5"	3XL

*Measurement taken 6" above mid-patella.



PROCARE®

ACTION™ Neoprene Brace/Wrap - 1/8"

Available in black 1/8" neoprene only. Front opening for easy on/off and custom fitting. Dual axis motion control hinges provide convenient dial-in adjustments with flexion/extension stops at 0°, 15°, 30°, 60° and 90°. Fully adjustable superior/inferior contact closure straps. Available in 13" or 17" length.

Indication: Ideal for mild to moderate ACL/PCL sprains; osteoarthritis, collateral ligament sprains and strains.

13"	17"	Circum.*	Size
79-94402	79-94412	13.5" – 15.5"	XS
79-94403	79-94413	15.5" – 18"	S
79-94405	79-94415	18" – 20.5"	M
79-94407	79-94417	20.5" – 23"	L
79-94408	79-94418	23" – 25.5"	XL
79-94409	79-94419	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PRO CARE®

ACTION™ Neoprene Brace/Sleeve - 1/8"

Available in choice of 13" or 17" lengths. Black 1/8" contoured neoprene sleeve with dual axis motion control hinges. 0°, 15°, 30°, 60° and 90° hinge adjustments. Molded condyle pads and four adjustable circumferential straps.

Indication: Ideal for mild to moderate ACL/PCL sprains; osteoarthritis, collateral ligament sprains, patellofemoral malalignment.

13"	17"	Circum.°	Size
79-94292	79-94302	13.5" – 15.5"	XS
79-94293	79-94303	15.5" – 18"	S
79-94295	79-94305	18" – 20.5"	M
79-94297	79-94307	20.5" – 23"	L
79-94298	79-94308	23" – 25.5"	XL
79-94299	79-94309	25.5" – 28"	2XL
79-94299-10		28" – 31"	3XL

[°]Measurement taken 6" above mid-patella.



PRO CARE®

WeekENDER™ Recreational Activity Brace

Lightweight, rigid, low profile. May be used for post-op, rehabilitation or functional application. Dual axis motion control hinges with foam/chamois condyle pads provide full and limited range-of-motion with flexion/extension stops at 0°, 15°, 30°, 60° and 90°.

Indication: Recommended for mild medial/lateral instability, including ACL/PCL laxity.

Right	Left	Circum.°	Size
79-94353	79-94363	14.5" – 18.5"	S
79-94355	79-94365	18.5" – 21"	M
79-94357	79-94367	21" – 23.5"	L
79-94358	79-94368	23.5" – 26.5"	XL
79-94359	79-94369	26.5" – 29"	2XL

Replacement Strap with Standard Foam Pads

79-94373	Thigh, Calf & Condyle Pads	S
79-94375	Thigh, Calf & Condyle Pads	M
79-94377	Thigh, Calf & Condyle Pads	L
79-94378	Thigh, Calf & Condyle Pads	XL
79-94379	Thigh, Calf & Condyle Pads	2XL

Chamois Replacement Pad Set

79-94323	Thigh, Calf & Condyle Pads	S
79-94325	Thigh, Calf & Condyle Pads	M
79-94327	Thigh, Calf & Condyle Pads	L
79-94328	Thigh, Calf & Condyle Pads	XL
79-94329	Thigh, Calf & Condyle Pads	2XL

[°]Measurement taken 6" above mid-patella.



PRO CARE®

EveryDAY™ Daily Activity Brace

Semi-rigid, lightweight combined instability brace affords maximum comfort and protection. Dual axis motion control hinges with durable shock absorbing condyle pads provide convenient dial-in adjustment with full and limited range-of-motion settings. Flexion/extension stops at 0°, 15°, 30°, 60° and 90°.

Indication: Recommended for mild medial/lateral instability, including ACL/PCL laxity.

	Circum.°	Size
79-94313	14.5" – 18.5"	S
79-94315	18.5" – 21"	M
79-94317	21" – 23.5"	L
79-94318	23.5" – 26.5"	XL

[°]Circumference measurement taken 6" above mid-patella.

ANKLE BRACES & SUPPORT – ANKLE STIRRUPS

AIRCAST®**Ankle Stirrups**

Aircast pioneered functional management of ankle sprains with the Air-Stirrup® Ankle Brace. Setting the standard of care for ankle sprains, the Air-Stirrup's design features the patented Duplex™ aircell technology that enhances circulation and reduces swelling.

Air-Stirrup® Ankle Brace

Anatomically designed semi-rigid shells lined with the patented Duplex™ aircell system for support and graduated compression during ambulation. The compression promotes efficient edema reduction in addition to helping accelerate rehabilitation from ankle sprains.

Indication: Acute injury, Post op; Ankle sprains grade I, II, III; Chronic instability.

Left	Right	Height	Size
02JL	02JR	6"	Ped.
02CL	02CR	8.75"	S
02BL	02BR	9"	M
02AL	02AR	10.5"	L
<i>Maximum Ankle Circumference 16.5"</i>			
0213	Heel assembly replacement		
0130	Replacement inflation tubing		





AIRCAST®

Ankle Sprain Care Kit™

Complete system for the functional management of an ankle sprain. Includes an Air-Stirrup Ankle Brace, ankle wrap, cold pack, exercise band, and a DVD plus booklet showing instructions for treatment and rehabilitation.

Indication: Acute injury; Post-op; Ankle sprains grade I, II; Chronic instability.

Left	Right	Height	Size
02CLK	02CRK	8.75"	S
02BLK	02BRK	9"	M
02ALK	02ARK	10.5"	L



AIRCAST®

Air-Stirrup Universe™

Anatomically designed semi-rigid shells lined with the patented Duplex™ aircell system for support and graduated compression during ambulation.

Indication: Ankle sprains; Chronic ankle instability.

	Height	Qty.	Size
02E	9.5"		Univ.
02E12	9.5"	12 pk.	Univ.
0212	Black heel assembly replacement		



AIRCAST®

Air-Stirrup® Universe™ Care Kit

Complete system for the functional management of an ankle sprain. Includes an Air-Stirrup Universe Ankle Brace, ankle wrap, cold pack, exercise band, and a DVD plus booklet showing instructions for treatment and rehabilitation.

Indication: Acute injury; Chronic instability.

	Height	Size
02EK	9.5"	Univ.



AIRCAST®

AirSport™

Incorporates clinically proven semi-rigid shell and aircells to provide comfort and support. Additional compression and stabilization is provided by anterior talofibular cross strap and integral forefoot and shin wraps. The unique "step-in" design (toes first inserted into the back of the brace) and automatic heel adjustment make the AirSport extremely fast and easy to put on.

Indication: Ankle sprains; Chronic instability; Prophylaxis.

Left	Right	Shoe Size		Size
		Men	Women	
02MXSL	02MXSR	Up to 5	Up to 5	XS
02MSL	02MSR	5.5 – 7	5.5 – 8.5	S
02MML	02MMR	7.5 – 11	9 – 12.5	M
02MLL	02MLR	11.5 – 13	13 – 14.5	L
02MXLL	02MXLR	13.5 +	15 +	XL



AIRCAS T

Air-Stirrup Plus

Combines the clinically proven Air-Stirrup Ankle Brace with an ankle wrap and a comprehensive instruction booklet providing information on treatment and rehabilitation — integral components of ankle sprain management.

Indication: Acute injury; Post-op; Ankle sprains grade I, II, III; Chronic instability.

Left	Right	Height	Size
02CLP	02CRP	8.75"	S
02BLP	02BRP	9"	M
02ALP	02ARP	10.5"	L
02K110	Ankle wrap replacement		



PROCARE

Surround® Gel Ankle

Rigid thermoplastic shells with adjustable heel strap to help prevent inversion/eversion. Air and gel bladder may be used for cold therapy to treat tender or swollen ankles. Fits right or left ankle.

Indication: Acute injury; Post-op; Ankle sprains grade I, II, III; Chronic instability.

	Height	Size
79-97863	8.5"	S
79-97865	9"	M
79-97867	10"	L



PROCARE

Surround® FLOAM™ Ankle

Provides stability to the ankle by minimizing ankle rotation. Comprised of two polypropylene outer shells connected with soft, adjustable foot piece. A unique FLOAM™ bladder provides uniform compression and a customized fit. FLOAM™ has no memory, and therefore prohibits pressure from being displaced to bony prominences. Fits right or left ankle.

Indication: Acute ankle sprains; chronic ankle instabilities; osteoarthritis.

	Height	Size
79-81193	6"	Ped.
79-81195	9"	M
79-81197	10"	L



PROCARE

Surround® Air Ankle

Injection molded shells help minimize ankle rotation and inversion/eversion without restricting plantar or dorsi-flexion. Two full-length pre-inflated air chambers lined with soft foam help reduce swelling and enhance circulation. A soft, plush liner covers exterior of bladders to help maximize patient comfort and prevent skin irritation. Air valve at proximal end of bladder allows for simple adjustment of air when desired by using filling tube. Soft foam heel pad provides easy adjustment for proper fit.

Indication: Acute ankle sprains and strains, chronic instabilities, and prophylactic use.

Right	Left	Height	Size
79-81705	79-81715	9"	M
79-81707	79-81717	10"	L



PRO CARE®

Excelerator™ Stirrup Ankle Splint

Constructed of a three-layer foam system for added comfort and bi-lateral rigid shells to help prevent inversion and eversion. Universal sizing fits right or left ankle. Contact closure straps for ease of application and removal.

Indication: Acute ankle sprain; chronic ankle instability; subtalar joint pain, osteoarthritis.

	Height	Size
79-81323	8"	S
79-81327	9.5"	L



PRO CARE®

Stirrup Ankle Support

Contoured, high impact molded shells help prevent inversion/eversion and provide protection for acute ankle sprains and chronic instabilities. Removable, padded nylon/neoprene liner is intended to provide superior fit and comfort around ankle joint. Contoured shells for left and right ankle.

Indication: Acute ankle sprains; ankle instability.

Right	Left	Height	Size
79-81250	79-81260	10"	Univ.
79-81251	79-81261	8.5"	Univ.



PRO CARE®

Universal Ankle Brace

Rigid thermoplastic one-piece outer shell helps prevent inversion/eversion. Cushioned liner incorporates neoprene forefoot closure. Elastic proximal strap provides additional compression.

Indication: Acute and chronic ankle instabilities, as well as, prophylactic use.

	Height	Size
79-81330	10"	Univ.



AIRCAST®

Leg Brace

Designed to provide functional management of stress fractures and the graduated management of certain stable fractures of the lower leg. Features anatomically designed shells lined with the patented Duplex aircell system. The Leg Brace can be ordered alone or with an optional Anterior Panel for additional tibial protection. A sock is included with each brace.

Indication: Stress fracture; Stable fracture (graduated fracture management).

Left	Right	Height	Size
03CL	03CR	13"	S
03DL	03DR	13"	S*
03AL	03AR	15.5"	Standard
03BL	03BR	15.5"	Standard*

*With Anterior Panel.

ANKLE BRACES & SUPPORT – SOFT ANKLE & SPECIALTY

PRO CARE®

Stabilized Ankle Support

Low profile design constructed with non-stretch nylon for comfort, performance and durability. Lace-up design and figure 8 strapping locks the calcaneus into place & controls inversion/eversion. A full circumferential elastic strap secures all strapping to provide superior tibia/fibula compression and stability. Removable medial/lateral stays for enhanced proprioception and stability. Fits easily into athletic or street shoe. Universal to fit right and left.

Indication: Ankle sprains; Chronic ankle instability.

	Circum.°	Size
79-81352	10" – 11"	XS
79-81353	11" – 12"	S
79-81355	12" – 13"	M
79-81357	13" – 14"	L
79-81358	14" – 15"	XL
79-81359	15" – 16"	2XL
79-81359-10	16" – 17"	3XL

°With athletic sock on, measure circumference of ankle from base of heel around ankle.





AIRCAST®

A60™ Ankle Support

Prophylactic support, protection, and comfort with simplified application. The design incorporates a stabilizer molded at a 60-degree angle to help guard against ankle sprains and prevent rollover. Manufactured with breathable material, for cool comfort. Light-weight anatomic design easily fits in athletic footwear. Easily adjustable application. Single strap replaces time-consuming lacing.

Indication: Chronic instability; Prophylaxis.

Left	Right	Shoe Size Male	Shoe Size Female	Size
02TSL	02TSR	Up to 7	Up to 8.5	S
02TML	02TMR	7.5 – 11.5	9 – 13	M
02TLL	02TLR	12+	13.5+	L



PROCARE®

Lace-Up Ankle Brace

Sturdy canvas stirrup design with removable, plastic medial/lateral stays provide support for acute ankle ligament sprains and chronic instabilities. Fits easily into any athletic shoe.

Indication: Acute ankle sprains; chronic ankle instability.

	Circum.*	Size
79-81312	5.5" – 7"	XS
79-81313	7" – 8.5"	S
79-81315	8.5" – 10"	M
79-81317	10" – 11.5"	L
79-81318	11.5" – 13.5"	XL

*Ankle circumference measurement.



PROCARE®

Kallassy Ankle Support®

Constructed of nylon lined neoprene which provides warmth and compression. Non-stretch lateral strapping system allows full range-of-motion. Counterbalance strap locks ankle in neutral position. Low profile design allows use in any athletic shoe. Black or white version available.

Indication: Ideal for chronic ankle instability, acute ankle sprain, prophylactic ankle support.

Right	Left	Color	Circum.*	Size
79-81411	79-81421	White	6" – 8.5"	2XS
79-81412	79-81422	White	8.5" – 9.5"	XS
79-81413	79-81423	White	9.5" – 10.5"	S/M
79-81417	79-81427	White	10.5" – 11.5"	L/XL
79-81418	79-81428	White	12.5" – 13.5"	2XL
79-81472	79-81482	Black	8.5" – 9.5"	XS
79-81473	79-81483	Black	9.5" – 10.5"	S/M
79-81477	79-81487	Black	10.5" – 11.5"	L/XL
79-81478	79-81488	Black	12.5" – 13.5"	2XL

*With athletic sock on measure circumference of ankle from base of heel around ankle.



PROCARE®

DS Ankle Wrap

Elastic/Lycra construction provides moderate compressive support without restriction. Closed heel design with low profile elastic contact closure figure-eight straps for medial/lateral and arch support. Lightweight material allows for ease of application. Fits right or left foot.

Indication: Mild ankle sprains and strains; edema control; tendonitis; osteoarthritis.

	Circum.*	Size
79-81362	6.5" – 7.5"	XS
79-81363	7.5" – 8.5"	S
79-81365	8.5" – 9.5"	M
79-81367	9.5" – 11"	L
79-81368	11" – 12.5"	XL
79-81369	12.5" – 14"	2XL

*Measurement taken 2" above medial/lateral malleoli.



PROCARE®

Double Strap Ankle Support

Durable mesh/elastic construction provides compressive support without restriction. Low profile elastic contact closure figure-eight straps provide medial/lateral and arch support without adding bulk. Fits right or left foot.

Indication: Mild ankle sprains and strains; edema control; tendonitis; medial/lateral instabilities.

	Circum.°	Size
79-81372	6.5" – 7.5"	XS
79-81373	7.5" – 8.5"	S
79-81375	8.5" – 9.5"	M
79-81377	9.5" – 11"	L
79-81378	11" – 12.5"	XL
79-81379	12.5" – 14"	2XL

[°]Measurement taken 2" above medial/lateral malleoli.



PROCARE®

Deluxe Ankle Support

Vinyl, foam, tricot construction with loop/lock contact closure and preformed aluminum medial/lateral stays. Fits right or left ankle.

Indication: Mild sprains and chronic ankle instabilities.

	Circum.°	Size
79-81063	7.75" – 8.5"	S
79-81065	8.5" – 10"	M
79-81067	10" – 11.5"	L
79-81068	11.5" – 13.5"	XL

[°]Ankle circumference measurement.



PROCARE®

Sport Ankle Splint

Quality, lace-up design with formed aluminum stays provide medial/lateral support. Durable perforated suede/vinyl construction.

Indication: Acute ankle sprains; chronic ankle instability.

	Circum.°	Size
79-81092	5.5" – 7"	XS
79-81093	7" – 8.5"	S
79-81095	8.5" – 10"	M
79-81097	10" – 11.5"	L
79-81098	11.5" – 13.5"	XL

[°]Ankle circumference measurement.



PROCARE®

Ankle Sleeve

1/8" neoprene sleeve provides warmth and compressive support for mild ankle sprains and soft tissue injuries. Easy on and off.

	Circum.°	Size
79-82302	7.25" – 7.75"	XS
79-82303	7.75" – 8.5"	S
79-82305	8.5" – 9.5"	M
79-82307	9.5" – 10.5"	L
79-82308	10.5" – 11.5"	XL

[°]Ankle circumference measurement.



AIRCAST®

AirLift™ PTTD Brace

Designed for the treatment of posterior tibial tendon dysfunction (PTTD), or for early signs and symptoms of adult acquired flat foot. Depending on the severity of the condition, it can be prescribed as part of a conservative treatment to stabilize the condition and help prevent its degeneration, or post-surgically and during rehabilitation.

Indication: Posterior Tibial Tendon Dysfunction in Stage I and early Stage II; Adult acquired flat foot; post-operative use; rehabilitation.

Left	Right	Shoe Size Male	Shoe Size Female	Size
02PSL	02PSR	Up to 7	Up to 8.5	S
02PML	02PMR	7 – 11	9 – 12.5	M
02PLL	02PLR	11.5+	13+	L
02PI70	Replacement inflation hand bulb			



PROCARE®

Elastic Ankle Support

Constructed of durable, breathable cotton elastic. Pull on design with open heel design for mild compression and support.

Indication: Mild ankle sprains and strains; edema; osteoarthritis.

	Circum.*	Size
79-81123	7.75" – 8.5"	S
79-81125	8.5" – 9.5"	M
79-81127	9.5" – 10.5"	L
79-81128	10.5" – 11.5"	XL
79-81129	11.5" – 12.5"	2XL

*Ankle circumference measurement.



PROCARE®

Bunion Splint

Flannel lined foam with brushed nylon backing for contact closure. Malleable stay and strap may be adjusted to accommodate indication.

Indication: Ideal for bunion alignment to help reduce pain.

Right	Left	Shoe Size Male	Shoe Size Female	Size
79-96133	79-96143	≤6	≤7.5	S
79-96135	79-96145	6.5 – 8.5	8 – 10	M
79-96137	79-96147	9+	10+	L



AIRCAST®

Bunion Aid®

Hinged bunion splint made with adjustable straps and a comfortable gel cushion to protect irritated tissue, distributes pressure, and provides corrective support. Ideal for bunion pain; mild or moderate malpositioning of the big toe; maintain post-surgical fixation of the big toe. Wear day or night, while walking or resting. One size fits most, left or right foot.

Indication: Bunion pain; mild or moderate malpositioning of the big toe; maintain post-surgical fixation of the big toe

	Description	Size
81-03H	Bunion Aid®	Universal
81-03E	Bunion Aid® Plus Replacement Set	Universal



PROCARE®

Silicone Heel Cups

Molded medical-grade soft silicone to help reduce shock and vibration loads to joints, to address Plantar Fasciitis, heel spurs and Achilles tendonitis. Softer center in heel region to help reduce strain to area sensitive to pain. Sold in pairs.

Indication: Plantar Fasciitis; heel spurs; Achilles tendonitis.

	Shoe Size Male	Shoe Size Female	Size
79-81102	Children/youth		XS
79-81103	5 – 9	5.5 – 9.5	S/M
79-81107	9.5+	10+	L/XL



PROCARE®

Tuli's® Heel Cups

Unique waffle design absorbs shock and pressure at heel strike. Standard style used for all activities. Pro style has wider, double-waffle design to absorb pressure from high impact sports. Sold in pairs.

Indication: Shock absorption; pressure relief.

Standard (Yellow)	Pro (Green)	Weight	Size
	79-72273	Youth	S
79-72280	79-72275	Adult (< 175 lbs.)	M
79-72281	79-72277	(175 lbs.+)	L

PLANTAR FASCIITIS



AIRCAST®

AirHeel/DNS Care Kit

For daytime, the unique Aircast AirHeel is a clinically-proven solution that provides comfort and massaging relief for Plantar Fasciitis. Two interconnected air cells located under the foot arch and in the back of the heel apply pulsating compression with every step, help reduce swelling and discomfort, and enhance circulation. The AirHeel is made of lightweight, breathable fabric, and is simple to apply and to adjust with one strap. For nighttime, the Aircast Dorsal Night Splint delivers both tailored flexion and optimal comfort for pain relief from Plantar Fasciitis. Thermoformed shell liner and anatomically designed footpad provide optimal comfort. The dorsiflexion strap allows for a custom stretch. Shipped with optional toe strap that secures fit within forefoot area. Designed for ease of use. Fits right or left.

Indication: Plantar Fasciitis.

	Shoe Size Male	Shoe Size Female	Size
09PKS	Up to 7	Up to 8	S
09PKM	7 - 10	9 - 11	M
09PKL	10+	11+	L



AIRCAST®

AirHeel™

Specifically designed to treat plantar fasciitis, Achilles tendonitis, and heel pain. Two interconnected aircells located under the foot arch and in back of the heel, apply pulsating compression with every step, help reduce swelling and discomfort, and enhance circulation. The AirHeel is made of lightweight breathable fabric and is simple to apply and adjust with one strap.

Indication: Plantar fasciitis; Achilles tendonitis; Heel pain.

	Shoe Size Male	Shoe Size Female	Size
09AS	Up to 7	Up to 8.5	S
09AM	7.5 – 11	8.5 – 12.5	M
09AL	11.5+	13+	L



AIRCAST®

Dorsal Night Splint

Provides tailored dorsiflexion and optimal comfort for pain relief from Plantar Fasciitis. Dorsiflexion adjustment strap provides custom stretch. Thermoformed liner and anatomically designed forefoot pad provide maximum comfort. Comes with optional toe assist strap to secure toes within forefoot area. Anterior shell design with flexible hinge allows for ease of application. Fits right or left foot.

Indication: Plantar Fasciitis.

	Shoe Size Male	Shoe Size Female	Size
09DS	5 – 9	6 – 10	S/M
09DL	10 – 14	11 – 15	L/XL



PROCARE®

Plantar Fasciitis Night Splint

Allows for controlled application of a prolonged stretch on the Plantar Fascia and Achilles Tendon. Lightweight, sturdy, one-piece flexible shell. Adjustable bilateral dorsi-flexion straps. Cool terry liner and padded calf and ankle straps for additional comfort. Shipped with toe wedge for additional stretch on plantar fascia. Fits right or left.

Indication: Plantar Fasciitis.

	Shoe Size Male	Shoe Size Female	Size
79-97753	≤ 6	≤ 7	S
79-97755	6.5 – 8.5	7.5 – 9.5	M
79-97757	9 – 11.5	10 – 12.5	L
79-97758	12+	13+	XL



PROCARE®

Dorsiwedge™ Night Splint

Soft padded, contoured posterior shell provides passive dorsi-flexion to address Plantar Fasciitis and Achilles Tendonitis. Removable liner may be washed for patient hygiene. Fits right or left foot. Shipped with full-length toe wedge for additional stretch on plantar fascia.

Indication: Plantar Fasciitis.

	Shoe Size Male	Shoe Size Female	Size
79-81403	≤ 6	≤ 6.5	S
79-81405	6.5 – 9.5	7 – 10	M
79-81407	10 – 12	10.5 – 12.5	L
79-81408	12+	13+	XL

POST-OP AND CAST SHOES

PROCARE®**Squared Toe Post-Op Shoe**

The square toe design acts as a bumper and provides additional room and comfort for the patient. Features high ankle strapping that maintains foot in proper location and a rigid rocker sole for ease of ambulation and pressure reduction. Design allows dressing to be kept clean and toes dry, providing more protection following post-operative procedures and forefoot trauma.

Indication: Post-op; forefoot trauma.


	Shoe Size Male	Shoe Size Female	Size
79-81232	3 – 5	4 – 6	XS
79-81233	5.5 – 7	6.5 – 8	S
79-81235	7.5 – 9	8.5 – 10	M
79-81237	9.5 – 12	10.5 – 13	L
79-81238	11.5+	13+	XL

**PROCARE®****IMPAX Grid Insole**

Designed for use with the ProCare Squared Toe Post-Op shoe, the dual-laminate insole is designed with Impax Grid technology to redistribute load and relieve pressure. Includes one Impax insole.

Indication: Redistribute load; relieve pressure.

	Shoe Size Male	Shoe Size Female	Size
79-81332	3 – 5	4 – 6	XS
79-81333	5.5 – 7	6.5 – 8	S
79-81335	7.5 – 9	8.5 – 10	M
79-81337	9.5 – 12	10.5 – 13	L
79-81338	12.5+	13.5+	XL





PROCARE®

Post-Op Shoe

Rigid E.V.A. rocker sole and contoured heel for added patient comfort. Low-profile foam/nylon mesh upper with tricot lining. Protective tongue provides maximum comfort. Available in loop-lock design.

Indication: Post foot surgery; metatarsal injuries; soft tissue procedures.

	Shoe Size	Size
Male Loop Lock Closure		
79-90182	5 – 7	XS
79-90183	7 – 9	S
79-90185	9 – 11	M
79-90187	11 – 13	L
79-90188	13+	XL
Female Loop Lock Closure		
79-90192	2 – 4	XS
79-90193	4 – 6	S
79-90195	6 – 8	M
79-90197	8 – 10	L



PROCARE®

Med/Surg Shoe

Designed for a variety of post-op and post-trauma applications. Soft foam mesh upper with padded tricot lining. Semi-flexible skid resistant sole and comfortable contoured heel. Large overlapping toe area accommodates for post-operative dressing and swelling. Loop lock contact closure straps. Fits right or left foot. Latex free.

Indication: Post-op; Post-Trauma applications.

Male	Shoe Size Male	Female	Shoe Size Female	Size
79-81132	5 – 7	79-81142	2 – 4	Ped./XS
79-81133	7 – 9	79-81143	4 – 6	S
79-81135	9 – 11	79-81145	6 – 8	M
79-81137	11 – 13	79-81147	8 – 10	L
79-83118	13+			XL



PROCARE®

Off-Loading Diabetic Shoe

Specifically designed to relieve pressure on ulcerated areas of the foot during the treatment of diabetic ulcers. Squared toe box and wide foot bed leave ample room for dressing. Removable reinforced toe cover protects from external aggression and allows easy access to the treatment area. Seamless design helps eliminating potential for skin breakdown. Impax Grid tri-laminate insole is designed to help redistribute load away from ulcers by simply removing Impax foam blocks. Fits Right or Left.

Indication: Pressure relief during treatment of diabetic ulcers.

	Shoe Size Male	Shoe Size Female	Size
79-81513	4 – 6	6.5 – 8.5	S
79-81515	6.5 – 8.5	9 – 11	M
79-81517	9 – 11	11.5 – 12.5	L
79-81518	11.5 – 14	13+	XL
Impax Diabetic Shoe Insole Accessory			
79-81523			S
79-81525			M
79-81527			L
79-81528			XL



PROCARE®

Rocker Cast Boot

Non-slip rocker sole to allow for more natural gait during ambulation. Sandal style with durable canvas upper. Forefoot and heel loop-lock contact closure for ease of application. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Color	Size
79-81110	4.5"	2.5"	Ped print	3XS
79-81111	5.5"	3.0"	Ped print	2XS
79-81112	7.25"	3.5"	Black	XS
79-81113	9.25"	4.25"	Black	S
79-81114	10.5"	4.5"	Black	S/M
79-81115	11.0"	4.75"	Black	M
79-81116	12.0"	5.0"	Black	M/L
79-81117	12.5"	5.25"	Black	L
79-81118	13.5"	6.0"	Black	XL



PROCARE®

Deluxe Cast Boot Sandal

Open toe/heel sandal with rubber rocker sole. Vinyl upper with forefoot and heel contact closure. Hook and loop closure for ease of application. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Size
79-81012	7.25"	3.5"	XS
79-81013	9.25"	4.25"	S
79-81015	11.0"	4.75"	M
79-81017	12.5"	5.25"	L
79-81018	13.5"	6.0"	XL



PROCARE®

Deluxe Cast Boot Weatherproof

Closed toe and heel protects cast from harsh weather. Vinyl upper with rubber rocker sole to help cushion cast during ambulation. Forefoot and heel contact closure. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Size
79-81022	7.25"	3.5"	XS
79-81023	9.25"	4.25"	S
79-81025	11.0"	4.75"	M
79-81027	12.5"	5.25"	L
79-81028	13.5"	6.0"	XL



PROCARE®

Deluxe Cast Boot Open Toe

Closed heel rocker sole with open toe design provides additional breathability. Constructed of durable vinyl materials that wipe clean with soap and water. Rocker sole for natural gait during ambulation. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Size
79-81042	7.25"	3.5"	XS
79-81043	9.25"	4.25"	S
79-81045	11.0"	4.75"	M
79-81047	12.5"	5.25"	L
79-81048	13.5"	6.0"	XL



PROCARE®

Beach Boy Cast Boot

Sandal style boot with rubber rocker sole. Sturdy elastic straps with contact closure. Open air design ideal for warmer climates. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Size
79-81072	7.5"	3.75"	XS
79-81073	9.0"	4.25"	S
79-81075	10.5"	4.5"	M
79-81077	12.0"	5.0"	L



PROCARE®

Rocker Cast Shoe w/Flex Sole

Constructed of durable canvas upper with flexible rocker sole for more natural gait during ambulation. Contact closure forefoot and heel straps for ease of application. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Size
79-81222	7.25"	3.5"	Ped
79-81223	9.25"	4.25"	S
79-81225	11.0"	4.75"	M
79-81227	12.5"	5.25"	L
79-81228	13.5"	6.0"	XL



PROCARE®

Orthowedge

Foam/mesh upper with contoured heel design. Wedge construction helps reduce pressure on the forefoot and keeps weight behind the MP joints to allow ulcerations under metatarsal heads and toes to heal faster.

Indication: Ideal for diabetic ulcers; foot surgery.

	Shoe Size Female	Shoe Size Male	Size
79-98442	2 – 4	3 – 5	XS
79-98443	4 – 6	5 – 7	S
79-98445	6 – 8	7 – 9	M
79-98447	8 – 10	9 – 11	L
79-98448	10+	11-13	XL



PROCARE®

Darco™ SlimLine Cast Boot

Rocker bottom with weatherproof, stretchable upper and contact closure.

Indication: For use with lower leg casted fractures.

	Shoe Size Male	Shoe Size Female	Size
79-98432		3 – 5	XS
79-98433		5.5 – 8	S
79-98435	6 – 8	8.5 – 10	M
79-98437	10.5 – 12	10.5 – 12	L
79-98438	12.5 – 14	12+	XL

WALKING BRACES

AIRCAST®

Elite

The Elite has a durable, semi-rigid shell that helps support the limb while providing full-shell protection. Overlapping Duplex aircells line the shell providing intermittent pneumatic compression for 3x faster edema reduction¹. The front panel is pre-inflated and combined with the overlapping aircells, provides a 'total contact' fit. SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities.

Indication: Stable fracture of foot/ankle/lower leg; Severe ankle sprain; Post-operative immobilization; Reduction of edema (swelling) 3x faster¹.

	Shoe Size		Size
	Male	Female	
01EP-XS	Up to 4	Up to 5	XS
01EP-S	4-7	5-8	S
01EP-M	7-10	8-11	M
01EP-L	10-13	11-15	L
01EP-XL	13+	15+	XL
Accessories			
01WXT-x	Weather Cover, Tall		
01WXS-x	Weather Cover, Short		
0130A-x	Hygiene Cover		
01T-x	Toe Cover		
01P414	Brace-Lok™		
01K-x	Heel Wedges		
0129P	Tube Stretch Sock		XS
0129A	Tube Stretch Sock		Univ.
Replacement Parts			
01E102x	Liner & Insole Replacement Kit, Elite		
01E104x	Strap Replacement Kit, Elite & Standard		
01E1050	Extension Strap Kit		Univ.

1. Conservative Therapy for Acute Lateral Ligament Lesions – Single Chamber vs. Two-Chamber Orthosis Systems; Schmidt, Mainers, Reintges, Lipke, Benesch, Gerngross - Surgery Dept of the Federal Arm Hospital, Ulm, Germany, 1999





AIRCAS[®]

Standard

The Standard has a durable, semi-rigid shell that supports the limb while providing protection. Housed within the shell are two adjustable distal aircells that provide compression and support the malleoli. SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities.

Indication: Stable fracture of foot/ankle/lower leg; Severe ankle sprain; Post-operative immobilization.

Shoe Size			
	Male	Female	Size
01EF-XS	Up to 4	Up to 5	XS
01EF-S	4-7	5-8	S
01EF-M	7-10	8-11	M
01EF-L	10-13	11-15	L
01EF-XL	13+	15+	XL
Accessories			
	Weather Cover, Tall		
01WXS-x	Weather Cover, Short		
0130A-x	Hygiene Cover		
01T-x	Toe Cover		
01P414	Brace-Lok™		
01K-x	Heel Wedges		
0129P	Tube Stretch Sock		XS
0129A	Tube Stretch Sock		Univ.
Replacement Parts			
01E102x	Liner & Insole Replacement Kit, Elite		
01E104x	Strap Replacement Kit, Elite & Standard		
01E1050	Extension Strap Kit		Univ.



AIRCAS[®]

Short

The Short has a durable, semi-rigid shell that helps support the limb while providing protection. The pre-inflated front panel and two adjustable distal aircells, housed within the shell, provide compression and support. SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities.

Indication: Metatarsal fracture; Stress and stable fracture of the foot; Bunionectomy, hallux valgus; Soft tissue injury; Severe ankle sprain.

Shoe Size			
	Male	Female	Size
01ES-XS	Up to 4	Up to 5	XS
01ES-S	4-7	5-8	S
01ES-M	7-10	8-11	M
01ES-L	10-13	11-15	L
01ES-XL	13+	15+	XL
Accessories			
01WXT-x	Weather Cover, Tall		
01WXS-x	Weather Cover, Short		
0130A-x	Hygiene Cover		
01T-x	Toe Cover		
01P414	Brace-Lok™		
01K-x	Heel Wedges		
0129P	Tube Stretch Sock		XS
0129A	Tube Stretch Sock		Univ.
Replacement Parts			
01E100x	Liner & Insole Replacement Kit, Elite		
01E100x	Strap Replacement Kit, Elite & Standard		
01E1050	Extension Strap Kit		Univ.



AIRCAST®

XP Walker™ (extra pneumatic)

The XP Walker has a lightweight, durable semi-rigid shell that helps support the limb while providing full-shell protection. Lining the shell are overlapping Duplex aircells providing intermittent pneumatic compression for edema reduction and callus formation. These aircells can be custom inflated for a “total contact” fit. Includes two socks and a hand bulb — Can be used with the Brace Lok™ system.

Indication: Stable foot and/or ankle fracture; severe ankle sprain; post-operative use.

	Shoe Size		Size
	Male	Female	
01P-P	Up to 4	Up to 5	Ped.
01P-S	4 – 7	5 – 8	S
01P-M	7 – 10	8 – 11	M
01P-L	10 – 13	11 – 15	L
01P-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
0193	Pneumatic Care Kit (2 socks, 1 hand bulb)		Univ.
0120	Replacement Hand Bulb		



AIRCAST®

FP Walker™ (foam pneumatic)

The FP Walker has a lightweight, durable semi-rigid shell that supports the limb while providing protection. Housed within the shell are two adjustable distal aircells that provide compression and support the malleoli. Includes a hand bulb.

Indication: Stable foot and/or ankle fracture; severe ankle sprain; post-operative use.

	Shoe Size		Size
	Male	Female	
01F-P	Up to 4	Up to 5	Ped.
01F-S	4 – 7	5 – 8	S
01F-M	7 – 10	8 – 11	M
01F-L	10 – 13	11 – 15	L
01F-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
1120SC	Replacement Hand Bulb		Univ.



AIRCAST®

SP Walker™ (short pneumatic)

The SP Walker has a lightweight, durable semi-rigid shell with a pre-inflated anterior aircell. Two adjustable aircells located at the malleoli may be individually inflated for “total contact fit”. The generous foot base has ample room for dressings without sacrificing comfort. Includes a hand bulb.

Indication: Metatarsal fracture; forefoot and midfoot injury; acute or post-operative use; bunionectomy; soft tissue injury; severe ankle sprain.

	Shoe Size		Size
	Male	Female	
01A-P	Up to 4	Up to 5	Ped.
01A-S	4 – 7	5 – 8	S
01A-M	7 – 10	8 – 11	M
01A-L	10 – 13	11 – 15	L
01A-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
1120SC	Replacement Hand Bulb		Univ.



AIRCAST®

XP Diabetic Walker Impax™ Replacement Kit

Designed for use with the XP Diabetic Walker System, the dual-laminate insole is designed with the Impax Grid technology to redistribute load away from ulcers. When ordering, Kit size should correspond to System size (i.e. medium Kit for medium System). Includes: one Impax insole and one stocking.

Indication: For use with the XP Diabetic Walker.

	Size
01PI-S	S
01PI-M	M
01PI-L	L
01PI-XL	XL



AIRCAST®

XP Diabetic Walker™ System

The XP Diabetic Walker incorporates an optimized rocker sole, semi-rigid shells, an innovative dual density insole and four aircells for maximum off-loading and “total contact” fit of the diabetic foot. Two extra insoles (including the dual-laminate Impax™ Grid insole), three stockings and one hand bulb with pressure gauge are included. Can be used with the Brace-Lok™ system.

Indication: Stable foot and/or ankle fracture; post-operative; hindfoot and heel foot protection, neuropathic ulcers, soft tissue injury.

	Shoe Size		Size
	Male	Female	
01PD-S	4 – 7	5 – 8	S
01PD-M	7 – 10	8 – 11	M
01PD-L	10 – 13	11 – 15	L
01PD-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
01PG	Hand Bulb with Pressure Gauge		Univ.



AIRCAST®

XP Diabetic Walker Replacement Kit

Designed for use with the XP Diabetic Walker System, the plastizote insole helps regulate shear stress. When ordering, Kit size should correspond to System size (i.e. medium Kit for medium System). Includes: one insole and one stocking.

Indication: For use with the XP Diabetic Walker.

	Size
01R-S	S
01R-M	M
01R-L	L
01R-XL	XL



AIRCAST®

Brace-Lok™

Compliance monitoring system for the XP Walker and XP Diabetic Walker.

Indication: The Brace-Lok is intended to monitor patient compliance.

	Qty.	Size
01P414	10 pk.	Univ.



AIRCAST®

Heel Wedges

The heel wedges can be placed on top of each other to achieve the desired heel rise (three heel wedges equal 22° of heel raise; two heel supports equal 16° of heel raise; one heel support equals 10° of heel raise).

Left	Right	Size
01K-L	01K-R	Standard (S, M, L)
01K-XLL	01K-XLR	XL



AIRCAST®

Tube Stretch Sock

The Tube Stretch Sock can be used with the XP Walker and the FP Walker. The Polyester Sock can be used with the XP Diabetic Walker System.

Indication: For use with the Aircast walking braces.

	Length	Size
0129P*	16"	Ped.
0129A*	20"	Univ.
Polyester Sock		
4003	30"	Univ.

*Tube Stretch Sock (Photo shown above)



AIRCAST®

Toe Cover

Made of a thin polycarbonate flexible plastic material and covered with a soft nylon foam, this optional adjustable accessory provides the ideal coverage for patients.

	Size
01T-S	Ped./S
01T-L	M/L/XL



AIRCAST®

Hygiene Cover

This easy-to-apply cover is used to provide a hygienic environment by preventing the spread of dirt, dust and germs brought in on a patients' walking brace.

	Size
0130A-PED	Ped./S
0130A-M/L	M/L
0130A-XL	XL



AIRCAST®

Weather Cover

The Aircast Walking Brace Weather Cover is designed to repel the elements and to keep the walking brace clean and dry. The Velcro closures on the back of the cover allow for easy application. Open bottom allows sure-footing for ambulatory patients. The Aircast Walking Brace Weather Cover comes in Small (fits: Ped. and Small), Large (fits: Medium and Large) and XLarge (fits: XLarge) and with an option for Tall or Short.

Tall	Short	Size
01WXT-S	01WXS-S	Ped./S
01WXT-L	01WXS-L	M/L
01WXT-XL	01WXS-XL	XL



PROCARE®

MiniTrax

MiniTrax™ is the only pediatric walking boot to feature kid-friendly design elements that make wearing a boot less intimidating. Beneath the fun exterior though, the MiniTrax is a serious boot with big benefits. Lightweight design includes anatomically shaped plastic uprights that help to facilitate a better fit. Easy to apply and adjust with a simple strapping and liner system. Low-profile rocker sole of the MiniTrax features a durable, non-slip tread tough enough to withstand the test of time—and of a kid.

Indication: Ideal following acute ankle sprains, metatarsal fractures stress and stable fractures of the lower leg.

MiniTrax Sizing Guide				
	Foot Length	Leg Length	Shoe Size	Size
79-95510	≤ 6"	6" - 7"	≤ 8.5	S (Ped)
79-95520	6" - 7"	≤ 11"	8.5 - 11.5	M (Ped)
79-955130	7" - 8"	≤ 14"	11.5 - 1.5	L (Youth)
Accessories				
79-95610	MiniTrax Replacement Liners			S (Ped)
79-95620	MiniTrax Replacement Liners			M (Ped)
79-95630	MiniTrax Replacement Liners			L (Youth)

PROCARE®

MaxTrax™ Air

Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures. Adjustable medial/lateral air cell system inflates both sides equally to accommodate changes in edema during rehab and maximize patient compliance. Cushioned inner/outer sole with low profile rocker bottom provides increased stability and comfort. Soft nylon/foam liner is constructed to ensure no seams are in critical post-surgical areas. Malleable uprights breakdown on standard height walking brace to accommodate shorter tibia/fibula.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size			Replacement Liner
	Male	Female	Junior	
MaxTrax Air				
79-95412		≤4	≤5	79-95362 Ped./XS
79-95413	≤5	4.5 - 6		79-95363 S
79-95415	5.5 - 10	6.5 - 11		79-95365 M
79-95417	10.5 - 13.5	11.5 - 14.5		79-95367 L
79-95418	14 - 17			79-95368 XL
MaxTrax Air Ankle				
79-95422		≤4	≤5	79-95372 Ped./XS
79-95423	≤5	4.5 - 6		79-95373 S
79-95425	5.5 - 10	6.5 - 11		79-95375 M
79-95427	10.5 - 13.5	11.5 - 14.5		79-95377 L
79-95428	14 - 17			79-95378 XL
Anterior/Posterior and Medial/Lateral Shells				
79-95402				XS/S
79-95403				M/XL



PROCARE®

XcelTrax™ Air

Designed for ideal support and comfort following trauma or post-operative procedure. The lightweight design and unique pre-shaped ergonomic frame help to modify the level of pressure on the injury at both the malleoli and lower leg. The soft nylon/foam liner has an integrated pump design making the pneumatic liner simple to inflate and adjust.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
XcelTrax Air				
79-95512	2 – 4	3.5 – 5	79-95532	XS
79-95513	4.5 – 7	6 – 8	79-95533	S
79-95515	7.5 – 10.5	8.5 – 11.5	79-95535	M
79-95517	10.5 – 12.5	11.5 – 13.5	79-95537	L
79-95518	12.5+	13.5+	79-95538	XL
XcelTrax Air Ankle				
79-95522	2 – 4	3.5 – 5	79-95342	XS
79-95523	4.5 – 7	6 – 8	79-95343	S
79-95525	7.5 – 10.5	8.5 – 11.5	79-95345	M
79-95527	10.5 – 12.5	11.5 – 13.5	79-95347	L
79-95528	12.5+	13.5+	79-95348	XL
Accessories				



PROCARE®

Nextep™ Contour Air

Provides protection and edema control following trauma or post-operative procedures. A medial/lateral dual air cell system inflates both sides of the liner equally and is easily adjusted with an integrated pump for proper compression, comfort, and custom fit. A secure, full wrap-around shell provides maximum support for ligament and tendon repair, osteotomy, pinning procedures, stress fractures or post-operative bunion surgery. Cushioned, non-slip rocker sole. Washable brush nylon foam liner with toe cover for maximum comfort.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
Nextep Contour Air				
79-95273	3.5 – 7.5	4.5 – 8.5	79-95363	S
79-95275	7 – 10	8 – 11.5	79-95365	M
79-95277	10 – 13	11 – 14	79-95367	L
Nextep Contour Air "Shortie"				
79-95443	3.5 – 7.5	4.5 – 8.5	79-95373	S
79-95445	7 – 10	8 – 11.5	79-95375	M
79-95447	10 – 13	11 – 14	79-95377	L
Accessories				
79-95440	Toe Cover 3/pk.			Univ.

PROCARE®

Nextep™ Contour 2 Air

Provides protection and edema control following trauma or post-operative procedures. A medial/lateral dual air cell system inflates both sides of the liner equally and is easily adjusted with an integrated pump for proper compression, comfort, and custom fit. Cushioned, non-slip rocker sole. Washable brush nylon foam liner for maximum comfort.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
Nextep Contour 2 Air				
79-95173	3.5 – 7.5	4.5 – 8.5	79-95363	S
79-95175	7 – 10	8 – 11.5	79-95365	M
79-95177	10 – 13	11 – 14	79-95367	L
Accessories				
79-95440	Toe Cover 3/pk.			Univ.



PROCARE®

MaxTrax™

Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures. A low profile rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Soft nylon/foam liner is constructed to ensure no seams are in critical post-surgical areas. Malleable uprights breakdown on standard height walking brace to accommodate shorter tibia/fibula.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size			Replacement Liner	Size
	Male	Female	Junior		
MaxTrax					
79-95322		≤4	≤5	79-95192	Ped./XS
79-95323	≤5	4.5 – 6		79-95193	S
79-95325	5.5 – 10	6.5 – 11		79-95195	M
79-95327	10.5 – 13.5	11.5 – 14.5		79-95197	L
79-95328	14 – 17			79-95198	XL
MaxTrax Ankle					
79-95342		≤4	≤5	79-95182	Ped./XS
79-95343	≤5	4.5 – 6		79-95183	S
79-95345	5.5 – 10	6.5 – 11		79-95185	M
79-95347	10.5 – 13.5	11.5 – 14.5		79-95187	L
79-95348	14 – 17			79-95188	XL
Anterior/Posterior and Medial/Lateral Shells					
79-95402					XS/S
79-95403					M/XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

XcelTrax™

Designed for ideal support and comfort following trauma or post-operative procedure. The lightweight design and unique pre-shaped ergonomic frame help to modify the level of pressure on the injury at both the malleoli and lower leg.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size			Replacement Liner	Size
	Male	Female	Junior		
XcelTrax					
79-95492	2 – 4	3.5 – 5		79-95552	XS
79-95493	4.5 – 7	6 – 8		79-95553	S
79-95495	7.5 – 10.5	8.5 – 11.5		79-95555	M
79-95497	10.5 – 12.5	11.5 – 13.5		79-95557	L
79-95498	12.5+	13.5+		79-95558	XL
XcelTrax Ankle					
79-95502	2 – 4	3.5 – 5		79-95562	XS
79-95503	4.5 – 7	6 – 8		79-95563	S
79-95505	7.5 – 10.5	8.5 – 11.5		79-95565	M
79-95507	10.5 – 12.5	11.5 – 13.5		79-95567	L
79-95508	12.5+	13.5+		79-95568	XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

Nextep™ Contour

Lightweight, durable, one-piece wrap-around shell, provides maximum immobilization and support following ligament and tendon repair, osteotomy, pinning procedures, stress fractures or post-operative bunion surgery. Unique design of circumferential shell also accommodates larger calf sizes. Low profile, cushioned rocker sole absorbs impact of heel strike to improve comfort during ambulation. Washable nylon foam liner with removable toe cover.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95063	3.5 – 7.5	4.5 – 8.5	79-95163	S
79-95065	7 – 10	8 – 11.5	79-95165	M
79-95067	10 – 13	11 – 14	79-95167	L
Nextep Contour Shortie				
79-95083	3.5 – 7.5	4.5 – 8.5		S
79-95085	7 – 10	8 – 11.5		M
79-95087	10 – 13	11 – 14		L
Accessories				
79-95210	Walker Extension Pad			Univ.
79-95440	Toe Cover 3/pk.			Univ.



PROCARE®

Nextep™ Contour 2

Lightweight, durable, one-piece wrap-around posterior shell, provides maximum immobilization and support. Brushed nylon/foam liner, low profile, cushioned non-skid sole. Unique design accommodates larger calf sizes. Rocker sole.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95073	3.5 – 7.5	4.5 – 8.5	79-95163	S
79-95075	7 – 10	8 – 11.5	79-95165	M
79-95077	10 – 13	11 – 14	79-95167	L
Accessories				
79-95210	Walker Extension Pad			Univ.
79-95440	Toe Cover 3/pk.			Univ.



PROCARE®

SideKICK™

The low profile design provides comfortable stability for acute ankle sprains, soft tissue and stress fractures to the lower leg and ankle, and stable fractures of the ankle. Other features include fixed, malleable uprights, a breathable fiber leg wrap, removable heel pads, and non-skid rocker sole.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95033	≤6	≤7	79-95123	S
79-95035	6 – 10	7 – 11	79-95125	M
79-95037	10+	11+	79-95127	L



PROCARE®

ProSTEP™

Lightweight, sturdy one-piece design with open heel. Non-skid rocker sole and washable fiber leg wrap with removable heel pads.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-98793	≤6	≤7	79-95123	S
79-98795	6 – 10	7 – 11	79-95125	M
79-98797	10+	11+	79-95127	L



PROCARE®

MaxTrax™ Air ROM

Protected range-of-motion and adjustable air cell system provides stability and helps decrease pain and edema following trauma or post-operative procedures. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion. Cushioned inner/outer sole with low profile rocker bottom provides increased stability and comfort. Soft nylon/foam liner.

Indication: Ideal following Achilles tendon repair for controlled rehab, and stable fractures of the lower leg.

	Shoe Size			Replacement Liner	Size
	Male	Female	Junior		
79-95432		≤4	≤5	79-95362	Ped./XS
79-95433	≤5	4.5 – 6		79-95363	S
79-95435	5.5 – 10	6.5 – 11		79-95365	M
79-95437	10.5 – 13.5	11.5 – 14.5		79-95367	L
79-95438	14 – 17			79-95368	XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

MaxTrax™ ROM

Provides protected range-of-motion following trauma or post-operative procedure. Ideal following Achilles tendon repair for controlled rehab. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion. Malleable uprights breakdown on standard height to accommodate shorter tibia/fibula. Cushioned inner/outer sole. Low profile rocker bottom.

Indication: Ideal following Achilles tendon repair for controlled rehab, and stable fractures of the lower leg.

	Shoe Size				
	Male	Female	Junior	Replacement Liner	Size
MaxTrax ROM					
79-95332		≤4	≤5	79-95192	Ped./XS
79-95333	≤5	4.5 – 6		79-95193	S
79-95335	5.5 – 10	6.5 – 11		79-95195	M
79-95337	10.5 – 13.5	11.5 – 14.5		79-95197	L
79-95338	14 – 17			79-95198	XL
MaxTrax ROM Ankle					
79-95352		≤4	≤5	79-95192	Ped./XS
79-95353	≤5	4.5 – 6		79-95193	S
79-95355	5.5 – 10	6.5 – 11		79-95195	M
79-95357	10.5 – 13.5	11.5 – 14.5		79-95197	L
79-95358	14 – 17			79-95198	XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

MaxTrax™ Diabetic Walker

Designed with a wider foot bed for increased stability and a soft nylon/foam liner with no seams in critical post-surgical areas and/or areas where neuropathy may occur. Tri-laminate insole is designed with Impax™ Grid technology to redistribute load away from ulcers. Low profile rocker bottom cushioned outer sole for shock absorption. Optional Toe Cover and a set of Upright Pads for additional patient comfort.

Indication: Ideal for replacement of total contact casting; treatment of pre-ulcerative or ulcerative conditions of the plantar surface of the foot; charcot foot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95453	≤5	4.5 – 6.5	79-95473	S
79-95455	5.5 – 10	6.5 – 11	79-95475	M
79-95457	10.5 – 13.5	11.5 – 14.5	79-95477	L
79-95458	14 – 17		79-95478	XL
79-95463	Impax Foam Insole			S
79-95465	Impax Foam Insole			M
79-95467	Impax Foam Insole			L
79-95468	Impax Foam Insole			XL
Accessories				
79-95450	MaxTrax Diabetic Toe Cover			Univ.
79-95460	MaxTrax Upright Pads (1 pair)			Univ.
79-95210	Walker Extension Pad			Univ.



PROCARE®

Achilles Wedge

4" neoprene padded foam wedge fits into most closed heel walkers. Each section is 3/4" in height and may be "peeled-away" during Achilles Tendon rehab.

	Size
79-95130	Univ.



PROCARE®

Evenup™

Simple, fast, easy to apply, the Evenup uses an innovative design that helps to equalize a patient's limb length and reduce body strain while walking in a cast or walker. Easily removable, the Evenup comes in three different sizes and two different thicknesses.

	Shoe Size		Size
	Male	Female	
79-95483	7 – 8.5	6.5 – 9.5	S
79-95485	9 – 11.5	9.5 – 12.5	M
79-95487	12+	12.5+	L



PROCARE®

Weather Cover

The ProCare Walking Brace Weather Cover is designed to repel the elements and to keep the walking brace clean and dry. The Velcro closures on the back of the cover allow for easy application. Open bottom allows sure-footing for ambulatory patients.

	Size
79-95093	S
79-95095	M
79-95097	L

PATIENT & PRESSURE CARE

PROCARE[®]**Body Wedge**

Stabilizes the patient securely in several positions without slipping or sliding. Supports arm or leg when lying on side or at foot of bed to prevent sliding in supine positions. Can also be used as a back support, footrest or blanket support. Comes compressed to reduce package sizing.

Indication: Pressure ulcer prevention, post-operative stabilization and positioning.

	Qty.	Length	Height	Width
79-90860	6/bx.	11.69"	8.375"	20"





PROCARE®

Hip Abduction Pillows

Designed for use following hip surgery when immobilization or post-operative positioning is required. All-foam construction for maximum support. Soft contact closure foam laminate straps help keep positioner in place and provide easy adjustment. Universal designs accommodate wide variety of patient sizes. 79-90180 has a 4" tear-off block on both ends of pillow for adjustment. Comes compressed to reduce package sizing.

Indication: Help keep hips in abduction (away from center of the body) following injury or surgery.

	Length	Height	Width Bottom	Size
79-90173	18"	6"	12"	S
79-90175	23"	6"	15"	M
79-90177	26"	6"	18"	L
79-90180	Max. 27" / Min. 19"	3"	6.5"	Univ.
79-90190*	18"	6"	18"	Univ.

*Foam wedge with removable, flannel lined washable cover and pile foam straps with contact closure.



PROCARE®

Leg Elevator

Full foam construction to elevate and maintain proper knee and lower leg position at 45° angle. Comes compressed to reduce package sizing.

Indication: Ideal for patients with leg casts, or those requiring stable elevation and internal/external rotation control, post-trauma or surgical procedures.

	Length	Height	Width
79-90191	31.5"	10"	10.25"




PROCARE®

Arm Elevation Pillow

Sculptured foam support aids in the reduction of pre/post-cast or post-surgical swelling. Provides comfortable, effective support with minimal restraint. Comes compressed to reduce package sizing. Excellent for use on wrist or shoulder surgery and mastectomy patients.

Indication: Intended for post-surgical or post-injury elevation of the hand and arm to promote venous return to minimize swelling.

	Height	Width	Length
79-90500	16.75"	11.5"	12.75"



PROCARE®


Donuts

Versatile positioners that provide anatomical support and comfort. May be used as a head rest; for protection of the ears and nerves of the head and face. Additional indications include protection of the elbows, knees and heels.

Indication: Cushions the head, ear or iliac crest. May be used to support elbows, heels and knees.

	Qty.	Height	Diameter
79-90880* Donuts w/Center (6", 4.5", 3" center)	36/bx.	2"	9"
79-90890 Donuts w/o (3" Center)	36/bx.	2"	7"
79-90891 Donuts w/o (4.5" Center)	36/bx.	2"	9"

*Donuts w/center (photo shown above).




PROCARE®

Foam Arm Boards

Helps maintain hemodynamic during surgery and enhances pressure redistribution. The Foam Arm Boards protect bony prominences and pressure areas such as elbows, knees, and heels. Comes compressed to reduce package sizing.

Indication: Supports entire arm during surgery. Can be used on an operating room table to enhance pressure redistribution.

	Qty.	Length	Height	Width
79-90820	12 pair/bx.	16"	1.75"	7.5"




PROCARE®

Head Cradles

Provide comfortable immobilization of head and neck. May be used for lateral, supine or prone positions. Designed to conform to sensitive facial features while maintaining head position. Comes compressed to reduce package sizing.

Indication: Provides comfortable immobilization of head and neck. May be used for lateral, supine or prone positions.

	Qty.	Length	Height	Width
79-90872 Head Cradle w/o Trach Slit	12/bx.	8"	4.5"	9"
79-90873 Head Cradle w/Trach Slit	36/bx.	8"	4"	9"



PROCARE®

Laminectomy Arm Cradles with Trough

Used to support the body off the chest while in a prone position. Offers cushioning for arms, shoulders, elbows and ulnar nerves. Design allows for varying arm positions. Comes compressed to reduce package sizing.

Indication: Offers cushioning for arms, shoulders, elbows, and ulnar nerves during laminectomy procedures or other prone positioned surgeries.

	Qty.	Length	Height	Width
79-90810	5 pair/bx.	23.5"	3"	5"



PROCARE®

Multi-Purpose Belt

Prevents patient movement on O.R. table or gurney.

Indication: Ideal for stretchers, treatment and procedure tables, this one piece design is used for patient positioning. The non-metal construction is safe for use in the MRI units.

	Qty.	Length	Width
79-90830	12/bx.	60"	2"




PROCARE®

O.R. Table Pad

A full-length convoluted foam table pad offers maximum protection and comfort for the patient. Used on top of an operating room table to enhance pressure redistribution and maintain hemodynamic during surgery. Comes compressed to reduce package sizing.

Indication: Comfortably supports patient during surgery.

	Qty.	Length	Height	Width
79-90850	6/bx.	72"	2"	20"



PROCARE®

Prone Head Positioners

Provides comfortable immobilization of head and neck. May be used for prone, lateral or supine positions. Designed to conform to sensitive facial features while maintaining head position. The Premium Prone Head Positioner maximizes comfort by using an extra soft comfortable support, reducing the pressure on the patients eyes and face. Comes compressed to reduce package sizing.

Indication: Provides comfortable immobilization of head and neck. May be used for lateral, supine or prone positions.

	Qty.	Length	Height	Width
79-90870	12/bx.	8.25"	4.125"	8.25"
79-90871*	12/bx.	10.5"	7"	9.5"

*Premium Prone Head Positioner (photo shown above).



PROCARE®

Ulnar Nerve Protectors

Protects against pinching and entrapment of the ulnar nerve. Provides cushioning on bony prominences of the arm. Available in standard and premium convoluted foam density.

Indication: Ideal for protecting against pinching and entrapment of the ulnar nerve. Provides cushioning on bony prominences of the arm.

	Qty.	Length	Height	Width
Standard 1.0 Density				
79-90800	36 pair/bx.	14.75"	2"	6.4"
Premium 1.5 Density				
79-90801*	36 pair/bx.	14.75"	1.75"	6.4"

*Premium 1.5 Density (photo shown above).



PROCARE®

Flex-Master® Non Sterile Clip Closure Bandage

Latex free, 100% woven cotton/elastic minimizes moisture retention and allergic reactions. Available in standard and double lengths.

Indication: Ideal for groin wraps, ice packs, knees, and ankles.

	Per Pk.	Length	Width
Standard Length			
79-98843	10	5.5 yds.	2"
79-98845	10	5.5 yds.	3"
79-98847	10	5.5 yds.	4"
79-98848	10	5.5 yds.	6"
Double Length			
79-98857	12	11 yds.	4"
79-98859	12	11 yds.	6"



PROCARE®

Flex-Master® Sterile Clip Closure Bandage

Standard and double lengths. Latex free, woven cotton/elastic construction with clip closure.

Indication: Ideal for post-surgical applications such as thoracic, knee, amputations, mastectomy, vascular, shoulder and burns.

	Per Pk.	Length	Width
Standard Length			
79-98723	20	5.5 yds.	2"
79-98725	20	5.5 yds.	3"
79-98727	20	5.5 yds.	4"
79-98728	20	5.5 yds.	6"
Double Length			
79-98717	20	11 yds.	4"
79-98718	20	11 yds.	6"



PROCARE®

Elastic Bandages - Self Closure

Constructed of breathable elastic. Each bandage is individually wrapped and may be secured by hook and loop closure or with two clip fasteners. Latex free.

Indication: Ideal for groin wraps, ice packs, knees, and ankles.

	Per Pk.	Length	Width
79-98233	10	5 yds.	2"
79-98235	10	5 yds.	3"
79-98237	10	5 yds.	4"
79-98238	10	5 yds.	6"



PROCARE®

Elastic Bandages w/Clip Closure

Cotton/Elastic wrap stretches to a full 4.5 yards and secures with two clips. Washable and auto clavable. Latex free.

Indication: Ideal for groin wraps, ice packs, knees, and ankles.

	Per Pk.	Length	Width
79-98583	10	4.5 yds.	2"
79-98585	10	4.5 yds.	3"
79-98587	10	4.5 yds.	4"
79-98588	10	4.5 yds.	6"



PROCARE®

Latex Rubber Bandage

Constructed of natural latex and designed to force blood out of an extremity by compressing the superficial vessels. Available in 3", 4", and 6" widths.

Indication: Used in conjunction with a tourniquet cuff to help accomplish exsanguinations (removal of blood) from a limb.

	Length	Width
79-98900	144"	3"
79-98910	144"	4"
79-98920	144"	6"



PROCARE®

Padded Hip Spica

Easy application facilitates easy dressing changes. Padding protects dressings and hip fracture. Compression helps prevent swelling and pain. Washable.

Indication: To provide protection to hip area and help prevent discomfort.

	Hip Measurement	Size
79-90003	24" – 32"	S
79-90005	30" – 38"	M
79-90007	36" – 42"	L
79-90008	40" – 50"	XL
79-90009	48" – 56"	2XL



PROCARE®

Catheter Strap

2" wide flannel lined elastic strap with contact closure. Special tab closure maintains tube position without occlusion.

Indication: Commonly used for catheter stabilization.

	Qty.	Length	Size
79-89200	12/pk.	24"	Univ.



PROCARE®

Heel and Elbow Protectors

Constructed of a breathable interfacing with fiberfill lining. Three-layer comfort cushion helps to eliminate pressure points. Anatomical design provides customized patient comfort. Hook and loop closure for ease of application and removal. For protection of the elbow and heel for injured and/or bedridden patients.

Indication: Ideal for pressure reduction of heel and/or elbow and pressure protection.

	Style	Qty.	Length
79-81001	Personal	Bulk Pairs, 18 Pair	Univ.
79-81002	Prevent	Indiv. Pairs, 18 Pair	Univ.



PROCARE®

Heel/Elbow Protector - Poly/Pile

Choice of polyester pile or vented foam with velfoam contact closure straps. Sold in pairs.

Indication: Ideal for pressure reduction of heel and/or elbow and pressure protection.

	Qty.	Size
79-81003	1 pair	S
79-81005	1 pair	M
79-81007	1 pair	L



PROCARE®

Mesh Heel/Elbow Protector

Knitted stretch stockinette covered comfort pad for use on heel, elbow or knee. Sold in pairs. One size fits all.

Indication: For protection of the elbow and heel for injured and/or bedridden patients.

	Description	Qty.	Size
79-81600	Foam Padded	1 Pair	Univ.



PROCARE®

Heel/Elbow Protectors - Universal

Choice of polyester pile or vented foam with velfoam contact closure straps. One size fits all. Sold in pairs.

Indication: For protection of the elbow and heel for injured and/or bedridden patients.

	Description	Qty.	Size
79-81040	Poly/Pile	1 Pair	Univ.
79-81050	Foam	1 Pair	Univ.
79-99740	Foam	3 Pair	Univ.



PROCARE®

Podous Boot

Lightweight, quality plastic shell with poly/pile liner help treat and prevent lower extremity disorders associated with trauma or immobility including pressure necrosis, ankle contractures, and post-operative hip abduction or rotation. Adjustable toe extension helps protect toes from bedding. Anti-rotation bar may be repositioned to assist with internal or external rotation. Removable transfer pad on dorsal side for limited patient ambulating. Fits right or left foot. Latex free.

Indication: Helps treat and prevent lower extremity disorders including pressure necrosis, ankle contractures, and post-operative hip abduction or rotation.

	Shoe Size			Length
	Male	Female	Calf Circ.	
79-90550	≤ 10	≤ 11	≤ 22"	Regular
79-90551	10.5 +	11.5 +	≤ 25"	XL
Replacement Liners				
79-90560	Liner Only	3/pk.		Regular
79-90561	Liner Only	3/pk.		XL



PROCARE®

Drop Foot Boot

90° upright frame with internal/external rotation in 10° increments helps control foot drop. Disposable poly/pile liner provides toe and heel protection.

Indication: Prevents Achilles contracture. Ideal for ankle contractures, muscle tightening, correct plantar flexion contractures.

	Description	Size
79-90160	Drop Foot Boot, Complete	Univ.
79-90170	Disposable Liner	Univ.
79-90161	Frame Only	Univ.



PROCARE®

Pre-Vent™ Ulnar Nerve Protector

Constructed of a breathable interfacing with fiberfill lining. Provides cushioned support for patient limb. Fluid repellent outer liner. Hook and loop closure for ease of application and removal. Can be used as post-operative support, insulated positioner for cryotherapy.

Indication: Traps fluids from patient, prevents pooling on O.R. surface area and expedites post-operative clean-up.

	Qty.	Size
79-91100	18 Pair	Univ.



PROCARE®

Pre-Vent™ O.R. Table Pad

Constructed of fluid-attracting face fabric with fiberfill lining. Impervious backing traps fluid, preventing pooling and expediting post-operative clean-up. Three section design can be customized to accommodate any surgical position. Does not interrupt heat transfer from warming blankets.

Indication: Traps fluids from patient, prevents pooling on O.R. surface area and expedites post-operative clean-up.

	Qty.	Size
79-91110	4/pk.	Univ.



PROCARE®

Pre-Vent™ Knee Crutch Pad

Constructed of a breathable interfacing with fiberfill lining. Provides cushioned support for patient limb. Fluid repellent outer liner. Hook and loop closure for ease of application and removal. Can be used as post-operative support, insulated positioner for cryotherapy applications.

Indication: Traps fluids from patient, prevents pooling on O.R. surface area and expedites post-operative clean-up.

	Qty.	Size
79-91120	12 Pair	Univ.



PROCARE®

Adjustable Anodized Aluminum Crutches

Packaged in pairs complete with tips, grips, and pads. Adjustable by using wing nut. 300 pound maximum weight. Latex free.

Indication: Can be used for many conditions including injuries, arthritis, or problems with balance.

	Approximate User Height	Size
79-91333	4'6" – 5'2"	Youth/S
79-91335	5'1" – 5'9"	Adult/M
79-91337	5'10" – 6'6"	Tall/L



PROCARE®

Adjustable Aluminum Cane

Anodized aluminum construction with adjustable range from 28" to 40" in height. 250 pound maximum weight.

Indication: Can be used for many conditions including injuries, arthritis, or problems with balance.

Adjustable Range	
79-91091	28" - 40"



PROCARE®

Anodized Aluminum Adjustable Fold-Up Walker

Includes hand grips and crutch tips. Folds flat to four-inch thickness for easy storage and car use, snap locking with a plated steel crossbar. Adjustable for proper height. 300 pound maximum weight.

Indication: Can be used for many conditions including injuries, arthritis, or problems with balance.

Adjustable Range	
79-91090	32" - 40"

TRACTION & PATIENT SAFETY

PROCARE[®]

Quick-Release Zipper Body Holder

Support in chairs, beds and wheelchairs for patients at risk for falls or those requiring a restraint to assist with medical treatment. Specifically designed with sleeves, an off-centered zipper and quick-release buckled straps for increased patient restraint

Indication: Provide comfort and security for wheelchair or in-bed patients.

	Color	Weight	Qty.	Size
79-91173	Red	80 – 120 lbs.	10/pk.	S
79-91175	Green	110 – 155 lbs.	10/pk.	M
79-91177	Yellow	145 – 190 lbs.	10/pk.	L
79-91178	Blue	180 – 230 lbs.	8/pk.	XL
79-91179	Purple	220 – 275 lbs.	8/pk.	2XL





PROCARE®

Vest Style Body Holders

Silklin body fabric with 65” straps. Cross points help to prevent excessive tightening. Arm holes help deter patient from working vest loose. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger him/herself or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

Regular	X-Strong Straps	Color	Weight	Qty.	Size
79-91643	79-91633	Pink	90 – 130 lbs.	10/pk.	S
79-91645	79-91635	Blue	120 – 180 lbs.	10/pk.	M
79-91647	79-91637	Yellow	170 – 220 lbs.	10/pk.	L
79-91648	79-91638	Green	210 – 300 lbs.	8/pk.	XL



PROCARE®

Secure-All™ Body Holder

Silklin body fabric with 65” straps. Can be used in hospital beds, wheelchair or geriatric chair. Cross points help prevent excessive tightening. Arm holes help deter patient from working vest loose. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger him/herself or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

	Color	Weight	Qty.	Size
79-91604	Lt. Blue	100 – 200 lbs.	10/pk.	S/M
79-91606	Navy	120 – 220 lbs.	20/pk.	M/L
79-91608	Green	210 – 300 lbs.	8/pk.	XL



PROCARE®

Vest Body Holder

Provides additional support and helps prevent unassisted bed or chair exits. Washable, light canvas lined with soft flannel. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger themselves or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

	Color	Weight	Size
79-91053	Red	80 lbs. – 120 lbs.	S
79-91055	Green	110 lbs. – 155 lbs.	M
79-91057	Yellow	145 lbs. – 190 lbs.	L
79-91058	Blue	180 lbs. – 230 lbs.	XL



PROCARE®

Mesh Vest Body Holder

Effective for patients in beds and wheelchairs. Lightweight, cool mesh construction. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger themselves or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

	Color	Weight	Size
79-91072	White	60 lbs. – 90 lbs.	XS
79-91073	Red	80 lbs. – 120 lbs.	S
79-91075	Green	110 lbs. – 155 lbs.	M
79-91077	Yellow	145 lbs. – 190 lbs.	L
79-91078	Blue	180 lbs. – 230 lbs.	XL



PROCARE®

Mesh Body Holder w/Sleeves

Support for patients in chairs, beds and wheelchairs. Sleeves provide increased patient restraint, cool mesh material. For patient comfort, the back has an off-center contact closure. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger themselves or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

	Color	Weight	Size
79-91063	Red	80 lbs. – 120 lbs.	S
79-91065	Green	110 lbs. – 155 lbs.	M
79-91067	Yellow	145 lbs. – 190 lbs.	L



PROCARE®

Vest Style Body Holder with Safety Sleeves

Silkin body fabric with 80” straps. Sleeve style helps to deter patient from working vest loose. Can be used in hospital bed, wheelchair or geriatric chair. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger him/herself or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

Regular	X-Strong Straps	Color	Weight	Qty.	Size
79-91623	79-91613	Pink	90 – 130 lbs.	10/pk.	S
79-91625	79-91615	Blue	120 – 180 lbs.	10/pk.	M
79-91627	79-91617	Yellow	170 – 220 lbs.	10/pk.	L
79-91628	79-91618	Green	210 – 300 lbs.	8/pk.	XL



PROCARE®

Safety Vest Body Holder

Versatile design for bed or chair use. Straps can secure outside of patient’s reach. Available in Cool-Mesh for maximum comfort with standard tie straps. Contraindicated for combative patients which exhibit high aggressiveness, who strike or endanger him/herself or others.

Indication: Provide comfort and security for wheelchair or in-bed patients.

	Color	Weight	Size
79-91033	Red	80 lbs. – 120 lbs.	S
79-91035	Green	110 lbs. – 155 lbs.	M
79-91037	Yellow	145 lbs. – 190 lbs.	L
79-91038	Blue	180 lbs. – 230 lbs.	XL



PROCARE®

Deluxe Quick-Release Limb Holder

Comfortable quilted cuffs adjust to most wrist or ankle sizes. Secures with hook and loop closure and quick-release buckles on both wrist and strap-end. 50” strap secures to movable position of bed.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91500	14 pair/pk.	Univ.



PROCARE®

Quick-Release Foam Limb Holder

Comfortable foam cuff design with single-strap attachment has convenient hook and loop closures and quick-release buckles. 58" strap secures to movable position of bed.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91510	24 pair/pk.	Univ.



PROCARE®

Quick Release Limb Holder

Constructed of soft quilted cuff with hook/loop closure for quick release as needed. 36" strap length.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91400	14 pair/pk.	Univ.
79-91420	14 pair/pk.	Univ./Extra Strong Straps
79-91490	50 pair/pk.	Ped.



PROCARE®

Personal Limb Holder

Soft quilted cuff design threads through tear drop opening and is secured with square knot for snug fit. 39" Sontara™ strap length.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Packaging
79-91460	40 pair/pk.	8 pair/box
79-91480	40 pair/pk.	Bulk Pairs



PROCARE®

Wraparound Limb Holder

Soft, quilted cuff, D-ring and hook/loop closures with 50" strap length.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91430	40 pair/pk.	Univ.
79-91440	40 pair/pk.	Univ./Extra Strong Straps
79-91450	40 pair/pk.	Univ./Double Extra Strong Straps



PROCARE®

Foam Limb Holder

Constructed of foam cuff and metal D-rings. 48" Sontara™ straps.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91410	24 pair/pk.	Univ.



PROCARE®

Secure-All™ Limb Holder

Soft quilted cuff design threads through tear drop opening and is secured with square knot for snug fit. 46" Sontara™ strap length.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91470	30 pair/pk.	Univ.



PROCARE®

Foam Limb Holders

Contact closure with double "D" ring cinch lock. Available in pairs. Sold in pediatric and universal sizes.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Color	Strap Length	Size
79-90100	10 pair	60"	Univ.
79-90110	1 pair	60"	Univ.
79-90113	1 pair	41"	Ped./XS



PROCARE®

Double Buckle Limb Holder

Foam restraint for ankle or elbow with double-buckle closure. Extra long tie. Sold 10 pairs per pack.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Color	Strap Length	Size
79-90120	10 pair	60"	Univ.
79-90122	10 pair	36"	Ped./XS



PROCARE®

Poly/Pile Limb Holder

Polyester pile wrist/ankle restraint with slide buckle closure. Available with standard buckle or quick release on strap. 46.5" strap length.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-90080	6 pair	Univ.
79-90410	1 pair	Univ.
Quick-Release Bundle		
79-90420	1 pair	Univ.



PROCARE®

Canvas Limb Holder

Sturdy canvas construction with web straps.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-90400	1 pair	Univ.



PROCARE®

Secure-All™ Finger Control Mitt

Constructed of silklin and fiberfill with hook/loop closures around wrist for secure fit. Option available with quick-release wrist strap.

Indication: Intended for patients who are at risk for disturbing medical treatment by pulling their IV line and/or catheter, or are prone to self-injury. This device is designed for use in a bed, on a stretcher or wheelchair.

	Qty.	Size
79-91220	6 Pair	Univ.
Quick-Release		
79-91221	6 pair	Univ.



PROCARE®

Elbow Immobilizer

Soft foam/nylon laminate construction for patient comfort. Aluminum stays help prevent patients from bending the elbow joint or accessing tubes. Contact closure straps for ease of application. Available in bulk pack of 50.

Indication: Limb positioning; Post-operative; IV therapy.

Each	50-Pack	Length	Width	Size
79-91512	79-91522	4"	7"	Infant/XS
79-91513	79-91523	5.75"	8.75"	Child/S
79-91515	79-91525	6.75"	11.25"	Youth/Adult/M
79-91517	79-91527	9.75"	15.5"	Adult/L



PROCARE®

Secure-All™ Body Belt

Silklin fabric covered foam belt with metal O-rings. 70" canvas straps.

Indication: Ideal for patient positioning and wheelchair support.

	Qty.	Size
79-91300	12/pk.	Univ.



PROCARE®

Roll Belt

Allows patient to sit and move from side to side within bed. Canvas fabric with web straps. One size fits most adults. Straps tie securely beneath bed, out of patient's reach.

Indication: Ideal for patients at risk for falls or those requiring a restraint for medical treatment.

	Qty.	Size
79-91280	1 Belt	Univ.



PROCARE®

Heavy Duty Gait Belt

Assists in safe patient transfer and ambulation, and helps prevent caregiver back injuries. Specially designed metal buckle offers maximum hold and constructed of poly-blend material for increased durability.

Indication: To assist cooperative patients that require ambulation and/or transfer assistance.

	Length	Size
79-91010	61"	Univ.



PROCARE®

Gait Belt

Used for supporting patient while walking or for securing into wheelchair. Sturdy, double back pile strap with quick-release buckle.

Indication: To assist cooperative patients that require ambulation and/or transfer assistance.

	Length	Size
79-91000	60"	Univ.



PROCARE®

Full Foam Bucks Traction

Convuluted foam helps minimize migration. Sewn-in medial and lateral stay, 6" spreader bar and loop/lock closure.

Indication: Ideal for hip or femur fractures.

	Length	Max. Calf Circum.	Size
79-83062	12"	11"	XS
79-83063	16"	13"	S
79-83065	20"	15.5"	M
79-83067	20"	18"	L
79-83068	20"	21"	XL



PROCARE®

Universal Bucks Traction with Foot Pad

Convuluted foam minimizes splint migration, increases air flow and reduces pressure on the limb. Metal medial/lateral stays provide stability. 6" spreader bar is incorporated into design. Foot pad with adjustable straps helps impede foot drop.

Indication: Ideal for hip or femur fractures.

	Length	Max Calf Circum.	Size
79-99890	20"	18.5"	Univ.



PROCARE®

Universal Bucks Traction

Two piece body is constructed of an inner liner with open-cell convuluted foam and foam laminate outer casing. Inner foam liner may be replaced for patient hygiene as needed. Universal sizing reduces stocking inventory. Three elastic lower leg straps secure with Velcro to either side of outer casing and can be adjusted for patient comfort. Adjustable elastic foot strap helps to impede foot drop. Removable aluminum medial/lateral stays. 6" spreader bar is incorporated into design.

Indication: Ideal for hip or femur fractures.

	Length	Max Calf Circum.	Size
79-99850	21"	21"	Univ.
79-99851	Replacement Liner		



PROCARE®

Padded Head Halter

Adjustable loop holder and flannel lined for patient comfort. Can be applied for hyperextension or flexion.

Indication: To be used in conjunction with traction treatment for cervical spine disorders including sprains and strains.

	Qty.	Size
79-83040	12/pk.	Univ.



PROCARE®

Deluxe Head Halter

Durable non-stretch design with adjustable straps. Fully padded. Halter easily adjusts with chin strap. Designed for overhead traction while providing occipital pressure. One size fits all.

Indication: Ideal for cervical strains and degenerative disc disease.

	Description	Size
79-83050	Deluxe Head Halter	Univ.



PROCARE®

Overdoor Cervical Traction Kit

Kit includes overdoor bracket, padded head halter, water weight bag, spreader bar and 86" traction cord to administer cervical traction at home or in physician's office.

Indication: Ideal for cervical strains and degenerative disc disease.

	Description	Size
79-90130	Cervical Traction Kit	Univ.



PROCARE®

Traction Cord

Flexible nylon cord, tensile-strength rated to 600 pounds. Available in spools of 100 and 33 yards. 5/32" diameter.

	Length	Size
79-98600	100 yard spool	Univ.
79-98601	33 yard spool	Univ.

FRACTURE BRACING & SPLINTS

PROCARE®

Colles Splint

Pre-molded, padded aluminum splint is designed for wrist support and immobilization prior to surgery or casting. Ventilation holes provide air circulation for patient comfort.

Indication: Immobilization of the wrist in the initial short and or most severe stage of injury.

Padded Left	Padded Right	Length x Width	Size
79-72112	79-72122	5.75" x 2.5"	Ped/XS
79-72113	79-72123	7" x 3.5"	S
79-72115	79-72125	8.5" x 3.5"	M
79-72117	79-72127	9" x 4.5"	L





PROCARE®

Forearm Splint

Padded aluminum splint is ideal for support of colles injuries prior to casting or surgery. Ventilation holes facilitate air circulation and ensure patient comfort.

Indication: Immobilization of the wrist and forearm in the initial short and/or most severe stage of injury.

Padded Left	Padded Right	Length x Width	Size
79-71972	79-71982	7" x 2.25"	XS
79-71973	79-71983	7.5" x 2.75"	S
79-71975	79-71985	9.5" x 3.75"	M
79-71977	79-71987	10.5" x 4.25"	L



PROCARE®

Mason Allen Splint

Constructed of aluminum and pre-molded for proper anatomical fit. Splint is designed to immobilize hand and wrist.

Indication: Immobilization of hand and wrist in the initial short and/or most severe stage of injury.

Unpadded	Padded	Length x Width	Size
79-72803	79-71043	9.5" x 5"	S
79-72807	79-71047	12.5" x 6"	L



PROCARE®

Cock-Up Splint

Padded aluminum splint designed for maximum use of fingers while maintaining wrist in cock-up position. Ideal for minor wrist injuries. Ventilation holes ensure proper circulation. Fits right or left hand.

Indication: Immobilization of the wrist in the initial short and/or most severe stage of injury.

Padded	Length x Width	Size
79-71943	6.5" x 2.5"	S
79-71945	8" x 3"	M
79-71947	9" x 3"	L



PROCARE®

Metacarpal Splints

Pre-formed, padded aluminum splint. Provides proper support for phalangeal, ulna, radius, and metacarpal fractures.

Indication: Immobilization of the wrist in the initial short and/or most severe stage of injury.

Padded Left	Padded Right	Length x Width	Size
79-72203	79-72213	8" x 3"	S
79-72205	79-72215	10.5" x 3.75"	M
79-72207	79-72217	11.5" x 4.5"	L



PROCARE®

Posterior Tibia/Fibula Splint

Aluminum pre-formed splint is designed to aid in treatment of tibia/fibula or ankle fractures. Ventilation holes facilitate air circulation to help ensure patient comfort.

Indication: For immobilization of the foot, ankle, and/or the distal 1/3 of the tibia and fibula in the acute stage.

Padded	Unpadded	Length x Width	Size
79-72283	79-72313	11" x 7"	S
79-72285	79-72315	13.5" x 8"	M
79-72287	79-72317	16" x 9.5"	L



PROCARE®

Femoral Leg Splint

Malleable coated aluminum is fully ventilated for patient comfort. Fits right or left leg.

Indication: For immobilization of the knee, tibia, fibula and/or the femur in the acute stage.

Padded	Length x Width	Size
79-72293	20" x 8"	S
79-72295	25" x 8"	M
79-72297	30" x 9.5"	L



PROCARE®

Baseball Splints

Padded aluminum splint is easily contoured for patient comfort.

Indication: Helps maintain stability by immobilizing phalangeal joints.

Padded	Length	Qty.	Size
79-72093	4"	12 pk.	S
79-72095	4.5"	12 pk.	M
79-72097	5"	12 pk.	L



PROCARE®

Finger Protectors

Lightweight aluminum, foam padded four prong splint protects injuries to fingertips. Easily molded to fit patient.

Indication: Ideal for use when total finger immobilization is required.

Padded	Length	Qty.	Size
79-71883	1.5"	12 pk.	S
79-71884	2.5"	12 pk.	S/M
79-71885	3"	12 pk.	M
79-71887	3.5"	12 pk.	L



PROCARE®

Finger Cots

Lightweight, padded aluminum protects injuries to fingertips and designed for ease of application and removal.

Indication: Ideal for use when total finger immobilization is required.

Padded	Length	Qty.	Size
79-71903	1.5"	12 pk.	S
79-71905	2.25"	12 pk.	M
79-71907	3.25"	12 pk.	L
79-71908	3.75"	12 pk.	XL



PROCARE®

Curved Fingers

Indicated for transverse fractures of the phalanges. Lightweight aluminum is curved to conform to finger in semi-flexed position.

Indication: Ideal for dislocations; phalangeal and metacarpal fractures.

Padded	Length x Width	Qty.	Size
79-71923	1.5" x .875"	12 pk.	S
79-71925	3" x .875"	12 pk.	M
79-71924	4" x .875"	12 pk.	M/Long
79-71926	5.5" x .875"	12 pk.	M/X-Long
79-71927	6" x .875"	12 pk.	L
79-71928	9" x .875"	12 pk.	XL



PROCARE®

Frog Splints

Easily contours for proper fit without bulky dressing. Padded aluminum splint helps maintain phalangeal joints in position.

Indication: Ideal for use when total finger immobilization is required.

Padded	Length x Width	Qty.	Size
79-71963	2.75" x 2.25"	12 pk.	S
79-71965	3.15" x 2.5"	12 pk.	M
79-71967	3.75" x 3"	12 pk.	L



PROCARE®

Plastic Finger Splints

Designed to aid in the repair of ruptured tendons. Lightweight clear or skin tone.

Indication: Fingertip and nail bed injuries; mallet finger.

Skin Tone	Clear	Length	Diameter	Qty.	Size
79-72241	79-72251	2.08"	.6"	12 pk.	#1
79-72242	79-72252	2.04"	.66"	12 pk.	#2
79-72243	79-72253	2.3"	.74"	12 pk.	#3
79-72244	79-72254	2.475"	.75"	12 pk.	#4
79-72245	79-72255	2.31"	.75"	12 pk.	#5
79-72246	79-72256	2.4"	.9"	12 pk.	#5.5
79-72247	79-72257	2.53"	.9"	12 pk.	#6
79-72248	79-72258	2.47"	1"	12 pk.	#7
79-72240	Finger Straps			25 pk.	



PROCARE®

Plastic Finger Splint Starter Set

30-piece assorted Plastic Finger Splints in choice of clear or skin tone colors.

Indication: Fingertip and nail bed injuries; mallet finger.

	Description	Size
79-72250	30 piece set, Skin tone	Univ.
79-72260	30 piece set, Clear	Univ.
79-72240	25 pk., Finger Straps	Univ.



PROCARE®

Aluminum Strips

Multi-purpose splints may be easily cut and contoured to desired shape. Closed cell hydroloc foam pads one side of each strip for patient comfort.

Indication: Phalangeal fractures.

	Length x Width	Qty.	Size
79-72163	18" x .5"	12 pk.	Univ.
79-72165	18" x .75"	12 pk.	Univ.
79-72167	18" x 1"	6 pk.	Univ.
79-72164	9" x .5"	24 pk.	Univ.
79-72166	9" x .75"	24 pk.	Univ.
79-72168	9" x 1"	12 pk.	Univ.
79-72010	18" x 2"	12 pk.	Univ.
79-72020	18" x 1"	12 pk.	Univ.
79-72030	9" x .5"	12 pk.	Univ.



PROCARE®

Plastalume™ Finger Splints and Kit

Coated aluminum with foam backing. 5/8" width. Kit contains six (6) of each size splint.

Indication: Applied to the extensor surface of the finger for stability.

	Type of End	Length	Qty.
79-73211	Splint Kit		48 pcs.
79-73212	Plain	1.25"	6 pk.
79-73213	Plain	2.25"	6 pk.
79-73214	Plain	3.25"	6 pk.
79-73215	Ball	3.25"	6 pk.
79-73216	Ball	4.25"	6 pk.
79-73217	Ball	5.25"	6 pk.
79-73218	Ball	6.25"	6 pk.
79-73219	Ball	7.25"	6 pk.



PROCARE®

Aluminum Finger Splint Set

Constructed of malleable aluminum foam. 48 total piece kit contains three each of each style; Sizes small, medium and large. Frog splint, Finger Cot, Curved Finger, Baseball Splint; Finger Protectors.

Indication: Ideal for finger phalange fracture, finger joint dislocation, finger tip injury, distal finger phalange fracture distal finger joint injury, and post partial finger amputation.

	Qty.	Size
79-71020	48 pcs.	Univ.



PROCARE®

Finger Trap

Smooth, tapered stainless steel, double mesh construction provides a safe, secure design. Spreader bar not included. Sold individually.

Indication: May be used to reduce forearm and colles fracture, prior to splint or cast application.

	Dimensions	Size
79-98513	4.25" x .625"	S
79-98515	6" x .75"	M
79-98517	4.5" x 1"	L



PROCARE®

Tibial Fracture Brace

A bilateral upright foot component to add medial/lateral support as well as to prevent distal slippage. Each brace comes complete with two (2) double thick cotton stockinettes. Measure length from superior border of patella to mid-level malleoli.

Indication: Foam lined preformed Sarmiento style patella tendon-bearing brace provides intimate contact needed for effective support of tibial fractures.

Right	Left	Calf Circ.	Length	Size
79-97703	79-97693	11" – 13"	15" – 16.5"	S
79-97704	79-97694	13" – 15"	18" – 20"	M/Long
79-97705	79-97695	13" – 15"	16.5" – 18.5"	M
79-97707	79-97697	15" – 18"	18" – 20"	L
79-97706	79-97696	15" – 18"	19.5" – 21.5"	L/Long



PROCARE®

Humeral Cuff

Provides total contact circumferential soft tissue compression for optimal control of humeral segments while allowing range of motion of both shoulder and elbow joints. Neutral design fits right or left. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Prefabricated foam lined orthosis designed for the management of mid-shaft and distal third humeral fractures.

	Bicep Circum.	Length	Size
79-99213	8" – 11"	6"	S
79-99215	10" – 13"	6.25"	M
79-99217	12" – 15"	7.25"	L
79-99218	14" – 17"	8"	XL



PROCARE®

Humeral Fracture Brace/Over the Shoulder

Provides total contact circumferential soft tissue compression for optimal control of humeral segments while allowing range of motion of both shoulder and elbow joints. Features a deltoid cap extension, which provides greater alignment control and prevents distal slippage. Neutral design fits right or left. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Prefabricated foam lined orthosis designed for the management of mid-shaft and distal third humeral fractures.

	Bicep Circum.	Length	Size
79-97953	8" – 11"	10"	S
79-97955	10" – 13"	10.25"	M
79-97957	12" – 15"	11.25"	L
79-97958	14" – 17"	12.75"	XL



PROCARE®

Ulnar Fracture Brace

Circumferential soft tissue compression provides cast-like rigidity for optimal control of ulnar segments. Full range of wrist and elbow motion is permitted to help prevent atrophy, increase vascularity and encourage bone growth. Neutral design fits right or left. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Lightweight, foam lined brace is designed for the management of mid-shaft and distal ulnar fractures.

	Distal Circ.	Proximal Circ.	Size
79-97943	6" – 7"	8" – 11"	S
79-97945	6.75" – 7.75"	9.5" – 12.5"	M
79-97947	7.5" – 8.5"	11" – 14"	L



PROCARE®

Wrist Hand Orthosis

Fabricated with a foam lined orthopedic grade thermoplastic, which provides durable performance and long lasting reliability. A foam tongue provides additional padding for the styloid process. Splint may be trimmed with cast scissors and formed using a heat gun. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Lightweight orthosis is designed to comfortably restrict motion of the wrist and hand while holding hand in neutral position.

Right	Left	Palm Width	Length	Size
79-97963	79-97973	2.5" – 3"	7.5"	S
79-97965	79-97975	3" – 3.5"	7.5"	M
79-97967	79-97977	3.5" – 4"	8"	L



PROCARE®

Wrist Hand Thumb Orthosis

Designed to restrict motion of the wrist and thumb while holding the hand in a neutral position. Fabricated of orthopedic grade thermoplastic and foam lined for comfortable, durable performance and long lasting reliability. Splint may be trimmed by using cast scissors and formed using a heat gun. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Orthosis provides rigid cast-like immobilization of the wrist, hand, and thumb without obstructing the web space between the thumb and fore finger.

Right	Left	Palm Width	Length	Size
79-97723	79-97713	2.5" – 3"	7.5"	S
79-97725	79-97715	3" – 3.5"	7.5"	M
79-97727	79-97717	3.5" – 4"	8"	L



PROCARE®

Super-Lite A.F.O.

Available in right or left, male and female sizes.

Indication: Low profile, trim-to-fit ankle orthosis for correction of drop foot.

	Right	Left	Size
Male	79-81150	79-81151	Univ.
Female	79-81160	79-81161	Univ.



PROCARE®

Posterior Leaf Splint (AFO)

Durable, orthopedic grade thermoplastic may be custom trimmed and reshaped. Will fit inside most lace-up shoes.

Indication: Designed to assist ambulation for those patients suffering mild to moderate neurological deficit resulting in drop foot.

Shoe Size		Men	Women	Height	Size
Righth	Left				
79-97743	79-97733	6-8	≤6	12.0"	S
79-97745	79-97735	8.5-10	6.5-8	13.5"	M
79-97747	79-97737	10.5-12	8.5-10	13.5"	L
79-97748	79-97738	12+	10+	14.0"	XL

COMPRESSION CARE

PROCARE®

Compression Pant

These body hugging pants are designed to be worn under clothing and feature a double tummy panel, 3" wide waistband, reinforced open crotch, side-separating cotton-lined zippers with a 1/8" thick, full length foam guard pad and full length contact closures for optimum compression and a comfortable, customizable fit. A strip of lace trim along the bottom of each pant leg helps prevent rollup. Latex free. Available in above and below the knee style in both standard and high waist configurations.



Indication: Provides excellent post-surgical support and compression for surgeries around the abdomen, buttocks, flank area, hips, thighs and knees.

Above Knee Compression Pant				
Above Knee Pant	Below Knee Pant	Hip Measurement	Color	Size
79-90613	79-90633	33" – 36"	White	S
79-90615	79-90635	37" – 40"	White	M
79-90617	79-90637	41" – 45"	White	L
79-90618	79-90638	46" – 50"	White	XL
79-90619	79-90639	51" – 55"	White	2XL
79-90623	79-90643	33" – 36"	Black	S
79-90625	79-90645	37" – 40"	Black	M
79-90627	79-90647	41" – 45"	Black	L
79-90628	79-90648	46" – 50"	Black	XL
79-90629	79-90649	51" – 55"	Black	2XL

High Waist Compression Pant				
Above Knee Pant	Below Knee Pant	Hip Measurement	Color	Size
79-90653	79-90673	33" – 36"	White	S
79-90655	79-90675	37" – 40"	White	M
79-90657	79-90677	41" – 45"	White	L
79-90658	79-90678	46" – 50"	White	XL
79-90659	79-90679	51" – 55"	White	2XL
79-90663	79-90683	33" – 36"	Black	S
79-90665	79-90685	37" – 40"	Black	M
79-90667	79-90687	41" – 45"	Black	L
79-90668	79-90688	46" – 50"	Black	XL
79-90669	79-90689	51" – 55"	Black	2XL



PROCARE®

Compression Bra

This non-binding, machine washable bra is so comfortable, your patients will want to wear it even after recovery is complete. The adjustable shoulder straps create the proper compression throughout the recovery process, and convenient front closures eliminate any twisting or straining to hook the bra. Latex free.

Indication: Delivers compression and support after breast augmentation surgery.

Beige	White	Black	Size
79-90932	79-90732	79-90832	32
79-90933	79-90733	79-90833	34
79-90935	79-90735	79-90835	36
79-90937	79-90737	79-90837	38
79-90938	79-90738	79-90838	40
79-90939	79-90739	79-90839	42



PROCARE®

Breast Wrap

Comfortably padded and the ideal 3" width for effective compression, this contoured band is easy to use thanks to its convenient front closure. One size fits most up to 45". Latex free.

Indication: Adds support and helps in implant stabilization and positioning.

	Color	Size
79-90690	White	Univ.



PROCARE®

Arm Compression Garment

Made from a soft spandex lined, body-hugging material for excellent compression and comfort. The Arm Compression Garment features two adjustable easy-to-reach front contact closure straps, an upper arm adjustment, and no internal seams for ultimate comfort. Latex free.

Indication: Provides post-surgical compression of the arms.

	Color	Bicep Circ.	Distance*	Size
79-90713	White	9" – 10"	10" – 12"	S
79-90715	White	10" – 11"	12" – 14"	M
79-90717	White	11" – 12"	14" – 16"	L
79-90718	White	12" – 13"	16" – 18"	XL
79-90719	White	13" – 14"	18" – 20"	2XL

*Measure distance between back of shoulders.



PROCARE®

Facial Wrap

This moisture-resistant wrap is made with breathable, latex-free fabric to keep your patients comfortable. Five separate contact closure tabs let you make adjustments as recovery progresses, all without irritating internal seams.

Indication: Gives compression and support after face and neck surgery.

	Color	Size
79-90700	White	Univ.



PROCARE®

Photo Modesty Panty

When discretion is paramount, these single-use panties made from stretchy polypropylene offer just the right amount of coverage to make your clients feel comfortable. The flexible elastic waistband makes application and removal easy.

Indication: Offers discretion and coverage during or post liposuction surgery.

	Color	Qty.	Size
79-90727	Blue	60/pk.	Univ.



PROCARE®

Imak Eye Pillow


Your patients will love the way this lightweight eye mask feels. Its breathable cotton Lycra fabric is soft to the touch, and smooth ergoBeads® gently massage the area around the eyes without the weight and clammy feeling of a cold mask - but with the same soothing result.

Indication: Helps reduce headaches, sinus pain, tired/puffy eyes and post surgical pain.

	Color	Size
79-90710	Black	Univ.



Visit
DJOGlobal.com



At DJO Global, we help people prevent injuries and support them on their journey to recovery. We call it 'Together in Motion,' and now all of our global brands are in one convenient location for you. DJOGlobal.com is home to comprehensive tools and information on our diverse portfolio of products and solutions to prevent injury, rehabilitate muscles, bones and joints and treat pain. Visit us online for detailed product information, specifications, education and more.

• 1,000+ DJO Global products

• Optimized for mobile devices

• Instructional Videos

Connect with us:



facebook.com/DJOGlobal



twitter.com/DJOGlobal

Take a look at other DJO Global products...



DonJoy pioneered functional knee bracing 30+ years ago and has led the profession ever since.



OA Nano

Revolutionary in weight and profile, OA Nano offers industry leading high-tech comfort and pain relief for mild to moderate knee OA. This safe, non-invasive treatment option is designed for people to stay active and move freely, making patient compliance a reality.



X-Act ROM Elbow

A versatile, easy to apply, comfortable telescoping brace that allows for ideal fixed or protected range of motion to help improve patient outcomes following surgery or injury. Easy shape-to-fit malleable forearm and bicep cuffs contour for a custom fit.



Retailed packaged orthopedic products catering to independent pharmacies and home healthcare dealers.



Fast Freeze

Fast Freeze delivers naturally cool pain therapy that provides rapid relief to sore muscles and joints. Available in convenient roll-on, gel and spray options, Fast Freeze is an invigorating, fast-absorbing pain reliever made with all natural ingredients.



Bell-Horn

Bell-Horn is your trusted single source for a competitive line of high-quality, innovative retail packaged orthopedic products for independent pharmacies and home healthcare dealers.



The Finest Quality Comfort Footwear **Period!**



Style

From boots to athletic shoes, and dress shoes to sandals, slippers and slip-ons, Dr. Comfort is the worldwide leading provider of diabetic footwear. Dr. Comfort was founded in 2002 and comfort and quality are exactly what the doctor ordered.



Comfort & Quality

Recognizing the need for style, comfort and quality in the therapeutic shoe industry, Dr. Comfort footwear line includes diabetic shoes, prescription and non-prescription inserts, diabetic socks, slippers, sandals and compression hosiery.

INDEX

Part No.	Description	Page No.	Part No.	Description	Page No.
2125	Gravity Cooler.....	12	11ASC	Self-Contained (SC) Knee Cryo/Cuff.....	13
2140	Cooler Hanger.....	12	11C01	Knee Cryo/Cuff.....	13
3008	VenaFlow Elite System Tubing.....	16	125-31-1111	Hand Based Ulnar Gutter Brace.....	5
3040	VenaFlow Elite Calf Cuff.....	16	126-31-1111	Hand Based Radial Gutter Brace.....	6
3045	VenaFlow Elite Thigh Cuff.....	16	12B01	Shoulder Cryo/Cuff.....	13
3046	VenaFlow Elite Foot Cuff.....	16	13A01	Thigh Cryo/Cuff.....	14
3225	ArterialFlow Upper Extremity Cuff.....	20	13C01	Calf Cryo/Cuff.....	13
3227	ArterialFlow Lower Extremity Cuff.....	20	14A01	Back/Hip/Rib Cryo/Cuff.....	14
0129A	Tube Stretch Sock.....	97	15A01	Elbow Cryo/Cuff.....	14
0130A	Hygiene Cover.....	98	16A01	Hand/Wrist Cryo/Cuff.....	14
01A	SP Walker™ (short pneumatic).....	96	200621-00	EXOS Form™ 621.....	9
01EF	Standard.....	94	200626-00	EXOS Form™ 626.....	9
01EP	Elite.....	93	200627-40	EXOS Form™ 627.....	10
01ES	Short.....	94	200631-40	EXOS Form™ 631.....	10
01F	FP Walker™ (foam pneumatic).....	95	210-21-XXXX	Short Thumb Spica.....	4
01K	Heel Wedges.....	97	211-31-1111	Extended Short Thumb Spica (ESTS).....	5
01P	XP Walker™ (extra pneumatic).....	95	220-31-1111	Wrist Brace.....	3
01PD	XP Diabetic Walker™ System.....	96	221-21-XXXX	Wrist Brace with Boa.....	3
01PI	XP Diabetic Walker Impax™ Replacement Kit.....	96	230-31-1111	Long Thumb Spica (LTS).....	4
01T	Toe Cover.....	98	231-21-XXXX	Long Thumb Spica with Boa.....	4
01WXT	Weather Cover.....	98	25-6755	X-Act ROM Knee Zip Ties.....	63
02E	Air-Stirrup Universe™.....	79	25-7525	Kit, Replacement Pad, X-Act™ ROM Hip.....	64
02EK	Air-Stirrup® Universe™ Care Kit.....	79	25-7528	Replacement Belt, X-Act™ ROM Hip.....	64
02MXXX	AirSport™.....	79	25-7529	Assembly, Hip Component, X-Act™ ROM Hip.....	64
02PXX	AirLift™ PTTD Brace.....	85	3007-PL	VenaFlow System Tubing.....	18
02TXX	A60™ Ankle Support.....	83	3010-PL	VenaFlow Calf Cuff.....	17
02XL	Air-Stirrup® Ankle Brace.....	78	3015-PL	VenaFlow Thigh Cuff.....	18
02XXK	Ankle Sprain Care Kit™.....	79	3016-PL	VenaFlow Foot Cuff.....	17
02XXX	Air-Stirrup Plus.....	80	307X	VenaFlow Elite System Repair.....	17
03XX	Leg Brace.....	81	30A	VenaFlow® System.....	17
05A-B	Pneumatic Armband.....	39	30B	VenaFlow Elite System.....	15
05EW	Mayo Clinic Elbow Brace.....	40	310-31-XXXX	Short Arm Fracture Brace with Boa.....	2
05FSX	StabilAir™ Wrist Brace.....	44	311-11-XXXX	Pediatric Short Arm Fracture Brace Open Thumb.....	3
05WX	A2™ Wrist Brace.....	44	312-31-XXXX	Short Arm Fracture Brace Open Thumb.....	3
06AB	Quick-Fit® Shoulder Immobilizer.....	31	315-31-1111	Thumb Spica Fracture Brace.....	4
06GMA	Arm Immobilizer.....	32	325-31-1111	Boxer's Fracture Brace.....	5
07X	Aircast Knee Immobilizer.....	62	326-31-1111	Radial Gutter Fracture Brace.....	5
08X	Infrapatellar Band.....	74	32A	ArterialFlow System*.....	19
09AX	AirHeel™.....	88	330-40-1111	Humeral Fracture Brace.....	6
09DX	Dorsal Night Splint.....	88	500R-X	VenaFlow Knee High Anti Embolism Stockings.....	21
09PKX	AirHeel/DNS Care Kit.....	87	501R-X	VenaFlow Thigh High Anti Embolism Stockings.....	22
100-001-EKIT	Exos Rapid Splint + Exothermix Kit.....	8	51A	IC Cooler.....	11
100-001-EXXX	Exos Rapid Splint + Exothermix.....	8	79-71020	Aluminum Finger Splint Set.....	127
100-001-PXXX	Exos Rapid Splint with Padding.....	8	79-7188X	Finger Protectors.....	125
100-001-XXX	Exos Rapid Splint.....	8	79-7190X	Finger Cots.....	126
100-002-PXXX	Exos Rapid Splint Extra with Padding.....	8	79-7192X	Curved Fingers.....	126
10A01	Ankle Cryo/Cuff.....	12	79-7194X	Cock-Up Splint.....	124
10C01	Foot Cryo/Cuff.....	12	79-7196X	Frog Splints.....	126
11-0201	VenaFlow System Repair Kits.....	18	79-7197X	Forearm Splint.....	124
11-0215-X-XXXXX	Reaction Knee® Brace.....	75	79-7209X	Baseball Splints.....	125
11-0541-X	Replacement Under Sleeve.....	75	79-7211X	Colles Splint.....	123
11-2151-9	X-Act ROM Knee Brace.....	63	79-7216X	Aluminum Strips.....	127
11-3233	X-Act™ ROM HIP.....	64	79-7220X	Metacarpal Splints.....	124
11-3233	X-Act™ ROM HIP.....	64	79-7224X	Plastic Finger Splints.....	126

Part No.	Description	Page No.	Part No.	Description	Page No.
79-7228X	Posterior Tibia/Fibula Splint	125	79-8231X	Elbow Sleeve	41
79-7229X	Femoral Leg Splint.....	125	79-8232X	Tennis Elbow Support	42
79-722XX	Tuli's® Heel Cups.....	86	79-8233X	Thigh Sleeve	70
79-722XX	Plastic Finger Splint Starter Set	127	79-8234X	Calf Sleeve	70
79-7280X	Mason Allen Splint.....	124	79-8239X	Reddie® Brace.....	75
79-7321X	Plastalume™ Finger Splints and Kit	127	79-8240X	Arthroscopy Knee Support	69
79-800XX	3-Panel Knee Splint	60	79-8246X	Universal Patella Knee with Buttress.....	73
79-800XX	Super Knee Splint.....	61	79-8247X	Universal Wrist-O-Prene™.....	48
79-8016X	Lites Visco Knee Support	73	79-824XX	Universal Knee Wrap with Closed Patella.....	68
79-8019X	Elastic Knee Support	69	79-8250X	Sacro-Lumbar Support w/Compression Straps.....	54
79-8023X	Elastic Spiral Knee Support.....	69	79-82570	Universal Surround® Elbow w/FLOAM™.....	40
79-8024X	Elastic Cartilage Knee Support	69	79-82600	Wrist/Thumb Wrap.....	51
79-8032X	Surround® Patella Strap.....	73	79-8261X	Clinic Knee Sleeve	68
79-806XX	Universal Foam Knee Splint	62	79-8262X	Knee Support with Closed Patella.....	68
79-808XX	Clinic 3-Panel Splint	61	79-8263X	Knee Support with Reinforced Patella.....	68
79-80XXX	Compression Knee Dressing.....	62	79-8270X	Universal Thumb-O-Prene™.....	51
79-8100X	Heel and Elbow Protectors	110	79-8270X	Sport Knee Reinforced	72
79-8100X	Heel/Elbow Protector - Poly/Pile.....	110	79-82710	ProCare® Thumb Splint.....	51
79-8101X	Deluxe Cast Boot Sandal.....	91	79-8272X	Stabilized Knee Support.....	71
79-8102X	Deluxe Cast Boot Weatherproof.....	91	79-83040	Padded Head Halter.....	121
79-81030	Deluxe Tennis Elbow.....	41	79-830XX	Universal Cervical Collar.....	24
79-81031	Clinic Tennis Elbow.....	41	79-830XX	Serpentine Cervical Collar.....	25
79-8104X	Deluxe Cast Boot Open Toe.....	91	79-830XX	One-Piece Flat Collar.....	26
79-8106X	Deluxe Ankle Support.....	84	79-83161	California Cervical Collar	26
79-8107X	Beach Boy Cast Boot	91	79-831XX	Patriot® Collar	27
79-8109X	Sport Ankle Splint.....	84	79-83242	Pediatric Collar Accessories.....	23
79-810X0	Heel/Elbow Protectors - Universal	111	79-8327X	Pediatric Collar Only	23
79-8110X	Silicone Heel Cups.....	86	79-8328X	ProCare® Transitional 172 Cervical Collar.....	23
79-8111X	Rocker Cast Boot	90	79-832XX	ProCare® Transitional 172 Cervical Collar.....	23
79-8112X	Elastic Ankle Support.....	85	79-832XX	Miami Jr.® Collar	28
79-8113X	Med/Surg Shoe	90	79-83303	Low Contour Cervical Colla.....	26
79-8118X	Tennis Elbow Support w/FLOAM™	41	79-8335X	Pediatric Collar w/ Repl. Pads.....	23
79-81190	IMAK Elbow Wrap.....	42	79-83370	Vista® Collar.....	24
79-8119X	Surround® FLOAM™ Ankle	80	79-83371	Vista Replacement Pad Set	24
79-811XX	Super-Lite A.F.O.	129	79-83372	Vista TX Repl. Pad Set.....	24
79-8121X	Elastic Elbow Support	42	79-8338X	ACC Transitional 172 Posterior Panel/Pad	23
79-8122X	Rocker Cast Shoe w/Flex Sole.....	92	79-83390	Vista® TX Collar.....	24
79-8123X	Squared Toe Post-Op Shoe	89	79-835XX	Universal Clinic Collar.....	25
79-8125X	Stirrup Ankle Support	81	79-835XX	Clinic Low Density Collar.....	25
79-8131X	Lace-Up Ankle Brace.....	83	79-839XX	911 First Response® Collar	28
79-8132X	Excelerator™ Stirrup Ankle Splint.....	81	79-8400X	Deluxe Arm Sling w/Pad.....	36
79-81330	Universal Ankle Brace.....	81	79-8401X	Clinic Shoulder Immobilizer.....	32
79-8133X	IMPAX Grid Insole	89	79-8402X	Clinic Arm Sling.....	36
79-8135X	Stabilized Ankle Support.....	82	79-8403X	Elastic Shoulder Immobilizer.....	33
79-8136X	DS Ankle Wrap	83	79-84052	Arm Elevator Sling with Pockets	38
79-8137X	Double Strap Ankle Support.....	84	79-8415X	Blue Vogue Arm Sling.....	35
79-8140X	Dorsiwedge™ Night Splint.....	88	79-8416X	Shoulder Immobilizer w/Foam Straps.....	32
79-814XX	Kallassy Ankle Support®.....	83	79-8417X	Chieftain™ Arm Sling	37
79-8151X	Off-Loading Diabetic Shoe	90	79-841XX	Deluxe Shoulder Immobilizer	32
79-81600	Mesh Heel/Elbow Protector	111	79-84230	Deluxe Sling and Swathe.....	33
79-8170X	Surround® Air Ankle.....	80	79-8424X	Sling and Swathe.....	33
79-8200X	Sport Knee Sleeve	67	79-8427X	Specialty Arm Sling.....	36
79-82050	Wrist Wrap	49	79-8429X	Quick Release Arm Sling	37
79-8215X	Hinged Knee Support.....	76	79-84300	Universal Arm Sling	37
79-821X0	Dynamic Pull Universal Knee Wrap	74	79-8434X	Quick Release Economy Shoulder Immobilizer	33
79-8222X	Safety Wrist.....	47	79-84350	Universal Mesh Shoulder Immobilizer	34
79-8230X	Ankle Sleeve	84	79-84500	Shoulder Abduction Kit.....	34

Part No.	Description	Page No.	Part No.	Description	Page No.
79-8500X	Clavicle Splint	29	79-906XX	Compression Pant.....	131
79-8502X	Heavy Padded Clavicle	30	79-90700	Facial Wrap.....	132
79-850XX	Clavicle Strap	30	79-90710	Imak Eye Pillow.....	133
79-85100	Universal Clavicle	30	79-9071X	Arm Compression Garment	132
79-8700X	Contoured Wrist Support.....	45	79-90727	Photo Modesty Panty.....	133
79-87021	Universal Cock-Up Splint.....	49	79-9080X	Ulnar Nerve Protectors.....	108
79-87030	Universal Foam Wrist and Forearm Splint.....	48	79-90810	Laminectomy Arm Cradles with Trough.....	107
79-8703X	Deluxe Wrist & Forearm Supports.....	46	79-90820	Foam Arm Boards.....	107
79-8707X	Elastic Wrist Brace.....	45	79-90830	Multi-Purpose Belt	107
79-8709X	Universal Elastic Wrist Brace	46	79-90850	O.R. Table Pad	108
79-870X0	Universal Wrist & Forearm Supports	47	79-90860	Body Wedge.....	105
79-87100	Abducted ThumbSPICA™.....	52	79-9087X	Head Cradles.....	107
79-8711X	ThumbSPICA™	52	79-9087X	Prone Head Positioners	108
79-8715X	CTS Wrist Support	45	79-908XX	Donuts.....	106
79-8720X	Lace-Up Wrist Support	46	79-9093X	Compression Bra	132
79-8726X	ThumbGUARD™	52	79-91000	Gait Belt	120
79-8728X	ComfortFORM™ Wrist.....	44	79-91010	Heavy Duty Gait Belt	120
79-8732X	Cock-Up Splint	48	79-9103X	Safety Vest Body Holder	116
79-8736X	Loop-Lock Cock-Up Splint.....	49	79-9105X	Vest Body Holder	115
79-873XX	ComfortFORM™ Wrist w/Abducted Thumb.....	51	79-9106X	Mesh Body Holder w/Sleeves.....	116
79-8742X	Foam Wrist Splint	48	79-9107X	Mesh Vest Body Holder	115
79-8744X	Canvas B.A.T.H.™ Wrist Splint.....	47	79-91090	Anodized Aluminum Adjustable Fold-Up Walker.....	113
79-8746X	Quick-Fit® Wrist	45	79-91091	Adjustable Aluminum Cane	113
79-874XX	Quick-Fit® W.T.O.	50	79-91100	Pre-Vent™ Ulnar Nerve Protector	112
79-8750X	Quick-Fit® Wrist & Forearm	44	79-91110	Pre-Vent™ O.R. Table Pad.....	112
79-8754X	IMAK Smart Glove®.....	49	79-91120	Pre-Vent™ Knee Crutch Pad	112
79-8756X	Quick-Fit® Wrist II	43	79-91173	Quick-Release Zipper Body Holder	114
79-875X0	IMAK Pil-O-Splint®	47	79-9122X	Secure-All™ Finger Control Mitt.....	119
79-8900X	10" Double-Pull Sacro-Lumbar Support.....	55	79-9133X	Adjustable Anodized Aluminum Crutches	112
79-8901X	Neoprene/Elastic Sacro-Lumbar w/Compression Straps..	55	79-914X0	Wraparound Limb Holder	117
79-8902X	Sacro-Lumbar Support	55	79-914X0	Personal Limb Holder	117
79-89040	10" Criss-Cross Support	54	79-91500	Deluxe Quick-Release Limb Holder.....	116
79-8906X	Deluxe Rib Belt	59	79-9151X	Elbow Immobilizer	119
79-89080	10" Universal Abdominal Support	58	79-9160X	Secure-All™ Body Holder	115
79-8908X	Industrial Back Support.....	56	79-916XX	Vest Style Body Holders	115
79-890X0	Universal Rib Belt	59	79-916XX	Vest Style Body Holder with Safety Sleeves	116
79-890XX	Abdominal Binder	58	79-9217X	Thumb Splint	52
79-8917X	Clinic Industrial Back Support.....	56	79-9285X	Patella Stabilizer with Buttress.....	72
79-8918X	Criss-Cross Support with Compression Straps	54	79-9400X	KneeRANGER® - Sized	64
79-891XX	Universal Deluxe Rib Belt.....	58	79-94040	KneeRANGER® Lite.....	65
79-891XX	8" Universal Rib Belt	59	79-94060	KneeRANGER® Range-of-Motion	65
79-89200	Catheter Strap.....	110	79-94120	KneeRANGER® Rehab Splint.....	65
79-8920X	Sacro-Lumbar Support with Mesh Back.....	55	79-94130	KneeRANGER® - Universal.....	64
79-8924X	Clinic Retention Support with Compression Straps.....	54	79-9413X	KneeRANGER® - Neoprene	66
79-8929X	Maternity Belt	56	79-94140	KneeRANGER® Fracture Hinges.....	66
79-8935X	ComfortFORM™ Back	53	79-9415X	KneeRANGER® Ankle/Shoe Insert	66
79-8936X	Personal Abdominal Binder	58	79-9425X	KneeRANGER® II - Sized.....	64
79-89XXX	Premium Panel Elastic Binder.....	57	79-9429X	ACTION™ Neoprene Brace/Sleeve - 1/8"	76
79-9000X	Padded Hip Spica.....	110	79-942XX	ElbowRANGER® Motion Control Splint	40
79-9012X	Double Buckle Limb Holder	118	79-942XX	KneeRANGER® II - Universal	63
79-9018X	Post-Op Shoe	90	79-9431X	EveryDAY™ Daily Activity Brace.....	77
79-90191	Leg Elevator	106	79-9435X	WeekENDER™ Recreational Activity Brace	77
79-901X0	Drop Foot Boot	111	79-9443X	Reinforced Patella Stabilizer	72
79-901XX	Foam Limb Holders	118	79-9448X	Hinged Patella Stabilizer.....	72
79-901XX	Hip Abduction Pillows.....	106	79-944XX	Dual Axis Hinge.....	66
79-90500	Arm Elevation Pillow.....	106	79-944XX	Lateral Patella Stabilizer	73
79-90690	Breast Wrap	132	79-944XX	ACTION™ Neoprene Brace/Wrap - 1/8.....	76

Part No.	Description	Page No.	Part No.	Description	Page No.
79-944XX	Dual Axis Hinge.....	77	800-SIZINGTAPE	Exos Sizing Tape	7
79-9503X	SideKICK™	102	800-XX-STRAP	Accessory Strap	7
79-9506X	Nextep™ Contour.....	102	81-03E	Bunion Aid® Plus Replacement Set.....	85
79-9507X	Nextep™ Contour 2	102	81-03H	Bunion Aid®.....	85
79-9508X	Nextep Contour Shortie.....	102	896-XX-1111	EXOS Wrist Undersleeve	7
79-9517X	Nextep™ Contour 2 Air	100			
79-9527X	Nextep™ Contour Air.....	100			
79-9532X	MaxTrax™.....	101			
79-9533X	MaxTrax™ ROM.....	103			
79-9534X	MaxTrax™ Ankle	101			
79-9543X	MaxTrax™ Air ROM	103			
79-9545X	MaxTrax™ Diabetic Walker.....	104			
79-9548X	Evenup™.....	104			
79-9549X	XcelTrax™.....	101			
79-954XX	MaxTrax™ Air.....	99			
79-9550X	XcelTrax™ Ankle.....	100			
79-955XX	MiniTrax.....	99			
79-955XX	XcelTrax™ Air	99			
79-960XX	Quick-Fit® Basic Knee Splint	61			
79-9613X	Bunion Splint.....	85			
79-97000	Arm Band w/Compression Pad	40			
79-9701X	Universal CTS Wrist Brace	46			
79-9703X	Industrial Back Support w/Compression Pad	56			
79-97210	Hot/Cold Therapy Wrap.....	34			
79-9770X	Tibial Fracture Brace	128			
79-9772X	Wrist Hand Thumb Orthosis.....	129			
79-9775X	Plantar Fasciitis Night Splint.....	88			
79-977XX	Posterior Leaf Splint (AFO).....	130			
79-9786X	Surround® Gel Ankle	80			
79-9794X	Ulnar Fracture Brace	129			
79-9795X	Humeral Fracture Brace/Over the Shoulder	128			
79-9796X	Wrist Hand Orthosis	129			
79-98233	Elastic Bandages - Self Closure, 5 yards x 2" 10 Pack.....	8			
79-98235	Elastic Bandages - Self Closure, 5 yards x 3" 10 Pack.....	8			
79-98237	Elastic Bandages - Self Closure, 5 yards x 4" 10 Pack.....	8			
79-98238	Elastic Bandages - Self Closure, 5 yards x 6" 10 Pack	8			
79-9823X	Elastic Bandages - Self Closure	109			
79-9823X	Elastic Bandages - Self Closure, 5 yards.....	8			
79-9843X	Darco™ SlimLine Cast Boot	92			
79-9851X	Finger Trap.....	128			
79-9858X	Elastic Bandages w/Clip Closure	109			
79-9860X	Traction Cord	122			
79-9867X	Contoured Clavicle	30			
79-9872X	Flex-Master® Sterile Clip Closure Bandage	109			
79-9879X	ProSTEP™.....	102			
79-9884X	Flex-Master® Non Sterile Clip Closure Bandage	108			
79-989XX	Latex Rubber Bandage.....	109			
79-9913X	Healthcare Bear Arm Sling	38			
79-9915X	Deep Pocket Economy Arm Sling.....	37			
79-9921X	Humeral Cuff	128			
79-99850	Universal Bucks Traction	121			
79-9844X	Orthowedge.....	92			
800-00	Exos Other Accessories.....	7			
800-40-STRAP	Accessory Strap.....	8			
800-HEATERBAG	Exos Rapid Splint Warmer	8			
800-LR	Boa Locking Ring Kit.....	6			
800-OVEN3P	Exos Top Loading Oven.....	6			

Aircast, Air-Stirrup, Sport-Stirrup, ArterialFlow, AutoChill and VenaFlow are registered trademarks of DJO, LLC.

Air-Stirrup Universe, Air-Stirrup Light, Ankle Sprain Care Kit, AirSport , A2, A6o, AirHeel, AirLift , ARC Brace-Lok, Cryo/Cuff, Cryo/Cuff IC, DJO, XP Walker, FP Walker, SP Walker and StabilAir are trademarks of DJO, LLC.

Surround, ElbowRANGER, KneeRANGER, Reddie, Flex-Master, Turtle Neck, First Response, Quick-Fit, Kallassy, and Procure are registered trademarks of DJO, LLC.

DJO, Impax, ComfortFORM, Wrist-O-Prene, ThumbSPICA, ThumbGUARD, Powerbelt, ACTION, WeekENDER, EveryDAY, SideKICK, ProSTEP, ProROM, Nextep, Thumb-O-Prene, MaxTrax, Chieftan, B.A.T.H., Excelerator, Pre-vent, Secure-All, and Never Stop Getting Better are trademarks of DJO, LLC.

Aspen and Vista are trademarks of Aspen Medical Products, Inc. Miami J and Miami Jr. are registered trademarks of Ossur, Lycra is a trademark of E.I. du Pont de Nemours, Floam is a trademark of Mattel Inc., Plastalume is a registered trademark of Brown Medical Industries, Dorsiwedge is a registered trademark of Walkwell Int., Tuli's is a registered trademark of International Comfort Products, Inc., and Darco is a trademark of Darco International, Velcro is a registered trademark of Velcro USA.

AIRCAST®

PROCARE®



Product Catalog



DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOglobal.com

Together in Motion™

