




Electrotherapy
Accessory Guide



Premium Electrodes

StimCare Premium electrodes provide even, comfortable stimulation eliminating uncomfortable spiked current provided by typical electrodes.

- Incorporates a highly conductive silver ink that controls the flow of electrical current over the entire surface of the electrode
- Better current flow results in more comfortable stimulation
- Features a thick, hypoallergenic gel formulated with pure aloe vera, known for its natural moisturizing effects
- Provides superior adherence to the skin and comfortable treatments over multiple uses
- A tinted, double-coated storage liner allows electrodes to be applied on either side for storage
- Latex Free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	199657-001	2"	Premium Electrode (Round)	4/pk
B	199660-001	2.75"	Premium Electrode (Round)	4/pk
C	199658-001	2" x 2"	Premium Electrode (Square)	4/pk
D	199659-001	2" x 4"	Premium Electrode (Rectangle)	4/pk

Carbon Foam Electrodes


- Flexible foam backing for maximum comfort
- Wide variety of shapes and sizes available for use with most TENS and NMES devices
- Resilient adhesive for longer life in extreme temperatures
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	86905215	2"	Reusable Electrode (Round)	4/pk
B	86905220	2" x 2"	Reusable Electrode (Square)	4/pk
C	86906320	1.5" x 2.5"	Reusable Electrode (Oval)	4/pk
D	86905240	2" x 3.5"	Reusable Electrode (Rectangle)	4/pk
E	86906240	2" x 4"	Reusable Electrode (Oval)	4/pk
F	86906350	2.75" x 5"	Reusable Electrode (Oval)	2/pk
G	86905250	2" x 5"	Reusable Electrode (Rectangle)	2/pk

Carbon Cloth Electrodes


- Low profile, cloth-backed electrodes provide the needed flexibility to conform to body surfaces
- Breathable fabric backing maximizes patient comfort, while anchored pigtail prevents pullout
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	199325-001	1.375"	Round	4/pk
B	198928-001	2"	Electrode 920 (Round)	4/pk
C	199327-001	2" x 2"	Electrode 922ST (Square)	4/pk
D	199024-001	2" x 3.5"	Electrode 924 (Rectangle)	4/pk

Blue Gel Sensitive Skin Electrodes



- Ideal for patients with sensitive skin
- Thicker layer gel increases patient comfort
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	199530-001	4" x 6"	Butterfly, Foam Paravertebral	1/pk
B	199173-001	1.75" x 3.75"	Rectangle	4/pk
C	199174-001	2" x 2"	Square	4/pk

Low Back Electrodes

- Unique design for treating the lower back
- Conforms well to the lower back
- 4" x 7" electrode features a 32" integrated lead wire for convenience
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	198816-001	4" x 6"	Butterfly, Paravertebral	1/pk
B	199530-001	4" x 6"	Butterfly, Skin Sensitive Foam Paravertebral	1/pk
C	86908600	4" x 7"	Pin w/32" Lead <i>*Only compatible with Empi right angle lead wires</i>	1/pk
D	199544-001	1.5" x 13"	Low Back Strip	1/pk

EasyWear™ Conductive Wrap

- Allows electrode placement without painful twists and turns
- Use at home, work or on-the-go
- Prescription required
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	199574-001	ONE SIZE FITS ALL (6" x 48")	Low Back Conductive Wrap	1/pk
	235643-001	4" x 10"	Extender (Not Shown)	1/pk
B	235641-001	4oz. Bottle	Conductive Spray	1/pk



A

SnapEase®

- Designed for localized stimulation with use of snap lead wires or snap adapter
- Unique gel leaves no residue
- Latex-free

	ORDER #	SIZE	PRODUCT	QUANTITY
A	199455-001	1.19"	SnapEase Electrode (Round)	48/pk



A

ComfortEase® Multi-Day Electrode

- Single-use electrode with durable adhesive is ideal for active patients
- May be worn in the shower if disconnected from the device
- Latex-free

	ORDER #	SIZE	PRODUCT	QUANTITY
A	87202104	2.25" x 2.5"	Pin	40/pk

NeuroAid®


- Carbon rubber electrode formulated for use in high heat and humidity
- Alternative gel formulated for excellent adhesion
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	56903799	1.5" x 1.75"	NeuroAid HH Reusable Electrode	4/pk
B	56907790	1.75" x 3.5"	NeuroAid HH Reusable Electrode	2/pk

Post-Operative Electrodes

- Sterile
- Designed for use with TENS devices for pain management after surgery
- Latex-free



	ORDER #	SIZE	PRODUCT	QUANTITY
A	199541-001	2"	Sterile Electrode (Square)	2/pk
B	199542-001	2" x 3.5"	Sterile Electrode (Rectangle)	2/pk



Active TENS™

Specialized wrap adds comfortable electrode support for predefined placement making it simple for patients.



A



B



C

	ORDER #	SIZE	PRODUCT	QUANTITY
A	199599-001	2"	Active Stimgel (Square)	4/pk
B	199648-001		Four Zone Lead Wire	1/pk
C	200049-001		Active Charger	1/pk

TENS/NMES Lead Wires

- Standard pin lead wires for TENS/ NMES devices are 40" in length (longer and shorter lengths also available)
- A coaxial bifurcated lead wire may be attached to a standard lead wire to allow for the use of four electrodes per lead wire or eight electrodes total
- Snap lead wires or a snap adapter may be used with snap electrodes



For use with:
Select, Epix VT



For use with:
Select, 300PV, Epix VT



For use with:
Select, 300PV, Epix VT



For use with:
Select, 300PV, Epix VT

	ORDER #	SIZE	PRODUCT	QUANTITY
A	193068-100	Black/Black 100 cm (40")	Touchproof Lead Wire (Select, Epx VT)	1/pk
B	193057-100	Red/Black 100 cm (40")	Touchproof Lead Wire (Select, 300 PV, Epix VT)	1/pk
C	193064-001	Red/Black 60 cm (24")	Pin Lead	1/pk
D	193057-150	Red/White 150 cm (60")	Touchproof Lead Wire (Select, 300 PV, Epix VT)	1/pk
E	700211-001	White/White 48 cm (19")	Coaxial Bifurcated Lead	1/pk
F	193052-100	Black/Black 100 cm (40")	Snap Touchproof Lead Wire	1/pk



For use with:
All TENS & NMES units

E



For use with:
IF 3Wave, ProMax, SportX

I



F



For use with:
IF 3Wave, ProMax, SportX

J



G



For use with:
Empi Active™

K



For use with:
Ortho DX

H



L

	ORDER #	SIZE	PRODUCT	QUANTITY
G	193014-001	Converts Pin to Snap	Snap Adapter	1/pk
H	630112	Coiled Cable	Touchproof Lead Wire (<i>OrthoDx</i>)	1/pk
I	620060-36	Red/Black 90 cm (36")	Touchproof Lead Wire (<i>IF 3Wave, ProMax, SportX</i>)	1/pk
J	620060-60	Red/Black 150 cm (60")	Touchproof Lead Wire (<i>IF 3Wave, ProMax, SportX</i>)	1/pk
K	199648-001	Black/Gray 90 cm (36")	Four Zone Lead Wire	1/pk
L	620049	Gray/Gray 35 cm (14")	Reverse Bifurcated Lead	1/pk

Skin Care Products

- Skin Prep Improves electrode adhesion, reduces irritation and helps extend the life of the electrode
- Empi Moisturizing Cream conditions irritated skin and contains vitamins A, D and E
- Empi Adhesive Remover easily removes adhesive residue



A



C



B



D



E



F



G



H



I



J



K

	ORDER #	SIZE	PRODUCT	QUANTITY
A	370190-01	12 oz. Bottle	Vitamin E Lotion	1/pk
B	370190-03	2 oz. Bottle	Vitamin E Lotion	1/pk
C	235509-001	—	Skin Prep Pads	54/bx
D	235510-001	—	Adhesive Remover Pads	100/bx
E	235292-001	3" x 30' Roll	Breathable Tape	1/pk
F	370203	12 oz. Jar	Sween Cream	1/pk
G	700200-001	50 grams	Tac Gel	1/pk
H	700201-001	60ml Tube	Res Off	1/pk
I	700737-002	2 oz. Bottle	Sween, Vitamin E Lotion	1/pk
J	700737-008	8 oz. Bottle	Sween, Vitamin E Lotion	1/pk
K	700705-001	2 oz. Bottle	Skin Prep Dab-o-matic	1/pk

Batteries and Chargers

- NiMH rechargeable batteries and battery charger are for use with TENS and 300PV devices



A



B



D



E



C



F



G

H

I

	ORDER #	PRODUCT	QUANTITY
A	200022-001	Battery Charger (TENS) <i>(batteries not included)</i>	1/pk
B	200030-001	NiMH 8.4V Rechargeable Battery	1/pk
C	200047-001	AA/AAA NiMH Battery Charger <i>(batteries not included)</i>	1/pk
D	200033-001	AA Battery NiMH	4/pk
E	200041-001	AAA Battery NiMH	6/pk
F	200049-001	Active Charger	1/pk
G	360102	Battery Pack (IF 3Wave)	1/pk
H	590029	Supply, Power, Cord	1/pk
I	201800	Docking Station (IF 3Wave)	1/pk

ELECTRODE USE AND CARE TIPS

For complete instructions, indications for use, contraindications, warnings, precautions, and adverse side effects, refer to the electrode instructions and to the user's manual of the electrotherapy device.

Reusable Electrodes

- Wash and dry treatment area prior to use
- Skin Prep improves adhesion, reduces irritation and helps extend the life of the electrode
- Attach electrodes to lead wires per device instructions
- Remove electrodes from liner. Do not pull by the wire.
- Attach electrodes per device instructions and clinician recommendations.
- Electrodes adhere better at room temperature
- If adhesive becomes over-saturated, allow electrode to air dry in a refrigerator with adhesive side up until gel regains tack.
- If the electrode gel appears dry, add a few drops of water to the electrode gel. Let rest to regain tack and apply to skin.
- Return electrode to the liner and seal in the provided package at room temperature
- Empi Moisturizing Cream may be applied after use to condition skin
- Electrode life varies due to skin condition and preparation, storage, climate and care
- Replace electrodes that do not stick or cause discomfort

Single-Use Electrodes

- Wash and dry treatment area prior to use
- Skin Prep may be applied to improve adhesion and protect skin
- Attach electrodes per device instructions and clinician recommendations
- Apply with care as the electrode may lose adhesion if removed and reapplied
- Use is acceptable in the shower if detached from device
- If adhesive remains after use, clean with Adhesive Remover
- Empi Moisturizing Cream may be applied after use to condition skin

Device Carriers

- Low Back Electrode Belt must be used with Velcro® strips for maximum effectiveness
- Device carriers, pouches and low back belts are available in black or white



	ORDER #	PRODUCT	QUANTITY
A	198935-001	White Device Carrier with 2" Wide Belt and Pouch	1/pk
B	198936-001	White Device Pouch	1/pk
C	199473-001	Black Device Carrier with 2" Wide Belt and Pouch	1/pk
D	199475-001	Black Device Pouch	1/pk
E	370163	Small Carrier (<i>Select, ProMax</i>)	1/pk
F	370164	Large Carrier (<i>OrthoDx, IF 3Wave</i>)	1/pk
G	199495-001	Carry Case with Empi Logo	1/pk
H	198805-001	Velcro Strips	32/pk
I	199474-001	Black Low Back Electrode Belt with Pouch	1/pk
J	198933-001	White Low Back Electrode Belt with Pouch (<i>Not Shown</i>)	1/pk

ORDERING INFORMATION

By Phone

Contact your Patient Services Representative at:
800-328-2536 or 651-415-9000

By Fax

Faxed orders may be sent to your Patient Services Representative at: 800-450-3593

By Email

Contact Empi using our web site at www.empi.com
Contact Your Empi representative to setup online ordering.

By Mail

Written orders may be sent to:

Empi
Attn: Order Department
P.O. Box 64640, St. Paul, Minnesota 55164-9856
For faster service, please include or have on hand:

Patient Name

Empi Patient Account Number

Patient Services Representative Name

Patient Services Representative Phone Number

To learn more about
Empi's Electrotherapy products visit
www.empi.com or call 800.328.2536



Empi, Inc. | A DJO Global Company

T 800.328.2536 **D** 651.415.9000 **F** 651.415.7414
205 Highway 22 East | Clear Lake, SD 57226 | U.S.A.
DJOglobal.com

Together in Motion™